TSH-IRMA

KIP1891 - KIP1894

History

Summary of change:

Previous Version: 170323-1	Current Version: 200224-1
No IVD logo	IVD logo added
LOT	Version :
No history	History added
PI number :	PI number removed
	Addition of the following sentence at the end of the English IFU: "Other translations of this Instruction for Use can be downloaded from our website: https://www.diasource-diagnostics.com/ "

Read entire protocol before use.

TSH-IRMA

I. INTENDED USE

Immunoradiometric assay kit for the in vitro quantitative measurement of human Thyroid Stimulating Hormone (TSH) in serum and plasma.

II. GENERAL INFORMATION

- A. Proprietary name :** DIAsource TSH-IRMA Kit
- B. Catalog number :** KIP1891 : 96 tests
KIP1894: 4 x 96 tests
- C. Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

Thyrotrope cells of the anterior pituitary synthesize and secrete human thyroid stimulating hormone (TSH), a glycoprotein of molecular weight 28,000 Da, comprising two subunits : α -TSH is very similar to a α subunit of LH, FSH and hCG, β -TSH differs from other hormone subunits and defines the immunological specificity.





TSH regulates the synthesis and release of thyroid hormones : thyroxin (T4) and triiodothyronine (T3). TSH secretion is stimulated by an hypothalamic peptide, TRH (TSH releasing hormone); a negative feedback on TSH secretion is exerted by T3 and T4.

Primary hyperthyroidism is now easily differentiated from euthyroidism by DIAsource ultrasensitive TSH-Irma, because of the high sensitivity (0.025 μ IU/ml) and high discrimination power.

IV. PRINCIPLES OF THE METHOD

The DIASource TSH-IRMA is an immunoradiometric assay based on coated tube separation. Mabs1, the capture antibodies, are attached to the lower and inner surface of the plastic tube. Calibrators or samples added to the tubes will at first show low affinity for Mabs1. Addition of Mab2, the signal antibody labelled with ¹²⁵I, will complete the system and trigger the immunological reaction. After washing, the remaining radioactivity bound to the tube reflects the antigen concentration. The use of several distinct Mabs avoids hyperspecificity.

V. REAGENTS PROVIDED

Reagents	96 tests Kit	4 x 96 tests Kit	Colour Code	Reconstitution
 Tubes coated with anti TSH (monoclonal antibodies)	2 x 48	8 x 48	yellow	Ready for use
 Anti-TSH- ¹²⁵ I (monoclonal antibodies) in TRIS maleate buffer with bovine serum albumin, azide (<0.1%), EDTA and inert red dye	1 vial 5.5 ml 700 kBq	4 vials 5.5 ml 4x700 kBq	red	Ready for use
 Calibrator 0 in bovine serum with thymol	1 vial lyophil.	2 vials lyophil.	yellow	Add 2.0 ml distilled water
 Calibrators 1 -7 in horse serum with thymol (see exact value on vial labels)	7 vials lyophil.	14 vials lyophil.	yellow	Add 2.0 ml distilled water
 Wash solution (TRIS-HCl)	1 vial 10 ml	4 vials 10 ml	brown	Dilute 70x with distilled water (use a magnetic stirrer).
 Controls 1 and 2 in human plasma and thymol	2 vials lyophil.	4 vials lyophil.	silver	Add 1 ml distilled water

Note: 1 µIU of the calibrator is equivalent to 1 µIU of the 2nd IRP 80/558

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- Distilled water
- Pipettes for delivery of: 50 µl, 200 µl, 1 ml and 2 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Vortex mixer
- Magnetic stirrer
- Tubes shaker (700 rpm)
- 5 ml automatic syringe (Cornwall type) for washing
- Aspiration system (optional)
- Any gamma counter capable of measuring ¹²⁵I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- Calibrators** : Reconstitute the calibrators 0-7 with 2.0 ml distilled water.
- Controls** : Reconstitute the controls with 1 ml distilled water.
- Working Wash solution** : Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.

- After reconstitution, calibrators and controls are stable for 8 days at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 3 months. Avoid subsequent freeze-thaw cycles.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum and plasma must be kept at 2-8°C.
- If the test is not run within 24 h., storage in aliquots at -20°C is recommended.
- Avoid subsequent freeze-thaw cycles.
- Serum or plasma (EDTA or heparine) provide similar results.

$$Y (\text{serum}) = 1.02x (\text{hep. plasma}) - 0.06 \quad r = 1 \quad n = 7$$

$$Y (\text{serum}) = 1.00x (\text{EDTA plasma}) + 0.05 \quad r = 1 \quad n = 7$$

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature prior to use. Thoroughly mix all reagents and samples by gentle agitation or swirling. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample. High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times. Prepare a calibration curve for each run, do not use data from previous runs.

B. Procedure

- Label coated tubes in duplicate for each calibrator, control and sample. For determination of total counts, label 2 normal tubes.
- Briefly vortex calibrators, controls, samples and dispense 200 µl of each into the respective tubes.
- Dispense 50 µl of anti-TSH-¹²⁵I tracer into each tube, including the uncoated tubes for total counts.
- Shake the rack containing the tubes gently by hand to liberate any trapped air bubbles.
- Incubate for 2 hours at room temperature on a shaker at 700 ± 100 rpm.
- Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
- Wash tubes with 2 ml Working Wash solution (except total counts). Avoid foaming during the addition of the Working Wash solution.
- Aspirate (or decant) the content of each tube (except total counts).
- Wash tubes again with 2 ml Wash solution (except total counts) and aspirate (or decant).
- After the last washing, let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
- Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

- Calculate the mean of duplicate determinations.
- On semi logarithmic or linear graph paper plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of TSH (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
- Read the concentration for each control and sample by interpolation on the calibration curve.
- Computer assisted data reduction will simplify these calculations. If automatic result processing is to be used, a 4-parameter logistic function curve fitting is recommended.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

TSH-IRMA		cpm	B/T (%)
Total count		359316	100
Calibrator	0.00 µIU/ml	123	0.03
	0.10 µIU/ml	572	0.16
	0.53 µIU/ml	1834	0.51
	1.54 µIU/ml	5284	1.47
	4.90 µIU/ml	16365	4.55
	14.00 µIU/ml	45387	12.63
	48.00 µIU/ml	122595	34.12
	90.00 µIU/ml	184397	51.32

XIII. PERFORMANCE AND LIMITATIONS

A. Detection limit

The Limit of Blank (LoB), Limit of Detection (LoD), and the Limit of Quantitation (LoQ), were determined in accordance with the CLSI guideline EP17-A.

The LoB was calculated by measuring the blank several times and calculating the 95th percentile of the distribution of the test values. The LoB was calculated to be 0.05 µIU/ml.

The LoD was calculated as described in the guideline. The LoD was calculated to be 0.09 µIU/ml.

The LoQ was calculated by testing 5 samples of low value 14 times in different test. The LoQ was calculated to be 0.12µIU/ml with CV of 20%.

B. Specificity

Cross-reacting hormones were added to a low and to a high TSH value calibrator. The apparent TSH response was measured.

added Hormone	TSH CAL 1 µIU/ml	TSH CAL 2 µIU/ml
-	0.09	49.00
LH 300 mIU/ml	0.80	47.86
FSH 300 mIU/ml	0.19	44.58
hCG 300000 mIU/ml	6.36	48.48

C. Precision

INTRA ASSAY				INTER ASSAY			
Serum	Replicate	<X> ± SD (µIU/ml)	CV (%)	Serum	Replicate	<X> ± SD (µIU/ml)	CV (%)
A	10	0.26 ± 0.02	6.0	A	20	1.34 ± 0.06	4.1
B	10	1.82 ± 0.03	1.4	B	20	13.69 ± 0.29	2.1
C	10	33.95 ± 0.20	0.6				

SD : Standard Deviation; CV: Coefficient of variation

D. Accuracy

RECOVERY TEST

Sample	Added TSH (µIU/ml)	Recovered TSH (µIU/ml)	Recovery (%)
Sample1			
Serum	105	107	102
Hep. plasma	105	105	100
EDTA plasma	105	107	102
Sample2			
Serum	0.62	0.65	105
Hep. plasma	0.62	0.63	102
EDTA plasma	0.62	0.65	105

DILUTION TEST			
Sample	Dilution	Theoretical Concent. (µIU/ml)	Measured Concent. (µIU/ml)
1	1/1	71.64	71.64
	1/2	35.82	35.66
	1/4	17.91	17.35
	1/8	8.95	8.42
	1/16	4.48	4.44
	1/32	2.24	2.17
	1/64	1.12	0.95
	1/128	0.56	0.51
	1/256	0.28	0.19
	1/512	0.14	0.08

Samples were diluted with zero calibrator.

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrator has been added to the coated tubes.

TIME DELAY				
	0'	10'	20'	30'
S 1 (µIU/ml)	0.17	0.16	0.15	0.16
S 2 (µIU/ml)	34	34	34	36

F. Hook-effect

A sample spiked with 2500 µIU/ml gives a result higher than the last calibration point.

XIV. LIMITATIONS

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays. Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed in case of the presence of heterophilic antibodies. Carefully evaluate the results of patients suspected of having these antibodies. If results are not consistent with other clinical observations, additional information should be required before diagnosis.

XV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises

XVI. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

Identification	Number of subjects	Range (µIU/ml)
Euthyroidism	216	0.2 - 3.2
Hyperthyroidism	59	< 0.01 - 0.09
Hypothyroidism	26	6.3- 158

XVII. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

This kit contains ¹²⁵I (half-life: 60 days) emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVIII. BIBLIOGRAPHY

1. BAYER M.F., KRIS J.P., McDOUGALL I.R.
Clinical experience with sensitive thyrotropin measurements: diagnostic and therapeutic Implication.
J. Nucl. Med. 1985, 26:1248-56.
2. CALDWELL G., KELLET H.A., GOW S.M., et al.
A new strategy for thyroid function testing.
Lancet 1985,i: 1117-9
3. FIELD J.B.
Pituitary thyrotropin: Mecanism of action.
In the thyroid, S. Werner and S.H. Ingbar. Ed. Harper Row, Hagestown, M.D., 1978.
4. MARDELL R.T., GAMBER T.R., WINTON M.R. J.
High sensitivity assay of thyroid stimulating hormone in patients receiving thyroxine for primary hypothyroidism and thyroid carcinoma.
Br. Med. J. 1985: 290:335-356.
5. MUSTO J.D., PIZZOLANTE J.M. and CHESARONE V.P.
A comment on thyrotropin measurement and evaluation.
Clin. Chem. 1984;30:329
6. PIERCE J.C.
Pituitary thyrotropin: Chemistry.
In the thyroid., S. Wemer and S.H. Ingbar. Ed. Harper Row, Hagerstown, M.D. 1978.
7. RODDIS M.J., BURRIN J.M., JOHANNSSSEN A., et al.
Serum thyrotropin: a first-line discriminatory test of thyroid function.
Lancet 1985; i: 277-8
8. ROSS D.S.
New Sensitive immunoradiometric assays for thyrotropin (Review).
Ann. Intern. Med. 1986; 104:718-21
9. PETER S.A. et al.
Elevated serum thyrotropin (TSH) levels in critically ill patients with acquired immunodeficiency syndrome (AIDS).
Exp. Clin. Endocrinol. 1993; 101(6);346-9

10. JAIMELA E. et al.
Ability of two new thyrotropin (YSH) assays to separate hyperthyroid patients from euthyroid patients with ow TSH.
Clin. Chem. 1994, Jan;40(1);101-105
11. KOMOROWSKI J. et al.
Stimulatory effect of thyrotropin (TSH) on interleukin 2 (IL2) release from human peripheral blood lymphocytes and dose-response study in vitro.
Horm. Metab. Res. 1993, Nov; 25(11); 598-9.
12. ADRIAANSE R. et al.
Pulsatile thyrotropin and prolactin secretion in a patient with an mixed thyrotropin and prolactin secreting pituitary adenoma.
Eur. J. Endocrinol. 1994, Feb; 130(2); 113-20.

XIX. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS ml	CALIBRATORS ml	SAMPLE(S) CONTROLS ml
Calibrators (0-7)	-	0.2	-
Samples, Controls	-	-	0.2
Tracer	0.05	0.05	0.05
Incubation	2 hours at RT with continous shaking		
Separation	-	aspirate (or decant)	
Working Wash solution	-	2.0	
Separation	-	aspirate (or decant)	
Working Wash solution	-	2.0	
Separation	-	aspirate (or decant)	
Counting	Count tubes for 60 seconds		

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

DIAsource Catalogue Nr : KIP1891 – KIP1894	Revision Number: 200224-1
-----------------------------------------------	------------------------------

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High sensitivity assay of thyroid stimulating hormone in patients receiving thyroxine for primary hypothyroidism and thyroid carcinoma.

Br. Med. J. 1985: 290:335-356.

5. MUSTO J.D., PIZZOLANTE J.M. and CHESARONE V.P.
A comment on thyrotropin measurement and evaluation.
Clin. Chem. 1984;30:329
6. PIERCE J.C.
Pituitary thyrotropin: Chemistry.
In the thyroid., S. Wemer and S.H. Ingbar. Ed. Harper Row, Hagerstown, M.D. 1978.
7. RODDIS M.J., BURRIN J.M., JOHANNSEN A., et al.
Serum thyrotropin: a first-line discriminatory test of thyroid function.
Lancet 1985; i: 277-8
8. ROSS D.S.
New Sensitive immunoradiometric assays for thyrotropin (Review).
Ann. Intern. Med. 1986; 104:718-21
9. PETER S.A. et al.
Elevated serum thyrotropin (TSH) levels in critically ill patients with acquired immunodeficiency syndrome (AIDS).
Exp. Clin. Endocrinol. 1993; 101(6);346-9
10. JAIMELA E. et al.
Ability of two new thyrotropin (YSH) assays to separate hyperthyroid patients from euthyroid patients with ow TSH.
Clin. Chem. 1994, Jan;40(1);101-105
11. KOMOROWSKI J. et al.
Stimulatory effect of thyrotropin (TSH) on interleukin 2 (IL2) release from human peripheral blood lymphocytes and dose-response study in vitro.
Horm. Metab. Res. 1993, Nov; 25(11); 598-9.
12. ADRIAANSE R. et al.
Pulsatile thyrotropin and prolactin secretion in a patient with an mixed thyrotropin and prolactin secreting pituitary adenoma.
Eur. J. Endocrinol. 1994, Feb; 130(2); 113-20.

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XVIII. LITERATUR

1. BAYER M.F., KRISS J.P., McDOUGALL I.R.
Clinical experience with sensitive thyrotropin measurements: diagnostic and therapeutic Implication.
 J. Nucl. Med. 1985, 26:1248-56.
2. CALDWELL G., KELLET H.A., GOW S.M., et al.
A new strategy for thyroid function testing.
 Lancet 1985,i: 1117-9
3. FIELD J.B.
Pituitary thyrotropin: Mecanism of action.
 In the thyroid, S. Werner and S.H. Ingbar. Ed. Harper Row, Hagestown, M.D., 1978.
4. MARDELL R.T., GAMBER T.R., WINTON M.R. J.
High sensitivity assay of thyroid stimulating hormone in patients receiving thyroxine for primary hypothyroidism and thyroid carcinoma.
 Br. Med. J. 1985: 290:335-356.
5. MUSTO J.D., PIZZOLANTE J.M. and CHESARONE V.P.
A comment on thyrotropin measurement and evaluation.
 Clin. Chem. 1984;30:329
6. PIERCE J.C.
Pituitary thyrotropin: Chemistry.
 In the thyroid., S. Werner and S.H. Ingbar. Ed. Harper Row, Hagerstown, M.D. 1978.
7. RODDIS M.J., BURRIN J.M., JOHANNSEN A., et al.
Serum thyrotropin: a first-line discriminatory test of thyroid function.
 Lancet 1985; i: 277-8

8. ROSS D.S.
New Sensitive immunoradiometric assays for thyrotropin (Review).
 Ann. Intern. Med. 1986; 104;718-21
9. PETER S.A. et al.
Elevated serum thyrotropin (TSH) levels in critically ill patients with acquired immunodeficiency syndrome (AIDS).
 Exp. Clin. Endocrinol. 1993; 101(6);346-9
10. JAIMELA E. et al.
Ability of two new thyrotropin (YSH) assays to separate hyperthyroid patients from euthyroid patients with ow TSH.
 Clin. Chem. 1994, Jan;40(1);101-105
11. KOMOROWSKI J. et al.
Stimulatory effect of thyrotropin (TSH) on interleukin 2 (IL2) release from human peripheral blood lymphocytes and dose-response study in vitro.
 Horm. Metab. Res. 1993, Nov; 25(11); 598-9.
12. ADRIAANSE R. et al.
Pulsatile thyrotropin and prolactin secretion in a patient with an mixed thyrotropin and prolactin secreting pituitary adenoma.
 Eur. J. Endocrinol. 1994, Feb; 130(2); 113-20.

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11. BAYER M.F., KRIS J.P., McDOUGALL I.R.
Clinical experience with sensitive thyrotropin measurements: diagnostic and therapeutic Implication.
J. Nucl. Med. 1985, 26:1248-56.
2. CALDWELL G., KELLET H.A., GOW S.M., et al.
A new strategy for thyroid function testing.
Lancet 1985,i: 1117-9
3. FIELD J.B.
Pituitary thyrotropin: Mecanism of action.
In the thyroid, S. Werner and S.H. Ingbar. Ed. Harper Row, Hagestown, M.D., 1978.
4. MARDELL R.T., GAMBER T.R., WINTON M.R. J.
High sensitivity assay of thyroid stimulating hormone in patients receiving thyroxine for primary hypothyroidism and thyroid carcinoma.
Br. Med. J. 1985: 290:335-356.
5. MUSTO J.D., PIZZOLANTE J.M. and CHESARONE V.P.
A comment on thyrotropin measurement and evaluation.
Clin. Chem. 1984;30:329
6. PIERCE J.C.
Pituitary thyrotropin: Chemistry.
In the thyroid., S. Wemer and S.H. Ingbar. Ed. Harper Row, Hagerstown, M.D. 1978.
7. RODDIS M.J., BURRIN J.M., JOHANNSEN A., et al.
Serum thyrotropin: a first-line discriminatory test of thyroid function.
Lancet 1985; i: 277-8
8. ROSS D.S.
New Sensitive immunoradiometric assays for thyrotropin (Review).
Ann. Intern. Med. 1986; 104:718-21
9. PETER S.A. et al.
Elevated serum thyrotropin (TSH) levels in critically ill patients with acquired immunodeficiency syndrome (AIDS).
Exp. Clin. Endocrinol. 1993; 101(6);346-9
10. JAIMELA E. et al.
Ability of two new thyrotropin (YSH) assays to separate hyperthyroid patients from euthyroid patients with ow TSH.
Clin. Chem. 1994, Jan;40(1);101-105

11. KOMOROWSKI J. et al.
Stimulatory effect of thyrotropin (TSH) on interleukin 2 (IL2) release from human peripheral blood lymphocytes and dose-response study in vitro.
Horm. Metab. Res. 1993, Nov; 25(11); 598-9.
12. ADRIAANSE R. et al.
Pulsatile thyrotropin and prolactin secretion in a patient with an mixed thyrotropin and prolactin secreting pituitary adenoma.
Eur. J. Endocrinol. 1994, Feb; 130(2); 113-20.

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- PETER S.A. et al.
Elevated serum thyrotropin (TSH) levels in critically ill patients with acquired immunodeficiency syndrome (AIDS).
Exp. Clin. Endocrinol. 1993; 101(6);346-9
- JAIMELA E. et al.
Ability of two new thyrotropin (YSH) assays to separate hyperthyroid patients from euthyroid patients with ow TSH.
Clin. Chem. 1994, Jan;40(1);101-105
- KOMOROWSKI J. et al.
Stimulatory effect of thyrotropin (TSH) on interleukin 2 (IL2) release from human peripheral blood lymphocytes and dose-response study in vitro.
Horm. Metab. Res. 1993, Nov; 25(11); 598-9.
- ADRIAANSE R. et al.
Pulsatile thyrotropin and prolactin secretion in a patient with an mixed thyrotropin and prolactin secreting pituitary adenoma.
Eur. J. Endocrinol. 1994, Feb; 130(2); 113-20.

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1. BAYER M.F., KRISS J.P., McDOUGALL I.R.
Clinical experience with sensitive thyrotropin measurements: diagnostic and therapeutic Implication.
J. Nucl. Med. 1985, 26:1248-56.
2. CALDWELL G., KELLET H.A., GOW S.M., et al.
A new strategy for thyroid function testing.
Lancet 1985,i: 1117-9
3. FIELD J.B.
Pituitary thyrotropin: Mecanism of action.
In the thyroid, S. Werner and S.H. Ingbar. Ed. Harper Row, Hagestown, M.D., 1978.
4. MARDELL R.T., GAMBER T.R., WINTON M.R. J.
High sensitivity assay of thyroid stimulating hormone in patients receiving thyroxine for primary hypothyroidism and thyroid carcinoma.
Br. Med. J. 1985: 290:335-356.
5. MUSTO J.D., PIZZOLANTE J.M. and CHESARONE V.P.
A comment on thyrotropin measurement and evaluation.
Clin. Chem. 1984;30:329
6. PIERCE J.C.
Pituitary thyrotropin: Chemistry.
In the thyroid., S. Wemer and S.H. Ingbar. Ed. Harper Row, Hagerstown, M.D. 1978.
7. RODDIS M.J., BURRIN J.M., JOHANNSEN A., et al.
Serum thyrotropin: a first-line discriminatory test of thyroid function.
Lancet 1985; i: 277-8
8. ROSS D.S.
New Sensitive immunoradiometric assays for thyrotropin (Review).
Ann. Intern. Med. 1986; 104;718-21
9. PETER S.A. et al.
Elevated serum thyrotropin (TSH) levels in critically ill patients with acquired immunodeficiency syndrome (AIDS).
Exp. Clin. Endocrinol. 1993; 101(6);346-9
10. JAIMELA E. et al.
Ability of two new thyrotropin (YSH) assays to separate hyperthyroid patients from euthyroid patients with ow TSH.
Clin. Chem. 1994, Jan;40(1);101-105
11. KOMOROWSKI J. et al.
Stimulatory effect of thyrotropin (TSH) on interleukin 2 (IL2) release from human peripheral blood lymphocytes and dose-response study in vitro.
Horm. Metab. Res. 1993, Nov; 25(11); 598-9.
12. ADRIAANSE R. et al.
Pulsatile thyrotropin and prolactin secretion in a patient with an mixed thyrotropin and prolactin secreting pituitary adenoma.
Eur. J. Endocrinol. 1994, Feb; 130(2); 113-20.

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XVIII. BIBLIOGRAFIA

1. BAYER M.F., KRIS J.P., McDOUGALL I.R.
Clinical experience with sensitive thyrotropin measurements: diagnostic and therapeutic Implication.
J. Nucl. Med. 1985, 26:1248-56.
2. CALDWELL G., KELLET H.A., GOW S.M., et al.
A new strategy for thyroid function testing.
Lancet 1985; i: 1117-9

3. FIELD J.B.
Pituitary thyrotropin: Mecanism of action.
In the thyroid, S. Werner and S.H. Ingbar. Ed. Harper Row, Hagestown, M.D., 1978.
4. MARDELL R.T., GAMBER T.R., WINTON M.R. J.
High sensitivity assay of thyroid stimulating hormone in patients receiving thyroxine for primary hypothyroidism and thyroid carcinoma.
Br. Med. J. 1985: 290:335-356.
5. MUSTO J.D., PIZZOLANTE J.M. and CHESARONE V.P.
A comment on thyrotropin measurement and evaluation.
Clin. Chem. 1984;30:329
6. PIERCE J.C.
Pituitary thyrotropin: Chemistry.
In the thyroid., S. Wemer and S.H. Ingbar. Ed. Harper Row, Hagerstown, M.D. 1978.
7. RODDIS M.J., BURRIN J.M., JOHANNSEN A., et al.
Serum thyrotropin: a first-line discriminatory test of thyroid function.
Lancet 1985; i: 277-8
8. ROSS D.S.
New Sensitive immunoradiometric assays for thyrotropin (Review).
Ann. Intern. Med. 1986; 104:718-21
9. PETER S.A. et al.
Elevated serum thyrotropin (TSH) levels in critically ill patients with acquired immunodeficiency syndrome (AIDS).
Exp. Clin. Endocrinol. 1993; 101(6):346-9
10. JAIMELA E. et al.
Ability of two new thyrotropin (YSH) assays to separate hyperthyroid patients from euthyroid patients with ow TSH.
Clin. Chem. 1994, Jan;40(1):101-105
11. KOMOROWSKI J. et al.
Stimulatory effect of thyrotropin (TSH) on interleukin 2 (IL2) release from human peripheral blood lymphocytes and dose-response study in vitro.
Horm. Metab. Res. 1993, Nov; 25(11); 598-9.
12. ADRIAANSE R. et al.
Pulsatile thyrotropin and prolactin secretion in a patient with an mixed thyrotropin and prolactin secreting pituitary adenoma.
Eur. J. Endocrinol. 1994, Feb; 130(2); 113-20.

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Br. Med. J. 1985: 290:335-356.

5. MUSTO J.D., PIZZOLANTE J.M. and CHESARONE V.P.
A comment on thyrotropin measurement and evaluation.
Clin. Chem. 1984;30:329
6. PIERCE J.C.
Pituitary thyrotropin: Chemistry.
In the thyroid., S. Wemer and S.H. Ingbar. Ed. Harper Row, Hagerstown, M.D. 1978.
7. RODDIS M.J., BURRIN J.M., JOHANNSEN A., et al.
Serum thyrotropin: a first-line discriminatory test of thyroid function.
Lancet 1985; i: 277-8
8. ROSS D.S.
New Sensitive immunoradiometric assays for thyrotropin (Review).
Ann. Intern. Med. 1986; 104;718-21
9. PETER S.A. et al.
Elevated serum thyrotropin (TSH) levels in critically ill patients with acquired immunodeficiency syndrome (AIDS).
Exp. Clin. Endocrinol. 1993; 101(6);346-9
10. JAIMELA E. et al.
Ability of two new thyrotropin (YSH) assays to separate hyperthyroid patients from euthyroid patients with ow TSH.
Clin. Chem. 1994, Jan;40(1);101-105
11. KOMOROWSKI J. et al.
Stimulatory effect of thyrotropin (TSH) on interleukin 2 (IL2) release from human peripheral blood lymphocytes and dose-response study in vitro.
Horm. Metab. Res. 1993, Nov; 25(11); 598-9.
12. ADRIAANSE R. et al.
Pulsatile thyrotropin and prolactin secretion in a patient with an mixed thyrotropin and prolactin secreting pituitary adenoma.
Eur. J. Endocrinol. 1994, Feb; 130(2); 113-20.

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