# **anti-GAD<sub>65</sub> RIA**

***KIPM2070***

***KIPM2071***

# History

---

## Summary of change :

<b>Previous Version :</b> 180921/2	<b>Current Version :</b> 200224/1
Multilanguage IFU	Addition of the following sentence at the end of the English IFU: "Other translations of this Instructions for Use can be downloaded from our website: <a href="https://www.diasource-diagnostics.com/">https://www.diasource-diagnostics.com/</a> "
/	Addition of Manufacturer symbol
PI number	PI number cleared

Read entire protocol before use.

## anti-GAD<sub>65</sub> RIA

### I. INTENDED USE

Radioligand assay for the determination of autoantibodies to Glutamic Acid Decarboxylase (GAD<sub>65</sub>Ab) in human serum.

### II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource anti-GAD<sub>65</sub> RIA Kit
- B. **Catalog number :** KIPM2070: 100 tests  
KIPM2071: 50 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.  
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

**For technical assistance or ordering information contact :**  
**Tel : +32 (0)10 84.99.11      Fax : +32 (0)10 84.99.91**

### III. CLINICAL BACKGROUND

Type 1 diabetes, also known as insulin-dependent diabetes mellitus (IDDM), results from a chronic autoimmune destruction of the insulin-secreting pancreatic beta cells, probably initiated by exposure of genetically susceptible host to an environmental agent. Autoimmune destruction of beta cells is thought to be completely asymptomatic until 80-90% of the cells are lost. This process may take years to complete and may occur at any time in all ages.

During the preclinical phase, this autoimmune process is marked by circulating autoantibodies to beta cell antigens. These autoantibodies are present years before the onset of type 1 diabetes and prior to clinical symptoms. Early studies utilized the immunofluorescence test for islet-cell antibodies (ICA), which has been difficult to standardize and is now replaced by a combination of several radioimmunoassays for antibodies against specific beta cell antigens, such as insulin (IAA), glutamic acid decarboxylase (GAD) and tyrosine phosphatase ICA 512 (IA2).

GAD, the enzyme that catalyzes the conversion of glutamate to GABA, has been identified in two isoforms, molecular weight 65.000 (GAD65) and 67.000 (GAD67). Although GAD autoantibodies are found in type 1 diabetes and in the rare neurological disorder Stiff man syndrome (SMS), the GAD autoantibodies profile in the two diseases differs. Auto-antibodies of SMS patients recognize predominantly linear epitopes on the N-terminus of the GAD protein. The use of a truncated antigen in the anti-GAD65 is limiting the use of this assay for the detection of neurological antibodies.

GAD65 autoantibodies are present in 70-80 % of newly diagnosed patients with type 1 diabetes and can be detected many years before clinical onset of the disease. The combination of the autoantibodies to GAD65 and IA2 is highly relevant for risk assessment of type 1 diabetes in children and adolescence. These tests in combination are more sensitive and predictive than ICA in risk groups, e.g. relatives of type 1 diabetes.

GAD65 autoantibodies also occur in a subset of adults with type 2 diabetes (non-insulin-dependent diabetes mellitus). These patients can have pronounced hyperglycemia, and after therapy with oral hypoglycemic agents for several months to years they may become insulin dependent. Therefore, these patients are thought to have a slowly progressive form of type 1 diabetes, sometimes called latent diabetes or latent autoimmune diabetes in adults (LADA). The presence of GAD65 autoantibodies in sera of such patients is a sensitive and specific marker for future insulin dependency.

#### IV. PRINCIPLES OF THE METHOD

The DIASource anti-GAD<sub>65</sub> RIA kit is a direct assay based on the principle of radioligand assays. Highly purified human recombinant GAD<sub>65</sub> is labeled with <sup>125</sup>Iodine. This tracer is used in excess and bound by the GAD<sub>65</sub> autoantibodies of the sample.

The DIASource anti-GAD<sub>65</sub> RIA tracer meets the highest requirements with regard to purity, enzymologic identity, fast reaction kinetics, cross reactivity at zero level and stability. These are the main prerequisites for the specific binding of the tracer and its exclusive recognition by the GAD<sub>65</sub> autoantibodies of the sample.

By adding dissolved Protein A which binds to the Fc moiety of the autoantibodies, sandwich-type complexes are formed. By adding a precipitation buffer the immune complex will be precipitated as a solid phase which facilitates the simple separation of the bound fraction (B) by centrifugation. After removing the supernatant which contains the non-bound tracer by aspiration or decantation, the radioactivity of the remaining precipitate is measured.

The concentration of GAD<sub>65</sub> autoantibodies (anti-GAD<sub>65</sub>) in the sample is reflected by the specifically bound tracer amount. The radioactive signal (cpm) of the bound fraction (B) is proportional to the autoantibody concentration.

No immune complex is formed if autoantibodies against GAD<sub>65</sub> are absent in the sample, as the tracer binds solely to GAD<sub>65</sub> autoantibodies, but not to Protein A.

A standard curve with a range of 0.1 - 120 U/ml is established by measuring cpm respectively the binding B/T % of the calibrators 1 - 6. The anti-GAD<sub>65</sub> concentration value of the patient's sample is directly read off against this curve.

#### V. REAGENTS PROVIDED

Reagents	50 Tests Kit	100 Tests Kit	Reconstitution		
<table border="1"> <tr> <td>Ag</td> <td>125I</td> </tr> </table> <p>TRACER: <sup>125</sup>Iodine labelled GAD<sub>65</sub> (human, recombinant)</p>	Ag	125I	1 vial Lyo <0.05 MBq	2 vials Lyo <0.05 MBq/vial	<b>Reconstitute</b> with 2.6ml Reconstitution buffer
Ag	125I				
<table border="1"> <tr> <td>CAL</td> <td>N</td> </tr> </table> <p>Calibrators - N = 1 to 6 Anti-GAD<sub>65</sub> calibrators in human serum (conc: see vial label)</p>	CAL	N	6 vials 0.25ml	6 vials 0.25ml	<b>Ready for use</b>
CAL	N				
<table border="1"> <tr> <td>PROTEIN A</td> </tr> </table> <p>Protein A component</p>	PROTEIN A	1 vial Lyo	2 vial Lyo	<b>Reconstitute</b> with 2.6ml Reconstitution buffer	
PROTEIN A					
<table border="1"> <tr> <td>REC</td> <td>SOLN</td> </tr> </table> <p>Reconstitution buffer</p>	REC	SOLN	1 vial 6 ml	1 vial 12 ml	<b>Ready for use</b>
REC	SOLN				
<table border="1"> <tr> <td>CONTROL</td> <td>N</td> </tr> </table> <p>Controls - N = 1 or 2 Anti- GAD<sub>65</sub> controls in human serum (conc: see vial label)</p>	CONTROL	N	2 vials 0.25ml	2 vials 0.25ml	<b>Ready for use</b>
CONTROL	N				
<table border="1"> <tr> <td>PREC</td> <td>AGENT</td> </tr> </table> <p>Precipitation buffer</p>	PREC	AGENT	1 vial 55 ml	1 vial 110 ml	<b>Ready for use</b>
PREC	AGENT				

#### VI. SUPPLIES NOT PROVIDED

- Precision pipettes 20 - 100 µl, 1000 µl
- Disposable pipette tips
- Distilled or de-ionized water
- Absorbent paper or paper towel
- Vortexer
- Centrifuge
- γ-counter

#### VII. REAGENT PREPARATION

Allow samples to reach room temperature prior to assay. Take care to agitate serum samples gently in order to ensure homogeneity.

##### Tracer:

Reconstitute with 2.6 ml Reconstitution buffer per vial. Reconstituted tracer remains stable for 2 weeks, stored at 2 - 8 °C.

##### Reconstitution Buffer:

Reconstitution buffer is ready for use and serves for the reconstitution of the tracer and the Protein A component.

##### Precipitation Buffer:

Precipitation buffer is ready for use and serves for the precipitation and washing of the radioactive immune complex. **Use cold.**

##### Protein A component:

Reconstitute with 2.6 ml Reconstitution buffer per vial. The reconstituted solution remains stable for 2 weeks stored at 2 - 8 °C.

**Calibrators 1 - 6:** Ready for use.

**Control serum 1-2:** Ready for use.

#### VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

The DIASource anti-GAD<sub>65</sub> RIA kit has been designed for 100 and 50 determinations, respectively. This is sufficient for the analysis of 41 or 16 unknown samples as well as for calibrators and control sera, assayed in duplicates.

The expiry date of each component is reported on its respective label, that of the complete kit on the box label.

Upon receipt, all components of the DIASource anti-GAD<sub>65</sub> RIA kit have to be kept at 2 - 8 °C, preferably in the original kit box.

#### IX. SPECIMEN COLLECTION AND PREPARATION

- Blood is taken by venipuncture. After clotting, the serum is separated by centrifugation. Lipaemic and hemolytic samples should not be employed.

- The samples may be kept at 2 - 8 °C up to three days. Long-term storage requires - 20 °C.

- Repeated freezing and thawing should be avoided. For multiple use, initially aliquot samples and keep at - 20 °C.

#### X. PROCEDURE

**Use of conical tubes for the DIASource anti-GAD<sub>65</sub> kit only.**

1. Label test tubes appropriately.
2. Pipette into the corresponding tubes according to assay scheme
  - 20 µl calibrators,
  - 20 µl control sera,
  - 20 µl patient's samples, each.
3. Add 50 µl tracer (prepared from tracer and reconstitution buffer), each, to **all tubes**, including those for total radioactivity T.  
*Tubes T are now separated until radioactivity is measured.*
4. Incubate for 2 hours (at room temperature).
5. Add 50 µl Protein A solution (prepared from Protein A component and reconstitution buffer), each.  
(Agitate the suspension gently prior to use - please cf. section VII. Reagent preparation before use).
6. Incubate for 1 hour (at room temperature).
7. Add 1 ml **cold** precipitation buffer, each.
8. Centrifuge the tubes for 30 minutes at 4-8°C at 2000 x g.
9. Aspirate supernatant completely or decant. For removal of any remaining liquid, turn tubes upside down (minimum 5 minutes) and absorb any droplets by carefully tapping on blotting paper.
10. Measure radioactivity of **all tubes including T**.  
Recommended counting time: 1 minute

#### XI. CALCULATION OF RESULTS

The standard curve is established by plotting the mean cpm values or preferentially the respective binding rates B related to the total radioactivity T (B/T %) of the calibrators 1 - 6 on the ordinate, y-axis, versus their respective anti-GAD<sub>65</sub>-concentrations on the abscissa, x-axis.

It is recommended to use both axes in logarithmic scale.

The anti-GAD<sub>65</sub> concentrations of the controls and the unknown samples are directly read off in U/ml against the respective cpm or B/T values.


DIAsource anti-GAD<sub>65</sub> M may be used also with Computer Assisted Analysis using software with spline smoothing calculation, such as for sandwich-type assays (IRMA).

## XII. TYPICAL DATA

**Do not use for evaluation!**

Test tubes	cpm (mean)	$\frac{B}{T}$ %	U/ml
Total radioactivity T	31977	100	---
Calibrator 1	742	2.3	<b>0.1</b>
Calibrator 2	2464	7.7	<b>1</b>
Calibrator 3	3838	12.0	<b>3</b>
Calibrator 4	6683	20.9	<b>10</b>
Calibrator 5	13172	41.2	<b>35</b>
Calibrator 6	25901	80.8	<b>120</b>
Control serum	---	---	---
Patient 1	10077	31.5	<b>20.7</b>

Typical standard curve:

## XIII. PERFORMANCE AND LIMITATIONS

### A. Detection limit

The most appropriate and statistically reasonable definition of the lower detection limit of any assay is at present the so-called **functional assay sensitivity**.

This functional assay sensitivity generally represents that concentration which corresponds to the 10 % (within-assay) and to the 20 % (between assay) coefficient of variation in the respective precision profiles of the assay in the lower concentration range. Upon correct and thorough performance of the DIAsource anti- GAD<sub>65</sub> RIA kit, this value is found at approx. 0.7 U/ml.

Anti-GAD<sub>65</sub> values below this defined level of functional assay sensitivity do not meet the statistical criteria for reliability according to GLP (Good Laboratory Practice) and therefore cannot be distinguished from zero due to the statistically necessary certainty.

Anti-GAD<sub>65</sub> concentrations above approx. 0.7 U/ml, however, fulfill these criteria and are consequently assessed as valid.

### B. Specificity

The DIAsource anti-GAD<sub>65</sub> RIA kit using <sup>125</sup>Iodine labelled GAD<sub>65</sub> has showing 100% specificity and 100% sensitivity when comparing with an established reference assay.

The high quality of the tracer (<sup>125</sup>I-GAD<sub>65</sub>) ensures the exclusive reaction of anti-GAD<sub>65</sub> autoantibodies in direct assay principle. Any detectable cross reaction with autoantibodies to thyroglobulin, thyroidal peroxidase, to the TSH receptor, to IAA and to IA2 does not exist.

### C. Calibration

The units in the DIAsource anti-GAD<sub>65</sub> RIA kit are arbitrary units.

### D. Parallelism of standards and serum samples

Dilutions of specimen in anti-GAD<sub>65</sub> free human serum are determined according to their expected theoretical values with the DIAsource anti-GAD<sub>65</sub> RIA kit. In some cases the expected theoretical values are not determined due to the heterogeneous nature of the autoantibody population in human serum and in view of epitope specificity and affinity of the autoantibodies.

Sample Nr.	Intraassay (n=20)		Interassay (n=5 x 10)	
	Mean (U/ml)	CV (%)	Mean (U/ml)	CV (%)
1	0.8	10.0	0.6	20.7
2	8.2	6.8	8.4	5.2
3	45.9	1.3	45.6	1.1

### E. Limitations

Healthy individuals should be tested negative by using the DIAsource anti-GAD<sub>65</sub> RIA kit.

However, GAD<sub>65</sub> autoantibodies may be also present in a rare neurological disorder, Stiff-man Syndrome (SMS). GAD autoantibodies from patients with SMS are directed mainly to linear epitopes of GAD<sub>65</sub> and GAD<sub>67</sub>, whereas type 1 diabetes antibodies recognize only conformational epitopes of GAD<sub>65</sub>. In SMS, the antibodies titers are found much higher. For further quantification such sera could be diluted with factor 1:50 and 1:100 with GAD<sub>65</sub> autoantibody negative sera.

Any clinical diagnosis should not be based on the results of in vitro diagnostic method alone. Physicians are supposed to consider all clinical and laboratory findings possible to state a diagnosis.

## XIV. REFERENCE INTERVALS

anti-GAD <sub>65</sub>	
GAD <sub>65</sub> antibodies negative	< 1.0 U/ml
GAD <sub>65</sub> antibodies Greyzone	≥ 1.0 – < 2.0 U/ml
GAD <sub>65</sub> antibodies positive	≥ 2.0 U/ml

Normal range was established by Receiver Operating Characteristics (ROC) analysis.

It is recommended that each laboratory establishes its own normal and pathological reference ranges for serum anti-GAD<sub>65</sub> antibodies levels as usually done for other diagnostic parameters, too. Therefore, the above mentioned reference values provide a guide only to values which might be expected.

## XV. PRECAUTIONS AND WARNINGS

### Safety

For *in vitro* diagnostic use only.

This kit contains <sup>125</sup>I (half-life: 59 days), emitting ionizing X (27 keV) and γ (35 keV) radiations.

**This kit is for in vitro use only.** Follow the working instructions carefully. This instruction manual is valid only for the present kit with the given composition. An exchange of single components is not in agreement with CE regulations.

The expiration dates stated on the respective labels are to be observed. The same relates to the stability stated for reconstituted reagents.

All reagents should be kept at 2 - 8 °C before use in the original shipping container.

Some of the reagents contain small amounts (< 0.1 %) of sodium azide as a preservative. They must not be swallowed or allowed to come into contact with skin or mucosa. The possible formation of heavy metal azides in the drainage has to be prevented by sufficient rinsing with water.

Source materials derived from human body fluids or organs used in the preparation of this kit were tested and found negative for both Hepatitis an HIV antibody. However, no known test guarantees the absence of such viral agents. Therefore, handle all components and all patient samples as if potentially hazardous.

Since the kit contains radioactive material the following precautions should be observed:

- Do not smoke, eat or drink while handling radioactive material in any room designated for working with radioactive material,
- Always use protective gloves,
- Never pipette radioactive material by mouth,
- Wipe up spills promptly, washing the affected surface thoroughly with a decontaminant,
- Place contaminated tissues, tubes, bench covers, gloves etc. in a specially marked container, discard liquid and solid radioactive waste only as permitted by federal, state or local authorities and regulations.

It is the responsibility of the user of this product to handle radioactive material in accordance to the national rules given by law or other statements of the local authorities.

In any case GLP with all general and individual regulation has to be applied to the use of this kit.

**XVI. SUMMARY OF THE PROTOCOL**

1	Label test tubes*	Cal	Ctrl	Pat.-Sera 1, 2 etc.	T
2	Pipette Cal 1 - 6 Ctrl 1 - 2 Patient sera	20 µl	20 µl	20 µl	
3	Pipette Reconstituted Tracer	50 µl	50 µl	50 µl	50 µl
4	Incubate**	2 hours (at room temperature)			
5	Pipette Reconstituted Protein A component	50 µl	50 µl	50 µl	
6	Incubate**	1 hours (at room temperature)			
7	Pipette Cold precipitation buffer	1 ml	1 ml	1 ml	
8	Centrifuge	30 minutes at 4-8°C at 2000 x g			
9	Decant supernatant or Aspirate supernatant	leave tubes upside down on absorbent paper for 5 minutes  quantitatively.			
10	Count radioactivity	Counting time: 1 minute			

\* use conical tubes

\*\* Prior to incubation, agitate the tubes briefly in order to ensure homogeneous reaction conditions.

DIAsource Catalogue Nr : KIPM2070/KIPM2071	Revision nr : 200224/1
---	---------------------------

Revision date : 2020-02-24

**Other translations of this Instructions for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>**

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번호	명칭	구성		비고
		50 Test Kit	100 Test Kit	
1	Tracer( <sup>125</sup> I labelled GAD <sub>65</sub> )	1 vial	2 vial	
2	Protein A	1 vial	2 vial	
3	Calibrator 1-6	6 vial, 0.25ml	6vial,0.25ml	
4	Control I, II	2 vial	2 vials,0.25ml	
5	Reconstitution Buffer	1 vial,6ml	1 vial, 12ml	
6	Precipitation Buffer	1 vial,55ml	1 vial, 110ml	

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