# **PIIIP-IRMA**

***OCFK07-PIIIP***

# History

## Summary of change:

Previous Version: 220408	Current Version: 220602																																																																
<p><b>E. Interferences</b> The presence of bilirubin at concentrations of up to 0.25 mg/mL, hemoglobin up to 20 g/L and triglycerides up to 20 g/L have no effect on the assay results. The immuno-assay is protected against heterophilic antibodies. However, we cannot guarantee that this protection is exhaustive.</p>	<p><b>E. Interferences</b> The presence of bilirubin at concentrations of up to 0.5 mg/mL, hemoglobin up to 20 mg/mL and triglycerides up to 25 mg/mL have no effect on the assay results. The immuno-assay is protected against heterophilic antibodies. However, we cannot guarantee that this protection is exhaustive.</p>																																																																
<p><b>C. Precision</b> This was evaluated using 2 samples assayed 20 times in the same series (Intra-Assay) and 3 samples assayed in 12 different runs (Inter-Assay).</p> <p style="text-align: center;">INTRA-ASSAY PRECISION                      INTER-ASSAY PRECISION</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Serum</th> <th>N</th> <th>Concent (U/ml)</th> <th>CV (%)</th> <th>Serum</th> <th>N</th> <th>Concent (U/ml)</th> <th>CV (%)</th> </tr> </thead> <tbody> <tr> <td>A</td> <td>20</td> <td>1.5</td> <td>3.0</td> <td>D</td> <td>12</td> <td>1.1</td> <td>2.6</td> </tr> <tr> <td>B</td> <td>20</td> <td>7.7</td> <td>1.7</td> <td>E</td> <td>12</td> <td>1.9</td> <td>5.7</td> </tr> <tr> <td></td> <td></td> <td></td> <td></td> <td>F</td> <td>12</td> <td>6.9</td> <td>6.5</td> </tr> </tbody> </table> <p>CV: Coefficient of variation</p>	Serum	N	Concent (U/ml)	CV (%)	Serum	N	Concent (U/ml)	CV (%)	A	20	1.5	3.0	D	12	1.1	2.6	B	20	7.7	1.7	E	12	1.9	5.7					F	12	6.9	6.5	<p><b>C. Precision</b> This was evaluated using 2 samples assayed 20 times in the same series (Intra-Assay) and 3 samples assayed in 12 different runs (Inter-Assay).</p> <p style="text-align: center;">INTRA-ASSAY PRECISION                      INTER-ASSAY PRECISION</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Serum</th> <th>N</th> <th>Concent (U/ml)</th> <th>CV (%)</th> <th>Serum</th> <th>N</th> <th>Concent (U/ml)</th> <th>CV (%)</th> </tr> </thead> <tbody> <tr> <td>A</td> <td>20</td> <td>1.5</td> <td>3.5</td> <td>D</td> <td>12</td> <td>1.1</td> <td>2.7</td> </tr> <tr> <td>B</td> <td>20</td> <td>7.7</td> <td>7.5</td> <td>E</td> <td>12</td> <td>1.9</td> <td>6.0</td> </tr> <tr> <td></td> <td></td> <td></td> <td></td> <td>F</td> <td>12</td> <td>6.9</td> <td>6.8</td> </tr> </tbody> </table> <p>CV: Coefficient of variation</p>	Serum	N	Concent (U/ml)	CV (%)	Serum	N	Concent (U/ml)	CV (%)	A	20	1.5	3.5	D	12	1.1	2.7	B	20	7.7	7.5	E	12	1.9	6.0					F	12	6.9	6.8
Serum	N	Concent (U/ml)	CV (%)	Serum	N	Concent (U/ml)	CV (%)																																																										
A	20	1.5	3.0	D	12	1.1	2.6																																																										
B	20	7.7	1.7	E	12	1.9	5.7																																																										
				F	12	6.9	6.5																																																										
Serum	N	Concent (U/ml)	CV (%)	Serum	N	Concent (U/ml)	CV (%)																																																										
A	20	1.5	3.5	D	12	1.1	2.7																																																										
B	20	7.7	7.5	E	12	1.9	6.0																																																										
				F	12	6.9	6.8																																																										

CE

en

Read entire protocol before use.

## PIIP-IRMA

### I. INTENDED USE

Radioimmunoassay for the *in vitro* quantitative measurement of procollagen III peptide (PIIP) in serum.

The kit is intended for professional use.

### II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource PIIP-IRMA Kit
- B. **Catalog number :** OCFK07-PIIP : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.  
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

**For technical assistance or ordering information contact :**  
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

### III. CLINICAL BACKGROUND

#### A. Biological activity

Procollagen type III is synthesized in fibroblasts as a biosynthetic precursor of collagen type III, and then released. The propeptides are split off in the extracellular space during the conversion into collagen. The N-terminal propeptide col 1-3 (PIIP) is formed during this process in equimolar proportions to collagen type III and enters the circulation. The serum levels of N-terminal propeptide can therefore be used as a measurement of collagen III synthesis.

The procollagen III peptide col 1-3 (MW 45000) can be broken down by proteolysis into fractions col 1 (MW 10 000) which are not detected in the serum with DIAsource PIIP-IRMA Kit.

#### B. Clinical applications

The main collagens found in the connective tissue of the liver are types I and III, in addition to small amounts of type IV. If, as a result of pathological conditions, there is active proliferation of connective tissue (fibrosis) in the liver, increasing amounts of procollagen III peptide col 1-3 are formed.

The transformation of functioning liver tissue into connective tissue is detectable by means of the raised serum procollagen III peptide level. This is, for instance, the case in alcoholic or virus induced forms of liver fibrosis and cirrhosis. Raised procollagen III peptide levels can also occur with some other disorders and for this reason this parameter is particularly suitable for monitoring the course of

liver disorders.

Pathological conditions affecting the liver, which are associated with active proliferation of connective tissue, give rise to raised serum procollagen III peptide values. The actual transformation of functioning liver tissue into connective tissue can be detected by measuring procollagen III peptide in the serum.

According to the degree of severity of the disease, procollagen III peptide in the serum is raised in chronic active hepatitis (CAH), fibrosis and cirrhosis of the liver. Levels in chronic persistent hepatitis (CPH) are generally in the normal range; procollagen III peptide may be positively raised in degeneration of the liver. In acute hepatitis procollagen III peptide in the serum is also raised. There are, however, other conditions in which procollagen III peptide is raised without a detectable change in the liver (e.g. pulmonary fibrosis, rheumatic disorders, myocardial infarction, acromegaly, multiple trauma).

The diagnostic importance is not in the initial diagnosis but in monitoring the course of the disease. There is a good correlation with histological findings in fibrosis and cirrhosis.

#### IV. PRINCIPLES OF THE METHOD

DIASource PIIP-IRMA Kit permits the in vitro determination of human procollagen III peptide in human serum (or plasma) using the principle of a 2-step sandwich assay. A complex of solid-phase anti-PIIP antibodies (monoclonal, mouse), procollagen III peptide in the sample and <sup>125</sup>I-labelled anti-PIIP antibodies (monoclonal, mouse) is formed during this process. At the end of the reaction the free tracer is removed by decanting (or aspiration) and subsequent washing.

The amount of tracer specifically bound to the coated test tubes is measured with a gamma scintillation counter.

Evaluation of the unknowns is carried out by reading off from a calibrator curve constructed under identical conditions.

#### V. REAGENTS PROVIDED

Reagents	96 Tests Kit	Reconstitution			
 Tubes coated with anti-PIIP mouse monoclonal antibody.	2 x 48	<b>Ready</b> for use			
<table border="1" data-bbox="124 694 263 739"> <tr> <td>Ab</td> <td><sup>125</sup>I</td> </tr> </table> TRACER: anti-PIIP <sup>125</sup> I mouse monoclonal antibody, bovine gelatine, sodium azide	Ab	<sup>125</sup> I	1 vial 42 ml 315 kBq	<b>Ready</b> for use	
Ab	<sup>125</sup> I				
<table border="1" data-bbox="124 840 263 884"> <tr> <td>CAL</td> <td>0</td> </tr> </table> Zero Calibrator in bovine albumine, PIIP of calf serum, sodium azide.	CAL	0	1 vial lyophil.	<b>Add 0.5 ml</b> distilled water	
CAL	0				
<table border="1" data-bbox="124 952 263 996"> <tr> <td>CAL</td> <td>N</td> </tr> </table> Calibrators - N = 1 to 6 (see exact values on vial labels) in bovine albumine, PIIP of calf serum, sodium azide.	CAL	N	6 vials lyophil.	<b>Add 0.5 ml</b> distilled water	
CAL	N				
<table border="1" data-bbox="79 1108 319 1153"> <tr> <td>INC</td> <td>BUF</td> </tr> </table> Incubation buffer : mouse monoclonal antibody, sodium azide	INC	BUF	1 vial 42 ml	<b>Ready</b> for use	
INC	BUF				
<table border="1" data-bbox="79 1232 319 1276"> <tr> <td>WASH</td> <td>SOLN</td> <td>CONC</td> </tr> </table> Wash soln conc (NaCl, tween 20)	WASH	SOLN	CONC	1 vial 40 ml	<b>Dilute</b> 20 x with distilled water (use a magnetic stirrer).
WASH	SOLN	CONC			
<table border="1" data-bbox="79 1366 295 1411"> <tr> <td>CONTROL</td> <td>N</td> </tr> </table> Controls - N = 1 or 2 In human serum, sodium azide.	CONTROL	N	2 vials lyophil.	<b>Add 0.5 ml</b> distilled water	
CONTROL	N				

**Note :** Use the incubation buffer for sera dilutions.

To obtain values expressed in weight quantities (µg/L), it is necessary to multiply the results with the DIASource PIIP-IRMA Kit (U/mL) by a factor of 8.

#### VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Pipettes for delivery of: 20 µl, 400 µl and 500 µl (the use of accurate pipettes with disposable plastic tips is recommended)
3. Vortex mixer
4. Magnetic stirrer
5. Tube shaker
6. 5 ml automatic syringe (Cornwall type) for washing
7. Aspiration system (optional)
8. Any gamma counter capable of measuring <sup>125</sup>I may be used (minimal yield 70%).

#### VII. REAGENT PREPARATION

- A. Calibrators:** Reconstitute the calibrators with 0.5 ml distilled water.
- B. Controls:** Reconstitute the controls with 0.5 ml distilled water.
- C. Working Wash solution:** Prepare an adequate volume of Working Wash solution by adding 19 volumes of distilled water to 1 volume of Washing Reagent (20x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

#### VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, calibrators and controls are stable at -20°C until the expiry date.
- Avoid subsequent freeze-thaw cycles.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

#### IX. SPECIMEN COLLECTION AND PREPARATION

When blood samples have been taken, serum or plasma is obtained by the usual methods. The serum or plasma is used directly in the assay or stored for up to 3 days at 2-8°C. If they are to be stored for a longer period, a temperature of -20°C is recommended. The samples should be carefully mixed after thawing. Repeated thawing should be avoided.

PIIP values may be outside the measuring range in some patients suffering for example from alcoholic hepatitis and multiple trauma. In such cases it is advisable to dilute the patient's sample 1/5 with the incubation buffer.

The DIASource PIIP-IRMA Kit enables concentrations of 0.1 to 14 U/mL of PIIP to be measured.

For concentrations above 8 U / mL, it is recommended to dilute with the incubation buffer available into the kit.

#### X. PROCEDURE

##### A. Handling notes

Do not use the kit or components beyond expiry date.

Do not mix materials from different kit lots.

Bring all the reagents to room temperature prior to use.

Thoroughly mix all reagents and samples by gentle agitation or swirling. Use a clean disposable pipette tip for addition of each different reagent and sample in order to avoid cross-contamination. High precision pipettes or automated pipetting equipment will improve the precision.

Respect the incubation times.

Prepare a calibration curve for each run, do not use data from previous runs.

##### B. Procedure

1. Label coated tubes in duplicate for each calibrator, sample and control. For determination of total counts, label 2 normal tubes.
2. Briefly vortex calibrators, samples, controls and dispense 20 µl of each into the respective tubes.
3. Dispense 400 µL incubation buffer into each test tube.
4. Shake the test tubes on a horizontal shaker for 2 hours at 18-25°C (300 rpm).
5. Aspirate
6. Dispense 1 mL wash buffer into each test tube, aspirate. Renew one time this operation.
7. Dispense 400 µL <sup>125</sup>I-anti-PIIP solution into each tube.
8. Shake on the horizontal shaker for 3 hours at 18-25°C (300 rpm).
9. Aspirate
10. Dispense 1 mL wash buffer into each test tube, aspirate. Renew one time this operation.
11. Measure the tubes for 1 minute on a gamma counter.

#### XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. Calculate the bound radioactivity as a percentage of the binding determined at the zero calibrator point (0) according to the following formula :


$$B/B_0 (\%) = \frac{\text{Counts (Calibrator or sample)}}{\text{(Zero Calibrator)}} \times 100 \text{ Counts}$$

3. Plot the (B/B<sub>0</sub>(%)) values for each calibrator point as a function of the PIIP concentration of each calibrator point. Reject obvious outliers.
4. Computer assisted methods can be used to construct the calibration curve. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.
5. By interpolation of the sample (B/B<sub>0</sub>(%)) values, determine the PIIP concentrations of the samples from the calibration curve.
6. For each assay, the percentage of total tracer bound in the absence of unlabelled PIIP (B<sub>0</sub>/T) must be checked.

### XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

PIIP	Cpm	Concentration (U/ml)
Calibrator 0	1010	0
Calibrator 1	5673	0.49
Calibrator 2	14379	1.07
Calibrator 3	30079	1.94
Calibrator 4	63531	4.13
Calibrator 5	85575	6.82
Calibrator 6	106620	16.6
Control 1	28030	1.82
Control 2	83616	6.53

### XIII. PERFORMANCE AND LIMITATIONS

#### A. Detection limit

The detection limit assessed by analytical method is defined as the smallest detectable concentration different from zero with a probability of 95%. It was measured below to 0.02 U / mL.

#### B. Specificity

The monoclonal antibodies used in the kit are specific for the N-terminal procollagen III peptide col 1-3. The possibility of a cross-reaction with other basal membrane proteins occurring in the concentration ranges that are physiologically relevant can be virtually ruled out.

#### C. Precision

This was evaluated using 2 samples assayed 20 times in the same series (Intra-Assay) and 3 samples assayed in 12 different runs (Inter-Assay).

##### INTRA-ASSAY PRECISION

##### INTER-ASSAY PRECISION

Serum	N	Concent (U/ml)	CV (%)	Serum	N	Concent (U/ml)	CV (%)
A	20	1.5	3.5	D	12	1.1	2.7
B	20	7.7	7.5	E	12	1.9	6.0
				F	12	6.9	6.8

CV: Coefficient of variation

#### D. Accuracy

Samples with high concentrations of PIIP were diluted. The recovery percentages obtained ranged from 69 to 103%.

Known quantities of PIIP were added to human sera. The recovery percentages from the samples ranged from 81 to 117%.

#### E. Interferences

The presence of bilirubin at concentrations of up to 0.5 mg/mL, hemoglobin up to 20 mg/mL and triglycerides up to 25 mg/mL have no effect on the assay results. The immuno-assay is protected against heterophilic antibodies. However, we cannot guarantee that this protection is exhaustive.

### XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises.

### XV. REFERENCE INTERVALS

The clinical trials were conducted in different centres on a collective of 158 healthy volunteers (both sex), over 20 years old. In addition of this adult study, a specific one for children and young people is done.

It is recommended that each laboratory establish its own normal values. The values given below are indicative.


Adults normal values (Normal values: 0.3 – 0.8 U/mL)

Normale values	n	Mean U/mL	5. Perc. U/mL	95. Perc. U/mL	Min. U/mL	Max. U/mL
	158	0.53	0.33	0.79	0.21	0.98


Children normal values

Age	n	Mean U/mL	5. Perc. U/mL	95. Perc. U/mL	Min. U/mL	Max. U/mL
0 - 10	19	0.98	0.34	6.1	0.29	9.4
10 - 15	30	1.0	0.64	1.7	0.52	2.4
15 - 20	34	0.76	0.37	1.8	0.29	1.95
> 20	158	0.53	0.33	0.79	0.21	0.98

#### NORMAL POPULATION

#### Concentrations of P III P in serum of groups of different age

## XVI. PRECAUTIONS AND WARNINGS

### Safety

For *in vitro* diagnostic use only.

This radioactive product may only be received, purchased, stored or used by persons so authorized, and by laboratories covered by such authorization. The solution should under no circumstances be administered to humans or to animals. The purchase, storage, use and exchange of radioactive products are subject to the laws in force in the user's country. Enforcement of the basic radiation protection rules will ensure adequate security.

A summary of these is given below:

Radioactive products must be stored in their original containers in a suitable area. A record of the reception and storage of radioactive products must be kept up to date.

Handling of radioactive products should take place in a suitably equipped area with restricted access (controlled zone).

Do not eat, drink, smoke or apply cosmetics in a controlled zone. Do not mouth pipette radioactive solutions. Avoid any direct contact with all radioactive products by using laboratory coats and protective gloves. Contaminated laboratory equipment and glassware must be disposed of immediately after contamination to prevent cross contamination of different isotopes. Any contamination or radioactive substance loss should be dealt with in accordance with the established procedures. All radioactive waste disposal must be carried out according to the regulations in force.

The raw materials of human origin contained in the reagents of this kit have been tested with licensed kits and have been found to be negative for anti-HIV 1, anti-HIV 2 and anti-HCV antibodies and the HBs antigen. However, as it is still impossible to strictly guarantee that such products are incapable of transmitting hepatitis, the HIV virus or any other viral infection, all raw materials of human origin, including the samples to be assayed, must be treated as potentially infectious.

Do not pipette by mouth.

Do not smoke, eat or drink in areas in which samples or kit reagents are handled. Wear disposable gloves while handling kit reagents or samples and wash hands thoroughly afterwards.

Avoid splashing. Decontaminate and dispose of samples and all potentially contaminated materials as if they contained infectious agents. The best decontamination method is autoclaving for a minimum of one hour at 121.5°C. Some reagents contain sodium azide as a preservative. Avoid any absorption of the reagents, as well as any contact with the skin or mucous membranes. Sodium azide may react with lead or copper piping to form highly explosive metal azides. Dilute well when disposing of waste.

## XVII. BIBLIOGRAPHY

Collazos J, Diaz F. Role of the measurement of serum procollagen type III N-terminal peptide in the evaluation of liver diseases. Clin Chim Acta. 1994;227:37-43.

Heickendorff L, Frost L, Madsen JK, et al. Serum propeptides of type I and III procollagens in renal transplant recipients. A comparison of cyclosporine and azathioprine treatment. Nephron. 1994;67:203-8.

Hiramatsu N, Hayashi N, Kasahara A et al. Improvement of liver fibrosis in chorionic hepatitis C patients treated with natural interferon alpha. J Hepatol. 1995;22:135-42.

Jeffers LJ, Coelho-Little ME, Cheinquer H, et al. Procollagen-III peptide and chronic viral C hepatitis. Am J Gastroenterol. 1995;90:1437-40.

Klappacher G, Franzen P, Haab D, et al. Measuring extracellular matrix turnover in the serum of patients with idiopathic or ischemic dilated cardiomyopathy and impact on diagnosis and prognosis. Am J Cardiol. 1995;75:913-8.

Plebani M, Burlina A. Biochemical markers of hepatic fibrosis. Clin Biochem. 1991;24:219-39.

Yamauchi M, Mizuhara Y, Maezawa Y, Toda G. Serum tenascin levels in chronic liver disease. Liver. 1994;14:148-53.

Yudoh K, Matsui H, Kanamori M, Ohmori K, Tsuji H, Tatezaki S. Serum levels of laminin, type IV collagen and type III procollagen peptide as markers for detection of metastasis. Jpn J Cancer Res. 1994;85:1263-9.

## XVIII. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS µl	CALIBRATORS µl	SAMPLE (S) CONTROLS µl
Calibrators (0 to 6)	-	20	-
Samples, Controls	-	-	20
Incubation buffer	-	400	400
Incubation	2 hours at 18-25°C with continuous shaking (300 rpm)		
Separation		Aspirate	
Working Wash solution		1.0 ml	
Separation		Aspirate	
Working Wash solution		1.0 ml	
Separation		Aspirate	
Tracer	400	400	400
Incubation	3 hours at 18-25°C with continuous shaking (300 rpm)		
Separation		Aspirate	
Working Wash solution		1.0 ml	
Separation		Aspirate	
Working Wash solution		1.0 ml	
Separation		Aspirate	
Counting	Count tubes for 60 seconds		

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

DIAsource Catalogue Nr : OCFK07-PIIP	Revision nr : 220602
---	-------------------------

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**XVII. LITERATUR**

Collazos J, Diaz F. Role of the measurement of serum procollagen type III N-terminal peptide in the evaluation of liver diseases. Clin Chim Acta. 1994;227:37-43.  
 Heckendorff L, Frost L, Madsen JK, et al. Serum propeptides of type I and III procollagens in renal transplant recipients. A comparison of cyclosporine and azathioprine treatment. Nephron. 1994;67:203-8.  
 Hiramatsu N, Hayashi N, Kasahara A et al. Improvement of liver fibrosis in

chorionic hepatitis C patients treated with natural interferon alpha. J Hepatol. 1995;22:135-42.  
 Jeffers LJ, Coelho-Little ME, Cheinquer H, et al. Procollagen-III peptide and chronic viral C hepatitis. Am J Gastroenterol. 1995;90:1437-40.  
 Klappacher G, Franzen P, Haab D, et al. Measuring extracellular matrix turnover in the serum of patients with idiopathic or ischemic dilated cardiomyopathy and impact on diagnosis and prognosis. Am J Cardiol. 1995;75:913-8.  
 Plebani M, Burlina A. Biochemical markers of hepatic fibrosis. Clin Biochem. 1991;24:219-39.  
 Yamauchi M, Mizuhara Y, Maezawa Y, Toda G. Serum tenascin levels in chronic liver disease. Liver. 1994;14:148-53.  
 Yudoh K, Matsui H, Kanamori M, Ohmori K, Tsuji H, Tatezaki S. Serum levels of laminin, type IV collagen and type III procollagen peptide as markers for detection of metastasis. Jpn J Cancer Res. 1994;85:1263-9.

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