IVD

CE

CA 50 IRMA

OCFM07-CA50

Version : 220524

History

Summary of change:

Previous Version:	Current Version:
/	220524
	First creation

en

Read entire protocol before use.

CA 50 IRMA

I. INTENDED USE

Immunoradiometric assay kit for the in vitro quantitative measurement of CA 50 antigen in human serum.

The kit is intended for professional use.

II. GENERAL INFORMATION

- A. Proprietary name : DIAsource CA 50 IRMA Kit
- B. Catalog number : OCFM07-CA50 : 96 tests
- C. Manufactured by : DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

Tumor cells express substances in the cell membrane which are not usually produced in healthy cell membranes. The determination of these tumor-associated structures is a valuable tool in the diagnosis of malignant disorders. Using the hybridoma technique of Köhler and Milstein, specific immunological reagents (monoclonal antibodies, MAB) can be obtained which recognize tumor-associated antigens. A monoclonal antibody of this type, C-50 MAB, was obtained after immunization using a colorectal adenocarcinoma cell line, Colo 205. The C-50 MAB recognizes two different carbohydrate chains, the sialylated Lewis-a and the hitherto unknown sialylated lactotetraose. Structures containing CA-50 are mainly found in gastrointestinal cancers (e.g. pancreatic, stomach, hepatic and colorectal cancers) but also sometimes in other malignant growths (endometrial cancers). The CA-50 antigens occur in the cell membrane in a lipid-bound form (as ganglioside) and in a form bound to a high molecular weight protein (as glycoprotein). The CA-50 antigens are released by the tumors into the blood stream where they can be specifically detected by means of immunological techniques based on C-50 MAB.

CA-50 is normally detectable only in low concentrations in the serum of healthy men and women. A slight elevation in CA-50 concentrations can sometimes be observed in patients with benign disorders.

Pathologically raised serum CA-50 levels are encountered in the presence of CA-50-producing tumors, e.g. tumors of the pancreas, gastrointestinal tract, endometrium and bladder.

If a malignant disorder is suspected and elevated CA-50 levels are observed, additional tests are recommended. The main indication for determination of CA-50 is the follow-up of tumor patients, i.e. the monitoring of treatment efficacy and the monitoring of the course of the disease over time and determination of its prognosis. Benign disorders can lead to CA-50 values above the normal range; this frequently occurs, for example, in cases of acute and chronic pancreatitis, ulcerative colitis, Crohn's disease, cirrhosis of the liver and hepatitis.

IV. PRINCIPLES OF THE METHOD

The DIAsource CA 50 IRMA in vitro assay of the CA-50 antigen in human serum (or plasma) using the 2-step "sandwich" test principle. A complex of anti-CA-50 antibodies (monoclonal, mouse) bound to the tube wall, CA-50 in the sample and ^{125}I -labeled anti-CA-50 antibodies (monoclonal, mouse) is formed. After the reaction, the free tracer fraction is removed by aspiration and washing.

The amount of tracer specifically bound to the coated test tubes is measured with a gamma scintillation counter.

Evaluation of results provided by unknown samples is performed by reading off from a standard curve constructed under identical conditions.

V. REAGENTS PROVIDED

Reagents	96 tests Kit	Reconstitution
Tubes coated with anti CA50 (mouse monoclonal antibodies)	2 x 48	Ready for use
Ab ^{125}I Anti-CA 50 ^{125}I (mouse monoclonal antibodies) in buffer with non-specific mouse immunoglobulins, bovine serum albumin, sodium azide, mouse and red dye .	1 vial 22 ml 344 kBq	Ready for use
CAL 0 Zero calibrator in buffer and sodium azide	1 vial 0.5 ml liquid	Ready for use
CAL N Calibrators 1-6:human CA 50 antigen in human serum, buffer and sodium azide. (see exact value on vial labels)	6 vials 0.5 ml liquid	Ready for use
CONTROL N Control: human CA 50 antigen in human serum and sodium azide. (see exact value on vial label)	2 vials 0.5 ml liquid	Ready for use
BUFFER Incubation buffer with bovine serum albumin, bovine immunoglobulins, blue dye and sodium azide .	1 vial 25 ml liquid	Ready for use
WASH SOLN CONC Wash solution (Phosphate Buffer)	1 vial 50 ml	Dilute 25x with distilled water (use a magnetic stirrer).

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Precision micropipettes or similar , for delivery of: 50 μl and 200 μl (the use of accurate pipettes with disposable plastic tips is recommended).
3. Disposable plastic tubes
4. Vortex mixer
5. Magnetic stirrer
6. Tubes shaker : circular horizontal (300 rpm)
7. 5 ml automatic syringe (Cornwall type) for washing
8. Aspiration system (optional)
9. Any gamma scintillation counter capable of measuring ^{125}I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

Working Wash solution : Prepare an adequate volume of Working Wash solution by adding 24 volumes of distilled water to 1 volume of Wash Solution (25x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening, all kit components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

The assay is performed directly on serum. If the test is to be carried out within 3 days after the sample is taken, the serum samples can be stored at 2-8°C. If not, they should be divided into aliquots which must be stored frozen (-20°C). Avoid re-freezing a sample after defrosting. After defrosting or removal from the refrigerator, shake the samples well. Dilutions: If high CA-50 levels are suspected, dilutions should be prepared with incubation buffer.

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature (18-25°C) at least 30 minutes prior to use. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample. Thoroughly mix all reagents and samples by gentle agitation or swirling. High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times. Prepare a calibration curve for each run, do not use data from previous runs.

B. Procedure

1. Label coated tubes in duplicate for each calibrator, sample and control. For determination of total counts, label 2 uncoated disposable plastic tubes.
2. Briefly vortex calibrators, samples, control.
3. Respect the order in which reagents are to be added:
4. Distribute 50 μl of calibrator, controls, or sample to be assayed at the bottom of the tubes prepared for this purpose. Use a new pipette tip for each sample.
5. Distribute 200 μl of buffer in each tube.
6. Shake the tubes on a horizontal agitation system (300 rpm) for 2 hours at a temperature of between 18 and 25°C.
7. Wash the coated tubes as follows:
 - Aspirate the contents of the tubes as completely as possible.
 - Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
 - Add 1.0 mL of Working Wash solution to each tube.
 - Avoid foaming during the addition of the Working Wash solution.
 Repeat this process twice.
8. Add 200 μl of anti-CA 50 ^{125}I to all the tubes.
9. Shake the tubes on a horizontal agitation system for 2 hours at a temperature of between 18 and 25°C.
10. Wash the coated tubes (except total counts) as previously described.
11. After the last washing, let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
12. Measure the remaining radioactivity bound to the coated tubes with a gamma scintillation counter for 60 seconds.


XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. Plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of CA50 (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
3. Read the concentration for each control and sample by interpolation on the calibration curve.
4. Computer assisted data reduction will simplify these calculations. If automatic result processing is to be used, a 4-parameter logistic function curve fitting is recommended.

XII. TYPICAL CALIBRATION CURVE

The following data are for illustration only and should never be used instead of the real time calibration curve.

CA 50 IRMA		cpm
Calibrator	0.0 U/ml	868
	6.3 U/ml	2267
	12.1 U/ml	3348
	19.6 U/ml	5592
	37.5 U/ml	11899
	72.7 U/ml	25262
	159.1 U/ml	51091
Control	15.3 U/ml	4313
	51.1 U/ml	17064

XIII. PERFORMANCE AND LIMITATIONS

A. Detection limit

The limit of detection is defined as the smallest detectable concentration above 0. It has been assessed as 0.9 U/ml.

B. Interferences

The presence of bilirubin at concentrations of up to 0.5 mg/mL, hemoglobin up to 20 mg/mL and triglycerides up to 25 mg/mL have no effect on the assay results. The immuno-assay is protected against heterophilic antibodies. However, we cannot guarantee that this protection is exhaustive.

C. Precision

INTRA ASSAY				INTER ASSAY			
Serum	Replicate	Concentration (U/ml)	CV (%)	Serum	Replicate	Concentration (U/ml)	CV (%)
A	20	13.1	6.5	C	13	50.9	6.5
B	20	51.0	4.1	D	13	13.0	6.7

CV: Coefficient of variation

D. Accuracy

RECOVERY TEST:

Known quantities of CA-50 were added to human sera. The recovery percentages from the samples ranged from 102 to 120%.

DILUTION TEST:

Samples with high concentrations of CA-50 were diluted. The recovery percentages obtained ranged from 91 to 108%.

E. Specificity

The monoclonal antibodies used in this assay guarantee specific measurement of the CA-50 antigen.

XIV. LIMITATIONS

- Do not extrapolate sample values beyond the last calibrator. Dilute the samples concerned and retest.

- The DIAsource CA-50 kit enables concentrations of 0.2 to 175 U/mL of CA-50 to be measured.

XV. INTERNAL QUALITY CONTROL

- If the results obtained for Control are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Do not freeze-thaw more than twice.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises

XVI. REFERENCE INTERVALS

The normal range for CA-50 was determined using 57 samples of serum from healthy men and women. Statistical evaluation demonstrated a concentration of 22.9 U/mL for the 90th percentile and of 23.5 U/mL for the 95th percentile.

XVII. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

This kit contains ^{125}I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVIII. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS µl	CALIBRATORS µl	SAMPLE(S) CONTROL µl
Calibrators (0-6) Control or Samples Buffer	- - -	50 - 200	- 50 200
Incubation	2 hours at 18-25°C , with continuous shaking (300 rpm)		
Separation Working Wash solution Separation Working Wash solution Separation	- - - - -	Aspirate 1.0 ml Aspirate 1.0 ml Aspirate	
Tracer	200	200	200
Incubation	2 hours at 18-25°C , with continuous shaking (300 rpm)		
Separation Working Wash solution Separation Working Wash solution Separation	- - - - -	Aspirate 1.0 ml Aspirate 1.0 ml Aspirate	
Counting	Count tubes for 60 seconds		

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

DIAsource Catalogue Nr : OCFM07-CA50	Revision Number : 220524
---	-----------------------------

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