25OH Vitamin D Total ELISA 90'

KAP1971-F1

Version: 220516

History

Summary of change:

Previous Version: 220405				Current Version: 220516																																																																																																											
Old DiaSource logo 				New DiaSource logo 																																																																																																											
V. REAGENTS PROVIDED				V. REAGENTS PROVIDED																																																																																																											
<table border="1"> <thead> <tr> <th>INC</th><th>BUF</th><th></th><th></th></tr> </thead> <tbody> <tr> <td colspan="2">Incubation Buffer with casein and proclin</td><td>1 vial 20 ml</td><td>green Ready for use</td></tr> </tbody> </table> <table border="1"> <thead> <tr> <th>CHROM</th><th>TMB</th><th></th><th></th></tr> </thead> <tbody> <tr> <td colspan="2">Chromogenic solution TMB (Tetramethylbenzydine)</td><td>1 vial 12 ml</td><td>orange Ready for use</td></tr> </tbody> </table> <table border="1"> <thead> <tr> <th>STOP</th><th>SOLN</th><th></th><th></th></tr> </thead> <tbody> <tr> <td colspan="2">Stop solution HCl 1M</td><td>1 vial 12 ml</td><td>Ready for use</td></tr> </tbody> </table>				INC	BUF			Incubation Buffer with casein and proclin		1 vial 20 ml	green Ready for use	CHROM	TMB			Chromogenic solution TMB (Tetramethylbenzydine)		1 vial 12 ml	orange Ready for use	STOP	SOLN			Stop solution HCl 1M		1 vial 12 ml	Ready for use	<table border="1"> <thead> <tr> <th>INC</th><th>BUF</th><th></th><th></th></tr> </thead> <tbody> <tr> <td colspan="2">Incubation Buffer with casein and proclin</td><td>1 vial 30 ml</td><td>Ready for use</td></tr> </tbody> </table> <table border="1"> <thead> <tr> <th>CHROM</th><th>TMB</th><th></th><th></th></tr> </thead> <tbody> <tr> <td colspan="2">Chromogenic solution TMB (Tetramethylbenzydine)</td><td>1 vial 13 ml</td><td>Ready for use</td></tr> </tbody> </table> <table border="1"> <thead> <tr> <th>STOP</th><th>SOLN</th><th></th><th></th></tr> </thead> <tbody> <tr> <td colspan="2">Stop solution 0.2 M H₂SO₄</td><td>1 vial 13 ml</td><td>Ready for use</td></tr> </tbody> </table>				INC	BUF			Incubation Buffer with casein and proclin		1 vial 30 ml	Ready for use	CHROM	TMB			Chromogenic solution TMB (Tetramethylbenzydine)		1 vial 13 ml	Ready for use	STOP	SOLN			Stop solution 0.2 M H ₂ SO ₄		1 vial 13 ml	Ready for use																																																								
INC	BUF																																																																																																														
Incubation Buffer with casein and proclin		1 vial 20 ml	green Ready for use																																																																																																												
CHROM	TMB																																																																																																														
Chromogenic solution TMB (Tetramethylbenzydine)		1 vial 12 ml	orange Ready for use																																																																																																												
STOP	SOLN																																																																																																														
Stop solution HCl 1M		1 vial 12 ml	Ready for use																																																																																																												
INC	BUF																																																																																																														
Incubation Buffer with casein and proclin		1 vial 30 ml	Ready for use																																																																																																												
CHROM	TMB																																																																																																														
Chromogenic solution TMB (Tetramethylbenzydine)		1 vial 13 ml	Ready for use																																																																																																												
STOP	SOLN																																																																																																														
Stop solution 0.2 M H ₂ SO ₄		1 vial 13 ml	Ready for use																																																																																																												
VII. REAGENT PREPARATION				VII. REAGENT PREPARATION																																																																																																											
D. Working HRP conjugate solution: Prepare the solution according to the number of used strips, as indicated in the below table: for example for 6 strips (48 wells): 60 µl of concentrated conjugate and 30 µl of concentrated HRP to 6 ml of conjugate buffer.				D. Working HRP conjugate solution: Prepare the solution according to the number of used strips, as indicated in the below table: for example for 6 strips (48 wells): 130 µl of concentrated conjugate and 65 µl of concentrated HRP to 13 ml of conjugate buffer.																																																																																																											
<table border="1"> <thead> <tr> <th>Nb of strips</th><th>Volume of Conjugate Buffer (ml)</th><th>Volume of Concentrated Conjugate (µl)</th><th>Volume of Concentrated HRP (µl)</th></tr> </thead> <tbody> <tr><td>1</td><td>1</td><td>10</td><td>5</td></tr> <tr><td>2</td><td>2</td><td>20</td><td>10</td></tr> <tr><td>3</td><td>3</td><td>30</td><td>15</td></tr> <tr><td>4</td><td>4</td><td>40</td><td>20</td></tr> <tr><td>5</td><td>5</td><td>50</td><td>25</td></tr> <tr><td>6</td><td>6</td><td>60</td><td>30</td></tr> <tr><td>7</td><td>7</td><td>70</td><td>35</td></tr> <tr><td>8</td><td>8</td><td>80</td><td>40</td></tr> <tr><td>9</td><td>9</td><td>90</td><td>45</td></tr> <tr><td>10</td><td>10</td><td>100</td><td>50</td></tr> <tr><td>11</td><td>11</td><td>110</td><td>55</td></tr> <tr><td>12</td><td>12</td><td>120</td><td>60</td></tr> </tbody> </table>				Nb of strips	Volume of Conjugate Buffer (ml)	Volume of Concentrated Conjugate (µl)	Volume of Concentrated HRP (µl)	1	1	10	5	2	2	20	10	3	3	30	15	4	4	40	20	5	5	50	25	6	6	60	30	7	7	70	35	8	8	80	40	9	9	90	45	10	10	100	50	11	11	110	55	12	12	120	60	<table border="1"> <thead> <tr> <th>Nb of strips</th><th>Volume of Conjugate Buffer (ml)</th><th>Volume of Concentrated Conjugate (µl)</th><th>Volume of Concentrated HRP (µl)</th></tr> </thead> <tbody> <tr><td>1</td><td>3</td><td>30</td><td>15</td></tr> <tr><td>2</td><td>5</td><td>50</td><td>25</td></tr> <tr><td>3</td><td>7</td><td>70</td><td>35</td></tr> <tr><td>4</td><td>9</td><td>90</td><td>45</td></tr> <tr><td>5</td><td>11</td><td>110</td><td>55</td></tr> <tr><td>6</td><td>13</td><td>130</td><td>65</td></tr> <tr><td>7</td><td>15</td><td>150</td><td>75</td></tr> <tr><td>8</td><td>17</td><td>170</td><td>85</td></tr> <tr><td>9</td><td>19</td><td>190</td><td>95</td></tr> <tr><td>10</td><td>21</td><td>210</td><td>105</td></tr> <tr><td>11</td><td>23</td><td>230</td><td>115</td></tr> <tr><td>12</td><td>25</td><td>250</td><td>125</td></tr> </tbody> </table>				Nb of strips	Volume of Conjugate Buffer (ml)	Volume of Concentrated Conjugate (µl)	Volume of Concentrated HRP (µl)	1	3	30	15	2	5	50	25	3	7	70	35	4	9	90	45	5	11	110	55	6	13	130	65	7	15	150	75	8	17	170	85	9	19	190	95	10	21	210	105	11	23	230	115	12	25	250	125
Nb of strips	Volume of Conjugate Buffer (ml)	Volume of Concentrated Conjugate (µl)	Volume of Concentrated HRP (µl)																																																																																																												
1	1	10	5																																																																																																												
2	2	20	10																																																																																																												
3	3	30	15																																																																																																												
4	4	40	20																																																																																																												
5	5	50	25																																																																																																												
6	6	60	30																																																																																																												
7	7	70	35																																																																																																												
8	8	80	40																																																																																																												
9	9	90	45																																																																																																												
10	10	100	50																																																																																																												
11	11	110	55																																																																																																												
12	12	120	60																																																																																																												
Nb of strips	Volume of Conjugate Buffer (ml)	Volume of Concentrated Conjugate (µl)	Volume of Concentrated HRP (µl)																																																																																																												
1	3	30	15																																																																																																												
2	5	50	25																																																																																																												
3	7	70	35																																																																																																												
4	9	90	45																																																																																																												
5	11	110	55																																																																																																												
6	13	130	65																																																																																																												
7	15	150	75																																																																																																												
8	17	170	85																																																																																																												
9	19	190	95																																																																																																												
10	21	210	105																																																																																																												
11	23	230	115																																																																																																												
12	25	250	125																																																																																																												
IX. SPECIMEN COLLECTION AND PREPARATION				IX. SPECIMEN COLLECTION AND PREPARATION																																																																																																											
<ul style="list-style-type: none"> - Serum and heparinized plasma provide similar results. Y (Heparin plasma) = 0.95 x (serum) + 2.8 ng/ml, $R^2 = 0.97$, $n = 19$ 				<ul style="list-style-type: none"> - Serum and heparinized plasma provide similar results. Y (Heparin plasma) = 1.0444 x (serum) – 0.7913 ng/ml, $R^2 = 0.9872$, $n = 10$ 																																																																																																											
B. Procedure <ol style="list-style-type: none"> Pipette 75 µl of Incubation Buffer into all the wells. Pipette 100 µl of the working HRP conjugate solution into each well. Incubate the microtiterplate for 15 minutes at room temperature (18-25°C), on a plate shaker (400 rpm) 				B. Procedure <ol style="list-style-type: none"> Pipette 200 µl of Incubation Buffer into all the wells. Pipette 250 µl of the working HRP conjugate solution into each well. Incubate the microtiterplate for 15 minutes at room temperature (18-25°C), on a plate shaker (400 rpm) 																																																																																																											

D. Accuracy

RECOVERY TEST	
Added 25OH-Vit.D ₃ (ng/ml)	Recovery (%)
0	100
12.8	102
Added 25OH-Vit.D ₂ (ng/ml)	Recovery (%)
0	100
11.2	97

D. Accuracy

RECOVERY TEST	
Added 25OH-Vit.D ₃ (ng/ml)	Recovery (%)
5	97.1
10	95.1
25	86.7
50	93.9
Added 25OH-Vit.D ₂ (ng/ml)	Recovery (%)
5	99.8
10	91.5
25	94.8
50	85.6

DILUTION TEST

Sample dilution	Theoretical concen. (ng/ml)	Measured concen. (ng/ml)	Recovery (%)
1/1	74.9	74.9	100
1/2	37.4	37.9	101
1/4	18.7	19.2	103
1/8	9.3	10.2	110
1/1	56.3	56.3	100
1/2	28.1	28.7	102
1/4	14.0	14.8	106
1/8	7.0	5.8	83

DILUTION TEST

Sample dilution	Theoretical concen. (ng/ml)	Measured concentration (ng/ml)	Recovery (%)
1/1	-	62.75	100
1/2	31.38	30.10	95.9
1/4	15.69	14.69	93.6
1/8	7.84	7.39	94.2

XVIII. SUMMARY OF THE PROTOCOL

	CALIBRATORS (μ l)	SAMPLE(S) CONTROLS (μ l)
Calibrators (0-5) Controls, Samples Incubation Buffer	25 - 75	- 25 75
Working HRP Conjugate	100	100

XVIII. SUMMARY OF THE PROTOCOL

	CALIBRATORS (μ l)	SAMPLE(S) CONTROLS (μ l)
Calibrators (0-5) Controls, Samples Incubation Buffer	25 - 200	- 25 200
Working HRP Conjugate	250	250

en

CE

Read entire protocol before use.

25OH Vitamin D Total ELISA 90'

I. INTENDED USE

Immunoenzymetric assay for the *in vitro* quantitative measurement of 25-hydroxyvitamin D2 and D3 (25OH-D2 and 25OH-D3) in serum and plasma.

II. GENERAL INFORMATION

- A. Proprietary name: DIAsource 25OH Vitamin D Total ELISA 90' Kit
- B. Catalog number: KAP1971-F1: 96 tests
- C. Manufactured by: DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2 B-1348 Louvain-la Neuve, Belgium.

For technical assistance or ordering information contact:
Tel: +32 (0)10 84 99 11 Fax: +32 (0)10 84 99 90

III. CLINICAL BACKGROUND

Vitamin D is the generic term used to designate Vitamin D2 or ergocalciferol and Vitamin D3 or cholecalciferol.

Humans naturally produce Vitamin D3 when the skin is exposed to ultraviolet sun rays.

In the liver mainly, Vitamin D3 is metabolised into 25-Hydroxyvitamin D3 (25OH D3) which is the main form of Vitamin D circulating in the body.

25OH D3 is a precursor for other Vitamin D metabolites and has also a limited activity by itself.

The most active derivative is 1,25-hydroxyvitamin D3, produced in the kidney (or placenta) by 1-hydroxylation of 25OH D3.

25OH Vitamin D stimulates the intestinal absorption of both calcium and phosphorus and also bone resorption and mineralisation.

25OH Vitamin D might also be active in other tissues responsible for calcium transport (placenta, kidney, mammary gland ...) and endocrine gland (parathyroid glands, beta cells...).

Vitamin D3 and Vitamin D2 are also available by ingestion through food or dietary supplementation.

As Vitamin D2 is metabolised in a similar way to Vitamin D3, both contribute to the overall Vitamin D status of an individual.

It is the reason why it is very important to measure both forms of 25OH Vitamin D equally for a correct diagnosis of Vitamin D deficiency, insufficiency or intoxication.

Vitamin D deficiency is an important risk factor for rickets, osteomalacia, senile osteoporosis, cancer and pregnancy outcomes.

The measurement of both 25OH Vitamin D forms is also required to determine the cause of abnormal serum calcium concentrations in patients.

Vitamin D intoxication has been shown to cause kidney and tissue damages.

IV. PRINCIPLES OF THE METHOD

The DiaSource 25OH Vitamin D Total ELISA 90' is a solid phase Enzyme Linked Immunosorbent Assay performed on microtiterplates. During a first 1 hour incubation step, at room temperature (18-25°C), total 25OH Vitamin D (D_2 and D_3) present in calibrators, controls and samples is dissociated from binding serum proteins to fix on binding sites of a specific monoclonal antibody. After 1 washing step, a fixed amount of 25OH Vitamin D-labelled with biotin in presence of horseradish peroxidase (HRP), compete with unlabelled 25OH Vitamin D_2 and 25OH Vitamin D_3 present on the binding sites of the specific monoclonal antibody. After a 15 minutes incubation at room temperature (18-25°C), the microtiterplate is washed to stop the competition reaction. The Chromogenic solution (TMB) is added and incubated for 15 minutes. The reaction is stopped with the addition of Stop Solution and the microtiterplate is then read at the appropriate wavelength. The amount of substrate turnover is determined colourimetrically by measuring the absorbance, which is inversely proportional to the total 25OH Vitamin D (D_2 and D_3) concentration.

A calibration curve is plotted and the total 25OH Vitamin D (D_2 and D_3) concentrations of the samples are determined by dose interpolation from the calibration curve.

V. REAGENTS PROVIDED

Reagents		KAP 1971-F1	Reconstitution
LJ	Microtiterplate (96 breakable wells) with anti 25OH Vit. D_2 and D_3 (monoclonal antibodies)	96 wells	Ready for use
CAL	0	1 vial lyophilised	Add 1 ml distilled water
CAL	N	5 vials lyophilised	Add 1 ml distilled water
CONTROL	N	2 vials lyophilised	Add 1 ml distilled water
INC	BUF	1 vial 30 ml	Ready for use
CONJ	CONC	1 vial 0.3 ml	Dilute 100 x with conjugate buffer
CONJ	BUF	1 vial 30 ml	Ready for use
HRP	CONC	1 vial 0.2 ml	Dilute 200 x with conjugate buffer
WASH	SOLN	1 vial 10 ml	Dilute 200 x with distilled water (use a magnetic stirrer).
CHROM	TMB	1 vial 13 ml	Ready for use
STOP	SOLN	1 vial 13 ml	Ready for use
Stop solution 0.2M H_2SO_4			

Note: Use Calibrator 0 for dilution of samples with values above the highest calibrator.

No international reference material is available.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Pipettes for delivery of: 25 μ l, 100 μ l, 200 μ l, 250 μ l and 1 ml (the use of accurate pipettes with disposable plastic tips is recommended)
3. Vortex mixer
4. Magnetic stirrer
5. Plate shaker (400 rpm)
6. Washer for microtiterplates
7. Microtiterplate reader capable of reading at 450 nm and 650 nm (bichromatic reading)

VII. REAGENT PREPARATION

- A. **Calibrator 0:** Reconstitute the calibrator 0 with 1 ml distilled water
- B. **Calibrators 1 - 5:** Reconstitute the calibrators 1-5 with 1 ml distilled water
- C. **Controls:** Reconstitute the controls with 1 ml distilled water.
- D. **Working HRP conjugate solution:**
The working HRP conjugate solution is to be prepared during the incubation and minimum 45 minutes before its use.
(cf X.B.5)

Prepare an adequate volume of working HRP conjugate solution by mixing the 3 reagents in the following sequence: (1) Conjugate buffer, (2) Concentrated Conjugate, (3) Vortex, (4) Concentrated HRP, (5) Vortex.

The order of addition of those 3 reagents is critical and should be rigorously respected to get reproducible Optical Densities.

Prepare the solution according to the number of used strips, as indicated in the below table: for example for 6 strips (48 wells): 130 μ l of concentrated conjugate and 65 μ l of concentrated HRP to 13 ml of conjugate buffer. Use a vortex to homogenize.

Until its use, keep the working HRP conjugate at room temperature (18-25°C) and avoid direct sunlight or use a brown glass vial for its preparation.

The preparation of working HRP conjugate is not stable and must be discarded if not used.

Nb of strips	Volume of Conjugate Buffer (ml)	Volume of Concentrated Conjugate (μ l)	Volume of Concentrated HRP (μ l)
1	3	30	15
2	5	50	25
3	7	70	35
4	9	90	45
5	11	110	55
6	13	130	65
7	15	150	75
8	17	170	85
9	19	190	95
10	21	210	105
11	23	230	115
12	25	250	125

E. Working Wash solution: Prepare an adequate volume of Working Wash solution by adding 199 volumes of distilled water to 1 volume of Wash Solution (200x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, calibrators and controls are stable for eight weeks at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 4 months. Avoid subsequent freeze-thaw cycles.
- Freshly prepared Working Wash solution should be used on the same day.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- This kit is suitable for serum and plasma heparinized samples.
- Serum and plasma heparinized samples samples must be kept at 2-8°C.
- If the test is not run within 24 hrs, **sampling and storage at -20°C is recommended.**
- Avoid subsequent freeze-thaw cycles.
- Serum and heparinized plasma provide similar results.
- Y (Heparin plasma) = $1.0444 \times (\text{serum}) - 0.7913 \text{ ng/ml}$, $R^2 = 0.9872$, $n = 10$

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date.
Do not mix materials from different kit lots.
Bring all the reagents to room temperature (18-25°C) prior to use.
Thoroughly mix all reagents and samples by gentle agitation or swirling.
Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.
Use a clean plastic container to prepare the Wash Solution.
In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
For the dispensing of the Chromogenic Solution and the Stop Solution avoid pipettes with metal parts.
High precision pipettes or automated pipetting equipment will improve the precision.
Respect the incubation times.
To avoid drift, the time between pipetting of the first calibrator and the last sample must be limited to the time mentioned in section XIII paragraph E (Time delay).
Prepare a calibration curve for each run, do not use data from previous runs.
Dispense the Chromogenic Solution within 15 minutes following the washing of the microtiterplate.
During incubation with Chromogenic Solution, avoid direct sunlight on the microtiterplate.
Each well can only be used once

B. Procedure

1. Select the required number of strips for the run. The unused strips should be resealed in the bag with a desiccant and stored at 2-8°C.
2. Secure the strips into the holding frame.
3. Pipette 25 µl of each Calibrator, Control and Sample into the appropriate wells.
4. Pipette 200 µl of Incubation Buffer into all the wells.
5. Incubate for 1 hour at room temperature (18-25°C), on a plate shaker (400 rpm)
6. Prepare the Working HRP conjugate solution during the incubation and minimum 45 minutes before its use.
7. Aspirate the liquid from each well.
8. Wash the plate 3 times by:
 - dispensing 0.35 ml of Wash Solution into each well
 - aspirating the content of each well
9. Pipette 250 µl of the working HRP conjugate solution into each well
10. Incubate the microtiterplate for 15 minutes at room temperature (18-25°C), on a plate shaker (400 rpm)
11. Aspirate the liquid from each well.
12. Wash the plate 3 times by:
 - dispensing 0.35 ml of Wash Solution into each well
 - aspirating the content of each well
13. Pipette 100 µl of the Chromogenic solution into each well within 15 minutes following the washing step.
14. Incubate the microtiterplate for 15 minutes at room temperature (18-25°C), on a plate shaker (400 rpm), avoid direct sunlight.
15. Pipette 100 µl of Stop Solution into each well.
16. Read the absorbances at 450 nm (reference filter 630 nm or 650 nm) within 1 hour and calculate the results as described in section XI

XI. CALCULATION OF RESULTS

1. Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).
2. Calculate the mean of duplicate determinations.
3. We recommend the use of computer assisted methods to construct the calibration curve. 4-parameter logistic function curve fitting is the preferred method. Reject obvious outliers.
4. By interpolation of the sample OD values, determine the 25OH Vitamin D concentrations of the samples from the calibration curve.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

25OH-ELISA		OD units
Calibrator	0 ng/ml 4.88 ng/ml 13.10 ng/ml 27.39 ng/ml 56.08 ng/ml 121.30 ng/ml	2.600 2.428 2.105 1.723 0.906 0.282

Note: 1 ng/ml = 2.5 pmol/ml

XIII. PERFORMANCE AND LIMITATIONS

A. Detection Limit

The LOB (Limit of blank) was calculated by measuring the blank several times and was calculated as the mean – 1.65 standard deviation of the distribution of the test values. The LOB was calculated to be 2.012 ng/ml.

The LOD (Limit of detection) was calculated as the LOB – 1.645 standard deviation of a low concentration sample tested in 10 different run. The LOD was calculated to be 3.035 ng/ml.

The LOQ (Limit of quantitation) was calculated by testing 8 samples of low values, 10 times. The LOQ was calculated to be 5.91 ng/ml.

B. Specificity

Cross reactivity of the 25OH Vitamin D Total ELISA 90' assay was determined by testing sera with spiked and unspiked cross reactants. The results are summarized in the following table:

Compound and Concentration	% Cross reaction
25OH-Vitamin D ₃ at 10 ng/mL	100
25OH-Vitamin D ₂ at 10 ng/mL	80.1
1,25(OH) ₂ -Vitamin D ₃ at 200 ng/mL	3.06
1,25(OH) ₂ -Vitamin D ₂ at 667 ng/mL	0.73
Vitamin D ₃ at 200 ng/mL	1.16
Vitamin D ₂ at 200 ng/mL	1.02
3-epi-25OH-Vitamin D ₃ at 20 µg/mL	0.07

The assay performance is not affected by hemolysis (2.5-5g/L hemoglobin tested), bilirubinemia (0.5-1g/L bilirubin conjugated and unconjugated tested) or triglycerides (0.625—2.5 g/L tested). Ascorbic acid (Vitamin C) (0.01-1g/L tested) and biotin (0.002-0.04g/L tested) don't interfere with this assay.

C. Precision

INTRA-ASSAY				INTER-ASSAY			
Sample	N	$\bar{X} \pm SD$ (ng/ml)	C.V. (%)	Sample	N	$\bar{X} \pm SD$ (ng/ml)	C.V. (%)
A	24	21.55 ± 0.85	3.92	A	10	18.72 ± 1.16	6.2
B	24	45.24 ± 1.45	3.20	B	10	67.90 ± 2.23	3.3

SD: Standard Deviation; CV: Coefficient of variation

D. Accuracy

RECOVERY TEST	
Added 25OH-Vit.D ₃ (ng/ml)	Recovery (%)
5	97.1
10	95.1
25	86.7
50	93.9
Added 25OH-Vit.D ₂ (ng/ml)	Recovery (%)
5	99.8
10	91.5
25	94.8
50	85.6

DILUTION TEST			
Sample dilution	Theoretical concentration (ng/ml)	Measured concentration (ng/ml)	Recovery (%)
1/1	-	62.75	100
1/2	31.38	30.10	95.9
1/4	15.69	14.69	93.6
1/8	7.84	7.39	94.2

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when incubation buffer is dispensed 10, 15 and 20 minutes after the sample has been added in the coated wells.

TIME DELAY				
	0' (ng/ml)	10' (ng/ml)	15' (ng/ml)	20' (ng/ml)
Sample 1	18.0	20.1	21.2	23.0
Sample 2	40.1	41.4	37.6	41.5

XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls which contain azide will interfere with the enzymatic reaction and cannot be used.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises
- It is recommended that Controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.
- It is good practise to check visually the curve fit selected by the computer.

XV. EXPECTED VALUES

Dietary intake, race, season and age are known to affect the normal levels of 25OH.Vit.D₃.

Each laboratory should establish its own range based on their local population.

Recent literature has suggested the following ranges for the classification of 25 OH Vitamin D status:

Level	ng/mL
Deficient	<10
Insufficient	10-29
Sufficient	30-100
Potential Toxicity	>100

XVI. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for

HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures. All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious. Avoid any skin contact with all reagents, Stop Solution contains H₂SO₄. In case of contact, wash thoroughly with water. Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves. For more information, see Material Safety Data Sheet (MSDS)

XVII. BIBLIOGRAPHY

- ZERWEKH J.E. (2008) **Blood biomarkers of Vitamin D status.** Am. J. Clin. Nutr., 87(suppl):1087S-91S.
- HOLICK M.F. (2006) **Resurrection of Vitamin D deficiency and rickets.** J. Clin. Invest., 116:2062-2072.
- HEANEY R.P. (2000) **Vitamin D: how much do we need, and how much is too much.** Osteoporos. Int., 11:553-555.
- DAWSON-HUGHES B., HEANEY R.P., HOLICK M.F., LIPS P., MEUNIER P.J. (1997) **Prevalence of Vitamin D insufficiency in an adult normal population.** Osteoporos. Int., 7:439-443.
- BISCHOFF-FERRARI H.A., GIOVANNUCCI E., WILLETT W.C., DIETRICH T., DAWSON-HUGHES B. (2006) **Estimation of optimal serum concentrations of 25-hydroxyvitamin D for multiple health outcomes.** Am. J. Clin. Nutr., 84:18-28.
- HOLICK M.F.(2004) **Sunlight and vitamin D for bone health and prevention of autoimmune diseases, cancers and cardiovascular disease.** Am. J. Clin. Nutr., 80:1678S-1688S.
- HEANEY R.P. (2010) **Defining deficiency of vitamin D.** Clinical Laboratory International, October 2010, vol.34 : 16-19.
- HOLICK M.F. (2007) **Vitamin D deficiency.** N. Engl. J. Med., 357:266-281.
- TAHA N. M , VIETH R.(2010) **The problem of an optimal target level for 25-Hydroxyvitamin D, the test for vitamin D nutritional status.** Clinical Laboratory International, November 2010, vol.34 : 28-30
- HOLICK M.F. (2009) **Vitamin D Status: Measurement, Interpretation, and Clinical Application** Ann. Epidemiol., 19:73-78.

XVIII. SUMMARY OF THE PROTOCOL

	CALIBRATORS (μ l)	SAMPLE(S) CONTROLS (μ l)
Calibrators (0-5)	25	-
Controls, Samples	-	25
Incubation Buffer	200	200
Incubate for 1 hour at room temperature (18-25°C) with continuous shaking at 400 rpm.		
Prepare the working HRP conjugate during the incubation and minimum 45 minutes before its use. The sequence of preparation is critical, see VII. Reagent Preparation.		
Aspirate the contents of each well.		
Wash 3 times with 350 μ l of Wash Solution and aspirate.		
Working HRP Conjugate	250	250
Incubate for 15 minutes at room temperature (18-25°C) with continuous shaking at 400 rpm.		
Aspirate the contents of each well.		
Wash 3 times with 350 μ l of Wash Solution and aspirate.		
Chromogenic Solution	100	100
Incubate for 15 min at room temperature (18-25°C) with continuous shaking at 400 rpm.		
Stop Solution	100	100
Read on a microtiterplate reader. Record the absorbance of each well at 450 nm (versus 630 or 650 nm).		

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

DIAsource Catalogue Nr : KAP1971-F1	Revision nr : 220516
--	-------------------------

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