



# IGFBP-2 Elisa

*KAPME05*



# History

## Summary of change :

Previous Version :	Current Version :
200224/1	210527
<b>INTENDED USE</b> Measurement of human IGFBP-2 in human serum, EDTA-plasma, cerebrospinal fluid, breast milk, amniotic fluid, saliva and in cell culture medium.	<b>INTENDED USE</b> Measurement of human IGFBP-2 in human serum and EDTA-plasma.
<b>Sample Preparation and Storage</b> Serum and plasma samples as well as cell culture medium, breast milk, amniotic fluid, cerebrospinal fluid and saliva are applicable. The blood sample for serum preparation should be gained according to standardized venipuncture procedure. The samples should be stored without anticoagulation reagents. Hemolytic reactions have to be avoided. The blood has to be allowed to clot and after complete clotting, serum is separated by centrifugation. Blood samples may be taken at any time of the day. Storage at RT      max. 2 days Storage at -20°C    max. 2 years Are not allowed to have more than 10 freeze/thaw cycles.	<b>Sample Preparation and Storage</b> Serum and plasma samples are applicable. The blood sample for serum preparation should be gained according to standardized venipuncture procedure. The samples should be stored without anticoagulation reagents. Hemolytic reactions have to be avoided. The blood has to be allowed to clot and after complete clotting, serum is separated by centrifugation. Blood samples may be taken at any time of the day and should be stored in firmly closed sample vials. Haemolytic reactions are to be avoided. Storage at RT      max. 24 hours Storage at -20°C    max. 2 years Are not allowed to have more than 10 freeze/thaw cycles.
<b>Materials not Provided</b> NA	<b>Materials not provided</b> Addition of "Polyethylene PE/Polypropylene PP tubes for dilution of samples.
<b>Technical hints</b> No information regarding the washing and shaking steps.	<b>Technical hints</b> Information added
<b>PROCEDURE</b> wash the wells 3 times with 250 µl of Washing Buffer / well respectively.	<b>Procedure</b> wash the wells 5 times with 300 µl of Washing Buffer / well respectively.
No information regarding the quality control	Information added
<b>PERFORMANCE CHARACTERISTICS</b> <b>Sensitivity</b> The analytical sensitivity of the assay yields 0.2 ng/ml (2x SD of zero calibrator)	<b>PERFORMANCE CHARACTERISTICS</b> <b>Sensitivity</b> Sensitivity was assed by measuring the blank and calculating the theoretical concentration of the blank + 2SD. The analytical sensitivity of the E05 is 0.2 ng/mL.
<b>PERFORMANCE CHARACTERISTICS</b> <b>Table 5:</b> Linearity of the sample dilution: Dilution test performed on Two serum samples used as matrix, Cerebrospinal fluid and Amniotic fluid.	<b>PERFORMANCE CHARACTERISTICS</b> <b>Table 5:</b> Linearity of the sample dilution: Dilution test performed on Two serum samples used as matrix
<b>Interference</b> No data for Hemoglobin interference	<b>Interference</b> Data for Hemoglobin interference added
No information regarding Instructions for use for scientific application	Information added



# IGFBP-2-ELISA

Enzyme Immunoassay for the Quantitative Determination of  
Insulin-like Growth Factor Binding Protein-2  
KAPME05  
*IN VITRO DIAGNOSTIC USE*

en

DIAsource ImmunoAssays S.A. - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 91

## INTENDED USE

Measurement of human IGFBP-2 in human serum and EDTA-plasma.

## CLINICAL IMPLICATIONS

The IGFBP-2 concentration is age-dependent in blood (3).

Normal values for healthy individuals (1.5 to > 70 years) were evaluated for this assay.

Supplementary parameter to IGFBP-3 in the diagnosis of growth disorders (IGFBP-2/IGFBP-3 ratio), IGFBP-2 is an inhibitor of growth hormone action (3,4).

Progression-dependent tumor marker in leukaemia (5), astrocytic CNS tumors (6,7), prostate- (8), suprarenal cortex-(9)-, hepatocellular (10) and other carcinomas.

Anti-aging parameter: IGFBP-2 as a marker of physiological functionality (20).

## INTRODUCTION


Insulin-like growth factors (IGFs) regulate the proliferation, differentiation, apoptose, cell adhesion and metabolism in various tissues and cell types. The IGF-I, which is produced mainly in liver under the influence of growth hormone (GH), regulates as hormone the linear growth of the bones and the process of sexual maturity, while IGF-II is mainly a growth factor of fetal tissue (11-13). The biological actions of IGF over the IGF-Type-I receptor are modulated variably through the IGF binding proteins (IGFBP-1 to-6) (14). IGFBP-2 is, after IGFBP-3, the second most frequent IGFBP in the human blood. IGFs, especially tumor typical pro-IGF-forms and hormones regulate the expression of IGFBP-2, GH effect is thereby inhibiting. At cellular level IGFBP-2 seems to stimulate the proliferation and dissemination of solid tumors via an IGF-independent mechanism (15,16).

## PHYSIOLOGICAL MEANING

IGFBP-2 is a unglycosylated polypeptide of 31.3 kDa, which forms binary IGF-complexes and shows no circadian rhythm in the circulation. The serum concentration of IGFBP-2 increases in fasting, after major surgery and after trauma, but the increasing of the concentration is most intensive in malignant diseases. The correlation of the IGFBP-2 level to the degree of progression is a striking feature in various tumor types as is the normalization of the IGFBP-serum levels after remission (5-10). During the GH-therapy, e.g. in short stature and in GH-abuse (doping) the IGFBP-2 level decreases. In Trisomy 18 IGFBP-2 in maternal serum is decreased and IGFBP-1 is increased; therefore the ratio IGFBP-2 /IGFBP-1 is a marker for this chromosome abnormality (17).

## MATERIALS

### Materials Provided

- 1)  **Microtiter plate**, ready for use: **Microtiter plate** with 96 wells, divided up in 12 strips with 8 wells separately breakable, coated with anti-IGFBP-2 antibody and packed in a laminate bag.
- 2) 

CAL	N
-----	---

**Calibrators 1-5**, lyophilized: contain recombinant IGFBP-2: Calibrator values are between 2 - 80 ng/ml (2, 10, 20, 40, 80 ng/ml) IGFBP-2 and have to be reconstituted with **750 µl Dilution Buffer** each.
- 3) 

CONTROL	N
---------	---

**2 Controls**, lyophilized: Contains human serum and has to be reconstituted with **100 µl Dilution Buffer**. The exact concentration of IGFBP-2 is given on the vial label.
- 4) 

DIL	BUF
-----	-----

**Dilution Buffer**, 50 ml, ready for use.
- 5) 

Ab	BIOT
SAV	HRP

**Antibody-Conjugate**, 12 ml, ready for use, contains a mixture of rabbit-anti-hIGFBP-2-antibody biotinylated + streptavidin horseradish peroxidase conjugate ready for use. Use 100 µl/well in the assay.
- 6) 

WASH	SOLN	CONC
------	------	------

**Washing Buffer**, 50 ml, 20-fold concentrated: The **Washing Buffer** has to be diluted 1:20 with distilled water or demineralised water before use (e.g. add the complete contents of the flask (50 ml) into a graduated flask and fill up with A.dest. to 1000 ml). Please dilute only according to requirements. The diluted washing buffer is stable for 4 weeks at 2-8°C. It has to be at room temperature for usage !

- 7) 

CHROM	TMB
-------	-----

**TMB-Substrate Solution** 12 ml, ready for use, horseradish-peroxidase-(HRP)-substrate, stabilised Tetramethylbencidine.
- 8) 

STOP	SOLN
------	------

**Stopping Solution**, 0.2 M sulphuric acid, 12 ml, ready for use. *Caution Acid !*
- 9) **Sealing tape** for covering of the microtiter plate, 2 x

#### Materials not Provided

- Precision pipettes (100 and 200µl) Micropipettes and multichannel pipettes with disposable plastic tips
- Distilled or Deionized water for dilution of the washing buffer
- Polyethylene PE/Polypropylene PP tubes for dilution of samples
- Micropipettes and multichannel pipettes with disposable plastic tips
- Vortex-mixer
- Plate washer and plate shaker ( $\geq 350$  rpm) (recommended)
- Calibrated Microplate reader ("ELISA-Reader") with filter for 450/620nm wavelength
- Timer (120 min. range)

#### WARNING AND PRECAUTIONS

The kit should not be used beyond the expiration date on the kit label.

For In Vitro Diagnostic Use only. For Professional use only.

The DIAsource kit is suitable only for in vitro diagnostics and not for internal use in humans and animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. DIAsource will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided.

Do not use obviously damaged or microbial contaminated or spilled material.

Caution: This kit contains material of human and/or animal origin. Therefore all components and patient's specimens should be treated as potentially infectious.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations. Material Safety Data Sheet is available on request.

#### Human Serum

Following components contain human serum: Controls.

Source human serum for the control sera provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV). No known methods can offer total security of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

#### Reagents Calibrators 1-5, Antibody-Conjugate, Dilution Buffer, Washing Buffer

Contain as preservative a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one ( $<0.015\%$ )

- |           |   |
|-----------|---|
| H317      | May cause an allergic skin reaction.  |
| P280      | Wear protective gloves/ protective clothing/ eye protection/ face protection.                           |
| P272      | Contaminated work clothing should not be allowed out of the workplace.                                  |
| P261      | Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.  |
| P333+P313 | If skin irritation or rash occurs: Get medical advice/ attention.                                       |
| P302+P352 | IF ON SKIN: Wash with plenty of soap and water.   |
| P501      | Dispose of contents/ container in accordance with local/ regional/ national/ international regulations. |

#### TMB-Substrate Solution

The TMB-Substrate (S) contains 3,3',5,5' Tetramethylbencidine ( $<0.05\%$ )

- |            |   |
|------------|---|
| H315       | Causes skin irritation.   |
| H319       | Causes serious eye irritation.                                      |
| H335       | May cause respiratory irritation.                                   |
| P261       | Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.              |
| P305+P351+ | IF IN EYES: Rinse cautiously with water for several minutes.        |
| P338       | Remove contact lenses, if present and easy to do. Continue rinsing. |

#### Stopping Solution

The Stopping Solution contains 0.2 M acid sulphur acid ( $H_2SO_4$ )

- |            |   |
|------------|---|
| H290       | May be corrosive to metals.   |
| H314       | Causes severe skin burns and eye damage.  |
| P280       | Wear protective gloves/ protective clothing/ eye protection/ face protection.           |
| P301+P330+ | IF SWALLOWED: rinse mouth.  |
| P331       | Do NOT induce vomiting.   |
| P305+P351+ | IF IN EYES: Rinse cautiously with water for several minutes.                            |
| P338       | Remove contact lenses, if present and easy to do. Continue rinsing.                     |
| P309+P310  | IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician. |

#### First aid procedures:

Skin contact: Wash affected area rinse immediately with plenty of water at least 15 minutes. Remove contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: After swallowing the product, if the affected person is conscious, rinse out the mouth with plenty of water: seek medical advice immediately.

## **METHODOLOGY**

### **Principle**

The DAsource IGFBP-2 Elisa kit is a so-called Sandwich-Assay using two specific and high affinity antibodies. The IGFBP-2 in the samples binds to the first antibody coated on the microtiter plate. In the following step the second specific anti-IGFBP-2-Antibody binds in turn to the immobilised IGFBP-2. The second antibody is biotinylated and will be applied in a mixture with a Streptavidin-Peroxidase-Enzyme Conjugate. In the closing substrate reaction the turn of the colour will be catalysed quantitatively depending on the IGFBP-2-level of the samples. It recognizes IGFBP-2 quantitatively and is not influenced by increased IGF-I or IGF-II values. Related molecules such as IGFBP-3 do not show any cross-reactions in the test.

### **Sample Preparation and Storage**

Serum and plasma EDTA samples are applicable.

The blood sample for serum preparation should be gained according to standardized venipuncture procedure. The samples should be stored without anticoagulation reagents. Hemolytic reactions have to be avoided. The blood has to be allowed to clot and after complete clotting, serum is separated by centrifugation.

Whole blood should ideally be separated after about 2 hours and processed after 12 hours at the latest and stored at -20 ° C until the measurement.

Blood samples may be taken at any time of the day and should be stored in firmly closed sample vials. Haemolytic reactions are to be avoided.

Storage at RT           max. 24 hours

Storage at -20°C       max. 2 years

Are not allowed to have more than 10 freeze/thaw cycles.

Triglycerides bilirubin and hemoglobin in the sample do not interfere up to a concentration of 100 mg/mL, 200 µg/mL and 1 mg/mL respectively. The use of hemolytic, lipemic or icteric samples should nevertheless be validated beforehand by the user.

#### Sample dilution

Dilution: 1:21 with Dilution Buffer

Example: 15 µL sample are added to 300 µL Dilution Buffer (dilution factor 21). After mixing this solution, 2 x 100 µL are used in the assay.

The serum samples must be diluted 1:10 – 1:30 with Dilution Buffer prior to measurement, depending on the expected IGFBP-2 values. In general, a dilution of 1:21 is suitable (and thus the minimum required sample volume for a duplicate determination is 15 µL serum). An extraction step is not required.

### **Technical hints**

The assay has to be conducted strictly according the test protocol herein.

Reagents with different lot numbers cannot be mixed.

Store the kit at 2-8°C after receipt until its expiry date. The lyophilized reagents should be stored at -20 °C after reconstitution. Avoid repeated thawing and freezing. Store the unused seal stripes of the microtiterplate together with the desiccant at 2-8°C.

#### Calibrators and Controls

For the reconstitution of the lyophilised components (calibrators 1-5 and Control Sera), the kit Dilution Buffer has to be used. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam !) with a Vortex mixer.

The reconstituted calibrator and controls can be stored for 4 weeks at -20°C. For further use, thaw quickly but gently (avoid temperature increase above room temperature and avoid excessive vortexing). Up to 3 of the freeze-thaw cycles did not influence the assay.

After reconstitution dilute the Controls with the Dilution Buffer in the same ratio (1:21) as the sample.

Bring all reagents to room temperature (20 – 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming.

The required volume of washing buffer is prepared by 1:20 dilution of the provided 20fold concentrate with deionised water. The diluted washing buffer is stable for 4 weeks at 2-8°C. It has to be at room temperature for usage!

Incubation at room temperature means: Incubation at 20 - 25°C. The TMB-Substrate Solution is photosensitive—store and incubation in the dark.

#### Shaking

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtiter plate shaker. We are recommending approx. 350 rpm. Due to certain technical differences deviations may occur, in case the rotation frequency must be adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/ or false values, excessive shaking may result in high optical densities and/ or false values.

### Washing

Proper washing is of basic importance for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

The danger of handling with potentially infectious material must be taken into account.

Manual washing should be used. Washing Buffer may be dispensed via a multistep device, a multichannel pipette, or a squirt bottle. Decant contents into a biohazard bin, then blot plate on absorbent tissue. Wash the plate by adding 300µL Washing Buffer WB/well, then decant and blot on absorbent tissue. Repeat this step 4 more times for total of 5 washes.

## **ASSAY PROCEDURE**

All determinations (Calibrators, Control Sera and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

When performing the assay, the **Calibrators 1-5, Controls and samples** should be pipetted as fast as possible (e.g., 15 Minutes). To avoid distortions due to differences in incubation times, antibody-Conjugate and the **TMB-Substrate solution** should be added to the plate in the same order and the same time interval as the samples. **Stopping Solution** should be added to the plate in the same order as the **TMB-Substrate Solution**.

- 1) After reconstitution dilute the Controls with the Dilution Buffer in the same ratio (1:21) as the sample.
- 2) Add **100 µl** Dilution Buffer to the first wells (blank). Subsequently, add **100µl** of each Calibrator or diluted **Control (1 & 2)** or diluted **Sample** to the following wells.
- 3) Cover the wells with sealing tape and incubate the plate for **1 hour** shaking at room temperature at **350 rpm**.
- 4) After incubation decant the contents of the wells into a disinfectant (risk of infection!) and wash the wells 5 times with **300 µl** of **Washing Buffer** / well respectively.
- 5) Following the last washing step pipette **100µl** of the **Antibody-Conjugate** in each well. Cover the wells with the sealing tape and incubate the plate for **30 Minutes** shaking at **room temperature at 350 rpm**.
- 6) After incubation wash the wells 5 times with **Washing Buffer** as described above.
- 7) Pipette **100 µl** of the **TMB-substrate** solution in each well.
- 8) Incubate the plate for **15 minutes** in the dark at **room temperature**.
- 9) Stop the reaction by adding **100 µl** of **Stopping Solution** to all the wells.
- 10) Measure the absorbance within **30 minutes** at **450 nm** with  $\geq 590$  nm as reference wavelength.

## **QUALITY CONTROL**

Good laboratory practice requires that controls are included in each assay. The test results are only valid if the test has been performed following the instructions. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state or local standards/laws.

### Quality criteria

For the evaluation of the assay it is required that the absorbance values of the blank should be below 0.25, and the absorbance of Calibrator 5 should be above 1.00. Samples, which yield higher absorbance values than Calibrator 5, should be re-tested with a higher dilution.

## **EVALUATION OF RESULTS**

### **Establishing the Calibration Curve:**

The calibrators provided contain the following concentrations of IGFBP-2:

<b>Calibrators</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>
ng/ml	2	10	20	40	80

- 1) Calculate the mean absorbance (MA) value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance (MA) of the blank from the mean absorbance of all other values.
- 3) Plot the calibrators concentrations 1-5 on the x-axis versus the mean value of the absorbance of the calibrators on the y-axis. Read the concentration of each control and sample by interpolation on the calibration curve.
- 4) Recommendation: Calculation of the calibration curve should be done by using a computer program, because the curve is in general (without respective transformation) not ideally described by linear regression. A four parametric logistic (4-PL) curve fit should be used for recalculation of IGFBP-2 concentrations.
- 5) The IGFBP-2 concentration of the diluted sample or the diluted control sera 1&2 in ng/ml is calculated in this way, the IGFBP-2 concentration of the **undiluted sample** and of 1 & 2 is calculated **by multiplication** with the respective dilution factor.

### **Example of a typical calibration curve**

The exemplary shown calibration curve in Fig.1 **cannot be used** for calculation of your test results. You have to establish a calibration curve for each test you conduct!

Exemplary calculation of the IGFBP-2 concentration of undiluted sample:

Measured extinction of your sample 0.37  
Measured extinction of the blank 0.06


Your **measurement program** will calculate the IGFBP-2 concentration of the diluted sample automatically by using the difference (0.31) of sample and blank for the calculation. You only have to determine the most suitable curve fit (here: polynomial 3<sup>rd</sup> degree).

In this exemplary case the following equation is solved by the program to calculate the IGFBP-2 concentration in the sample:

$$0.31 = -0.0012048 + 0.039581x + 5.1788 \cdot 10^{-0.005} \cdot x^2 - 1.8929x \cdot 10^{-0.006} \cdot x^3$$

$$7.93 = x$$

if the dilution factor (1:21) is taken into account the IGFBP-2 concentration of the undiluted sample is  $7.93 \times 21 = 166.55$  ng/ml



**Fig. 1: Exemplary Calibration Curve** with a polynomial 3<sup>rd</sup> degree as curve fit.

## INTERPRETATION OF RESULTS

The test results should not be the only base for therapeutic decisions. The results should be interpreted in regard to anamnesis, further clinical observations and results of other diagnostic investigations. Further, it is recommended to establish reference and cut-off values corresponding to the relevant group of patients for each laboratory.

## PERFORMANCE CHARACTERISTICS

### Sensitivity

Sensitivity was assed by measuring the blank and calculating the theoretical concentration of the blank + 2SD. The analytical sensitivity of the KAPME05 is 0.2 ng/mL.

### Specificity

The specificity of this ELISA was investigated. Neither 100 ng/mL IGFBP-1 nor 20 µg/mL IGFBP-3 containing samples resulted in a signal significantly different of the blank.

### Interference

Interference of bilirubin, triglycerides and hemoglobin was tested by adding different amounts of these substances to human serum containing IGFBP-2. For comparison the same amount of buffer without any substance was also added to the serum. Table 1 demonstrates that neither bilirubin nor triglycerides exert any influence on the measurement of IGFBP-2 in human serum.

**Table 1:** Interference

Bilirubin		Triglycerides		Hemoglobin	
[µg/ml]	% of control	[mg/ml]	% of control	[mg/mL]	% of control
25	95.07	12.5	100.79	0.125	99
50	92.80	25	101.01	0.25	105
100	93.83	50	103.65	0.5	100
200	88.15	100	101.34	1	100



### Recovery

Recombinant IGFBP-2 was added in three different concentration to human serum. The IGFBP-2 concentration was measured and the mean relative recovery in comparison to buffer was 108%. Some exemplary data are shown in table 2.

**Table 2:** Recovery of recombinant human IGFBP-2 in Serum

IGFBP-2 [ng/ml]	+1000 ng/ml	+500 ng/ml	+100 ng/ml	Mean [ng/ml]
% Recovery	100,00	112,00	114,00	108,67

### Reproducibility and Precision

The inter- and intra assay coefficients of variability are below 10%. Exemplary determinations are shown in table 3 and 4.

**Table 3 :** Inter-Assay-Variation

Sample1 (ng/ml)	137	159	152	8% CV
Sample 2 (ng/ml)	672	697	688	2% CV
Sample 3 (ng/ml)	928	929	956	2% CV

**Table 4:** Intra-Assay-Variation

Sample 1 ng/ml	322	375	298	305	318	311	320	325	302	301	305	317	6.48% CV
Sample 2 ng/ml	612	609	616	648	594	597	620	613	617	611	636	698	4.49% CV

### Linearity

The test can be diluted over a very wide range. The linearity of serum dilutions is shown in table 5.

**Table 5:** Linearity of the sample dilution:

Dilution	Serum 1 (ng/ml)	Serum 2 (ng/ml)
1:10	938	582
1:20	1061	673
1:40	1055	719
1:80	1004	691
1:160	952	668

## EXPECTED VALUES

The IGFBP-2 concentration in serum is depended on age (Table 6) and on Body Mass Index (BMI; Table 7). For data collection of these reference values IGFBP-2 levels were determined in serum of over 400 normal children and adults (see table 7 figure 2); (3).

**Table 6:** BMI dependent reference values, adults between 20 and 80 years


BMI [kg/m <sup>2</sup> ]	N	IGFBP-2 [ng/ml]		percentiles		
		mean	SD	5th	50th	95th
15	12	612	110	431	612	793
17,5	14	568	126	361	568	775
20	76	509	144	271	509	746
22,5	124	449	162	182	449	716
25	101	398	165	127	398	670
27,5	52	348	147	106	348	590
30	25	315	118	120	315	510
32,5	15	282	90	135	282	430
35	4	251	80	119	251	383
37,5	4	220	71	104	220	336

**Table 7 :** IGFBP-2 serum levels (in ng/ml) of > 400 healthy individuals. The normal range is given by the 5., 50. and 95. percentile for age classes.

### Age-dependent normal range of serum IGFBP-2

age (years)	5. percentile (ng/ml)	50. percentile (ng/ml)	95. percentile (ng/ml)
1	408	545	728
2	359	500	696
3	317	460	668
4	277	421	640
5	243	388	617
6	217	361	602
7	194	336	583
8	173	312	562
9	154	289	542
10	138	268	522
11	123	249	503
12	111	232	486
13	101	219	477
14	94	212	470
15	89	207	465
16	86	207	460
17	84	214	466
18	84	223	483
19	84	232	500
25	99	280	580
35	110	381	686
45	130	403	702
55	140	410	715

65	151	418	727
75	153	427	740
80	156	430	744



Assembled by Dr. R. Schweizer, Tübingen, Germany

**Figure 2** : IGFBP-2 serum levels (in ng/mL) of > 400 healthy individuals. The normal range is given by the 5th, 50th and 95th percentile.

#### LIMITATION

Deviation from the reference range can be expected especially in hypothyroidism, after major surgery, in polytrauma, in Diabetes mellitus (due to insulin therapy), in fasting and in malignant diseases.

The DIAsource IGFBP-2 ELISA is based on mono- and polyclonal antibodies. Generally the result of any immunological test system can be influenced by heterophilic antibodies, anti-species antibodies or rheumatic factors. The assay design reduces these potential influences to a minimum but they cannot be excluded in any case.

## Instructions for use for scientific application

#### SCIENTIFIC APPLICATION

IGFBP-2 is present in different concentrations in various body fluids and in conditioned cell culture media of many cell lines.

#### Samples suitable for scientific application

Cerebrospinal fluid, breast milk, amniotic fluid, saliva and in cell culture medium.

**Table 8** : Linearity of the sample dilution.

Dilution	Cerebrospinal fluid (ng/mL)	Amniotic fluid (ng/mL)
1:10	426	Not determined
1:20	428	460
1:40	379	483
1:80	318	431
1:160	426	415

**Table 9 :** Expected values of IGFBP-2 in body fluids of human origin and in cell culture supernatants:

Sample	Expected Value [ng/mL]
Serum	[ 100 - 1000 ]
Cerebrospinal fluid	[ 100 - 300 ]
Amniotic fluid	[ 200 - 10000 ]
Seminal plasma	[ 5000 - 15000 ]
Breast milk	[ 1500 - 3000 ]
Cell culture supernatants	[ 5 - 300 ]

**Table 10:** Results of saliva-sample matrix tests:

Recombinant IGFBP-2: 38,2 [ng/mL] % Recovery		
Saliva + DIL	2.02	n.a.
Saliva + rec. IGFBP-2	35.7	89
Saliva + rec. IGFBP-2	37.36	93
Saliva + rec. IGFBP-2	35.67	89
Saliva + rec. IGFBP-2	36.17	90

**Cross reactions with animal samples**

This assay is specific for human IGFBP-2, low degree of cross reactions was found with commercial dog, horse, donkey, cat and goat sera. No cross-reactivity was found with pig, bovine, rabbit, mouse, chicken, rat, guinea pig, sheep sera.

**REFERENCES**

- Elmlinger MW, Wimmer K, Biemer E, Blum WF, Ranke MB, Dannecker GE, (1996) Insulin-like growth factor binding protein 2 is differentially expressed in leukaemic T- and B-cell lines. *Growth Regulation* 6, 152-157
- Hoeflich A, Nedbal S, Blum WF, Erhard M, Lahm H, Brem G, Kolb HJ, Wanke R, Wolf E. (2001) Growth inhibition in giant growth hormone transgenic mice by overexpression of insulin-like growth factor-binding protein-2. *Endocrinology* 142:1889-98
- Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2000) Significance of basal IGF-I, IGFBP-3 and IGFBP-2 measurements in the diagnostics of short stature in children. *Hormone Research* 54: 60-68
- Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2001). Relevance of IGF-I, IGFBP-3, and IGFBP-2 measurements during GH treatment of GH-deficient and non-GH-deficient children and adolescents. *Hormone Research* 55: 115-124
- Crofton PM, Ahmed SF, Wade JC, Elmlinger MW, Ranke MB, Kelnar CJH, Wallace WHB (1999) Effects of a third intensification block of chemotherapy on bone and collagen turnover, insulin-like growth factor I, its binding proteins and short term growth in children with acute lymphoblastic leukemia. *European Journal of Cancer* 35: 960-967
- Muller HL, Oh Y, Lehnbecher T, Blum WF, Rosenfeld RG. (1994) Insulin-like growth factor-binding protein-2 concentrations in cerebrospinal fluid and serum of children with malignant solid tumors or acute leukemia. *J Clin Endocrinol Metab.* 79: 428-34
- Martin W Elmlinger, Martin H Deininger, Burkhardt S Schuett, Richard Meyermann, Frank Duffner, Ernst H Grote, and Michael B Ranke (2001) In vivo expression of the insulin-like growth factor binding protein-2 in human gliomas increases with the tumor grade. *Endocrinology* 142: 1652-1658
- Cohen P, Peehl DM, Stamey TA, Wilson KF, Clemmons DR, Rosenfeld RG (1993) Elevated levels of insulin-like growth factor-binding protein-2 in the serum of prostate cancer patients. *J Clin Endocrinol Metab.* 76: 1031-1035
- Boulle N, Baudin E, Gicquel C, Logie A, Bertherat J, Penfornis A, Bertagna X, Luton JP, Schlumberger M, Le Bouc Y. (2001) Evaluation of plasma insulin-like growth factor binding protein-2 as a marker for adrenocortical tumors. *Eur J Endocrinol.* 144: 29-36
- Ranke MB, Maier KP, Schweizer R, Stadler B, Schleicher S, Elmlinger MW, Flehmig B (2003), Pilot study of elevated levels of insulin-like growth factor-binding protein-2 as indicators of hepatocellular carcinoma. *Horm Res* 60:174-180
- Rosenfeld RG, Roberts CT Jr. (eds.) (1999) *The IGF system: Contemporary Endocrinology Series*; Humana Press
- Jones JI, Clemmons DR. (1995) Insulin-like growth factors and their binding proteins: biological actions. *Endocr Rev.* 16: 3-34
- Chard T (1994) Insulin-like growth factors and their binding proteins in normal and abnormal fetal growth. *Growth Reg* 4: 91-100
- Ranke MB, Elmlinger MW (1997) Functional role of insulin-like growth factor binding proteins. *Hormone Research* 48 (suppl 4): 9-15
- Elmlinger MW, Bell M, Schütt, BS, Langkamp M, Kutoh E, Ranke MB (2001) Transactivation of the IGFBP-2 promoter in human tumor cell lines. *Molecular and Cellular Endocrinology* 175: 211-218
- Hoeflich A, Reisinger R, Lahm H, Kiess W, Blum WF, Kolb HJ, Weber MM, Wolf E. (2001) Insulin-like growth factor-binding protein 2 in tumorigenesis: protector or promoter? *Cancer Res.* 61: 8601-8610
- Miell JP, Langford KS, Jones JS, Noble P, Westwood M, White A, Nicolaides KH. (1997) The maternal insulin-like growth factor (IGF) and IGF-binding protein response to trisomic pregnancy during the first trimester: a possible diagnostic tool for trisomy 18 pregnancies. *J Clin Endocrinol Metab.* 82: 287-292
- Elmlinger MW., Regine Grund, Michael Buck, Hartmut A Wollmann, Norman Feist, Matthias M. Weber, Christian C. Speer and Michael B. Ranke (1999). Limited proteolysis of the insulin-like growth factor binding protein (IGFBP)-2 by a specific serine protease activity in early breast milk. *Pediatric Research* 46: 76-81

19. Elmlinger MW, Michael S Sanatani, Mathias Bell, Günther E Dannecker, Michael B Ranke (1998) Elevated insulin-like growth factor (IGF) binding protein (IGFBP)-2 and IGFBP-4 expression of leukemic T-cells is affected by autocrine/paracrine IGF-II action but not by IGF type I receptor expression. European Journal of Endocrinology 138: 337-343
20. van den Beld, Blum W, Pols H, Grobbee D, Lamberts S (2003) Serum insulin-like growth factor binding protein-2 levels as an indicator of functional ability in elderly men. European Journal of Endocrinology 148: 627-634

#### SUMMARY OF THE ASSAY

Reagent:	Reconstitution:	dilution:
Calibrators 1-5	in 750 µl Dilution Buffer	
Controls Serum	in 100 µl Dilution Buffer	1:21 with Dilution Buffer
Washing Buffer		1:20 with Distilled water
Dilute samples 1:21 in Dilution Buffer, mix immediately, Use 100 µL for each well in the assay.		
Before assay procedure bring all reagents to room temperature +20-25°C.		

#### Assay Procedure for Double Determination

Pipette	Reagents	Position
100 µl	Dilution Buffer	A1/2
100 µl	Calibrator 1 ( 2 ng/ml)	B1/2
100 µl	Calibrator 2 (10 ng/ml)	C1/2
100 µl	Calibrator 3 (20 ng/ml)	D1/2
100 µl	Calibrator 4 (40 ng/ml)	E1/2
100 µl	Calibrator 5 (80 ng/ml)	F1/2
100 µl	Control Serum 1 (1:21 diluted)	G1/2
100 µl	Control Serum 2 (1:21 diluted)	H1/2
100 µl	Sample dilution (1:21 diluted)	following wells
Cover the wells with the sealing tape.		
<b>Incubation: 1 h at RT, 350 rpm</b>		
5x 300 µl	Decant the contents of the wells and wash 5x with 300 µl Wash Buffer	each well
100 µl	Antibody-Conjugate	each well
Cover the wells with the sealing tape.		
<b>Incubation: 30 min at RT, 350 rpm</b>		
5 x 300 µl	Decant the contents of the wells and wash 5 x with 300 µl Wash Buffer	each well
100 µl	TMB-Substrate Solution	each well
<b>Incubation: 15 min in the Dark at RT</b>		
100 µl	Stopping Solution	each well
Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength.		

Revision Date: 27-05-2021