IVD

CE

CA 72-4 IRMA

ELSA-CA72-4

Version : 230123

Date of issue : 23/01/2023

Revision date: 23/01/2023

History

Summary of change:

Current Version:
230123
New logo

CE

en

Read entire protocol before use.

CA 72-4 IRMA

I. INTENDED USE

Immunoradiometric assay kit for the in vitro quantitative measurement of TAG 72 antigen in human serum.

II. GENERAL INFORMATION

- A. Proprietary name : DIAsource CA 72-4 IRMA Kit
- B. Catalog number : ELSA-CA72-4 : 96 tests
- C. Manufactured by : DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

TAG 72 is a high molecular weight glycoprotein associated with tumors, and recognized by B 72-3* and CC 49 monoclonal antibodies*.

The B 72-3 antibody was obtained by immunizing a mouse with a metastatic human breast cancer-enriched membrane fraction. In immunohistochemistry, it shows good affinity for gastro-intestinal and human mammary carcinomas, as opposed to the corresponding benign or normal tissues. The CC 49 antibody was generated after immunizing a mouse with TAG 72 previously purified by affinity chromatography. It recognizes a different epitope from that recognized by B 72-3.

This assay system was used in numerous clinical studies, which demonstrated its good level of sensitivity to gastric cancers(70% for the metastatic stages and 20% for the non-metastatic) and its outstanding specificity. The assay also enables the follow-up of gastric cancer cases under treatment or in remission.

IV. PRINCIPLES OF THE METHOD

The DIAsource CA 72-4 IRMA is a solid-phase sandwich immunoradiometric assay. Two monoclonal antibodies were prepared against sterically remote sites on the TAG 72 molecule, the first being coated on the lower and inner surface of the plastic tubes, while the second radiolabelled with iodine 125, is used as a tracer.

TAG 72 molecules present in the calibrators or the samples to be tested are sandwiched between the two antibodies. Excess unbound tracer is easily removed during the procedure's washing step, and the coated tubes retain only the adsorbed antibody/antigen/tracer antibody combination.

The amount of radioactivity bound to the coated tube is proportional to the concentration of TAG 72 present in the sample.

V. REAGENTS PROVIDED

Reagents	96 tests Kit	Reconstitution
Tubes coated with anti CA72-4 (monoclonal antibodies)	2 x 48	Ready for use
Anti-CA 72-4 ^{125}I (monoclonal antibodies) in phosphate buffer with bovine serum albumin, Proclin, mouse immunoglobulins and inert red dye .	1 vial 31 ml 388 kBq	Ready for use
Calibrators 1:5:human TAG 72 antigen in human serum and azide. (see exact value on vial labels)	5 vials 0.8 ml liquid	Ready for use
Control: human TAG 72 antigen in human serum and azide. (see exact value on vial label)	1 vial 0.8 ml liquid	Ready for use
Buffer with bovine serum albumin and azide .	1 vial 20 ml liquid	Ready for use
Buffer with bovine serum albumin and azide .	1 vial 10 ml liquid	Ready for use
Wash solution (Tween 20 - NaCl)	1 vial 50 ml	Dilute 28x with distilled water (use a magnetic stirrer).

Note:

Use the diluent as calibrator 0.
Use the diluent for samples dilutions.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Precision micropipettes or similar , for delivery of: 100 μl , 200 μl and 300 μl (the use of accurate pipettes with disposable plastic tips is recommended).
3. Disposable plastic tubes
4. Vortex mixer
5. Magnetic stirrer
6. Tubes shaker : circular horizontal (400 rpm)
7. 5 ml automatic syringe (Cornwall type) for washing
8. Aspiration system (optional)
9. Any gamma scintillation counter capable of measuring 125 Iodine may be used (minimal yield 70%).

VII. REAGENT PREPARATION

Working Wash solution : Prepare an adequate volume of Working Wash solution by adding 27 volumes of distilled water to 1 volume of Wash Solution (28x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening, all kit components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

The assay is performed directly on serum. If the test is to be carried out within 24 hours, the samples may be stored at 2-8°C. Otherwise, they should be divided into aliquots and deep frozen (-20°C) until needed. Avoid subsequent freeze-thaw cycles.

Dilutions

Should elevated CA 72-4 levels be suspected, dilution is performed with the diluent found in the kit.

It is recommended that disposable plastic tubes be used when carrying out dilutions.

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature (18-25°C) at least 30 minutes prior to use.

In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample. Thoroughly mix all reagents and samples by gentle agitation or swirling.

High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times.

Prepare a calibration curve for each run, do not use data from previous runs.

B. Procedure

1. Label coated tubes in duplicate for each calibrator, sample and control. For determination of total counts, label 2 uncoated disposable plastic tubes.
2. Briefly vortex calibrators, samples, control.
3. Respect the order in which reagents are to be added:
4. Dispense 200 μl of Buffer into each coated tube .
4. Add 100 μL of calibrators, control or samples to the corresponding groups of tubes.
5. Gently mix each tube with a Vortex-type mixer.
6. Incubate 4 hours \pm 5 min at room temperature (18-25°C) under continuous shaking (400 rpm).
7. Wash the coated tubes as follows:
 - Aspirate the contents of the tubes as completely as possible.
 - Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
 - Add 2.0 mL of Working Wash solution to each tube.
 - Avoid foaming during the addition of the Working Wash solution.
 Repeat this process twice.
8. Add 300 μL of anti-CA 72-4 ^{125}I monoclonal antibody to all the tubes.
9. Gently mix each tube with a Vortex-type mixer.
10. Incubate overnight (16-20 h) at 2-6°C.
11. The incubation temperature should not exceed 6°C.
12. Wash the coated tubes (except total counts) as previously described.
13. After the last washing, let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
13. Measure the remaining radioactivity bound to the coated tubes with a gamma scintillation counter for 60 seconds.

XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. Plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of CA72-4 (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
3. Read the concentration for each control and sample by interpolation on the calibration curve.

4. Computer assisted data reduction will simplify these calculations. If automatic result processing is to be used, a 4-parameter logistic function curve fitting is recommended.

XII. TYPICAL CALIBRATION CURVE

The following data are for illustration only and should never be used instead of the real time calibration curve.

CA72-4 IRMA		cpm	B/T (%)
Total count		170106	100
Calibrator			
	0.0 U/ml	763	0.45
	4.3 U/ml	1449	0.85
	11.1 U/ml	2913	1.71
	24.1 U/ml	4935	2.90
	47.9 U/ml	8729	5.13
	102.7 U/ml	16803	9.88

Since no international reference material is available for CA72-4 antigen, DIAsource CA 72-4 IRMA calibrator values are assigned against a set of in-house reference standards.

XIII. PERFORMANCE AND LIMITATIONS

A. Detection limit

The LoB (Limit of blank) was calculated by measuring the blank several times and was calculated as the mean + 1.645 standard deviation of the distribution of the test values . The LoB was calculated to be 0.88 U/ml.

The LoD (Limit of Detection) was calculated as the LoB + 1.645 standard deviation of a low concentration sample tested in 10 different runs .The LoD was calculated to be 1.9 U/ml.

The LoQ (Limit of Quantification) was calculated by testing 5 low values samples, 10 times. The LoQ was calculated to be 3.4 U/ml.

B. Interferences

No interference has been observed :

with haemoglobin up to 2.5 g/l and 5 g/l ,

with triglycerides up to 0.5 and 2.5 g/l ,

with bilirubin non conjugated up to 20 and 150 mg/l

with conjugated bilirubin up to 4 and 20 mg/l

The immunoassay is protected against any human anti-mouse antibody (HAMA) interference by the addition of a protector to the tracer (non-specific mouse immunoglobulins). However, we cannot guarantee that this protection is exhaustive.

C. Precision

INTRA ASSAY				INTER ASSAY			
Serum	Replicate	$\text{} \pm \text{SD}$ (U/ml)	CV (%)	Serum	Replicate	$\text{} \pm \text{SD}$ (U/ml)	CV (%)
A	30	6.2 ± 0.6	8.9	A	10	5.4 ± 0.8	14.3
B	30	17.2 ± 1.2	7.2	B	10	16.2 ± 1.7	10.5

SD : Standard Deviation; CV: Coefficient of variation

D. Accuracy

RECOVERY TEST:

Known quantities of CA 72-4 were added to low human serum samples. The recovery percentages of CA 72-4 in serum samples ranged from 80.3 to 117.1%.

Added CA 72-4 (U/ml)	Recovered CA 72-4 (U/ml)	Recovery (%)
Sample 1		
2.05	2.13	103.9
4.30	4.80	111.6
8.70	7.73	88.9
16.59	14.80	89.2
25.50	24.64	96.6
48.20	45.65	94.7
91.85	78.54	85.5
99.06	109.79	110.8
Sample 2		
2.05	2.40	117.1
4.30	3.98	92.6
8.70	6.99	80.3
16.59	15.47	93.3
25.50	23.19	90.9
48.20	45.67	94.8
91.85	82.14	89.4
99.06	105.87	106.9

DILUTION TEST:

An elevated sample was diluted with the Diluent, with the recovery percentages ranging from 81.7 to 112.5 %.

Dilution factor	Theoretical concentration (U/ml)	Measured concentration (U/ml)	recovery (%)
1	152.1		
0.6	91.2	92.9	101.9
0.5	76.0	68.0	89.4
0.4	60.8	57.5	94.6
0.3	45.6	40.4	88.6
0.2	30.4	24.9	81.7
0.1	15.2	12.5	81.9
0.05	7.6	7.3	96.1
0.01	1.5	1.7	108.6
0.005	0.8	0.9	112.5

Samples were diluted with the Diluent.

XIV. LIMITATIONS

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays. Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed in case of the presence of heterophilic antibodies. Carefully evaluate the results of patients suspected of having these antibodies.
- If results are not consistent with other clinical observations, additional information should be required before diagnosis.

XV. INTERNAL QUALITY CONTROL

- If the results obtained for Control are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Do not freeze-thaw more than twice.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises

XVI. REFERENCE INTERVALS

Among 120 apparently healthy individuals, 94% of the results were below 4 U/ml, with no difference observed between sex, age and smoking habits.

These results are given as an indication and it is recommended that each laboratory establish its own range of clinical values.

For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history and other findings.

XVII. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

This kit contains ^{125}I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVIII. BIBLIOGRAPHY

Brandt B, Liffers E, Sasse W, Assmann G. CA 72-4 : the tumor marker for gastric carcinoma. Dans : Klapdor R, ed. Tumor associated antigens, oncogenes, receptors, cytokines in tumor diagnosis and therapy at the beginning of the nineties. Cancer of the breast. State and trends in diagnosis and therapy. München, Bern, Wien, New-York : W. Zuckschwerdt Verlag. 1992:10-2.

Byrne DJ, Browning MCK, Cuschieri A. CA 72-4 : a new tumor marker for gastric cancer. Br J Surg. 1990;77:1010-3.

Colcher D, Paterson A, Sears H, Schlom J. Use of monoclonal antibody b72.3 to detect a circulating tumor-associated glycoprotein (TAG 72) in the serum of colon carcinoma patients. Hybridoma. 1985;4(1):71.

Gero EG, Colcher D, Ferroni P et al. CA 72-4 radioimmunoassay for the detection of the TAG 72 carcinoma-associated antigen in serum of patients. J Clin Lab Anal. 1989;3:360-9.

Guadagni F, Roselli M, Amato T et al. CA 72-4 measurement of tumor-associated glycoprotein 72 (TAG 72) as a serum marker in the management of gastric carcinoma. Cancer Res. 1992;52:1222-7.

Johnson VG, Schlom J, Paterson AJ, Bennet J, Magnani JL, Colcher D. Analysis of human tumor-associated glycoprotein (TAG 72) identified by monoclonal antibody B72.3. Cancer Res. 1986;46:850-7.

Ohuchi N, Gero E, Mori S et al. Clinical evaluation of CA 72-4 immunoradiometric assay for serum TAG 72 antigen in patients with carcinoma. J of Tumor Marker Oncol. 1990;5(1):1-10.

XIX. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS ml	CALIBRATORS ml	SAMPLE(S) CONTROL ml
Buffer Calibrators (0-5) Control or Samples	-	0.2 0.1 -	0.2 - 0.1
Incubation	4 hours at 18-25°C , with continuous shaking(400 rpm)		
Separation Working Wash solution Separation Working Wash solution Separation	- - - -	Aspirate (or decant) 2.0 Aspirate (or decant) 2.0 Aspirate (or decant)	
Tracer	0.3	0.3	0.3
Incubation	20 ± 2 hours at 2-6°C		
Separation Working Wash solution Separation Working Wash solution Separation	- - - -	Aspirate 2.0 Aspirate 2.0 Aspirate	
Counting	Count tubes for 60 seconds		

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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Colcher D, Paterson A, Sears H, Schlom J. Use of monoclonal antibody b72.3 to detect a circulating tumor-associated glycoprotein (TAG 72) in the serum of colon carcinoma patients. Hybridoma.1985;4(1):71.

Gero EG, Colcher D, Ferroni P et al. CA 72-4 radioimmunoassay for the detection of the TAG 72 carcinoma-associated antigen in serum of patients. J Clin Lab Anal.1989;3:360-9.

Guadagni F, Roselli M, Amato T et al. CA 72-4 measurement of tumor-associated glycoprotein 72 (TAG 72) as a serum marker in the management of gastric carcinoma. Cancer Res. 1992 52 1222 7

Johnson VG, Schlom J, Paterson AJ, Bennet J, Magnani JL, Colcher D. Analysis of human tumor-associated glycoprotein (TAG 72) identified by monoclonal antibody B72.3. Cancer Res. 1986 46 850 7

Ohuchi N, Gero E, Mori S et al. Clinical evaluation of CA 72-4 immunoradiometric assay for serum TAG 72 antigen in patients with carcinoma. J of Tumor Marker Oncol. 1990;5(1):1-10.

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