# **CYFRA 21-1 IRMA**

***ELSA-CYFRA***

---

Version : 230123

Date of issue : 23/01/2023

Revision date: 23/01/2023

# History

---

## Summary of change:

<b>Current Version:</b> 230123
New logo

Read entire protocol before use.

## CYFRA 21-1 IRMA

### *I. INTENDED USE*

Immunoradiometric assay kit for the in vitro quantitative measurement of CYFRA 21-1 in human serum.

### *II. GENERAL INFORMATION*

- A. **Proprietary name :** DIAsource CYFRA 21-1 IRMA Kit
- B. **Catalog number :** ELSA-CYFRA : 48 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.  
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

**For technical assistance or ordering information contact :**  
**Tel : +32 (0)10 84.99.11                      Fax : +32 (0)10 84.99.91**

### *III. CLINICAL BACKGROUND*

CYFRA 21-1 is a cytokeratin 19 fragment found in serum of cancer patients. Precise recognition of this fragment is made with two monoclonal antibodies (BM 19-21 and KS 19-1)\* which were obtained after immunisation of mice with MCF-7 cells.

Cytokeratin 19 (CK19) is a member of the intermediate filament group of proteins, whose physiological role remains unclear. It is an acid-type cytoplasmic protein, with a molecular weight of 40 000 D, expressed in simple epithelium. On the death of the cell, it is released into the serum in the form of soluble fragments.

In immunohistochemistry, CK19 is found in the cytoplasm of the epithelial tumor cell, including that of bronchial cancers. Preliminary clinical studies of bronchial cancer patients sera have shown that CYFRA 21-1 assay is useful in the diagnosis and follow-up of non-small cell lung carcinoma and particularly of squamous cell carcinoma of the lung.




#### IV. PRINCIPLES OF THE METHOD

The DIAsource CYFRA 21-1 IRMA is a solid-phase sandwich immunoradiometric assay. Two monoclonal antibodies were prepared against sterically remote sites on the CYFRA 21-1 molecule, the first being coated on the lower and inner surface of the plastic tubes, while the second radiolabelled with iodine 125, is used as a tracer.

CYFRA 21-1 molecules present in the calibrators or the samples to be tested are sandwiched between the two antibodies. Excess unbound tracer is easily removed during the procedure's washing step, and the coated tubes retain only the adsorbed antibody/antigen/tracer antibody combination.

The amount of radioactivity bound to the coated tube is proportional to the amount of CYFRA 21-1 present at the beginning of the assay.

#### V. REAGENTS PROVIDED

Reagents	48 tests Kit	Reconstitution
 Tubes coated with anti CYFRA21-1 (monoclonal antibodies)	1 x 48	Ready for use
 Ab <sup>125</sup> I Anti-CYFRA21-1 <sup>125</sup> I (monoclonal antibodies) in phosphate buffer with bovine serum albumin, azide, mouse immunoglobulins and inert red dye .	1 vial 16 ml 311 kBq	Ready for use
 CAL 0 Calibrator 0 in bovine serum	1 vial lyophil.	Add 4 ml distilled water
 CAL N Calibrators 1-4 :human CYFRA 21-1 in bovine serum (see exact value on vial labels)	4 vials lyophil.	Add 1 ml distilled water
 WASH SOLN CONC Wash solution (TRIS-HCl)	1 vial 10 ml	Dilute 70x with distilled water (use a magnetic stirrer).
 CONTROL N Control :human CYFRA 21-1 in bovine serum( see exact value on vial label)	1 vial lyophil.	Add 1 ml distilled water

**Note:** Use the zero calibrator for sera dilutions.

#### VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Precision micropipettes or similar , for delivery of: 100 µl, 300 µl ,1 ml and 4 ml (the use of accurate pipettes with disposable plastic tips is recommended).
3. Disposable plastic tubes
4. Vortex mixer
5. Magnetic stirrer
6. 5 ml automatic syringe (Cornwall type) for washing
7. Aspiration system (optional)
8. Any gamma scintillation counter capable of measuring 125 Iodine may be used (minimal yield 70%).

#### VII. REAGENT PREPARATION

- A. Calibrators :** Reconstitute the calibrator 0 with 4 ml distilled water and the calibrators 1-4 with 1 ml.
- B. Control :** Reconstitute the control with 1 ml distilled water.
- C. Working Wash solution :** Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash

Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

#### VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.
- After reconstitution , the calibrator 0, the calibrators 1 - 4 and the control, should be kept at 2-8°C for maximum 5 days or at -20°C for maximum 2 months.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

#### IX. SPECIMEN COLLECTION AND PREPARATION

- The assay is performed directly on serum.
- If the test is to be carried out within 24 hours, the samples should be refrigerated at 2-8°C. Otherwise, they should be divided into aliquots and deep frozen (-20°C) until needed.
- Avoid subsequent freeze-thaw cycles.
- Do not use plasma samples.

#### X. PROCEDURE

##### A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature (18-25°C) at least 30 minutes prior to use.

In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample. Thoroughly mix all reagents and samples by gentle agitation or swirling.

High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times.

Prepare a calibration curve for each run, do not use data from previous runs.

##### B. Procedure

1. Label coated tubes in duplicate for each calibrator, sample and control. For determination of total counts, label 2 uncoated disposable plastic tubes.
2. Briefly vortex calibrators, samples, control and dispense 100 µl of each into the respective tubes.
3. Dispense 300 µl of anti-CYFRA 21-1 <sup>125</sup>I tracer into each tube, including the uncoated tubes for total counts.
4. Mix each tube gently with a Vortex mixer.
5. Incubate for 20 h ± 2 h at 2-8°C.
6. Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
7. Wash tubes with 3 ml Working Wash solution (except total counts). Avoid foaming during the addition of the Working Wash solution.
8. Aspirate (or decant) the content of each tube (except total counts).
9. Repeat step 7 and step 8 , twice.
10. After the last washing, let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
11. Count tubes in a gamma counter for 60 seconds.

#### XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. Plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of CYFRA 21-1 (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
3. Read the concentration for each control and sample by interpolation on the calibration curve.
4. Computer assisted data reduction will simplify these calculations. If automatic result processing is to be used, a 4-parameter logistic function curve fitting is recommended.

### XII. TYPICAL CALIBRATION CURVE

The following data are for illustration only and should never be used instead of the real time calibration curve.

CYFRA 21-1 IRMA		cpm	B/T (%)
Total count		241728	100
Calibrator	0.0 ng/ml	54	0.02
	1.4 ng/ml	1047	0.41
	8.4 ng/ml	5497	2.25
	24.6 ng/ml	19168	7.91
	48.7 ng/ml	38754	16.01

Since no international reference material is available for CYFRA antigen, DIASource CYFRA 21-1 IRMA calibrator values are assigned against a set of in-house reference standards.

### XIII. PERFORMANCE AND LIMITATIONS

#### A. Detection limit

The LoB (Limit of blank) was calculated by measuring the blank several times and was calculated as the mean + 1.645 standard deviation of the distribution of the test values. The LoB was calculated to be 0.05 ng/ml.

The LoD (Limit of Detection) was calculated as the LoB + 1.645 standard deviation of a low concentration sample tested in 10 different run. The LoD was calculated to be 0.25 ng/ml.

The LoQ (Limit of Quantification) was calculated by testing 5 low values samples, 10 times. The LoQ was calculated to be 0.85 ng/ml.

#### B. Interferences

No interference has been observed :

- with haemoglobin up to 2.5 g/l and 5 g/l ,
- with triglycerides up to 0.5 and 2.5 g/l ,
- with bilirubin non conjugated up to 20 and 150 mg/l
- with conjugated bilirubin up to 4 and 20 mg/l

The immunoassay is protected against any human anti-mouse antibody (HAMA) interference by the addition of a protector to the tracer (non-specific mouse immunoglobulins). However, we cannot guarantee that this protection is exhaustive.

#### C. Precision

INTRA ASSAY				INTER ASSAY			
Serum	Replicate	<X> ± SD (ng/ml)	CV (%)	Serum	Replicate	<X> ± SD (ng/ml)	CV (%)
A	30	3.1 ± 0.2	6.7	A	10	2.6 ± 0.1	5.6
B	30	27.2 ± 1.6	6.0	B	10	22.8 ± 2.0	9.0

SD : Standard Deviation; CV: Coefficient of variation

#### D. Accuracy

##### RECOVERY TEST:

Known quantities of CYFRA 21-1 were added to a low human serum. The recovery percentages of CYFRA 21-1 in the samples ranged from 85.2 to 102.6%.

Added CYFRA 21-1 (ng/ml)	Recovered CYFRA 21-1 (ng/ml)	Recovery (%)
1.15	0.98	85.2
3.09	3.17	102.6
7.26	6.39	88.0
15.07	13.04	86.5
21.43	19.71	92.0
30.08	29.05	96.6
36.3	36.16	99.6
45.39	40.62	89.5

### DILUTION TEST:

An elevated sample was diluted with zero calibrator, with the recovery percentages ranging from 90.8 to 115.8 %.

Dilution factor	Theoretical concentration (ng/ml)	Measured concentration (ng/ml)	recovery (%)
1	176.9		100
0.4	70.8	70.7	112.9
0.3	53.1	59.9	115.3
0.2	35.4	40.8	115.8
0.1	17.7	20.5	101.8
0.05	8.8	9.0	94.0
0.01	1.8	1.7	90.8
0.005	0.9	0.8	

### XIV. LIMITATIONS

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays. Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed in case of the presence of heterophilic antibodies. Carefully evaluate the results of patients suspected of having these antibodies. If results are not consistent with other clinical observations, additional information should be required before diagnosis.

### XV. INTERNAL QUALITY CONTROL

- If the results obtained for Control are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Do not freeze-thaw more than twice.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises

### XVI. REFERENCE INTERVALS

Among 196 apparently healthy individuals, 97.4% of the results were below 3.3 ng/ml, with no difference observed between sex, age and smoking habits.

96.1 % of patients with benign pulmonary disease (n=51) had a CYFRA 21-1 level below 3.3 ng/ml.

These results are given as an indication and it is recommended that each laboratory establish its own range of clinical values.

**For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history and other findings.**

### XVII. PRECAUTIONS AND WARNINGS

#### Safety

For *in vitro* diagnostic use only.

This kit contains <sup>125</sup>I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding

jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

### XVIII. BIBLIOGRAPHY

Broers JLV, Ramaekers FCS, Klein Rot M, et al. Cytokeratins in different types of human lung cancer as monitored by chain-specific monoclonal antibodies. *Cancer Res.* 1988;48:3221-9.

Fuchs E, Tyner AL, Giudice GJ, Marchuk D, Chaudhury AR, Rosenberg M. The human keratin genes and their differential expression. Dans : Sawyer RH, ed. *The molecular and developmental biology of keratins.* In : Moscona AA, Monroy A, eds. *Current topics in developmental biology.* San Diego: Academic Press;1987:5-34 (Vol 22).

Kasper M, Stosiek P, Typlt H, Karsten U. Histological evaluation of three monoclonal antibodies. 1. Normal tissues. *Eur J Cancer Clin Oncol.* 1987;23(2):137-47.

Karsten U, Papsdorf G, Roloff G et al. Monoclonal anti-cytokeratin antibody from a hybridoma clone generated by electrofusion. *Eur J Cancer Clin Oncol.* 1985;21(6):733-40.

Mc Cormack P. Current surgical approach to non-small-cell-lung cancer. *Oncol.* 1991;5(11):39-46.

Moll R, Franke WW, Shiller DL, Geiger B, Krepler R. The catalog of human cytokeratins: patterns of expression in normal epithelia, tumors and cultured cells. *Cell.* 1982;31:11-24.

Pujol JL, Grenier J, Daures JP, Daver A, Pujol H, Michel FB. Serum fragments of cytokeratin subunit 19 measured by CYFRA 21-1 immunoradiometric assay as a marker of lung cancer. *Cancer Res.* 1993;53:1-6.

Schegel R, Banks-Shegel S, Pinkus GS. Immunohistochemical localization of keratin in normal human tissues. *Lab Invest.* 1980;42(1):91-6.

Stasiak PC, Purkis PE, Leigh IM, Lane EB. Keratin 19 : predicted amino acid sequence and broad tissue distribution suggest it evolved from keratinocyte keratins. *J Invest Dermatol.* 1989;92(5):707-16.

Sun TT, Eichner R, Nelson WG, Tseng SCG, Weiss RA, Jarvinen M, Woodcock-Mitchell J. Keratin classes : molecular markers for different types of epithelial differentiation. *J Invest Dermatol.* 1983;81:109s-15s.

### XIX. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS ml	CALIBRATORS ml	SAMPLE(S) CONTROL ml
Calibrators (0-4) Samples, Control Tracer	- - 0.3	0.1 - 0.3	- 0.1 0.3
Incubation	20 ± 2 hours at 2-8°C		
Separation Working Wash solution(ml)	- - -	Aspirate (or decant)  3.0	
Separation Working Wash solution(ml)	- - -	Aspirate (or decant)  3.0	
Separation Working Wash solution	- - -	Aspirate (or decant)  3.0	
Separation	- -	Aspirate (or decant)	
Counting	Count tubes for 60 seconds		

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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## **.XVIII. BIBLIOGRAPHIE**

Broers JLV, Ramaekers FCS, Klein Rot M, et al. Cytokeratins in different types of human lung cancer as monitored by chain-specific monoclonal antibodies. *Cancer Res.* 1988 48 3221 9

Fuchs E, Tyner AL, Giudice GJ, Marchuk D, Chaudhury AR, Rosenberg M. The human keratin genes and their differential expression. Dans : Sawyer RH, ed. *The molecular and developmental biology of keratins.* In : Moscona AA, Monroy A, eds. *Current topics in developmental biology.* San Diego : Academic Press;1987:5-34 (Vol 22).

Kasper M, Stosiek P, Typl H, Karsten U. Histological evaluation of three monoclonal antibodies. 1. Normal tissues. *Eur J Cancer Clin Oncol.* 1987 (23)2/137.47

Karsten U, Papsdorf G, Roloff G et al. Monoclonal anti-cytokeratin antibody from a hybridoma clone generated by electrofusion. *Eur J Cancer Clin Oncol.* 1985 (21)6/733.40

Mc Cormack P. Current surgical approach to non-small-cell-lung cancer. *Oncol.* 1991 (5)11/39.46

Moll R, Franke WW, Shiller DL, Geiger B, Krepler R. The catalog of human cytokeratins: patterns of expression in normal epithelia, tumors and cultured cells. *Cell.* 1982 31 11 24

Pujol JL, Grenier J, Daures JP, Daver A, Pujol H, Michel FB. Serum fragments of cytokeratin subunit 19 measured by CYFRA 21-1 immunoradiometric assay as a marker of lung cancer. *Cancer Res.* 1993 53 1 6

Schegel R, Banks-Schegel S, Pinkus GS. Immunohistochemical localization of keratin in normal human tissues. *Lab Invest.* 1980 (42)1/91.6

Stasiak PC, Purkis PE, Leigh IM, Lane EB. Keratin 19 : predicted amino acid sequence and broad tissue distribution suggest it evolved from keratinocyte keratins. *J Invest Dermatol.* 1989 (92)5/707.16

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Fuchs E, Tyner AL, Giudice GJ, Marchuk D, Chaudhury AR, Rosenberg M. The human keratin genes and their differential expression. In: Sawyer RH, ed. *The molecular and developmental biology of keratins.*

In : Moscona AA, Monroy A, eds. *Current topics in developmental biology.* San Diego: Academic Press; 1987:5-34 (Vol 22).

Kasper M, Stosiek P, Typl H, Karsten U. Histological evaluation of three monoclonal antibodies. 1. Normal tissues. *Eur J Cancer Clin Oncol.* 1987;23(2):137-47.

Karsten U, Papsdorf G, Roloff G et al. Monoclonal anti-cytokeratin antibody from a hybridoma clone generated by electrofusion. *Eur J Cancer Clin Oncol.* 1985;21(6):733-40.

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Schegel R, Banks-Shegel S, Pinkus GS. Immunohistochemical localization of keratin in normal human tissues.

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Fuchs E, Tyner AL, Giudice GJ, Marchuk D, Chaudhury AR, Rosenberg M. The human keratin genes and their differential expression. Dans: Sawyer RH, ed. *The molecular and developmental biology of keratins.* In: Moscona AA, Monroy A, eds. *Current topics in developmental biology.* San Diego: Academic Press;1987:5-34 (Vol 22).

Kasper M, Stosiek P, Typl H, Karsten U. Histological evaluation of three monoclonal antibodies. I. Normal tissues. *Eur J Cancer Clin Oncol.* 1987;23(2):137-47.

Karsten U, Papsdorf G, Roloff G et al. Monoclonal anti-cytokeratin antibody from a hybridoma clone generated by electrofusion. *Eur J Cancer Clin Oncol.* 1985;21(6):733-40.

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Schegel R, Banks-Schegel S, Pinkus GS. Immunohistochemical localization of keratin in normal human tissues. *Lab Invest.* 1980;42(1):91-6.

Stasiak PC, Purkis PE, Leigh IM, Lane EB. Keratin 19: predicted amino acid sequence and broad tissue distribution suggest it evolved from keratinocyte keratins. *J Invest Dermatol.* 1989;92(5):707-16.

Sun TT, Eichner R, Nelson WG, Tseng SCG, Weiss RA, Jarvinen M, Woodcock-Mitchell J. Keratin classes: molecular markers for different types of epithelial differentiation. *J Invest Dermatol.* 1983;81:109s-15s.

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