IVD

CE

IGFBP-3-ELISA

KAP1171

Version : 230123

Date of issue : 23/01/2023

Revision date: 23/01/2023

History

Summary of change:

Current Version:
230123
New logo

CE

en

Read entire protocol before use.

IGFBP-3-ELISA

I. INTENDED USE

Immunoenzymatic assay for the *in vitro* quantitative measurement of the human Insulin-like Growth Factor Binding Protein-3 (IGFBP-3) in serum.

II. GENERAL INFORMATION

A. Proprietary name : DIAsource IGFBP-3-ELISA Kit

B. Catalogue number : KAP1171 : 96 tests

C. Manufactured by : DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :

Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

The insulin-like growth factor (IGF) system is the primary regulator of normal body growth and regeneration, affecting cell proliferation, differentiation and apoptosis. In addition, the IGF-system appears to modify insulin sensitivity and long-term glucose metabolism. Finally, numerous epidemiological, experimental and clinical data indicate that the IGF-system is also involved in the development of several common cancers as well as frequent diseases such as atherosclerosis and type 2 diabetes mellitus.

The IGF-system consists of a family of closely related peptides, which includes the two primary growth promoting peptides, IGF-I and IGF-II, 6 specific high-affinity IGF-binding proteins (IGFBP-1 to -6), and a large non-IGF-binding glycoprotein, the "acid labile subunit" (ALS).

IGFBP-3 is the most abundant IGF-binding protein, accounting for as much as 75% or more of the circulating IGF-binding capacity in healthy subjects. IGFBP-3 shares functional properties with IGFBP-5 in that both peptides are able to form high molecular weight ternary complexes of ~150 kilo Dalton with ALS and either IGF-I or -II. However, IGFBP-5 circulates in much lower concentrations than IGFBP-3, and in healthy subjects the ternary complexes carry as much as 90% of IGFBP-3 but only about 50% of IGFBP-5.

Originally, the IGFBPs were thought to serve as IGF-carrier proteins, stabilizing plasma IGF levels and controlling the egress of IGF from the circulation to the extra-vascular compartment. Furthermore, it was assumed that IGFBP-complexed IGF was biologically more or less inactive, being deprived its ability to interact with the IGF-I receptor. However, it soon became apparent that in some experimental settings the IGFBPs stimulated rather than inhibited IGF-I mediated actions, and accordingly, the IGFBPs are now often referred to as *modulators* of IGF-I bioactivity. In addition, the majority of the IGFBPs, and in particular IGFBP-3, exerts IGF-I and IGF-II receptor independent effects, possibly involving interactions with specific receptors located at the cell surface and intracellular. For example, IGFBP-3 is nowadays considered to serve as an anti-cancer molecule, apparently protecting against several common cancers, and effects of IGFBP-3 on insulin signaling in cultured adipocytes have also been suggested.

The turnover of the ternary complexes is very slow, and the plasma concentration of IGFBP-3 remains stable throughout the day, being unaffected by short-term nutritional changes. Thus, the level of IGFBP-3 may be determined by one single measurement. GH is the primary regulator of IGFBP-3 as well as of IGF-I and ALS and therefore, all three peptides increase during the pubertal growth spurt, where after levels gradually decline with increasing age. In children, IGFBP-3 has been shown to correlate with the 24-h integrated GH secretion and in particular in children IGFBP-3 may be helpful in the diagnosis of GH deficiency.

IV. PRINCIPLES OF THE METHOD

The DIAsource IGFBP-3-ELISA is a solid phase Enzyme Amplified Sensitivity Immunoassay performed on microtiterplates. Calibrators and samples react with the capture monoclonal antibody (MAb 1) coated on microtiter well and with a monoclonal antibody (MAb 2) labelled with horseradish peroxidase (HRP). After an incubation period allowing the formation of a sandwich: coated MAb 1 – human IGFBP-3– MAb 2 – HRP, the microtiterplate is washed to remove unbound enzyme labelled antibody. Bound enzyme-labelled antibody is measured through a chromogenic reaction. Chromogenic solution (TMB) is added and incubated. The reaction is stopped with the addition of Stop Solution and the microtiterplate is then read at the appropriate wavelength. The amount of substrate turnover is determined colourimetrically by measuring the absorbance, which is proportional to the human IGFBP-3 concentration.

A calibration curve is plotted and IGFBP-3 concentration in samples is determined by interpolation from the calibration curve.

V. REAGENTS PROVIDED

Reagents	96 tests Kit	Color Code	Reconstitution
 Microtiterplate with 96 anti IGFBP3 (monoclonal antibodies) coated breakable wells	96 wells	blue	Ready for use
Ab HRP CONC Conjugate: HRP labelled anti-IGFBP-3 (monoclonal antibodies) in citrate buffer with bovine serum albumin and Proclin	1 vial 0.5 ml	yellow	Dilute 20 x with conjugate buffer
CONJ BUF Conjugate buffer: Tris buffer with bovine serum albumin and thymol	1 vial 10 ml	red	Ready for use
CAL N Calibrators - N = 1 to 5 , in phosphate buffer with bovine serum and thymol. (see exact values on vials labels) Calibrators are prediluted. ! Use dilution buffer as zero calibrator.	5 vials lyophilized	yellow	Add 1 ml distilled water
DIL BUF Dilution buffer: Phosphate buffer with bovine albumin, bovine serum and thymol.	1 vial 100 ml	black	Ready for use
CONTROL N Controls - N = 1 or 2 in human plasma with thymol. (see exact values on vials labels) Controls are prediluted.	2 vials lyophilized	silver	Add 1 ml distilled water
WASH SOLN CONC Wash Solution (Tris-HCl)	1 vial 10 ml	brown	Dilute 200 x with distilled water (use a magnetic stirrer).
CHROM TMB Chromogenic TMB Solution (Tetramethylbenzidine)	1 vial 12 ml	white	Ready for use
STOP SOLN Stop solution: HCl 1.0 N	1 vial 12 ml	black	Ready for use

Note: use the dilution buffer as zero calibrator.

The calibrators are standardized against the NIBSC/WHO recombinant IGFBP-3, reference reagent coded 93/560.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- High quality distilled water
- Pipettes for delivery of: 50 µl, 100 µl and 1 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Plastic tubes for dilution of samples
- Vortex mixer
- Magnetic stirrer
- Washer for microtiterplates
- Microtiterplate reader capable of reading at 450 nm and 650 nm (monochromatic reading)

VII. REAGENT PREPARATION

A. Calibrators: Reconstitute calibrators with 1 ml distilled water.

! Use dilution buffer as zero calibrator

B. Controls : Reconstitute the controls with 1 ml distilled water.

C. Working IGFBP-3-HRP conjugate : Prepare an adequate volume of conjugate solution by adding for example :100 µl of the 20 x concentrated IGFBP-3-HRP conjugate to 2 ml of conjugate buffer. Use a vortex to homogenize. Extemporaneous preparation is recommended.

D. Working Wash solution : Prepare an adequate volume of Working Wash solution by adding 199 volumes of distilled water to 1 volume of Wash Solution (200x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.
- Unused strips must be stored, at 2-8°C, in a sealed bag containing a desiccant until expiration date.
- After reconstitution, calibrators and controls are stable for one week at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 3 months. Avoid successive freezing and thawing
- The concentrated Wash Solution is stable at room temperature (18-25°C) until expiration date.
- The Working IGFBP-3-HRP conjugate is stable for 4 hours at room temperature (18-25°C), avoid direct sunlight.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

▪ Serum must be kept at 2 - 8°C.

- If the test is not run within 24 hours, storage in aliquots at -20°C is recommended. Avoid subsequent freeze thaw cycles.
- Prior to use, all samples should be at room temperature (18-25°C). It is recommended to vortex the samples before use.
- Do not use haemolysed samples.

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date.

Do not mix materials from different kit lots.

Bring all the reagents to room temperature (18-25°C) prior to use.

Thoroughly mix all reagents and samples by gentle agitation or swirling. Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.

Use a clean plastic container to prepare the Wash Solution.

In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.

For the dispensing of the Chromogenic Solution and the Stop Solution avoid pipettes with metal parts.

High precision pipettes or automated pipetting equipment will improve the precision.

Respect the incubation times.

To avoid drift, the time between pipetting of the first calibrator and the last sample must be limited to the time mentioned in section XIII paragraph E (Time delay).

Prepare a calibration curve for each run, do not use data from previous runs.

Dispense the Chromogenic Solution within 15 minutes following the washing of the microtiterplate.

During incubation with Chromogenic Solution, avoid direct sunlight on the microtiterplate.

Each well of the microplate can only be used once.

B. Procedure

1. Label one plain plastic tube for each sample.
- 2.. Dispense 1 ml of Dilution Buffer into each tube.
3. Add 10 µl of sample into these tubes.
4. Vortex pre-diluted samples ,reconstituted calibrators and controls.
5. Select the required number of wells for the run. The unused wells should be resealed in the bag with a desiccant and stored at 2-8°C.
6. Secure the wells into the holding frame.
7. Pipette 100 µl of dilution buffer as zero calibrator. Pipette 100 µl of each calibrator, control and diluted sample into the appropriate wells.
8. Pipette 50 µl of IGFBP-3-HRP conjugate solution into all the wells.
9. Incubate for 2 hours at room temperature (18-25°C).
10. Aspirate the liquid from each well.
11. Wash the plate 3 times with 400 µl wash solution and aspirate
12. Pipette 100 µl of the Chromogenic solution into each well within 15 minutes following the washing step.
13. Incubate the microtiterplate for 30 minutes at room temperature (18-25°C), avoid direct sunlight
14. Pipette 100 µl of Stop Solution into each well.
15. Read at 450 nm (reference filter 630 nm or 650 nm) within 1 hour and calculate the results as described in section XI

XI. CALCULATION OF RESULTS

1. Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).
2. Calculate the mean of duplicate determinations.
3. Plot the OD values (ordinate) for each calibrator against the corresponding concentration of Hu IGFBP-3 (abscissa) and draw a calibration curve through the calibrator points by connecting the plotted points with straight lines.
4. Read the concentration for each control and sample by interpolation on the calibration curve.
5. Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

IGFBP-3-ELISA		OD units
Calibrator	0 ng/ml	0.028
	460 ng/ml	0.114
	1270 ng/ml	0.311
	3020 ng/ml	0.778
	6710 ng/ml	1.403
	16070 ng/ml	2.588

XIII. PERFORMANCE AND LIMITATIONS

A. Detection Limit

Sixteen zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations above the average OD at zero binding, was 10 ng/ml

B. Specificity

Some potentially interfering hormones have been tested. At concentrations up to 10 µg/ml, none of the following hormones showed significant interference :

- rhIGF-BP1
- rhIGF-BP2
- rhIGF-BP4
- rhIGF-BP5
- rhIGF-BP6
- rhIGF-I
- rhIGF-II

C. Precision

INTRA ASSAY				INTER ASSAY			
Serum	N	\timesSD (ng/ml)	CV (%)	Serum	N	\timesSD (ng/ml)	CV (%)
A	22	827.3 ± 41.99	5.1	A	10	3074 ± 198.67	6.4
B	24	2081.7 ± 104.5	5.0	B	10	4951 ± 296.4	6

SD : Standard Deviation; CV: Coefficient of variation

D. Accuracy

RECOVERY TEST

Sample	Added IGFBP-3 (ng/ml)	Recovered IGFBP-3 (ng/ml)	Recovery (%)
Serum 1	3700	3880	104.8%
	5940	6680	112.0%
	10680	11700	109.5%
Serum 2	3700	3760	101.6%
	5940	6620	111.0%
	10680	11620	108.8%

DILUTION TEST

Sample	Dilution	Theoretical Concent. (ng/ml)	Measured Concent. (ng/ml)
Serum A	1/1		6720
	1/2	3360	3510
	1/4	1680	1540
	1/8	840	760
	1/16	420	440
Serum B	1/1		6510
	1/2	3250	3190
	1/4	1620	1410
	1/8	810	710
	1/16	400	440

Serial dilutions were made after the initial dilution with dilution buffer as described in the procedure section X. B. 2-3 .

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrators have been added to the coated wells.

TIME DELAY		
	0 min	30 min
S1	(ng/ml) 4060	(ng/ml) 4800
S2	5930	5900

XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls which contain azide will interfere with the enzymatic reaction and cannot be used.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises
- It is recommended that Controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.
- It is good practise to check visually the curve fit selected by the computer.

XV REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

Age Group	MALES (ng/ml)			FEMALES (ng/ml)		
	Mean	Range	N	Mean	Range	N
0 - 2 years	2639	1481-4481	15	2348	1398-3485	12
3 - 5 years	2405	1479-3053	12	2752	2059-3325	13
6 - 8 years	3186	2174-5128	17	3282	2469-4495	13
9 - 11 years	3263	2020-4705	21	3298	2343-4640	11
12 - 14 years	3672	2239-5971	19	4241	3000-7022	14
15 - 17 years	4031	2710-5235	21	4181	2539-6607	20
18 - 20 years	3826	2304-5537	10	3709	2272-6102	9
21 - 30 years	3372	2093-4553	11	3766	2704-5595	10
31 - 40 years	2704	1190-4140	14	3372	2660-4533	12
41 - 50 years	3886	2318-6897	18	3247	2323-4046	16
51 - 60 years	3176	2113-4625	16	3830	1603-5998	15
> 60 years	2827	1155-3877	23	3621	1995-6505	21

XVI. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with all reagents. Stop Solution contains HCl. In case of contact, wash thoroughly with water.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

For more information, see Material Safety Data Sheet (MSDS)

XVII. BIBLIOGRAPHY

1. LEROITH D., BONDY C., YAKAR S., LIU JL., BUTLER A. **The somatomedin hypothesis : 2001.** Endocr Rev 2001; 22:53-74.
2. POLLAK MN., SCHERNHAMMER ES., HANKINSON SE. **Insulin-like growth factors and neoplasia.** Nat Rev Cancer 2004; 4:505-518.
3. YUEN K., FRYSTYK J., UMPLERBY M., FRYKLUND L., DUNGER D. **Changes in free rather than total insulin-like growth factor-I enhance insulin sensitivity and suppress endogenous peak growth hormone (GH) release following short-term low dose GH administration in young healthy adults.** J Clin Endocrinol Metab 2004; 89:3956-3964.
4. KHANDWALA HM., McCUTCHEON IE., FLYVBJERG A., FRIEND KE. **The effects of insulin-like growth factors on tumorigenesis and neoplastic growth.** Endocr Rev 2000; 21:215-244.
5. RENEHAN AG., ZWAHLEN M., MINDER PC., O'DWYSER ST., SHALET PS., EGGER PM. **Insulin-like growth factor (IGF)-I, IGF binding protein-3, and cancer risk : systematic review and meta-regression analysis.** The Lancet 2004; 363:1346-1353.
6. JUUL A., SCHEIKE T., DAVIDSEN M., GYLLENborg J., JORGENSEN T. **Low serum insulin-like growth factor I is associated with increased risk of ischemic heart disease : a population-based case-control study.** Circulation 2002; 106:939-944.
7. SANDHU MS., HEALD AH., GIBSON JM., CRUICKSHANK JK., DUNGER DB., WAREHAM NJ. **Circulating concentrations of insulin-like growth factor-I and development of glucose intolerance : a prospective observational study.** The Lancet 2002; 359:1740-1745.
8. VAESEN N., HEUTINK P., JANSSEN JA., WITTEMAN JC., TESTERS L., HOFMAN A., LAMBERTS SW., OOSTRA BA., POLS HA., VAN DUIJN CM. **A polymorphism in the gene for IGF-I : functional properties and risk for type 2 diabetes and myocardial infarction.** Diabetes 2001; 50:637-642.
9. JUUL A. **Serum levels of insulin-like growth factor-I and its binding proteins in health and disease.** Growth Horm IGF Res 2003; 13:113-170
10. FIRTH SM., BAXTER RC. **Cellular actions of the insulin-like growth factor binding proteins** Endocr Rev 2002; 23:824-854.
11. BAXTER RC., MEKA S., FIRTH SM. **Molecular distribution of IGF binding protein-5 in human serum.** J Clin Endocrinol Metab 2002; 87:271-276.
12. RICORT JM. **Insulin-like growth factor binding protein (IGFBP) signaling.** Growth Horm IGF Res 2004; 14:277-286.
13. JONES JI., CLEMMONS DR. **Insulin-like growth factors and their binding proteins : biological actions.** Endocr Rev 1995; 16:3-34.
14. ALI O., COHEN P., LEE KW. **Epidemiology and biology of insulin-like growth factor binding protein-3 (IGFBP-3) as an anti-cancer molecule.** Horm Metab Res 2003; 35:726-733.
15. CHAN SS., TWIGG SM., FIRTH SM., BAXTER RC. **Insulin-like growth factor binding protein-3 leads to insulin resistance in adipocytes.** J Clin Endocrinol Metab 2005; 90:6588-6595.
16. JUUL A., MAIN K., BLUM WF., LINDHOLM J., RANKE MB., SKAKKEBAEK NE. **The ration between serum levels of insulin-like growth factor (IGF)-I and the IGF binding proteins (IGFBP-1, 2 and 3) decreases with age in healthy adults and is increased in acromegalic patients.** Clin Endocrinol (Oxf) 1994; 41:85-93.
17. BLUM WF., ALBERTSSON-WIKLAND K., ROSBERG S., RANKE MB. **Serum levels of insulin-like growth factor I (IGF-I) and IGF binding protein 3 reflect spontaneous growth hormone secretion.** J Clin Endocrinol Metab 1994; 76:1610-1616.
18. FRYSTYK J., IVARSEN P., SKJAERBAEK C., FLYVBJERG A., PEDERSEN EB., ORSKOV H. **Serum-free insulin-like growth factor I correlates with clearance in patients with chronic renal failure.** Kidney Int 1999; 56:2076-2084.
19. FRYSTYK J. **In Endocrinology and Metabolism – Clinics of North America 2005 : Endocrinology of aging, Chapter XI : Aging somatotropic axis mechanisms and implications of IGFBP adaptation.**

XVIII.**SUMMARY OF THE PROTOCOL**

	CALIBRATORS CONTROLS µl	SAMPLE(S) µl
DILUTION OF SAMPLES		
Dilution buffer	-	1000
Sample	-	10
Shaking	Vortex	
INCUBATION		
Calibrators (0 to 5), controls	100	-
Diluted Samples,	-	100
Diluted Conjugate	50	50
Incubate for 2 hours at room temperature (18-25°C). Aspirate the contents of each well. Wash 3 times with 400 µl of Wash Solution and aspirate.		
Chromogenic Solution	100	100
Incubate for 30 min at room temperature (18-25°C)		
Stop Solution	100	100
Read on a microtiterplate reader and record the absorbance of each well at 450 nm (versus 630 or 650 nm)		

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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XVII. BIBLIOGRAPHIE

1. LEROITH D., BONDY C., YAKAR S., LIU JL., BUTLER A. **The somatomedin hypothesis : 2001.** Endocr Rev 2001; 22:53-74.
2. POLLAK MN., SCHERNHAMMER ES., HANKINSON SE. **Insulin-like growth factors and neoplasia.** Nat Rev Cancer 2004; 4:505-518.
3. YUEN K., FRYSTYK J., UMPLEBY M., FRYKLUND L., DUNGER D. **Changes in free rather than total insulin-like growth factor-I enhance insulin sensitivity and suppress endogenous peak growth hormone (GH) release following short-term low dose GH administration in young healthy adults.** J Clin Endocrinol Metab 2004; 89:3956-3964.
4. KHANDWALA HM., McCUTCHEON IE., FLYVBJERG A., FRIEND KE. **The effects of insulin-like growth factors on tumorigenesis and neoplastic growth.** Endocr Rev 2000; 21:215-244.
5. RENEHAN AG., ZWAHLEN M., MINDER PC., O'DWYSER ST., SHALET PS., EGGER PM. **Insulin-like growth factor (IGF)-I, IGF binding protein-3, and cancer risk : systematic review and meta-regression analysis.** The Lancet 2004; 363:1346-1353.
6. JUUL A., SCHEIKE T., DAVIDSEN M., GYLLENborg J., JORGENSEN T. **Low serum insulin-like growth factor I is associated with increased risk of ischemic heart disease : a population-based case-control study.** Circulation 2002; 106:939-944.
7. SANDHU MS., HEALD AH., GIBSON JM., CRUCKSHANK JK., DUNGER DB., WAREHAM NJ. **Circulating concentrations of insulin-like growth factor-I and development of glucose intolerance : a prospective observational study.** The Lancet 2002; 359:1740-1745.
8. VAESEN N., HEUTINK P., JANSSEN JA., WITTEMAN JC., TESTERS L., HOFMAN A., LAMBERTS SW., OOSTRA BA., POLS HA., VAN DUIJN CM. **A polymorphism in the gene for IGF-I : functional properties and risk for type 2 diabetes and myocardial infarction.** Diabetes 2001; 50:637-642.
9. JUUL A. **Serum levels of insulin-like growth factor-I and its binding proteins in health and disease.** Growth Horm IGF Res 2003; 13:113-170
10. FIRTH SM., BAXTER RC. **Cellular actions of the insulin-like growth factor binding proteins** Endocr Rev 2002; 23:824-854.
11. BAXTER RC., MEKA S., FIRTH SM. **Molecular distribution of IGF binding protein-5 in human serum.** J Clin Endocrinol Metab 2002; 87:271-276.
12. RICORT JM. **Insulin-like growth factor binding protein (IGFBP) signaling.** Growth Horm IGF Res 2004; 14:277-286.
13. JONES JI., CLEMMONS DR. **Insulin-like growth factors and their binding proteins : biological actions.** Endocr Rev 1995; 16:3-34.
14. ALI O., COHEN P., LEE KW. **Epidemiology and biology of insulin-like growth factor binding protein-3 (IGFBP-3) as an anti-cancer molecule.** Horm Metab Res 2003; 35:726-733.
15. CHAN SS., TWIGG SM., FIRTH SM., BAXTER RC. **Insulin-like growth factor binding protein-3 leads to insulin resistance in adipocytes.** J Clin Endocrinol Metab 2005; 90:6588-6595.
16. JUUL A., MAIN K., BLUM WF., LINDHOLM J., RANKE MB., SKAKKEBAEK NE. **The ration between serum levels of insulin-like growth factor (IGF)-I and the IGF binding proteins (IGFBP-1, 2 and 3) decreases with age in healthy adults and is increased in acromegalic patients.** Clin Endocrinol (Oxf) 1994; 41:85-93.
17. BLUM WF., ALBERTSSON-WIKLAND K., ROSBERG S., RANKE MB. **Serum levels of insulin-like growth factor I (IGF-I) and IGF binding protein 3 reflect spontaneous growth hormone secretion.** J Clin Endocrinol Metab 1994; 76:1610-1616.
18. FRYSTYK J., IVARSEN P., SKJAERBAEK C., FLYVBJERG A., PEDERSEN EB., ORSKOV H. **Serum-free insulin-like growth factor I correlates with clearance in patients with chronic renal failure.** Kidney Int 1999; 56:2076-2084.
19. FRYSTYK J. In Endocrinology and Metabolism – Clinics of North America 2005 : Endocrinology of aging, Chapter XI : Aging somatotropic axis mechanisms and implications of IGFBP adaptation.

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1. LEROITH D., BONDY C., YAKAR S., LIU JL., BUTLER A. **The somatomedin hypothesis : 2001.** Endocr Rev 2001; 22:53-74.
2. POLLAK MN., SCHERNHAMMER ES., HANKINSON SE. **Insulin-like growth factors and neoplasia.** Nat Rev Cancer 2004; 4:505-518.
3. YUEN K., FRYSTYK J., UMPLEBY M., FRYKLUND L., DUNGER D. **Changes in free rather than total insulin-like growth factor-I enhance insulin sensitivity and suppress endogenous peak growth hormone (GH) release following short-term low dose GH administration in young healthy adults.** J Clin Endocrinol Metab 2004; 89:3956-3964.
4. KHANDWALA HM., McCUTCHEON IE., FLYVBJERG A., FRIEND KE. **The effects of insulin-like growth factors on tumorigenesis and neoplastic growth.** Endocr Rev 2000; 21:215-244.
5. RENEHAN AG., ZWAHLEN M., MINDER PC., O'DWYSER ST., SHALET PS., EGGER PM. **Insulin-like growth factor (IGF)-I, IGF binding protein-3, and cancer risk : systematic review and meta-regression analysis.** The Lancet 2004; 363:1346-1353.
6. JUUL A., SCHEIKE T., DAVIDSEN M., GYLLENborg J., JORGENSEN T. **Low serum insulin-like growth factor I is associated with increased risk of ischemic heart disease : a population-based case-control study.** Circulation 2002; 106:939-944.
7. SANDHU MS., HEALD AH., GIBSON JM., CRUICKSHANK JK., DUNGER DB., WAREHAM NJ. **Circulating concentrations of insulin-like growth factor-I and development of glucose intolerance : a prospective observational study.** The Lancet 2002; 359:1740-1745.
8. VAESEN N., HEUTINK P., JANSSEN JA., WITTEMAN JC., TESTERS L., HOFMAN A., LAMBERTS SW., OOSTRA BA., POLS HA., VAN DUIJN CM. **A polymorphism in the gene for IGF-I : functional properties and risk for type 2 diabetes and myocardial infarction.** Diabetes 2001; 50:637-642.
9. JUUL A. **Serum levels of insulin-like growth factor-I and its binding proteins in health and disease.** Growth Horm IGF Res 2003; 13:113-170
10. FIRTH SM., BAXTER RC. **Cellular actions of the insulin-like growth factor binding proteins** Endocr Rev 2002; 23:824-854.
11. BAXTER RC., MEKA S., FIRTH SM. **Molecular distribution of IGF binding protein-5 in human serum.** J Clin Endocrinol Metab 2002; 87:271-276.
12. RICORT JM. **Insulin-like growth factor binding protein (IGFBP) signaling.** Growth Horm IGF Res 2004; 14:277-286.
13. JONES JI., CLEMMONS DR. **Insulin-like growth factors and their binding proteins : biological actions.** Endocr Rev 1995; 16:3-34.
14. ALI O., COHEN P., LEE KW. **Epidemiology and biology of insulin-like growth factor binding protein-3 (IGFBP-3) as an anti-cancer molecule.** Horm Metab Res 2003; 35:726-733.
15. CHAN SS., TWIGG SM., FIRTH SM., BAXTER RC. **Insulin-like growth factor binding protein-3 leads to insulin resistance in adipocytes.** J Clin Endocrinol Metab 2005; 90:6588-6595.
16. JUUL A., MAIN K., BLUM WF., LINDHOLM J., RANKE MB., SKAKKEBAEK NE. **The ration between serum levels of insulin-like growth factor (IGF)-I and the IGF binding proteins (IGFBP-1, 2 and 3) decreases with age in healthy adults and is increased in acromegalic patients.** Clin Endocrinol (Oxf) 1994; 41:85-93.
17. BLUM WF., ALBERTSSON-WIKLAND K., ROSBERG S., RANKE MB. **Serum levels of insulin-like growth factor I (IGF-I) and IGF binding protein 3 reflect spontaneous growth hormone secretion.** J Clin Endocrinol Metab 1994; 76:1610-1616.
18. FRYSTYK J., IVARSEN P., SKJAERBAEK C., FLYVBJERG A., PEDERSEN EB., ORSKOV H. **Serum-free insulin-like growth factor I correlates with clearance in patients with chronic renal failure.** Kidney Int 1999; 56:2076-2084.
19. FRYSTYK J. In Endocrinology and Metabolism – Clinics of North America 2005 : Endocrinology of aging, Chapter XI : Aging somatotropic axis mechanisms and implications of IGFBP adaptation.

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XVII. BIBLIOGRAFIA

1. LEROITH D., BONDY C., YAKAR S., LIU JL., BUTLER A. **The somatomedin hypothesis : 2001.** Endocr Rev 2001; 22:53-74.
2. POLLAK MN., SCHERNHAMMER ES., HANKINSON SE. **Insulin-like growth factors and neoplasia.** Nat Rev Cancer 2004; 4:505-518.
3. YUEN K., FRYSTYK J., UMPLEBY M., FRYKLUND L., DUNGER D. **Changes in free rather than total insulin-like growth factor-I enhance insulin sensitivity and suppress endogenous peak growth hormone (GH) release following short-term low dose GH administration in young healthy adults.** J Clin Endocrinol Metab 2004; 89:3956-3964.
4. KHANDWALA HM., McCUTCHEON IE., FLYVBJERG A., FRIEND KE. **The effects of insulin-like growth factors on tumorigenesis and neoplastic growth.** Endocr Rev 2000; 21:215-244.
5. RENEHAN AG., ZWAHLEN M., MINDER PC., O'DWYSER ST., SHALET PS., EGGER PM. **Insulin-like growth factor (IGF)-I, IGF binding protein-3, and cancer risk : systematic review and meta-regression analysis.** The Lancet 2004; 363:1346-1353.
6. JUUL A., SCHEIKE T., DAVIDSEN M., GYLLENborg J., JORGENSEN T. **Low serum insulin-like growth factor I is associated with increased risk of ischemic heart disease : a population-based case-control study.** Circulation 2002; 106:939-944.
7. SANDHU MS., HEALD AH., GIBSON JM., CRUICKSHANK JK., DUNGER DB., WAREHAM NJ. **Circulating concentrations of insulin-like growth factor-I and development of glucose intolerance : a prospective observational study.** The Lancet 2002; 359:1740-1745.
8. VAESEN N., HEUTINK P., JANSSEN JA., WITTEMAN JC., TESTERS L., HOFMAN A., LAMBERTS SW., OOSTRA BA., POLS HA., VAN DUIJN CM. **A polymorphism in the gene for IGF-I : functional properties and risk for type 2 diabetes and myocardial infarction.** Diabetes 2001; 50:637-642.
9. JUUL A. **Serum levels of insulin-like growth factor-I and its binding proteins in health and disease.** Growth Horm IGF Res 2003; 13:113-170
10. FIRTH SM., BAXTER RC. **Cellular actions of the insulin-like growth factor binding proteins** Endocr Rev 2002; 23:824-854.
11. BAXTER RC., MEKA S., FIRTH SM. **Molecular distribution of IGF binding protein-5 in human serum.** J Clin Endocrinol Metab 2002; 87:271-276.
12. RICORT JM. **Insulin-like growth factor binding protein (IGFBP) signaling.** Growth Horm IGF Res 2004; 14:277-286.
13. JONES JI., CLEMMONS DR. **Insulin-like growth factors and their binding proteins : biological actions.** Endocr Rev 1995; 16:3-34.
14. ALI O., COHEN P., LEE KW. **Epidemiology and biology of insulin-like growth factor binding protein-3 (IGFBP-3) as an anti-cancer molecule.** Horm Metab Res 2003; 35:726-733.
15. CHAN SS., TWIGG SM., FIRTH SM., BAXTER RC. **Insulin-like growth factor binding protein-3 leads to insulin resistance in adipocytes.** J Clin Endocrinol Metab 2005; 90:6588-6595.
16. JUUL A., MAIN K., BLUM WF., LINDHOLM J., RANKE MB., SKAKKEBAEK NE. **The ration between serum levels of insulin-like growth factor (IGF)-I and the IGF binding proteins (IGFBP-1, 2 and 3) decreases with age in healthy adults and is increased in acromegalic patients.** Clin Endocrinol (Oxf) 1994; 41:85-93.
17. BLUM WF., ALBERTSSON-WIKLAND K., ROSBERG S., RANKE MB. **Serum levels of insulin-like growth factor I (IGF-I) and IGF binding protein 3 reflect spontaneous growth hormone secretion.** J Clin Endocrinol Metab 1994; 76:1610-1616.
18. FRYSTYK J., IVARSEN P., SKJAERBAEK C., FLYVBJERG A., PEDERSEN EB., ORSKOV H. **Serum-free insulin-like growth factor I correlates with clearance in patients with chronic renal failure.** Kidney Int 1999; 56:2076-2084.
19. FRYSTYK J. In Endocrinology and Metabolism – Clinics of North America 2005 : Endocrinology of aging, Chapter XI : Aging somatotrophic axis mechanisms and implications of IGFBP adaptation.

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- LEROITH D., BONDY C., YAKAR S., LIU JL., BUTLER A. **The somatomedin hypothesis : 2001.** Endocr Rev 2001; 22:53-74.
- POLLAK MN., SCHERNHAMMER ES., HANKINSON SE. Insulin-like growth factors and neoplasia. Nat Rev Cancer 2004; 4:505-518.
- YUEN K., FRYSTYK J., UMPLEBY M., FRYKLUND L., DUNGER D. Changes in free rather than total insulin-like growth factor-I enhance insulin sensitivity and suppress endogenous peak growth hormone (GH) release following short-term low dose GH administration in young healthy adults. J Clin Endocrinol Metab 2004; 89:3956-3964.
- KHANDWALA HM., McCUTCHEON IE., FLYVBJERG A., FRIEND KE. The effects of insulin-like growth factors on tumorigenesis and neoplastic growth. Endocr Rev 2000; 21:215-244.
- RENEHAN AG., ZWAHLEN M., MINDER PC., O'DWYSER ST., SHALET PS., EGGER PM. Insulin-like growth factor (IGF)-I, IGF binding protein-3, and cancer risk : systematic review and meta-regression analysis. The Lancet 2004; 363:1346-1353.
- JUUL A., SCHEIKE T., DAVIDSEN M., GYLLENborg J., JORGENSEN T. Low serum insulin-like growth factor I is associated with increased risk of ischemic heart disease : a population-based case-control study. Circulation 2002; 106:939-944.
- SANDHU MS., HEALD AH., GIBSON JM., CRUICKSHANK JK., DUNGER DB., WAREHAM NJ. Circulating concentrations of insulin-like growth factor-I and development of glucose intolerance : a prospective observational study. The Lancet 2002; 359:1740-1745.
- VAESSEN N., HEUTINK P., JANSEN JA., WITTEMAN JC., TESTERS L., HOFMAN A., LAMBERTS SW., OOSTRA BA., POLS HA., VAN DUIJN CM. A polymorphism in the gene for IGF-I : functional properties and risk for type 2 diabetes and myocardial infarction. Diabetes 2001; 50:637-642.
- JUUL A. Serum levels of insulin-like growth factor-I and its binding proteins in health and disease. Growth Horm IGF Res 2003; 13:113-170
- FIRTH SM., BAXTER RC. Cellular actions of the insulin-like growth factor binding proteins Endocr Rev 2002; 23:824-854.
- BAXTER RC., MEKA S., FIRTH SM. Molecular distribution of IGF binding protein-5 in human serum. J Clin Endocrinol Metab 2002; 87:271-276.
- RICORT JM. Insulin-like growth factor binding protein (IGFBP) signaling. Growth Horm IGF Res 2004; 14:277-286.
- JONES JI., CLEMMONS DR. Insulin-like growth factors and their binding proteins : biological actions. Endocr Rev 1995; 16:3-34.
- ALI O., COHEN P., LEE KW. Epidemiology and biology of insulin-like growth factor binding protein-3 (IGFBP-3) as an anti-cancer molecule. Horm Metab Res 2003; 35:726-733.
- CHAN SS., TWIGG SM., FIRTH SM., BAXTER RC. Insulin-like growth factor binding protein-3 leads to insulin resistance in adipocytes. J Clin Endocrinol Metab 2005; 90:6588-6595.
- JUUL A., MAIN K., BLUM WF., LINDHOLM J., RANKE MB., SKAKKEBAEK NE. The ration between serum levels of insulin-like growth factor (IGF)-I and the IGF binding proteins (IGFBP-1, 2 and 3) decreases with age in healthy adults and is increased in acromegalic patients. Clin Endocrinol (Oxf) 1994; 41:85-93.
- BLUM WF., ALBERTSSON-WIKLAND K., ROSBERG S., RANKE MB. Serum levels of insulin-like growth factor I (IGF-I) and IGF binding protein 3 reflect spontaneous growth hormone secretion. J Clin Endocrinol Metab 1994; 76:1610-1616.
- FRYSTYK J., IVARSEN P., SKJAERBAEK C., FLYVBJERG A., PEDERSEN EB., ORSKOV H. Serum-free insulin-like growth factor I correlates with clearance in patients with chronic renal failure. Kidney Int 1999; 56:2076-2084.
- FRYSTYK J. In Endocrinology and Metabolism – Clinics of North America 2005 : Endocrinology of aging, Chapter XI : Aging somatotrophic axis mechanisms and implications of IGFBP adaptation.

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XVII. BIBLIOGRAFIA

- LEROITH D., BONDY C., YAKAR S., LIU JL., BUTLER A. **The somatomedin hypothesis : 2001.** Endocr Rev 2001; 22:53-74.
- POLLAK MN., SCHERNHAMMER ES., HANKINSON SE. **Insulin-like growth factors and neoplasia.** Nat Rev Cancer 2004; 4:505-518.
- YUEN K., FRYSTYK J., UMPLEBY M., FRYKLUND L., DUNGER D. **Changes in free rather than total insulin-like growth factor-I enhance insulin sensitivity and suppress endogenous peak growth hormone (GH) release following short-term low dose GH administration in young healthy adults.** J Clin Endocrinol Metab 2004; 89:3956-3964.
- KHANDWALA HM., McCUTCHEON IE., FLYVBJERG A., FRIEND KE. **The effects of insulin-like growth factors on tumorigenesis and neoplastic growth.** Endocr Rev 2000; 21:215-244.
- RENEHAN AG., ZWAHLEN M., MINDER PC., O'DWYSER ST., SHALET PS., EGGER PM. **Insulin-like growth factor (IGF)-I, IGF binding protein-3, and cancer risk : systematic review and meta-regression analysis.** The Lancet 2004; 363:1346-1353.
- JUUL A., SCHEIKE T., DAVIDSEN M., GYLLENborg J., JORGENSEN T. **Low serum insulin-like growth factor I is associated with increased risk of ischemic heart disease : a population-based case-control study.** Circulation 2002; 106:939-944.
- SANDHU MS., HEALD AH., GIBSON JM., CRUICKSHANK JK., DUNGER DB., WAREHAM NJ. **Circulating concentrations of insulin-like growth factor-I and development of glucose intolerance : a prospective observational study.** The Lancet 2002; 359:1740-1745.
- VAESSEN N., HEUTINK P., JANSSEN JA., WITTEMAN JC., TESTERS L., HOFMAN A., LAMBERTS SW., OOSTRA BA., POLS HA., VAN DUIJN CM. **A polymorphism in the gene for IGF-I : functional properties and risk for type 2 diabetes and myocardial infarction.** Diabetes 2001; 50:637-642.
- JUUL A. **Serum levels of insulin-like growth factor-I and its binding proteins in health and disease.** Growth Horm IGF Res 2003; 13:113-170
- FIRTH SM., BAXTER RC. **Cellular actions of the insulin-like growth factor binding proteins** Endocr Rev 2002; 23:824-854.
- BAXTER RC., MEKA S., FIRTH SM. **Molecular distribution of IGF binding protein-5 in human serum.** J Clin Endocrinol Metab 2002; 87:271-276.
- RICORT JM. **Insulin-like growth factor binding protein (IGFBP) signaling.** Growth Horm IGF Res 2004; 14:277-286.
- JONES JI., CLEMMONS DR. **Insulin-like growth factors and their binding proteins : biological actions.** Endocr Rev 1995; 16:3-34.
- ALI O., COHEN P., LEE KW. **Epidemiology and biology of insulin-like growth factor binding protein-3 (IGFBP-3) as an anti-cancer molecule.** Horm Metab Res 2003; 35:726-733.
- CHAN SS., TWIGG SM., FIRTH SM., BAXTER RC. **Insulin-like growth factor binding protein-3 leads to insulin resistance in adipocytes.** J Clin Endocrinol Metab 2005; 90:6588-6595.

16. JUUL A., MAIN K., BLUM WF., LINDHOLM J., RANKE MB., SKAKKEBAEK NE.
The ration between serum levels of insulin-like growth factor (IGF)-I and the IGF binding proteins (IGFBP-1, 2 and 3) decreases with age in healthy adults and is increased in acromegalic patients.
 Clin Endocrinol (Oxf) 1994; 41:85-93.
17. BLUM WF., ALBERTSSON-WIKLAND K., ROSBERG S., RANKE MB.
Serum levels of insulin-like growth factor I (IGF-I) and IGF binding protein 3 reflect spontaneous growth hormone secretion.
 J Clin Endocrinol Metab 1994; 76:1610-1616.
18. FRYSTYK J., IVARSEN P., SKJAERBAEK C., FLYVBJERG A., PEDERSEN EB., ORSKOV H.
Serum-free insulin-like growth factor I correlates with clearance in patients with chronic renal failure.
 Kidney Int 1999; 56:2076-2084.
19. FRYSTYK J.
 In Endocrinology and Metabolism – Clinics of North America 2005 : Endocrinology of aging, Chapter XI : Aging somatotropic axis mechanisms and implications of IGFBP adaptation.

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