IVD

CE

# 17OH-PROGESTERONE ELISA

KAP1401

# History

## Summary of change:

Previous Version:	Current Version:																																																
200224-1	220516																																																
Old DiaSource logo 	New DiaSource logo 																																																
<b>IX SPECIMEN COLLECTION AND PREPARATION</b>	<b>IX SPECIMEN COLLECTION AND PREPARATION</b>																																																
<ul style="list-style-type: none"> <li>This kit is suitable for serum, heparinized plasma or EDTA plasma samples.</li> <li>Plasma samples provide similar results than serum:  <math>Y(\text{EDTA plasma}) = 0.97 X(\text{serum}) - 0.01</math> r = 0.99 n=20  <math>Y(\text{heparinized plasma}) = 0.95 X(\text{serum}) - 0.07</math> r = 0.99 n =20</li> <li>Serum or plasma samples must be kept at 2-8°C.</li> <li>If the test is not run within 24 hrs, <b>sampling and storage at -20°C is recommended.</b></li> <li>Avoid subsequent freeze-thaw cycles.</li> </ul>	<ul style="list-style-type: none"> <li>This kit is suitable for serum, heparinized plasma or EDTA plasma samples.</li> <li>Plasma samples provide similar results than serum:  <math>Y(\text{EDTA plasma}) = 0.93 X(\text{serum}) - 0.03</math> r = 0.98 n=20  <math>Y(\text{heparinized plasma}) = 0.96 X(\text{serum}) - 0.01</math> r = 0.99 n=20</li> <li>Serum or plasma samples must be kept at 2-8°C.</li> <li>If the test is not run within 24 hrs, <b>sampling and storage at -20°C is recommended.</b></li> <li>Avoid subsequent freeze-thaw cycles.</li> </ul>																																																
<b>V. REAGENTS PROVIDED</b>	<b>V. REAGENTS PROVIDED</b>																																																
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 25%;">CHROM</th><th style="width: 25%;">TMB</th><th style="width: 25%;"></th><th style="width: 25%;"></th></tr> </thead> <tbody> <tr> <td>Chromogenic solution TMB (Tetramethylbenzidine)</td><td>1 vial 12 ml</td><td>orange</td><td>Ready for use</td></tr> <tr> <td>STOP</td><td>SOLN</td><td>1 vial 12 ml</td><td>Ready for use</td></tr> <tr> <td colspan="4">Stop solution HCl 1M</td></tr> </tbody> </table>	CHROM	TMB			Chromogenic solution TMB (Tetramethylbenzidine)	1 vial 12 ml	orange	Ready for use	STOP	SOLN	1 vial 12 ml	Ready for use	Stop solution HCl 1M				<p>Removal of the column relative to color code.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 25%;">CHROM</th><th style="width: 25%;">TMB</th><th style="width: 25%;"></th><th style="width: 25%;"></th></tr> </thead> <tbody> <tr> <td>Chromogenic solution TMB (Tetramethylbenzidine)</td><td>1 vial 13 ml</td><td></td><td>Ready for use</td></tr> <tr> <td>STOP</td><td>SOLN</td><td>1 vial 13 ml</td><td>Ready for use</td></tr> <tr> <td colspan="4">Stop solution 0.2M H<sub>2</sub>SO<sub>4</sub></td></tr> </tbody> </table>	CHROM	TMB			Chromogenic solution TMB (Tetramethylbenzidine)	1 vial 13 ml		Ready for use	STOP	SOLN	1 vial 13 ml	Ready for use	Stop solution 0.2M H <sub>2</sub> SO <sub>4</sub>																			
CHROM	TMB																																																
Chromogenic solution TMB (Tetramethylbenzidine)	1 vial 12 ml	orange	Ready for use																																														
STOP	SOLN	1 vial 12 ml	Ready for use																																														
Stop solution HCl 1M																																																	
CHROM	TMB																																																
Chromogenic solution TMB (Tetramethylbenzidine)	1 vial 13 ml		Ready for use																																														
STOP	SOLN	1 vial 13 ml	Ready for use																																														
Stop solution 0.2M H <sub>2</sub> SO <sub>4</sub>																																																	
<b>XII. TYPICAL DATA</b>	<b>XII. TYPICAL DATA</b>																																																
The following data are for illustration only and should never be used instead of the real time calibration curve.	The following data are for illustration only and should never be used instead of the real time calibration curve.																																																
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">17OH-Progesterone ELISA</th><th>mod units</th></tr> </thead> <tbody> <tr> <td>Calibrator</td><td>0 ng/ml</td><td>2922</td></tr> <tr> <td></td><td>0.08 ng/ml</td><td>2768</td></tr> <tr> <td></td><td>0.4 ng/ml</td><td>1904</td></tr> <tr> <td></td><td>1.3 ng/ml</td><td>1088</td></tr> <tr> <td></td><td>2.5 ng/ml</td><td>655</td></tr> <tr> <td></td><td>7.5 ng/ml</td><td>241</td></tr> <tr> <td></td><td>15 ng/ml</td><td>164</td></tr> </tbody> </table>	17OH-Progesterone ELISA		mod units	Calibrator	0 ng/ml	2922		0.08 ng/ml	2768		0.4 ng/ml	1904		1.3 ng/ml	1088		2.5 ng/ml	655		7.5 ng/ml	241		15 ng/ml	164	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">17OH-Progesterone ELISA</th><th>OD units</th></tr> </thead> <tbody> <tr> <td>Calibrator</td><td>0 ng/ml</td><td>2.306</td></tr> <tr> <td></td><td>0.09 ng/ml</td><td>2.205</td></tr> <tr> <td></td><td>0.4 ng/ml</td><td>1.744</td></tr> <tr> <td></td><td>1.4 ng/ml</td><td>1.075</td></tr> <tr> <td></td><td>2.7 ng/ml</td><td>0.737</td></tr> <tr> <td></td><td>7.9 ng/ml</td><td>0.355</td></tr> <tr> <td></td><td>15.6 ng/ml</td><td>0.235</td></tr> </tbody> </table>	17OH-Progesterone ELISA		OD units	Calibrator	0 ng/ml	2.306		0.09 ng/ml	2.205		0.4 ng/ml	1.744		1.4 ng/ml	1.075		2.7 ng/ml	0.737		7.9 ng/ml	0.355		15.6 ng/ml	0.235
17OH-Progesterone ELISA		mod units																																															
Calibrator	0 ng/ml	2922																																															
	0.08 ng/ml	2768																																															
	0.4 ng/ml	1904																																															
	1.3 ng/ml	1088																																															
	2.5 ng/ml	655																																															
	7.5 ng/ml	241																																															
	15 ng/ml	164																																															
17OH-Progesterone ELISA		OD units																																															
Calibrator	0 ng/ml	2.306																																															
	0.09 ng/ml	2.205																																															
	0.4 ng/ml	1.744																																															
	1.4 ng/ml	1.075																																															
	2.7 ng/ml	0.737																																															
	7.9 ng/ml	0.355																																															
	15.6 ng/ml	0.235																																															
<b>XIII. PERFORMANCE AND LIMITATIONS</b>	<b>XIII. PERFORMANCE AND LIMITATIONS</b>																																																
<b>A. Detection Limit</b>	<b>A. Detection Limit</b>																																																
The LOB was calculated to be 0.04 ng/ml. The LOD was calculated to be 0.07 ng/ml. The LOQ was calculated to be 0.1 ng/ml with CV of 20%.	The LOB was calculated to be 0.03 ng/ml. The LOD was calculated to be 0.06 ng/ml. The LOQ was calculated to be 0.1 ng/ml with CV of 20%.																																																

## B. Specificity

Compound	Added amount (ng/ml)	Cross-Reactivity (%)
17OH-Progesterone	-	100
Progesterone	1000	0.58
17- $\alpha$ -hydroxypregnenolone	1000	0.5
21-deoxycortisol	1000	4.3
Pregnenolone	1000	0.03
11-deoxycortisol	1000	0.29
Corticosterone	10000	0.01
11-deoxycorticosterone	1000	0.28
Testosterone	5000	0.01
Androstenedione	5000	0.01
Estradiol	10000	ND
17- $\alpha$ -hydroxypregnenolone sulfate	5000	0.44

## B. Specificity

Compound	Added amount (ng/ml)	Cross-Reactivity (%)
17OH-Progesterone	-	100
Progesterone	1000	0.45
17- $\alpha$ -hydroxypregnenolone	1000	0.22
21-deoxycortisol	1000	0.96
Pregnenolone	1000	0.01
11-deoxycortisol	1000	0.13
Corticosterone	10000	0.00
11-deoxycorticosterone	1000	0.06
Testosterone	5000	0.01
Androstenedione	5000	0.01
Estradiol	10000	0.00
17- $\alpha$ -hydroxypregnenolone sulfate	5000	ND

Substance	17OH-P (ng/mL)	Concentration of Interferent (mg/dL)	Sample concentration after addition of the 17OH-P (ng/ml)	Mean % Variation
Hemoglobin	0.81	250	0.76	-10.1%
		500	0.74	
	1.83	250	1.69	
		500	1.5	
Bilirubin	0.81	50	0.78	-2.2%
		100	0.78	
	1.83	50	1.92	
		100	1.73	
Triglyceride	0.72	62.5	0.66	-4.6%
		125	0.66	
		250	0.67	
	1.56	62.5	1.65	
		125	1.55	
		250	1.41	

Substance	17OH-P (ng/mL)	Concentration of Interferent (mg/dL)	Sample concentration after addition of the 17OH-P (ng/ml)	Mean % Variation
Hemoglobin	0.55	250	0.60	5.3%
		500	0.60	
	1.24	250	1.25	
		500	1.26	
Bilirubin	0.55	50	0.56	4.5%
		100	0.58	
	1.24	50	1.23	
		100	1.32	
Triglyceride	0.55	62.5	0.56	2.7%
		125	0.54	
		250	0.55	
	1.24	62.5	1.19	
		125	1.23	
		250	1.24	

## C. Precision

INTRA-ASSAY			
Sample	N	$<X> \pm SD$ (ng/ml)	C.V. (%)
A	24	$1.16 \pm 0.05$	4.3
B	24	$2.93 \pm 0.12$	4.1

## C. Precision

INTRA-ASSAY			
Sample	N	$<X> \pm SD$ (ng/ml)	C.V. (%)
A	24	$0.87 \pm 0.04$	4.2
B	24	$1.64 \pm 0.08$	4.8

INTER-ASSAY			
Sample	N	$<X> \pm SD$ (ng/ml)	C.V. (%)
A	14	$0.99 \pm 0.06$	6.1
B	14	$2.32 \pm 0.14$	6.0

INTER-ASSAY			
Sample	N	$<X> \pm SD$ (ng/ml)	C.V. (%)
A	14	$0.35 \pm 0.03$	8.5
B	14	$1.35 \pm 0.11$	8.4

## D. Accuracy

RECOVERY TEST:			
Sample	Added 17-OHP (ng/ml)	Recovered 17-OHP (ng/ml)	Recovery (%)
1	10	9.61	96
	5	4.84	97
	1	0.94	94
	0.5	0.43	86
2	10	9.18	92
	5	4.5	90
	1	0.9	90
	0.5	0.5	100

## D. Accuracy

RECOVERY TEST			
Sample	Expected concentration (ng/ml)	Measured concentration (ng/ml)	Recovery (%)
1	-	0.59	-
	10.59	10.17	96.0
	5.59	4.53	81.0
2	1.59	1.85	116.1
	-	0.82	-
	10.82	12.86	118.9
	5.82	6.08	104.6
	1.82	2.03	111.8

DILUTION TEST:				DILUTION TEST			
Sample dilution	Theoretical concentration (ng/ml)	Measured concentration (ng/ml)	Recovery (%)	Sample dilution	Theoretical concentration (ng/ml)	Measured concentration (ng/ml)	Recovery (%)
1/1	12.01	12.01	100	1/1	-	8.98	-
1/2	6.01	5.51	91.7	1/2	4.49	4.19	93.2
1/4	3.00	2.80	93.2	1/4	2.25	2.17	96.7
1/8	1.50	1.58	105.2	1/8	1.12	1.22	108.4
1/1	10.73	10.73	100	1/16	0.56	0.65	115.0
1/2	5.37	4.77	88.9	1/1	-	8.16	-
1/4	2.68	2.68	99.9	1/2	4.08	3.28	80.3
1/8	1.34	1.43	106.6	1/4	2.04	1.65	80.6
				1/8	1.02	1.03	101.1
				1/16	0.51	0.60	118.4

en

Read entire protocol before use.

## 17OH-PROGESTERONE ELISA

### I. INTENDED USE

Immunoenzymetric assay for the *in vitro* quantitative measurement of human 17- $\alpha$ -hydroxyprogesterone (17-OHP) in serum and plasma.

### II. GENERAL INFORMATION

- A. Proprietary name: DIAsource 17OH-PROGESTERONE ELISA Kit
- B. Catalog number: KAP1401: 96 tests
- C. Manufactured by: DIAsource ImmunoAssays S.A.  
Rue du Bosquet, 2 B-1348 Louvain-la Neuve, Belgium.

For technical assistance or ordering information contact:  
Tel: +32 (0)10 84 99 11    Fax: +32 (0)10 84 99 90

### III. CLINICAL BACKGROUND

#### A. Biological activities of 17- $\alpha$ -hydroxyprogesterone

17- $\alpha$ -hydroxyprogesterone (17-OHP) is a C-21 steroid hormone (molecular weight 330.3) which is produced from 17- $\alpha$ -hydroxy pregnenolone in the adrenals and also in the ovaries, testes and placenta.

17-OHP is hydroxylated at positions 11 and 21 to produce cortisol via 11-deoxycortisol.

#### B. Clinical applications of 17- $\alpha$ -hydroxyprogesterone determination

As a rule, serum or amniotic fluid 17-OHP dosages are relevant to diagnose congenital adrenal hyperplasia (CAH). This CAH is due to a specific enzyme defect (six distinct enzyme deficiencies have been described). As a result of these deficiencies, ACTH increases and produces adrenal hyperplasia and the raise of many steroid precursors. But it is also very interesting to know the value of 17-OHP in patients with varicocele (17-OHP and testosterone represent markers of Leydig cell function) and in ageing male patients to detect Benign Prostatic Hypertrophy (BPH) and carcinoma of the prostate (PCA) (plasma17-OHP is significantly lower in PCA and BPH groups than in normal men).

There are other domains for 17-OHP investigations as: male infertility, girls with peripubertal virilization, children with premature adrenache (in these cases, the values of 17-OHP are increased without or after ACTH stimulation).

#### IV. PRINCIPLES OF THE METHOD

The DIAsource 17OH-PROGESTERONE ELISA kit is a solid phase Enzyme Linked Immunosorbent Assay performed on microtiterplates. During a 1 hour incubation step, at room temperature (18-25°C), with shaking, 17-OHP present in calibrators, controls and samples competes with 17-OHP-HRP conjugate for binding sites of a specific antibody immobilized on the wells of the microtiterplate. After this incubation, microtiterplate is washed to stop the competition reaction. The Chromogenic solution (TMB) is added and incubated for 30 minutes. The reaction is stopped with the addition of Stop Solution and the microtiterplate is then read at the appropriate wavelength. The amount of substrate turnover is determined colourimetrically by measuring the absorbance, which is inversely proportional to the 17-OHP concentration.

A calibration curve is plotted and the 17-OHP concentrations of the samples are determined by dose interpolation from the calibration curve.

#### V. REAGENTS PROVIDED

Reagents			KAP1401	Reconstitution
 Microtiterplate with 96 breakable wells coated with anti17- $\alpha$ -OH- progesterone			1 x 96 wells	Ready for use
CAL   0 Calibrator 0: human serum with gentamycin and Proclin			1 vial 1 ml	Ready for use
CAL   N 17- $\alpha$ -OH-progesterone N = 1 to 6 (see exact values on vial labels) in human serum with gentamycin and Proclin			6 vials 1 ml	Ready for use
CONTROL   N Controls N = 2 in human serum with gentamycin and Proclin			2 vials lyophilised	Add 1 ml distilled water
Ag   HRP   CONC Conjugate: HRP labelled 17- $\alpha$ -OH-progesterone			1 vial 0.3 ml	Dilute 100 x with conjugate buffer
CONJ   BUF Conjugate Buffer with bovine serum albumin			1 vial 30 ml	Ready for use
WASH   SOLN   CONC Wash solution (TRIS-HCl)			1 vial 10 ml	Dilute 200 x with distilled water (use a magnetic stirrer).
CHROM   TMB Chromogenic solution TMB (Tetramethylbenzidine)			1 vial 13 ml	Ready for use
STOP   SOLN Stop solution 0.2M H <sub>2</sub> SO <sub>4</sub>			1 vial 13 ml	Ready for use

**Note:** Use Calibrator 0 for dilution of samples with values above the highest calibrator.

#### VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Pipettes for delivery of: 25  $\mu$ l, 100  $\mu$ l, 200  $\mu$ l and 1000  $\mu$ l (the use of accurate pipettes with disposable plastic tips is recommended)
3. Vortex mixer

4. Magnetic stirrer
5. Plate shaker (300 to 700 rpm)
6. Washer for microtiterplates
7. Microtiterplate reader capable of reading at 450 nm and 650 nm (bichromatic reading)

#### Optional (for extraction of newborns serum samples)

1. Diethyl ether (analytical grade; purity > 98%)
2. Disposable glass tubes (12 x 75 m), with stoppers.
3. DIAsource Reconstitution Solution for 17-OHP -ELISA ref: 4214024

#### VII. REAGENT PREPARATION

- Controls: Reconstitute the controls with 1 ml distilled water.
- Working 17-OHP-HRP conjugate: Prepare an adequate volume of conjugate solution by adding for example: 20  $\mu$ l of the 100x concentrated 17-OHP -HRP conjugate to 2 ml of conjugate buffer. Use a vortex to homogenize.  
**Extemporaneous preparation is recommended.**
- Working Wash solution: Prepare an adequate volume of Working Wash solution by adding 199 volumes of distilled water to 1 volume of Wash Solution (200x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

#### VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, controls are stable for 8 weeks at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 4 months. Avoid subsequent freeze-thaw cycles.
- Freshly prepared Working Wash solution should be used on the same day.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.
- The DIAsource reconstitution solution must be kept at 2 to 8°C for 1 month after the first opening.

#### IX. SPECIMEN COLLECTION AND PREPARATION

- This kit is suitable for serum, heparinized plasma or EDTA plasma samples.
- Plasma samples provide similar results than serum:  

$$Y(\text{EDTA plasma}) = 0.93 X(\text{serum}) - 0.03 \quad r = 0.98 \quad n=20$$

$$Y(\text{heparinized plasma}) = 0.96 X(\text{serum}) - 0.01 \quad r = 0.99 \quad n=20$$
- Serum or plasma samples must be kept at 2-8°C.
- If the test is not run within 24 hrs, **sampling and storage at -20°C is recommended.**
- Avoid subsequent freeze-thaw cycles.

#### X. PROCEDURE

##### A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature (18-25°C) prior to use. Thoroughly mix all reagents and samples by gentle agitation or swirling. Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended. Use a clean plastic container to prepare the Wash Solution. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample. For the dispensing of the Chromogenic Solution and the Stop Solution avoid pipettes with metal parts. High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times.

**To avoid drift, the time between pipetting of the first calibrator and the last sample must be limited to the time mentioned in section XIII paragraph E (Time delay).**

Prepare a calibration curve for each run, do not use data from previous runs. Dispense the Chromogenic Solution within 15 minutes following the washing of the microtiterplate.

During incubation with Chromogenic Solution, avoid direct sunlight on the microtiterplate.

## B. Extraction for NEWBORN SERUM SAMPLES (optional)

1. Label one glass tube for each sample (do not extract calibrators or controls).
2. Pipette 100 µl of serum, followed by 1.5 ml of diethyl ether (analytical grade; purity > 98%).
3. Vortex all the tubes vigorously (2 x 1 minute).
4. Let stand the tubes for 15 minutes to separate well aqueous phase (lower phase) and organic phase (upper phase).
5. Then place the tubes at -20°C in order to freeze the aqueous phase.
6. Prepare a second series of glass tubes and, for each sample, transfer the organic phase (upper phase) into these new tubes. Avoid contamination by the aqueous phase.
7. Evaporate the organic phase (diethyl ether) completely under a stream of air by placing the tubes at 37°C (water bath). Manipulate under a hood.
8. Dissolve the dry organic extract with 100 µl of Reconstitution solution (not provided, see VI.). Vortex vigorously for 1 minute.
9. Let stand for 10-15 minutes and vortex again for 1 minute. These volumes allow to perform the determination in duplicate.

## C. Procedure

1. Select the required number of strips for the run. The unused strips should be resealed in the bag with a desiccant and stored at 2-8°C.
2. Secure the strips into the holding frame.
3. Pipette 25 µl of each Calibrator, Control and Sample into the appropriate wells.
4. Pipette 200 µl of the working 17-OHP-HRP conjugate solution into all the wells.
5. Incubate for 1 hour at room temperature (18-25°C), on a plate shaker (400 rpm).
6. Aspirate the liquid from each well.
7. Wash the plate 3 times by:
  - dispensing 0.4 ml of Wash Solution into each well
  - aspirating the content of each well
8. Pipette 100 µl of the Chromogenic solution into each well within 15 minutes following the washing step.
9. Incubate the microtiterplate for 30 minutes at room temperature (18-25°C), on a plate shaker (400 rpm), avoid direct sunlight.
10. Pipette 100 µl of Stop Solution into each well.
11. Read the absorbances at 450 nm (reference filter 630 nm or 650 nm) within 1 hour and calculate the results as described in section XI

## XI. CALCULATION OF RESULTS

1. Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).
  2. Calculate the mean of duplicate determinations.
  3. Calculate for each calibrator, control and sample:
- $$\text{B/B0 (\%)} = \frac{\text{OD (Calibrator, Control or Sample)}}{\text{OD (Zero Calibrator)}} \times 100$$
4. Plot the (B/B0(%)) values for each calibrator point as a function of the 17OH Progesterone concentration of each calibrator point. Reject obvious outliers.
  5. Computer assisted methods can also be used to construct the calibration curve. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.
  6. By interpolation of the sample (B/B0 (%)) values, determine the 17OH-Progesterone concentrations of the samples from the calibration curve

## XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

17OH-Progesterone ELISA		mOD units
Calibrator	0 ng/ml	2.306
	0.09 ng/ml	2.205
	0.4 ng/ml	1.744
	1.4 ng/ml	1.075
	2.7 ng/ml	0.737
	7.9 ng/ml	0.355
	15.6 ng/ml	0.235

### Conversion factor:

From ng/ml to nmol/L: x 3.03

From nmol/L to ng/ml: x 0.33

To the best of our knowledge, no international reference material exists for this parameter.

## XIII. PERFORMANCE AND LIMITATIONS

### A. Detection Limit

The Limit of Blank (LoB), Limit of Detection (LoD), and the Limit of Quantitation (LoQ), were determined in accordance with the CLSI guideline EP17-A.

The LOB (Limit of Blank) was calculated by measuring the blank several times and was calculated as the mean – 1.65 standard deviations of the distribution of the test values. The LOB was calculated to be 0.03 ng/ml.

The LOD (Limit of Detection) was calculated as the LOB – 1.65 standard deviations of a low concentration sample tested in 10 different runs. The LOD was calculated to be 0.06 ng/ml.

The LOQ (Limit of Quantification) was calculated by testing 7 samples of low value 10 times. The LOQ was calculated to be 0.1 ng/ml with CV of 20%.

### B. Specificity

The specificity was estimated by spiking a pool of 17-OHP samples ( $\pm 0.6$  ng/ml) with steroids which might be present in patient samples.

Compound	Added amount (ng/ml)	Cross-Reactivity (%)
17OH-Progesterone	-	100
Progesterone	1000	0.45
17- $\alpha$ -hydroxypregnenolone	1000	0.22
21-deoxycortisol	1000	0.96
Pregnenolone	1000	0.01
11-deoxycortisol	1000	0.13
Corticosterone	10000	0.00
11-deoxycorticosterone	1000	0.06
Testosterone	5000	0.01
Androstanedione	5000	0.01
Estradiol	10000	0.00
17- $\alpha$ -hydroxypregnenolone sulfate	5000	ND

ND: non detectable

The effect of potential interfering substances on samples using the DIAsource 17-OHP ELISA test was evaluated. Different levels of Hemoglobin, Bilirubin and Triglyceride in serum samples were tested on samples with different 17-OHP concentrations. Our acceptance criteria was to have interference of less than 15%. The tested substances did not affect the performance of the DIAsoure 17-OHP ELISA.

Substance	17OH-P (ng/mL)	Concentration of Interferent (mg/dL)	Sample concentration after addition of the 17OH-P (ng/ml)	Mean % Variation
Hemoglobin	0.55	250	0.60	5.3%
		500	0.60	
	1.24	250	1.25	
		500	1.26	
Bilirubin	0.55	50	0.56	4.5%
		100	0.58	
	1.24	50	1.23	
		100	1.32	
Triglyceride	0.55	62.5	0.56	2.7%
		125	0.54	
		250	0.55	
	1.24	62.5	1.19	
		125	1.23	
		250	1.24	

### Conversion factor:

From ng/ml to nmol/L: x 3.03

From nmol/L to ng/ml: x 0.33

## C. Precision

INTRA-ASSAY			
Sample	N	$\text{\bar{x}} \pm \text{SD}$ (ng/ml)	C.V. (%)
A	24	$0.87 \pm 0.04$	4.2
B	24	$1.64 \pm 0.08$	4.8

INTER-ASSAY			
Sample	N	$\text{\bar{x}} \pm \text{SD}$ (ng/ml)	C.V. (%)
A	14	$0.35 \pm 0.03$	8.5
B	14	$1.35 \pm 0.11$	8.4

SD: Standard Deviation; CV: Coefficient of variation

## D. Accuracy

RECOVERY TEST			
Sample	Expected concentration (ng/ml)	Measured concentration (ng/ml)	Recovery (%)
1	-	0.59	
	10.59	10.17	96.0
	5.59	4.53	81.0
	1.59	1.85	116.1
2	-	0.82	
	10.82	12.86	118.9
	5.82	6.08	104.6
	1.82	2.03	111.8

DILUTION TEST			
Sample dilution	Theoretical concentration (ng/ml)	Measured concentration (ng/ml)	Recovery (%)
1/1	-	8.98	-
1/2	4.49	4.19	93.2
1/4	2.25	2.17	96.7
1/8	1.12	1.22	108.4
1/16	0.56	0.65	115.0
1/1	-	8.16	-
1/2	4.08	3.28	80.3
1/4	2.04	1.65	80.6
1/8	1.02	1.03	101.1
1/16	0.51	0.60	118.4

Samples were diluted with the zero calibrator.

## E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when Working 17-OHP-HRP Conjugate is dispensed 10 and 15 minutes after the calibrator has been added in the coated wells.

TIME DELAY			
	0' (ng/ml)	10' (ng/ml)	15' (ng/ml)
Sample 1	0.68	0.71	0.74
Sample 2	2.81	3.30	3.14

## F. Limitations for newborn samples

Very young children (< 3 months old) present high levels of 17OH pregnenolone sulphate. Despite the fact that there is almost no cross reaction with this molecule, we recommend to apply the extraction procedure (X. B.)

## XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls which contain azide will interfere with the enzymatic reaction and cannot be used.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises.
- It is recommended that Controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.
- It is good practise to check visually the curve fit selected by the computer.

## XV. EXPECTED VALUES

These values are given only for guidance; each laboratory should establish its own normal range of values.

	Concentration range (2.5 to 97.5 percentiles) (ng/ml)	Number of subjects
Normal males	0.37 – 2.38	38
Normal Females		
Follicular phase	0.32 – 0.82	28
Luteal phase	0.79 - 3.29	34
Pregnancy		
First and second trimester	1.3 – 3.9	25
Third trimester	2.2 – 12.4	14

Children:	Nb of subjects:	KAP1401 without extraction Concentration range (2.5 to 97.5 percentiles) ng/ml	KAP1401 with extraction Concentration range (2.5 to 97.5 percentiles) ng/ml
Newborns:			
0 to 1 month	20	1.41 – 16.81	0.36 - 2.5
1 to 2 months	44	1.44 – 8.27	0.37 - 1.68
2 to 3 months	27	1.11 – 6.01	0.25 - 1.83
3 to 4 months	4	1.58 – 6.39	0.56 - 2.14
1 to 2 years	20	0.11 – 0.50	NA
2 to 10 years	40	0.11 – 1.04	NA
11 to 14 years	12	0.40 – 2.00	NA
14 to 18 years	11	0.69 – 2.95	NA

## XVI. PRECAUTIONS AND WARNINGS

### Safety

For *in vitro* diagnostic use only.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with all reagents, Stop Solution contains H<sub>2</sub>SO<sub>4</sub>. In case of contact, wash thoroughly with water.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

## XVII. BIBLIOGRAPHY

- ANDO S. et al (1983) Plasma levels of 17-OH-progesterone and testosterone in patients with varicoceles. Acta Endocrinologica, 102,463-469.

2. BONACCORSI A. et al (1987)  
**Male infertility due to congenital adrenal hyperplasia: testicular biopsy findings, hormonal evaluation, and therapeutic results in three patients.**  
Fertility and Sterility, 47, 4, 664-670.

3. CACCIARE et al (1987)  
**Neonatal screening program for congenital adrenal hyperplasia in a homogeneous Caucasian population.**

Annals New York Academy of sciences 85-88.

4. DRAFTA D. et al (1982)  
**Plasma steroids in benign prostatic hypertrophy and carcinoma of the prostate.**  
J. Steroid Biochem. 17, 689-693.

5. GRANOFF A. et al (1985)  
**17-hydroxyprogesterone responses to adrenocorticotropin in children with premature adrenarche.**  
Journal of Clinical Endocrinology and Metabolism, 3, 409-415.

6. NEW M. et al (1983)  
**Genotyping steroid 21-hydroxylase deficiency hormonal reference data.**  
Journal of Clinical Endocrinology and Metabolism 57, 320-326.

7. SOLYOM J. et al (1987)  
**Detection of late-onset adrenal hyperplasia in girls with peripubertal virilization.**  
Acta Endocrinologica (Copenh.) 115, 413-418.

#### XVIII. SUMMARY OF THE PROTOCOL

	CALIBRATORS ( $\mu$ l)	SAMPLE(S) CONTROLS ( $\mu$ l)
Calibrators (0-6)	25	-
Controls, Samples or Extracted Samples	-	25
Working 17-OHP HRP Conjugate	200	200
Incubate for 1 hour at room temperature (18-25°C) with continuous shaking at 400 rpm. Aspirate the contents of each well. Wash 3 times with 400 $\mu$ l of Wash Solution and aspirate.		
Chromogenic Solution	100	100
Incubate for 30 min at room temperature (18-25°C) with continuous shaking at 400 rpm		
Stop Solution	100	100
Read on a microtiterplate reader. Record the absorbance of each well at 450 nm (versus 630 or 650 nm).		

DIAsource's Instrumentation Service can offer you the possibility of acquiring a protocol adapted for this kit to be used on the Stratec Gemini 2PS + Combo platform including protocol assay file, reagent file and tips for the good use of the kit on the instrument.

If you need any additional information please contact [Instrumentation@diасource.be](mailto:Instrumentation@diасource.be)

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diасource-diagnostics.com/>

DIAsource Catalogue Nr: KAP1401	Revision nr: 220516
------------------------------------	------------------------

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1. ANDO S. et al (1983)

**Plasma levels of 17-OH-progesterone and testosterone in patients with varicoceles.**

Acta Endocrinologica, 102, 463-469.

2. BONACCORSI A. et al (1987)

**Male infertility due to congenital adrenal hyperplasia: testicular biopsy findings, hormonal evaluation, and therapeutic results in three patients.**

Fertility and Sterility, 47, 4, 664-670.

3. CACCIARE et al (1987)

**Neonatal screening program for congenital adrenal hyperplasia in a homogeneous Caucasian population.**

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4. DRAFTA D. et al (1982)

**Plasma steroids in benign prostatic hypertrophy and carcinoma of the prostate.**

J. Steroid Biochem. 17, 689-693.

5. GRANOFF A. et al (1985)

**17-hydroxyprogesterone responses to adrenocorticotropin in children with premature adrenarche.**

Journal of Clinical Endocrinology and Metabolism, 3, 409-415.

6. NEW M. et al (1983)

**Genotyping steroid 21-hydroxylase deficiency hormonal reference data.**

Journal of Clinical Endocrinology and Metabolism 57, 320-326.

7. SOLYOM J. et al (1987)

**Detection of late-onset adrenal hyperplasia in girls with peripubertal virilization.**

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Other translations of this Instructions for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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## XVII. BIBLIOGRAFIA

1 ANDO S. et al (1983)

**Plasma levels of 17OHprogesterone and testosterone in patients with varicoceles.**

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Annals New York Academy of sciences 8588.

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J. Steroid Biochem. 17/689693

5 ANDO S. et al (1985)

**17hydroxyprogesterone responses to adrenocorticotropin in children with premature adrenarche.**

Journal of Clinical Endocrinology and Metabolism, 3, 409415.

6 ANDO S. et al (1983)

**Genotyping steroid 21hydroxylase deficiency hormonal reference data.**

Journal of Clinical Endocrinology and Metabolism, 57, 320326.

7 SOLYOM J. et al (1987)

**Detection of lateonset adrenal hyperplasia in girls with peripubertal virilization.**

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## XVII. BIBLIOGRAPHY

1. ANDO S. et al (1983)  
**Plasma levels of 17-OH-progesterone and testosterone in patients with varicoceles.**  
Acta Endocrinologica, 102, 463-469.
2. BONACCORSI A. et al (1987)  
**Male infertility due to congenital adrenal hyperplasia: testicular biopsy findings, hormonal evaluation, and therapeutic results in three patients.**  
Fertility and Sterility, 47, 4, 664-670.
3. CACCIARE et al (1987)  
**Neonatal screening program for congenital adrenal hyperplasia in a homogeneous Caucasian population.**  
Annals New York Academy of sciences 85-88.
4. DRAFTA D. et al (1982)  
**Plasma steroids in benign prostatic hypertrophy and carcinoma of the prostate.**  
J. Steroid Biochem. 17, 689-693.
5. GRANOFF A. et al (1985)  
**17-hydroxyprogesterone responses to adrenocorticotropin in children with premature adrenarche.**  
Journal of Clinical Endocrinology and Metabolism, 3, 409-415.
6. NEW M. et al (1983)  
**Genotyping steroid 21-hydroxylase deficiency hormonal reference data.**  
Journal of Clinical Endocrinology and Metabolism 57, 320-326.
7. SOLYOM J. et al (1987)  
**Detection of late-onset adrenal hyperplasia in girls with peripubertal virilization.**  
Acta Endocrinologica (Copenh.) 115, 413-418.

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