CA 125-IRMA

KIP0301

Version : 230123

Date of issue : 23/01/2023

Revision date: 23/01/2023

History

Summary of change:

Current Version:
230123
New logo

CE

en

Read entire protocol before use.

CA 125-IRMA

I. INTENDED USE

Immunoradiometric assay kit for the in vitro quantitative determination of the cancer associated antigen CA 125 in serum.

II. GENERAL INFORMATION

- A. Proprietary name : DIAsource CA 125-IRMA Kit
- B. Catalog number : KIP0301 : 96 tests
- C. Manufactured by : DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

A. Biological activities

CA 125 is a high molecular weight mucin type glycoprotein, originally defined by the Oc 125 monoclonal antibody (Mab) established by Bast et al. Different epitopes, co-expressed with the Oc 125 epitope on the CA 125 antigen, have been used for the development of heterologous assays for determination of the CA 125 antigen.

B. Clinical applications

Approximately 90% of ovarian cancers are celomic epithelial carcinomas and contain a celomic epithelium-related glycoprotein designated CA 125. CA 125 can be localized in most serous, endometrioid, and clear cell ovarian carcinomas; mucinous tumors express this antigen less frequently. CA 125 is also found in the epithelium of the fallopian tubes, endometrium, and uterine cervix.

Elevated CA 125 levels can result from abdominal diseases other than ovarian cancer. Although CA 125 is most consistently elevated in patients with epithelial ovarian cancer, it can be expressed in a number of gynaecologic (eg, endometrium, fallopian tube) and non-gynaecologic (eg, pancreas, breast, colon, lung) cancers. CA 125 levels are frequently elevated with tumor spread beyond the uterus.

IV. PRINCIPLES OF THE METHOD

The DIAsource CA 125-IRMA is a two-step immunoradiometric assay based on coated-tube separation. Mab1, the capture antibody, is attached to the lower and inner surface of the plastic tube. Add calibrators or samples to the tubes. After incubation, washing removes the occasional excess of antigen. Addition of Mab2, the signal antibody labelled with ^{125}I , will complete the system and trigger the immunological reaction. After washing, the remaining radioactivity bound to the tube reflects the antigen concentration.

V. REAGENTS PROVIDED

Reagents	96 Test Kit	Colour Code	Reconstitution
Tubes coated with anti CA 125(monoclonal antibody).	2 x 48	orange	Ready for use
Ab ^{125}I	1 vial 10,5 ml 760 kBq	red	Ready for use
Anti-CA 125- ^{125}I (monoclonal antibody) in phosphate buffer with bovine serum albumin, azide (<0.1%), bovine serum and inert red dye.			
CAL 0	1 vial lyophilised	yellow	Add 3.0 ml distilled water
Calibrators 0 in bovine serum with thymol .			
CAL N	5 vials lyophilised	yellow	Add 1.0 ml distilled water
Calibrators 1-5 in bovine serum with thymol. See exact value on vial labels.			
INC BUF	2 vials 3 ml	black	Ready for use
Incubation buffer: Tris Maleate buffer with bovine serum albumin, bovine serum and azide.			
WASH SOLN CONC	2 vials 10 ml	brown	Dilute 70 x with distilled water (use a magnetic stirrer).
Wash solution (TRIS-HCl).			
CONTROL N	2 vials lyophilised	silver	Add 1.0 ml distilled water
Controls - N = 1 or 2 in human plasma and thymol.			

Note: use the zero calibrator for sera dilutions.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Pipettes for delivery of: 50 μl , 100 μl , 1 ml and 3 ml (the use of accurate pipettes with disposable plastic tips is recommended)
3. Pipette for delivery of 1 to 5 ml of distilled water
4. Vortex mixer
5. Magnetic stirrer
6. Tube shaker (400rpm)
7. 5 ml automatic syringe (Cornwall type) for washing
8. Aspiration system (optional)
9. Any gamma counter capable of measuring ^{125}I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- A. **Calibrators:** Reconstitute the calibrator 0 with 3 ml distilled water and the calibrators 1-5 with 1.0 ml.
- B. **Controls:** Reconstitute the controls with 1.0 ml distilled water.

C. **Working Wash solution:** Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- The calibrators and controls are very unstable, use them immediately after reconstitution, freeze immediately in aliquots and keep them at -20°C for maximum 3 months. Avoid subsequent freeze-thaw cycles.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum samples must be kept at 2-8°C.
- If the test is not run within 24 hrs, storage at -20°C is recommended.
- Avoid subsequent freeze-thaw cycles.

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date.
Do not mix materials from different kit lots.
Bring all the reagents to room temperature prior to use.
Thoroughly mix all reagents and samples by gentle agitation or swirling.
Use a clean disposable pipette tip for addition of each different reagent and sample in order to avoid cross-contamination. High precision pipettes or automated pipetting equipment will improve the precision.
Respect the incubation times.
Prepare a calibration curve for each run, do not use data from previous runs.

B. Procedure.

1. Label coated tubes in duplicate for each calibrator, sample and control. For the determination of total counts, label 2 normal tubes.
2. Briefly vortex calibrators, samples and controls and dispense 100 μl of each into the respective tubes.
3. Dispense 50 μl Incubation buffer in each tube except those for total counts.
4. Incubate for 120 minutes at room temperature on a tube shaker (400 rpm).
5. Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
6. Wash tubes with 2 ml Working Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
7. Wash tubes again with 2 ml Wash solution (except total counts) and aspirate (or decant).
8. Let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
9. Dispense 100 μl of ^{125}I odine labelled anti CA 125 into each tube, including the uncoated tubes for total counts.
10. Incubate for 60 minutes at room temperature on a tube shaker(400 rpm).
11. Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
12. Wash tubes with 2 ml Working Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
13. Wash tubes again with 2 ml Wash solution (except total counts) and aspirate (or decant).
14. Let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
15. Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. On semi logarithmic or linear graph paper plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of CA 125 (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
3. Read the concentration for each control and sample by interpolation on the calibration curve.

4. Computer assisted data reduction will simplify these calculations. If automatic result processing is to be used, a 4-parameter logistic function curve fitting is recommended

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

CA 125-IRMA		cpm	B/T (%)
Total count		299443	100
Calibrator	0.0 U/ml	349	0.12
	6.4 U/ml	1261	0.42
	16.2 U/ml	3752	1.25
	49.0 U/ml	12488	4.17
	157.0 U/ml	39441	13.17
	568.0 U/ml	124992	41.74

Since no international reference material is available for CA 125 antigen, DIAsource CA 125-IRMA calibrator values are assigned against a set of in-house reference standards.

XIII. PERFORMANCE AND LIMITATIONS

A. Detection limit

Twelve zero calibrators were assayed along with a set of the other calibrators. The detection limit, defined as the apparent concentration of the average count at zero binding plus two standard deviations, was 0.7 U/ml.

B. Precision

INTRA-ASSAY PRECISION

INTER-ASSAY PRECISION

Serum	N	$\text{\langle X \rangle \pm SD (U/ml)}$	CV (%)	Serum	N	$\text{\langle X \rangle \pm SD (U/ml)}$	CV (%)
A	20	31.5 ± 1.2	3.9	A	20	38.0 ± 2.0	5.3
B	20	120.8 ± 2.8	2.3	B	20	114.5 ± 5.7	5.0

SD: Standard Deviation; CV: Coefficient of variation

C. Accuracy

DILUTION TEST

Sample	Dilution	Theoretical Concent. (U/ml)	Measured Concent. (U/ml)
A	1/1	-	618.8
	1/2	309.4	310.7
	1/4	154.7	161.4
	1/8	77.4	80.9
	1/16	38.7	43.2
	1/32	19.3	21.6
	1/64	9.7	10.3

Samples were diluted with calibrator 0.

RECOVERY TEST

Added CA 125 (U/ml)	Recovered CA 125 (U/ml)	Recovery (%)
31.1	35.5	114 %
65.8	67.7	103 %
124.3	135.7	109 %
257.9	268.9	104 %
531.8	552.4	104 %

D. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrator has been added to the coated tubes.

TIME DELAY				
	0'	10'	20'	30'
S 1 (U/ml)	35.9	37.2	39.3	38.5
S 2 (U/ml)	104.5	113.5	117.7	116.6

E. Hook-effect

Samples containing 340.000 U/ml CA 125 give a result higher than the last calibration point.

F. Specificity

The DIAsource CA 125 IRMA is based on two mouse monoclonal antibodies, OvK95 and Ov185, directed against two independent epitopes of the protein core of the CA 125 antigen.

XIV. LIMITATIONS

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays. Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed in case of the presence of heterophilic antibodies. Carefully evaluate the results of patients suspected of having these antibodies. If results are not consistent with other clinical observations, additional information should be required before diagnosis.

XV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises.

XVI. REFERENCE INTERVALS

Among 154 apparently healthy individuals, 98% of the results were below 35.0 U/ml.

For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history and other findings.

XVII. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

This kit contains ^{125}I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations. This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where bovine spongiform encephalopathy has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious. Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up. Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVIII. BIBLIOGRAPHY

1. BAST R., FEENEY M., LAZARUS H., NADLER L., COLVIN R., KNAPP R.
Reactivity of a monoclonal antibody with human ovarian carcinoma.
J Clin Invest 1981;68:1331-1337.
2. O'BRIEN T., RAYMOND L., BANNON G., FORD D., HARDARDOTTIR H., MILLER F., QUIRK J..
New monoclonal antibodies identify the glycoprotein carrying the CA 125 epitope.
Am J Obstet Gynecol 1991;165:1857-1864.
3. BAST RC Jr., XU FJ., YU YH., BARNHILL S., ZHANG Z., MILSS GB.
CA 125 : The past and the future.
Int J Biol Markers 1998;13:179-187.
4. BONFRER JMG., DUFFY MJ., RADTKE M., SEGURADO O., TORRE GC., VAN DALEN A., et al.
Tumor markers in gynecological cancers : EGTM recommendations.
Anticancer Res 1999;19:2807-2810.
5. NIH CONSENSUS DEVELOPMENTAL PANEL ON OVARIAN CANCER :
NIH CONSENSUS CONFERENCE.
Ovarian cancer : screening, treatment and follow-up.
JAMA 1995;273:491-497.
6. PETIGNAT P., JORIS F., OBRIST R.
How CA 125 is used in routine clinical practice.
Eur J Cancer 2000;36:1933-1937.
7. VON SCHLIPPE M., RUSTIN GJ.
Circulating tumor markers in ovarian tumors.
Forum (Genova) 2000;10:383-392.
8. COLAKOVIC S., LUKIC V., MITROVIC L., JELIC S., SUSNjar S., MARINKOVIC J.
Prognostic value of CA 125 kinetics and half-life in advanced ovarian cancer.
Int J Biol Markers 2000;15:147-152.

XIX. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS µl	CALIBRATORS µl	SAMPLE(S) CONTROLS µl
Calibrators (0 to 5)	-	100	-
Sample(s), controls	-	-	100
Incubation buffer	-	50	50
Incubation	120 minutes at room temperature with shaking at 400 rpm		
Separation	-	Aspirate (or decant)	
Working Wash solution		2 ml	
Separation		Aspirate (or decant)	
Working Wash solution		2 ml	
Separation		Aspirate (or decant)	
Tracer	100	100	100
Incubation	60 minutes at room temperature with shaking at 400 rpm		
Separation	-	Aspirate (or decant)	
Working Wash solution		2 ml	
Separation		Aspirate (or decant)	
Working Wash solution		2 ml	
Separation		Aspirate (or decant)	
Counting	Count tubes for 60 seconds		

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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XVIII. BIBLIOGRAPHIE

1. BAST R., FEENEY M., LAZARUS H., NADLER L., COLVIN R., KNAPP R.
Reactivity of a monoclonal antibody with human ovarian carcinoma.
J Clin Invest 1981;68:1331-1337.
2. O'BRIEN T., RAYMOND L., BANNON G., FORD D., HARDARDOTTIR H., MILLER F., QUIRK J.
New monoclonal antibodies identify the glycoprotein carrying the CA 125 epitope.
Am J Obstet Gynecol 1991;165:1857-1864.
3. BAST RC Jr., XU FJ., YU YH., BARNHILL S., ZHANG Z., MILSS GB.
CA 125 : The past and the future.
Int J Biol Markers 1998;13:179-187.
4. BONFRER JMG., DUFFY MJ., RADTKE M., SEGURADO O., TORRE GC., VAN DALEN A., et al.
Tumor markers in gynaecological cancers : EGTM recommendations.
Anticancer Res 1999;19:2807-2810.
5. NIH CONSENSUS DEVELOPMENTAL PANEL ON OVARIAN CANCER : NIH CONSENSUS CONFERENCE.
Ovarian cancer : screening, treatment and follow-up.
JAMA 1995;273:491-497.
6. PETIGNAT P., JORIS F., OBRIST R.
How CA 125 is used in routine clinical practice.
Eur J Cancer 2000;36:1933-1937.
7. VON SCHLIPPE M., RUSTIN GJ.
Circulating tumor markers in ovarian tumors.
Forum (Genova) 2000;10:383-392.

8. COLAKOVIC S., LUKIC V., MITROVIC L., JELIC S., SUSNjar S., MARINKOVIC J.

Prognostic value of CA 125 kinetics and half-life in advanced ovarian cancer.

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XVIII. BIBLIOGRAFIA

1. BAST R., FEENEY M., LAZARUS H., NADLER L., COLVIN R., KNAPP R.

Reactivity of a monoclonal antibody with human ovarian carcinoma.

J Clin Invest 1981;68:1331-1337.

2. O'BRIEN T., RAYMOND L., BANNON G., FORD D., HARDARDOTTIR H., MILLER F., QUIRK J..

New monoclonal antibodies identify the glycoprotein carrying the CA 125 epitope.

Am J Obstet Gynecol 1991;165:1857-1864.

3. BAST RC Jr., XU FJ., YU YH., BARNHILL S., ZHANG Z., MILSS GB.

CA 125 : The past and the future.

Int J Biol Markers 1998;13:179-187.

4. BONFRER JMG., DUFFY MJ., RADTKE M., SEGURADO O., TORRE GC., VAN DALEN A., et al.

Tumor markers in gynecological cancers : EGT recommendations.

Anticancer Res 1999;19:2807-2810.

5. NIH CONSENSUS DEVELOPMENTAL PANEL ON OVARIAN CANCER : NIH CONSENSUS CONFERENCE.

Ovarian cancer : screening, treatment and follow-up.

JAMA 1995;273:491-497.

6. PETIGNAT P., JORIS F., OBRIST R.

How CA 125 is used in routine clinical practice.

Eur J Cancer 2000;36:1933-1937.

7. VON SCHLIPPE M., RUSTIN GJ.

Circulating tumor markers in ovarian tumors.

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