# DHEA-S-RIA-CT

*KIP0481*

---

Version : 230123

Date of issue : 23/01/2023

Revision date: 23/01/2023

# History

---

## Summary of change:

<b>Current Version:</b>
<b>230123</b>
New logo

en

Read entire protocol before use.

## DHEA-S-RIA-CT

### I. INTENDED USE

Radioimmunoassay for the *in vitro* quantitative measurement of human 5-Androsten-3 $\beta$ -OL-17-Dione Sulphate (Dehydroepiandrosterone Sulphate or DHEA Sulphate) in serum.

### II. GENERAL INFORMATION

- A. Proprietary name : DIAsource DHEA-S-RIA-CT Kit
- B. Catalog number : KIP0481 : 96 tests
- C. Manufactured by : DIAsource ImmunoAssays S.A.  
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :  
Tel : +32 (0)10 84.99.11                  Fax : +32 (0)10 84.99.91

### III. CLINICAL BACKGROUND

#### A. Biological Activity

Dehydroepiandrosterone (DHEA) and its sulphated form, DHEA-S, are the most abundant steroid products of the adrenal glands. The adrenal is the sole source of these steroids in women, but in men there is also a small contribution from the testes. Formation of DHEA is stimulated by ACTH producing a significant diurnal variation in serum concentrations. DHEA is rapidly converted to DHEA-S by an enzyme present in the adrenals, liver and small intestine. DHEA-S is present at concentrations greater than 200 times that of DHEA and has a longer half-life, which largely removes the diurnal variation.

The aging influences the endocrine temporal structure, including DHEA-S which can be considered as a biomarker of aging, since its levels gradually decrease in older subjects. Plasma levels of DHEA-S increase steadily from about the seventh year of life, than gradually decline after the third decade.

DHEA is a 19-carbon steroid with a molecular weight of 288 Kd and a half-life in plasma of about 1-3 hours. The molecular weight of DHEA-S is 368.5 Kd and the half-life is about 10-20 hours. DHEA is formed from pregnenolone by the enzyme 17,20 desmolase and metabolised to androstenedione or testosterone by 3 beta- or 17 beta-hydrosteroid dehydrogenase repectively. Hydrosteroid sulphatase converts DHEA to DHEA-S and sulphohydrolase reverses this reaction.

#### B. Clinical Applications

Raised levels of DHEA-S are found in the plasma of patients with adrenal tumors or with congenital adrenal hyperplasia. DHEA-S may also be slightly elevated in patients with polycystic ovaries, supporting an adrenal component to the virilization seen in this condition. HCG-production tumors in men may lead to increased testicular DHEA production.

DHEA-S is often assayed in conjunction with free testosterone as an initial screen for hyperandrogenism in hirsutism. Sometimes DHEA-S is the only hormone circulating at a level above normal, and is apparently more likely to be elevated during the early stages of hirsutism than most other androgens.

DHEA-S is usually undetectable with adrenal insufficiency or panhypopituitarism. Concentrations are slightly decreased in pregnancy and with oral contraceptive use and markedly decreased following glucocorticoid administration.

Low circulating concentrations are seen with severe illness and in patients with AIDS.

#### IV. PRINCIPLES OF THE METHOD

A fixed amount of  $^{125}\text{I}$  labelled DHEA-S competes with the DHEA-S to be measured present in the sample or in the calibrator for a fixed amount of antibody sites being immobilized to the wall of a polystyrene tube. After 1-hour incubation at 37°C in water bath, an aspiration step terminates the competition reaction. The tubes are then washed twice with 2 ml of Working Wash Solution and aspirated. A calibration curve is plotted and the DHEA-S concentrations of the samples are determined by dose interpolation from the calibration curve.

#### V. REAGENTS PROVIDED

Reagents	96 Test Kit	Colour Code	Reconstitution
Tubes coated with anti DHEA-S	2 x 48	ivory	<b>Ready for use</b>
Ag $^{125}\text{I}$	1 vial 55 ml 450 kBq	red	<b>Ready for use</b>
TRACER: $^{125}\text{I}$ odine labelled DHEA-S (HPLC grade) in buffer with bovine serum albumin and azide (<0.1%)			
CAL      0	1 vial lyophilised	yellow	<b>Add 0.5 ml distilled water</b>
Zero calibrator in human serum, thymol			
CAL      N	5 vials lyophilised	yellow	<b>Add 0.5 ml distilled water</b>
Calibrators DHEA-S - N = 1 to 5 (see exact values on vial labels) in human plasma, thymol			
WASH      SOLN      CONC	1 vial 10 ml	brown	<b>Dilute 70 x with distilled water (use a magnetic stirrer).</b>
Wash solution (TRIS-HCl)			
CONTROL      N	2 vials Lyophilised	silver	<b>Add 0.5 ml distilled water</b>
Controls - N = 1 or 2 in human plasma, thymol (see exact values on vial labels)			

Note : Use the zero calibrator for sample dilutions.

#### VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Pipettes for delivery of: 10  $\mu\text{l}$  and 500  $\mu\text{l}$  (the use of accurate pipettes with disposable plastic tips is recommended)
3. Vortex mixer
4. Magnetic stirrer
5. Water bath at 37°C
6. 5 ml automatic syringe (Cornwall type) for washing
7. Aspiration system (optional)
8. Any gamma counter capable of measuring  $^{125}\text{I}$  may be used (minimal yield 70%).

#### VII. REAGENT PREPARATION

- A. **Calibrators** : Reconstitute the calibrators with 0.5 ml distilled water
- B. **Controls**: Reconstitute the controls with 0.5 ml distilled water.
- C. **Working Wash solution**: Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

#### VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, calibrators and controls are stable for one week at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 3 months. Avoid successive freezing and thawing.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, the tracer is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.

- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

#### IX. SPECIMEN COLLECTION AND PREPARATION

- Serum samples must be kept at 2-8°C.
- If the test is not run within 24 hours, storage in aliquots at -20°C is recommended.
- Avoid successive freezing and thawing.

#### X. PROCEDURE

##### A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature (18-25°C) prior to use. Thoroughly mix all reagents and samples by gentle agitation or swirling. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample. High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times. Prepare a calibration curve for each run, do not use data from previous runs. Each tube can only be used once.

##### B. Procedure

1. Label coated tubes in duplicate for each calibrator, control and sample. For the determination of total counts, label 2 normal tubes.
2. Briefly vortex calibrators, controls and samples and dispense 10  $\mu\text{l}$  of each into respective tubes.
3. Dispense 500  $\mu\text{l}$  of  $^{125}\text{I}$ odine labelled DHEA-S into each tube, including the uncoated tubes for total counts.
4. Shake the tube rack gently by hand to liberate any trapped air bubbles.
5. Incubate for 1 hour at 37°C in water bath.
6. Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
7. Wash tubes with 2 ml Working Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
8. Wash tubes with 2 ml Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
9. Let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
10. Count tubes in a gamma counter for 60 seconds.

#### XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. Calculate the bound radioactivity as a percentage of the binding determined at the zero calibrator point (0) according to the following formula:

$$\text{B/B}_0 (\%) = \frac{\text{Counts (Calibrator or sample)}}{\text{Counts (Zero Calibrator)}} \times 100$$

3. Plot the (B/B<sub>0</sub>(%)) values for each calibrator point as a function of the DHEA-S concentration of each calibrator point. Reject obvious outliers.
4. Computer assisted methods can also be used to construct the calibration curve. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.
5. By interpolation of the sample (B/B<sub>0</sub> (%)) values, determine the DHEA-S concentrations of the samples from the calibration curve.
6. For each assay, the percentage of total tracer bound in the absence of unlabelled DHEA-S (B<sub>0</sub>/T) must be checked.

#### XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

DHEA-S-RIA-CT	cpm	B/B <sub>0</sub> (%)
Total count	164710	
Calibrator	0.0 $\mu\text{g}/\text{dl}$	100.0
	1.7 $\mu\text{g}/\text{dl}$	91.6
	9.7 $\mu\text{g}/\text{dl}$	62.5
	41.0 $\mu\text{g}/\text{dl}$	32.3
	219.0 $\mu\text{g}/\text{dl}$	9.9
	1090.0 $\mu\text{g}/\text{dl}$	2.6
	877	

**Conversion factor:**

From µg/dl to µmol/l : x 0.0271

From µmol/l to µg/dl : x 36.85

**XIII. PERFORMANCE AND LIMITATIONS****A. Detection limit**

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations below the average counts at zero binding, was 1.23 µg/dl.

**B. Specificity**

The specificity of the polyclonal antibody used in this assay was evaluated by RIA, as the ratio of the amount of DHEA-S, which yields 50% inhibition of binding of labelled DHEA-S and the amount of cross-reacting compounds (analog) giving the same inhibition.

Compound	Cross-Reactivity (%)
DHEA-S	100.00
DHEA	< 0.05
DHEA-Glucuronide	< 0.05
Androstenediol	< 0.05
Estrone-Sulfate	0.08
5α-Dihydrotestosterone	< 0.05
Testosterone	< 0.05

Note : This table shows the cross-reactivity for the anti DHEA-S.

**C. Precision****INTRA-ASSAY PRECISION****INTER-ASSAY PRECISION**

Serum	N	$\text{\bar{X}} \pm \text{SD}$ (µg/dl)	CV (%)	Serum	N	$\text{\bar{X}} \pm \text{SD}$ (µg/dl)	CV (%)
A	20	$32.1 \pm 1.9$	5.9	A	12	$30.8 \pm 2.0$	6.5
B	20	$107.4 \pm 5.1$	3.6	B	12	$105.5 \pm 6.5$	6.5

SD: Standard Deviation; CV: Coefficient of variation

**D. Accuracy****DILUTION TEST**

Sample	Dilution	Theoretical Concent. (µg/dl)	Measured Concent. (µg/dl)
Serum 1	1/1	389.6	168.1
	1/2	194.8	
	1/4	97.4	
	1/8	48.7	
	1/16	24.3	
Serum 2	1/1	85.0	38.3
	1/2	42.5	
	1/4	21.2	
	1/8	10.6	
	1/16	5.3	

Samples were diluted with the zero calibrator.

**RECOVERY TEST**

Sample (µg/dl)	added DHEA-S (µg/dl)	Theoretical DHEA-S (µg/dl)	Measured DHEA-S (µg/dl)	Recovered %
Serum	0	NA	26.0	NA
	25	51	59.7	117
	50	76	80.6	106
	100	126	128.1	101.6

To the best of our knowledge, no international reference material exists for this parameter.

**E. Time delay between last calibrator and sample dispensing**

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrator has been added to coated tubes.

**TIME DELAY**

Serum (µg/dl)	0'	10'	20'	30'
Serum 1	30.8	33.8	31.5	31.2
Serum 2	98.5	113.5	118.2	102.5

**XIV. INTERNAL QUALITY CONTROL**

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises.

**XV. REFERENCE INTERVALS**

These values are given only for guidance; each laboratory should establish its own normal range of values.

The ranges are expressed as 2.5% to 97.5% percentiles.

Subjects (age group)	N	Mean (µg/dl)	Range *
Women	15	129.7	37.8 - 264.5
	19	276.7	117.3 - 531.0
	16	244.0	118.7 - 428.9
	33	197.5	70.2 - 413.4
	29	175.1	75.2 - 398.7
	23	162.0	63.3 - 291.4
	19	144.4	37.3 - 255.8
	13	109.8	20.9 - 198.0
	16	117.0	26.8 - 375.4
	13	102.9	34.5 - 178.1
Men	10	92.8	24.6 - 182.3
	15	114.3	6.3 - 377.0
	14	138.9	34.9 - 284.1
	42	262.2	80.5 - 479.4
	37	177.6	51.8 - 470.7

**XVI. PRECAUTIONS AND WARNINGS****Safety**

For *in vitro* diagnostic use only.

This kit contains  $^{125}\text{I}$  (half-life: 60 days), emitting ionizing X (28 keV) and  $\gamma$  (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections.

Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

For more information, see Material Safety Data Sheet (MSDS)

#### XVII. BIBLIOGRAPHY

1. Orentreich N, Brind J, Rizer R, Vogelman J : **Age changes and sex differences in serum dehydroepiandrosterone sulphate concentrations throughout adulthood.** J Clin Endocrinol Metab 59:551-555, 1984.
2. Orentreich N, Brind J, Vogelman J, Andres R, Baldwin H. : **Long-term longitudinal measurements of plasma dehydroepiandrosterone sulphate in normal men.** J Clin Endocrinol Metab 75:1002-1004, 1992.
3. Abbasi A, Drinka P, Mattson D, Rudman D : **Low circulating levels of insulin-like growth factors and testosterone in chronically institutionalized elderly men.** J Am Geriatr Soc 41:975-981, 1993.
4. Belanger A, Candas B, Dupont A et al. : **Changes in serum concentrations of conjugated and unconjugated steroid in 40 to 80 year old men.** J Clin Endocrinol Metab 79:1086-1090, 1994.
5. Borst S, Millard W, Lowenthal D : **Growth hormone, exercise, and aging : the future of therapy for the frail elderly.** J Am Geriatr Soc 42:528-535, 1994.
6. Morales A, Nolan J, Nelson J, Yen S : **Effects of replacement dose of dehydroepiandrosterone in men and women of advancing age.** J Clin Endocrinol Metab 78:1360-1367, 1994.
7. Tenover J : **Androgen administration to aging men.** Clin Androl 23:877-892, 1994.
8. Rosenfeld R, Rosenberg B, Fukushima D, Heuman L : **24 hour secretory pattern of dehydroepiandrosterone and dehydroandrosterone sulphate.** J Clin Endocrinol 40:850-855, 1975.
9. Liu C, Laughlin G, Fischer U, Yen S : **Marked attenuation of ultradian and circadian rhythms of dehydroepiandrosterone in postmenopausal women : evidence for a reduced 17,20-desmolase enzymatic activity.** J Clin Endocrinol Metab 71:900-906, 1990.
10. Thomas G, Frenoy N, Legrain S, Sebag-Lanoe R, Baulieu E, Debuire B : **Serum dehydroepiandrosterone sulphate levels as an individual marker.** J Clin Endocrinol Metab 79:1273-1276, 1994.
11. Tissandier O, Peres G, Fiet J, Piette F : **Testosterone, dehydroepiandrosterone, insulin-like growth factor 1, and insulin in sedentary and physically trained aged men.** Eur J Appl Physiol 85:177-184; 2001.
12. Guagnano MT, Del Ponte A, Manigrasso MR, Merlitti D, Pace-Palitti V, Sensis S : **Age-related circadian rhythm of DHEA-S plasma levels in male subjects.** Biological Rhythm Research 32,3:323-331, 2001.
13. Derkken J, Nagesser S.K, Meinders A.E, Haak H.R, Van De Velde C.J.H : **Identification of virilizing adrenal tumors in hirsute women.** New Engl J Med 331(15):968-973, 1994.
14. Lobo R.A, Paul W.L, Goebelmann U : **dehydroepiandrosterone sulphate as an indicator of adrenal androgen function.** Obstet Gynecol 57:69-73, 1981.
15. Cacciari E : **Gonadal and adrenal secretion of dehydroepiandrosterone sulphate in prepubertal and pubertal subjects.** J Endocrinol Invest 4:197-202, 1981.
16. Ebeling P, Koivisto V.A : **Physiological importance of dehydroepiandrosterone.** Lancet 343:1479-1481, 1994.

#### XVIII. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS µl	CALIBRATORS µl	SAMPLE(S) CONTROLS µl
Calibrators (0-5) Samples, controls Tracer	- - 500	10 - 500	- 10 500
Incubation	1 hour at 37°C in water bath		
Separation Working Wash solution Separation Working Wash solution Separation	-	Aspirate 2.0 ml Aspirate 2.0 ml Aspirate carefully	
Counting	Count tubes for 60 seconds		

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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## XVII. BIBLIOGRAPHIE

- Orentreich N, Brind J, Rizer R, Vogelman J : **Age changes and sex differences in serum dehydroepiandrosterone sulphate concentrations throughout adulthood.** J Clin Endocrinol Metab 59:551-555, 1984.
- Orentreich N, Brind J, Vogelman J, Andres R, Baldwin H. : **Long-term longitudinal measurements of plasma dehydroepiandrosterone sulphate in normal men.** J Clin Endocrinol Metab 75:1002-1004, 1992.
- Abbasi A, Drinka P, Mattson D, Rudman D : **Low circulating levels of insulin-like growth factors and testosterone in chronically institutionalized elderly men.** J Am Geriatr Soc 41:975-981, 1993.
- Belanger A, Candas B, Dupont A et al. : **Changes in serum concentrations of conjugated and unconjugated steroid in 40 to 80 year old men.** J Clin Endocrinol Metab 79:1086-1090, 1994.
- Borst S, Millard W, Lowenthal D : **Growth hormone, exercise, and aging : the future of therapy for the frail elderly.** J Am Geriatr Soc 42:528-535, 1994.
- Morales A, Nolan J, Nelson J, Yen S : **Effects of replacement dose of dehydroepiandrosterone in men and women of advancing age.** J Clin Endocrinol Metab 78:1360-1367, 1994.
- Tenover J : **Androgen administration to aging men.** Clin Androl 23:877-892, 1994.
- Rosenfeld R, Rosenberg B, Fukushima D, Heuman L : **24 hour secretory pattern of dehydroepiandrosterone and dehydroandrosterone sulphate.** J Clin Endocrinol 40:850-855, 1975.
- Liu C, Laughlin G, Fischer U, Yen S : **Marked attenuation of ultradian and circadian rhythms of dehydroepiandrosterone in postmenopausal women : evidence for a reduced 17,20-desmolase enzymatic activity.** J Clin Endocrinol Metab 71:900-906, 1990.

- Thomas G, Frenoy N, Legrain S, Sebag-Lanoe R, Baulieu E, Debuire B : **Serum dehydroepiandrosterone sulphate levels as an individual marker.** J Clin Endocrinol Metab 79:1273-1276, 1994.
- Tissandier O, Peres G, Fiet J, Piette F : **Testosterone, dehydroepiandrosterone, insulin-like growth factor 1, and insulin in sedentary and physically trained aged men.** Eur J Appl Physiol 85:177-184; 2001.
- Guagnano MT, Del Ponte A, Manigrasso MR, Merlitti D, Pace-Palitti V, Sensis S : **Age-related circadian rhythm of DHEA-S plasma levels in male subjects.** Biological Rhythm Research 32,3:323-331, 2001.
- DerkSEN J, Nagesser S.K, Meinders A.E, Haak H.R, Van De Velde C.J.H : **Identification of virilizing adrenal tumors in hirsute women.** New Engl J Med 331(15):968-973, 1994.
- Lobo R.A, Paul W.L, Goebelmann U : **dehydroepiandrosterone sulphate as an indicator of adrenal androgen function.** Obstet Gynecol 57:69-73, 1981.
- Cacciari E : **Gonadal and adrenal secretion of dehydroepiandrosterone sulphate in prepubertal and pubertal subjects.** J Endocrinol Invest 4:197-202, 1981.
- Ebeling P, Koivisto V.A : **Physiological importance of dehydroepiandrosterone.** Lancet 343:1479-1481, 1994.

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## XVII. BIBLIOGRAFIE

- Orentreich N, Brind J, Rizer R, Vogelman J : **Age changes and sex differences in serum dehydroepiandrosterone sulphate concentrations throughout adulthood.** J Clin Endocrinol Metab 59:551-555, 1984.
- Orentreich N, Brind J, Vogelman J, Andres R, Baldwin H. : **Long-term longitudinal measurements of plasma dehydroepiandrosterone sulphate in normal men.** J Clin Endocrinol Metab 75:1002-1004, 1992.
- Abbasi A, Drinka P, Mattson D, Rudman D : **Low circulating levels of insulin-like growth factors and testosterone in chronically institutionalized elderly men.** J Am Geriatr Soc 41:975-981, 1993.
- Belanger A, Candas B, Dupont A et al. : **Changes in serum concentrations of conjugated and unconjugated steroid in 40 to 80 year old men.** J Clin Endocrinol Metab 79:1086-1090, 1994.
- Borst S, Millard W, Lowenthal D : **Growth hormone, exercise, and aging : the future of therapy for the frail elderly.** J Am Geriatr Soc 42:528-535, 1994.
- Morales A, Nolan J, Nelson J, Yen S : **Effects of replacement dose of dehydroepiandrosterone in men and women of advancing age.** J Clin Endocrinol Metab 78:1360-1367, 1994.
- Tenover J : **Androgen administration to aging men.** Clin Androl 23:877-892, 1994.
- Rosenfeld R, Rosenberg B, Fukushima D, Heuman L : **24 hour secretory pattern of dehydroepiandrosterone and dehydroandrosterone sulphate.** J Clin Endocrinol 40:850-855, 1975.
- Liu C, Laughlin G, Fischer U, Yen S : **Marked attenuation of ultradian and circadian rhythms of dehydroepiandrosterone in postmenopausal women : evidence for a reduced 17,20-desmolase enzymatic activity.** J Clin Endocrinol Metab 71:900-906, 1990.
- Thomas G, Frenoy N, Legrain S, Sebag-Lanoe R, Baulieu E, Debuire B : **Serum dehydroepiandrosterone sulphate levels as an individual marker.** J Clin Endocrinol Metab 79:1273-1276, 1994.
- Tissandier O, Peres G, Fiet J, Piette F : **Testosterone, dehydroepiandrosterone, insulin-like growth factor 1, and insulin in sedentary and physically trained aged men.** Eur J Appl Physiol 85:177-184; 2001.
- Guagnano MT, Del Ponte A, Manigrasso MR, Merlitti D, Pace-Palitti V, Sensis S : **Age-related circadian rhythm of DHEA-S plasma levels in male subjects.** Biological Rhythm Research 32,3:323-331, 2001.
- Derkzen J, Nagesser S.K, Meinders A.E, Haak H.R, Van De Velde C.J.H : **Identification of virilizing adrenal tumors in hirsute women.** New Engl J Med 331(15):968-973, 1994.
- Lobo R.A, Paul W.L, Goebelsmann U : **dehydroepiandrosterone sulphate as an indicator of adrenal androgen function.** Obstet Gynecol 57:69-73, 1981.
- Cacciari E : **Gonadal and adrenal secretion of dehydroepiandrosterone sulphate in prepubertal and pubertal subjects.** J Endocrinol Invest 4:197-202, 1981.
- Ebeling P, Koivisto V.A : **Physiological importance of dehydroepiandrosterone.** Lancet 343:1479-1481, 1994.

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- Abbas A, Drinka P, Mattson D, Rudman D : **Low circulating levels of insulin-like growth factors and testosterone in chronically institutionalized elderly men.** J Am Geriatr Soc 41:975-981, 1993.
- Belanger A, Candas B, Dupont A et al. : **Changes in serum concentrations of conjugated and unconjugated steroid in 40 to 80 year old men.** J Clin Endocrinol Metab 79:1086-1090, 1994.
- Borst S, Millard W, Lowenthal D : **Growth hormone, exercise, and aging : the future of therapy for the frail elderly.** J Am Geriatr Soc 42:528-535, 1994.
- Morales A, Nolan J, Nelson J, Yen S : **Effects of replacement dose of dehydroepiandrosterone in men and women of advancing age.** J Clin Endocrinol Metab 78:1360-1367, 1994.
- Tenover J : **Androgen administration to aging men.** Clin Androl 23:877-892, 1994.
- Rosenfeld R, Rosenberg B, Fukushima D, Heuman L : **24 hour secretory pattern of dehydroepiandrosterone and dehydroandrosterone sulphate.** J Clin Endocrinol 40:850-855, 1975.
- Liu C, Laughlin G, Fischer U, Yen S : **Marked attenuation of ultradian and circadian rhythms of dehydroepiandrosterone in postmenopausal women : evidence for a reduced 17,20-desmolase enzymatic activity.** J Clin Endocrinol Metab 71:900-906, 1990.
- Thomas G, Frenoy N, Legrain S, Sebag-Lanoe R, Baulieu E, Debuire B : **Serum dehydroepiandrosterone sulphate levels as an individual marker.** J Clin Endocrinol Metab 79:1273-1276, 1994.
- Tissandier O, Peres G, Fiet J, Piette F : **Testosterone, dehydroepiandrosterone, insulin-like growth factor 1, and insulin in sedentary and physically trained aged men.** Eur J Appl Physiol 85:177-184; 2001.
- Guagnano MT, Del Ponte A, Manigrasso MR, Merlitti D, Pace-Palitti V, Sensis S : **Age-related circadian rhythm of DHEA-S plasma levels in male subjects.** Biological Rhythm Research 32,3:323-331, 2001.
- Derksen J, Nagesser S.K, Meinders A.E, Haak H.R, Van De Velde C.J.H : **Identification of virilizing adrenal tumors in hirsute women.** New Engl J Med 331(15):968-973, 1994.
- Lobo R.A, Paul W.L, Goebelsmann U : **dehydroepiandrosterone sulphate as an indicator of adrenal androgen function.** Obstet Gynecol 57:69-73, 1981.
- Cacciari E : **Gonadal and adrenal secretion of dehydroepiandrosterone sulphate in prepubertal and pubertal subjects.** J Endocrinol Invest 4:197-202, 1981.
- Ebeling P, Koivisto V.A : **Physiological importance of dehydroepiandrosterone.** Lancet 343:1479-1481, 1994.

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1. Orentreich N, Brind J, Rizer R, Vogelman J : Age changes and sex differences in serum dehydroepiandrosterone sulphate concentrations throughout adulthood. *J Clin Endocrinol Metab* 59:551-555, 1984.
2. Orentreich N, Brind J, Vogelman J, Andres R, Baldwin H. : Long-term longitudinal measurements of plasma dehydroepiandrosterone sulphate in normal men. *J Clin Endocrinol Metab* 75:1002-1004, 1992.
3. Abbasi A, Drinka P, Mattson D, Rudman D : Low circulating levels of insulin-like growth factors and testosterone in chronically institutionalized elderly men. *J Am Geriatr Soc* 41:975-981, 1993.
4. Belanger A, Candas B, Dupont A et al. : Changes in serum concentrations of conjugated and unconjugated steroid in 40 to 80 year old men. *J Clin Endocrinol Metab* 79:1086-1090, 1994.
5. Borst S, Millard W, Lowenthal D : Growth hormone, exercise, and aging : the future of therapy for the frail elderly. *J Am Geriatr Soc* 42:528-535, 1994.
6. Morales A, Nolan J, Nelson J, Yen S : Effects of replacement dose of dehydroepiandrosterone in men and women of advancing age. *J Clin Endocrinol Metab* 78:1360-1367, 1994.
7. Tenover J : Androgen administration to aging men. *Clin Androl* 23:877-892, 1994.
8. Rosenfeld R, Rosenberg B, Fukushima D, Heuman L : 24 hour secretory pattern of dehydroepiandrosterone and dehydroandrosterone sulphate. *J Clin Endocrinol* 40:850-855, 1975.

9. Liu C, Laughlin G, Fischer U, Yen S : Marked attenuation of ultradian and circadian rhythms of dehydroepiandrosterone in postmenopausal women : evidence for a reduced 17,20-desmolase enzymatic activity. *J Clin Endocrinol Metab* 71:900-906, 1990.
10. Thomas G, Frenoy N, Legrain S, Sebag-Lanoë R, Baulieu E, Debuire B : Serum dehydroepiandrosterone sulphate levels as an individual marker. *J Clin Endocrinol Metab* 79:1273-1276, 1994.
11. Tissandier O, Peres G, Fiet J, Piette F : Testosterone, dehydroepiandrosterone, insulin-like growth factor 1, and insulin in sedentary and physically trained aged men. *Eur J Appl Physiol* 85:177-184; 2001.
12. Guagnano MT, Del Ponte A, Manigrasso MR, Merlitti D, Pace-Palitti V, Sensis S : Age-related circadian rhythm of DHEA-S plasma levels in male subjects. *Biological Rhythm Research* 32,3:323-331, 2001.
13. Derksen J, Nagesser S.K, Meinders A.E, Haak H.R, Van De Velde C.J.H : Identification of virilizing adrenal tumors in hirsute women. *New Engl J Med* 331(15):968-973, 1994.
14. Lobo R.A, Paul W.L, Goebelsmann U : dehydroepiandrosterone sulphate as an indicator of adrenal androgen function. *Obstet Gynecol* 57:69-73, 1981.
15. Cacciari E : Gonadal and adrenal secretion of dehydroepiandrosterone sulphate in prepubertal and pubertal subjects. *J Endocrinol Invest* 4:197-202, 1981.
16. Ebeling P, Koivisto V.A : Physiological importance of dehydroepiandrosterone. *Lancet* 343:1479-1481, 1994.

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## XVII. BIBLIOGRAFIA

- Orentreich N, Brind J, Rizer R, Vogelman J : **Age changes and sex differences in serum dehydroepiandrosterone sulphate concentrations throughout adulthood.** J Clin Endocrinol Metab 59:551-555, 1984.
- Orentreich N, Brind J, Vogelman J, Andres R, Baldwin H. : **Long-term longitudinal measurements of plasma dehydroepiandrosterone sulphate in normal men.** J Clin Endocrinol Metab 75:1002-1004, 1992.
- Abbas A, Drinka P, Mattson D, Rudman D : **Low circulating levels of insulin-like growth factors and testosterone in chronically institutionalized elderly men.** J Am Geriatr Soc 41:975-981, 1993.
- Belanger A, Candas B, Dupont A et al. : **Changes in serum concentrations of conjugated and unconjugated steroid in 40 to 80 year old men.** J Clin Endocrinol Metab 79:1086-1090, 1994.
- Borst S, Millard W, Lowenthal D : **Growth hormone, exercise, and aging : the future of therapy for the frail elderly.** J Am Geriatr Soc 42:528-535, 1994.
- Morales A, Nolan J, Nelson J, Yen S : **Effects of replacement dose of dehydroepiandrosterone in men and women of advancing age.** J Clin Endocrinol Metab 78:1360-1367, 1994.
- Tenover J : **Androgen administration to aging men.** Clin Andol 23:877-892, 1994.

- Rosenfeld R, Rosenberg B, Fukushima D, Heuman L : **24 hour secretory pattern of dehydroepiandrosterone and dehydroandrosterone sulphate.** J Clin Endocrinol 40:850-855, 1975.
- Liu C, Laughlin G, Fischer U, Yen S : **Marked attenuation of ultradian and circadian rhythms of dehydroepiandrosterone in postmenopausal women : evidence for a reduced 17,20-desmolase enzymatic activity.** J Clin Endocrinol Metab 71:900-906, 1990.
- Thomas G, Frenoy N, Legrain S, Sebag-Lanoe R, Baulieu E, Debuire B : **Serum dehydroepiandrosterone sulphate levels as an individual marker.** J Clin Endocrinol Metab 79:1273-1276, 1994.
- Tissandier O, Peres G, Fiet J, Piette F : **Testosterone, dehydroepiandrosterone, insulin-like growth factor 1, and insulin in sedentary and physically trained aged men.** Eur J Appl Physiol 85:177-184; 2001.
- Guagnano MT, Del Ponte A, Manigrasso MR, Merlitti D, Pace-Palitti V, Sensis S : **Age-related circadian rhythm of DHEA-S plasma levels in male subjects.** Biological Rhythm Research 32,3:323-331, 2001.
- Derkzen J, Nagesser S.K, Meinders A.E, Haak H.R, Van De Velde C.J.H : **Identification of virilizing adrenal tumors in hirsute women.** New Engl J Med 331(15):968-973, 1994.
- Lobo R.A, Paul W.L, Goebelmann U : **dehydroepiandrosterone sulphate as an indicator of adrenal androgen function.** Obstet Gynecol 57:69-73, 1981.
- Cacciari E : **Gonadal and adrenal secretion of dehydroepiandrosterone sulphate in prepubertal and pubertal subjects.** J Endocrinol Invest 4:197-202, 1981.
- Ebeling P, Koivisto V.A : **Physiological importance of dehydroepiandrosterone.** Lancet 343:1479-1481, 1994.

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## XVII. BIBLIOGRAFIA

1. Orentreich N, Brind J, Rizer R, Vogelman J : **Age changes and sex differences in serum dehydroepiandrosterone sulphate concentrations throughout adulthoold.** J Clin Endocrinol Metab 59:551-555, 1984.
2. Orentreich N, Brind J, Vogelman J, Andres R, Baldwin H. : **Long-term longitudinal measurements of plasma dehydroepiandrosterone sulphate in normal men.** J Clin Endocrinol Metab 75:1002-1004, 1992.
3. Abbasi A, Drinka P, Mattson D, Rudman D : **Low circulating levels of insulin-like growth factors and testosterone in chronically institutionalized elderly men.** J Am Geriatr Soc 41:975-981, 1993.
4. Belanger A, Candas B, Dupont A et al. : **Changes in serum concentrations of conjugated and unconjugated steroid in 40 to 80 year old men.** J Clin Endocrinol Metab 79:1086-1090, 1994.
5. Borst S, Millard W, Lowenthal D : **Growth hormone, exercise, and aging : the future of therapy for the frail elderly.** J Am Geriatr Soc 42:528-535, 1994.
6. Morales A, Nolan J, Nelson J, Yen S : **Effects of replacement dose of dehydroepiandrosterone in men and women of advancing age.** J Clin Endocrinol Metab 78:1360-1367, 1994.
7. Tenover J : **Androgen administration to aging men.** Clin Androl 23:877-892, 1994.
8. Rosenfeld R, Rosenberg B, Fukushima D, Heuman L : **24 hour secretory pattern of dehydroepiandrosterone and dehydroandrosterone sulphate.** J Clin Endocrinol 40:850-855, 1975.
9. Liu C, Laughlin G, Fischer U, Yen S : **Marked attenuation of ultradian and circadian rhythms of dehydroepiandrosterone in postmenopausal women : evidence for a reduced 17,20-desmolase enzymatic activity.** J Clin Endocrinol Metab 71:900-906, 1990.
10. Thomas G, Frenoy N, Legrain S, Sebag-Lanoe R, Baulieu E, Debuire B : **Serum dehydroepiandrosterone sulphate levels as an individual marker.** J Clin Endocrinol Metab 79:1273-1276, 1994.

11. Tissandier O, Peres G, Fiet J, Piette F : **Testosterone, dehydroepiandrosterone, insulin-like growth factor 1, and insulin in sedentary and physically trained aged men.** Eur J Appl Physiol 85:177-184; 2001.
12. Guagnano MT, Del Ponte A, Manigrasso MR, Merlitti D, Pace-Palitti V, Sensis S : **Age-related circadian rhythm of DHEA-S plasma levels in male subjects.** Biological Rhythm Research 32,3:323-331, 2001.
13. Derkzen J, Nagesser S.K, Meinders A.E, Haak H.R, Van De Velde C.J.H : **Identification of virilizing adrenal tumors in hirsute women.** New Engl J Med 331(15):968-973, 1994.
14. Lobo R.A, Paul W.L, Goebelsmann U : **dehydroepiandrosterone sulphate as an indicator of adrenal androgen function.** Obstet Gynecol 57:69-73, 1981.
15. Cacciari E : **Gonadal and adrenal secretion of dehydroepiandrosterone sulphate in prepubertal and pubertal subjects.** J Endocrinol Invest 4:197-202, 1981.
16. Ebeling P, Koivisto V.A : **Physiological importance of dehydroepiandrosterone.** Lancet 343:1479-1481, 1994.

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- Orentreich N, Brind J, Rizer R, Vogelman J : Age changes and sex differences in serum dehydroepiandrosterone sulphate concentrations throughout adulthood. J Clin Endocrinol Metab 59:551-555, 1984.

- Orentreich N, Brind J, Vogelman J, Andres R, Baldwin H : Long-term longitudinal measurements of plasma dehydroepiandrosterone sulphate in normal men. J Clin Endocrinol Metab 75:1002-1004, 1992.
- Abbas A, Drinka P, Mattson D, Rudman D : Low circulating levels of insulin-like growth factors and testosterone in chronically institutionalized elderly men. J Am Geriatr Soc 41:975-981, 1993.
- Belanger A, Candas B, Dupont A et al. : Changes in serum concentrations of conjugated and unconjugated steroid in 40 to 80 year old men. J Clin Endocrinol Metab 79:1086-1090, 1994.
- Borst S, Millard W, Lowenthal D : Growth hormone, exercise, and aging : the future of therapy for the frail elderly. J Am Geriatr Soc 42:528-535, 1994.
- Morales A, Nolan J, Nelson J, Yen S : Effects of replacement dose of dehydroepiandrosterone in men and women of advancing age. J Clin Endocrinol Metab 78:1360-1367, 1994.
- Tenover J : Androgen administration to aging men. Clin Androl 23:877-892, 1994.
- Rosenfeld R, Rosenberg B, Fukushima D, Heuman L : 24 hour secretory pattern of dehydroepiandrosterone and dehydroandrosterone sulphate. J Clin Endocrinol 40:850-855, 1975.
- Liu C, Laughlin G, Fischer U, Yen S : Marked attenuation of ultradian and circadian rhythms of dehydroepiandrosterone in postmenopausal women : evidence for a reduced 17,20-desmolase enzymatic activity. J Clin Endocrinol Metab 71:900-906, 1990.
- Thomas G, Frenoy N, Legrain S, Sebag-Lanoe R, Baulieu E, Debuire B : Serum dehydroepiandrosterone sulphate levels as an individual marker. J Clin Endocrinol Metab 79:1273-1276, 1994.
- Tissandier O, Peres G, Fiet J, Piette F : Testosterone, dehydroepiandrosterone, insulin-like growth factor 1, and insulin in sedentary and physically trained aged men. Eur J Appl Physiol 85:177-184; 2001.
- Guagnano MT, Del Ponte A, Manigrasso MR, Merlitti D, Pace-Palitti V, Sensis S : Age-related circadian rhythm of DHEA-S plasma levels in male subjects. Biological Rhythm Research 32,3:323-331, 2001.
- Derkzen J, Nagesser S.K, Meinders A.E, Haak H.R, Van De Velde C.J.H : Identification of virilizing adrenal tumors in hirsute women. New Engl J Med 331(15):968-973, 1994.
- Lobo R.A, Paul W.L, Goebelmann U : dehydroepiandrosterone sulphate as an indicator of adrenal androgen function. Obstet Gynecol 57:69-73, 1981.
- Cacciari E : Gonadal and adrenal secretion of dehydroepiandrosterone sulphate in prepubertal and pubertal subjects. J Endocrinol Invest 4:197-202, 1981.
- Ebeling P, Koivisto V.A : Physiological importance of dehydroepiandrosterone. Lancet 343:1479-1481, 1994.

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2. Orentreich N, Brind J, Vogelman J, Andres R, Baldwin H. : **Long-term longitudinal measurements of plasma dehydroepiandrosterone sulphate in normal men.** J Clin Endocrinol Metab 75:1002-1004, 1992.
3. Abbasi A, Drinka P, Mattson D, Rudman D : **Low circulating levels of insulin-like growth factors and testosterone in chronically institutionalized elderly men.** J Am Geriatr Soc 41:975-981, 1993.
4. Belanger A, Candas B, Dupont A et al. : **Changes in serum concentrations of conjugated and unconjugated steroid in 40 to 80 year old men.** J Clin Endocrinol Metab 79:1086-1090, 1994.
5. Borst S, Millard W, Lowenthal D : **Growth hormone, exercise, and aging : the future of therapy for the frail elderly.** J Am Geriatr Soc 42:528-535, 1994.
6. Morales A, Nolan J, Nelson J, Yen S : **Effects of replacement dose of dehydroepiandrosterone in men and women of advancing age.** J Clin Endocrinol Metab 78:1360-1367, 1994.
7. Tenover J : **Androgen administration to aging men.** Clin Androl 23:877-892, 1994.
8. Rosenfeld R, Rosenberg B, Fukushima D, Heuman L : **24 hour secretory pattern of dehydroepiandrosterone and dehydroandrosterone sulphate.** J Clin Endocrinol 40:850-855, 1975.
9. Liu C, Laughlin G, Fischer U, Yen S : **Marked attenuation of ultradian and circadian rhythms of dehydroepiandrosterone in postmenopausal women : evidence for a reduced 17,20-desmolase enzymatic activity.** J Clin Endocrinol Metab 71:900-906, 1990.
10. Thomas G, Frenoy N, Legrain S, Sebag-Lanoe R, Baulieu E, Debuire B : **Serum dehydroepiandrosterone sulphate levels as an individual marker.** J Clin Endocrinol Metab 79:1273-1276, 1994.
11. Tissandier O, Peres G, Fiet J, Piette F : **Testosterone, dehydroepiandrosterone, insulin-like growth factor 1, and insulin in sedentary and physically trained aged men.** Eur J Appl Physiol 85:177-184; 2001.
12. Guagnano MT, Del Ponte A, Manigrasso MR, Merlitti D, Pace-Palitti V, Sensis S : **Age-related circadian rhythm of DHEA-S plasma levels in male subjects.** Biological Rhythm Research 32,3:323-331, 2001.
13. Derkken J, Nagesser S.K, Meinders A.E, Haak H.R, Van De Velde C.J.H : **Identification of virilizing adrenal tumors in hirsute women.** New Engl J Med 331(15):968-973, 1994.
14. Lobo R.A, Paul W.L, Goebelsmann U : **dehydroepiandrosterone sulphate as an indicator of adrenal androgen function.** Obstet Gynecol 57:69-73, 1981.
15. Cacciari E : **Gonadal and adrenal secretion of dehydroepiandrosterone sulphate in prepubertal and pubertal subjects.** J Endocrinol Invest 4:197-202, 1981.
16. Ebeling P, Koivisto V.A : **Physiological importance of dehydroepiandrosterone.** Lancet 343:1479-1481, 1994.

## XVII. BIBLIOGRAFIA

1. Orentreich N, Brind J, Rizer R, Vogelman J : **Age changes and sex differences in serum dehydroepiandrosterone sulphate concentrations throughout adulthood.** J Clin Endocrinol Metab 59:551-555, 1984.

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- Orentreich N, Brind J, Rizer R, Vogelman J : Age changes and sex differences in serum dehydroepiandrosterone sulphate concentrations throughout adulthood. *J Clin Endocrinol Metab* 59:551-555, 1984.

- Orentreich N, Brind J, Vogelman J, Andres R, Baldwin H. : Long-term longitudinal measurements of plasma dehydroepiandrosterone sulphate in normal men. *J Clin Endocrinol Metab* 75:1002-1004, 1992.
- Abbasi A, Drinka P, Mattson D, Rudman D : Low circulating levels of insulin-like growth factors and testosterone in chronically institutionalized elderly men. *J Am Geriatr Soc* 41:975-981, 1993.
- Belanger A, Candas B, Dupont A et al. : Changes in serum concentrations of conjugated and unconjugated steroid in 40 to 80 year old men. *J Clin Endocrinol Metab* 79:1086-1090, 1994.
- Borst S, Millard W, Lowenthal D : Growth hormone, exercise, and aging : the future of therapy for the frail elderly. *J Am Geriatr Soc* 42:528-535, 1994.
- Morales A, Nolan J, Nelson J, Yen S : Effects of replacement dose of dehydroepiandrosterone in men and women of advancing age. *J Clin Endocrinol Metab* 78:1360-1367, 1994.
- Tenover J : Androgen administration to aging men. *Clin Androl* 23:877-892, 1994.
- Rosenfeld R, Rosenberg B, Fukushima D, Heuman L : 24 hour secretory pattern of dehydroepiandrosterone and dehydroandrosterone sulphate. *J Clin Endocrinol* 40:850-855, 1975.
- Liu C, Laughlin G, Fischer U, Yen S : Marked attenuation of ultradian and circadian rhythms of dehydroepiandrosterone in postmenopausal women : evidence for a reduced 17,20-desmolase enzymatic activity. *J Clin Endocrinol Metab* 71:900-906, 1990.
- Thomas G, Frenoy N, Legrain S, Sebag-Lanoe R, Baulieu E, Debuire B : Serum dehydroepiandrosterone sulphate levels as an individual marker. *J Clin Endocrinol Metab* 79:1273-1276, 1994.
- Tissandier O, Peres G, Fiet J, Piette F : Testosterone, dehydroepiandrosterone, insulin-like growth factor 1, and insulin in sedentary and physically trained aged men. *Eur J Appl Physiol* 85:177-184; 2001.
- Guagnano MT, Del Ponte A, Manigrasso MR, Merlitti D, Pace-Palitti V, Sensis S : Age-related circadian rhythm of DHEA-S plasma levels in male subjects. *Biological Rhythm Research* 32,3:323-331, 2001.
- Derkzen J, Nagesser S.K, Meinders A.E, Haak H.R, Van De Velde C.J.H : Identification of virilizing adrenal tumors in hirsute women. *New Engl J Med* 331(15):968-973, 1994.
- Lobo R.A, Paul W.L, Goebelmann U : dehydroepiandrosterone sulphate as an indicator of adrenal androgen function. *Obstet Gynecol* 57:69-73, 1981.
- Cacciari E : Gonadal and adrenal secretion of dehydroepiandrosterone sulphate in prepubertal and pubertal subjects. *J Endocrinol Invest* 4:197-202, 1981.
- Ebeling P, Koivisto V.A : Physiological importance of dehydroepiandrosterone. *Lancet* 343:1479-1481, 1994.

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## XVII. IRODALOM

- Orentreich N, Brind J, Rizer R, Vogelman J : Age changes and sex differences in serum dehydroepiandrosterone sulphate concentrations throughout adulthood. *J Clin Endocrinol Metab* 59:551-555, 1984.
- Orentreich N, Brind J, Vogelman J, Andres R, Baldwin H : Long-term longitudinal measurements of plasma dehydroepiandrosterone sulphate in normal men. *J Clin Endocrinol Metab* 75:1002-1004, 1992.

3. Abbasi A, Drinka P, Mattson D, Rudman D : Low circulating levels of insulin-like growth factors and testosterone in chronically institutionalized elderly men. *J Am Geriatr Soc* 41:975-981, 1993.

- Belanger A, Candas B, Dupont A et al. : Changes in serum concentrations of conjugated and unconjugated steroid in 40 to 80 year old men. *J Clin Endocrinol Metab* 79:1086-1090, 1994.
- Borst S, Millard W, Lowenthal D : Growth hormone, exercise, and aging : the future of therapy for the frail elderly. *J Am Geriatr Soc* 42:528-535, 1994.
- Morales A, Nolan J, Nelson J, Yen S : Effects of replacement dose of dehydroepiandrosterone in men and women of advancing age. *J Clin Endocrinol Metab* 78:1360-1367, 1994.
- Tenover J : Androgen administration to aging men. *Clin Andol* 23:877-892, 1994.
- Rosenfeld R, Rosenberg B, Fukushima D, Heuman L : 24 hour secretory pattern of dehydroepiandrosterone and dehydroandrosterone sulphate. *J Clin Endocrinol* 40:850-855, 1975.
- Liu C, Laughlin G, Fischer U, Yen S : Marked attenuation of ultradian and circadian rhythms of dehydroepiandrosterone in postmenopausal women : evidence for a reduced 17,20-desmolase enzymatic activity. *J Clin Endocrinol Metab* 71:900-906, 1990.
- Thomas G, Frenoy N, Legrain S, Sebag-Lanoe R, Baulieu E, Debuire B : Serum dehydroepiandrosterone sulphate levels as an individual marker. *J Clin Endocrinol Metab* 79:1273-1276, 1994.
- Tissandier O, Peres G, Fiet J, Piette F : Testosterone, dehydroepiandrosterone, insulin-like growth factor 1, and insulin in sedentary and physically trained aged men. *Eur J Appl Physiol* 85:177-184; 2001.
- Guagnano MT, Del Ponte A, Manigrasso MR, Merlitti D, Pace-Palitti V, Sensis S : Age-related circadian rhythm of DHEA-S plasma levels in male subjects. *Biological Rhythm Research* 32,3:323-331, 2001.
- DerkSEN J, Nagesser S.K, Meinders A.E, Haak H.R, Van De Velde C.J.H : Identification of virilizing adrenal tumors in hirsute women. *New Engl J Med* 331(15):968-973, 1994.
- Lobo R.A, Paul W.L, Goebelmann U : dehydroepiandrosterone sulphate as an indicator of adrenal androgen function. *Obstet Gynecol* 57:69-73, 1981.
- Cacciari E : Gonadal and adrenal secretion of dehydroepiandrosterone sulphate in prepubertal and pubertal subjects. *J Endocrinol Invest* 4:197-202, 1981.
- Ebeling P, Koivisto V.A : Physiological importance of dehydroepiandrosterone. *Lancet* 343:1479-1481, 1994.

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