FSH-IRMA

KIP0841

Version : 230123

Date of issue : 23/01/2023

Revision date: 23/01/2023

History

Summary of change:

Current Version:
230123
New logo

en

Read entire protocol before use.

FSH-IRMA

I. INTENDED USE

Immunoradiometric assay kit for the *in vitro* quantitative measurement of human Follicle Stimulating Hormone (FSH) in serum and plasma.

II. GENERAL INFORMATION

A. Proprietary name : DIAsource FSH-IRMA Kit

B. Catalog number : KIP0841: 96 tests

C. Manufactured by : DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

The measurement of LH and FSH concentrations in serum is essential for investigating fertility and especially disorders of the hypothalamic/pituitary/gonadal axis.

Both LH and FSH are secreted by the basophil cells of the anterior pituitary as a result of gonadotropin releasing hormone (GnRH) secretion from hypothalamic cells.

In adults, LH and FSH hormones control gonadal functions; mainly gametogenesis and steroid secretion.

Circulating levels of FSH are controlled by a negative feedback effect on the hypothalamus by steroid hormones and gonadal peptides.

IV. PRINCIPLES OF THE METHOD

The DiaSource FSH-IRMA is an immunoradiometric assay based on coated tubes separation. Mabs1, the capture antibodies, are attached to the lower and inner surface of the plastic tube. Calibrators or samples added to the tubes will at first show low affinity for Mabs1. Addition of Mab2, the signal antibody labelled with ^{125}I , will complete the system and trigger the immunological reaction. After washing, the remaining radioactivity bound to the tube reflects the antigen concentration. The use of several distinct Mabs avoids hyperspecificity.

V. REAGENTS PROVIDED

Reagents	Quantity 96 tests	Quantity 4x96 tests	Colour Code	Reconstitution
Tubes coated with anti FSH (monoclonal antibodies)	2 x 48	8 x 48	red	Ready for use
Anti-FSH- ^{125}I (monoclonal antibodies) in borate-PO ₄ buffer with bovine serum albumin, azide (<0.1%) and inert red dye	1 vial 5.5 ml 700 kBq	4 vials 5.5 ml 4x700 kBq	red	Ready for use
Zero Calibrator in bovine serum with thymol	1 vial lyophil.	2 vials lyophil.	yellow	Add 2 ml distilled water
Calibrators 1-6 in bovine serum with thymol (see exact values on vial labels)	6 vials lyophil.	12 vials lyophil.	yellow	Add 1 ml distilled water
Wash solution (TRIS-HCl)	1 vial 10 ml	4 vials 10 ml	brown	Dilute 70x with distilled water (use a magnetic stirrer).
Controls 1 and 2 in human plasma with thymol	2 vials lyophil.	4 vials lyophil.	silver	Add 0.5 ml distilled water

Note: 1. Use the zero calibrator for sera dilutions.
 2. 1 mIU of the calibrator is equivalent to 1 mIU of the 1st IS 92/510.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Pipettes for delivery of: 50 µl, 100 µl, 500 µl, 1 ml and 2 ml. (the use of accurate pipettes with disposable plastic tips is recommended)
3. Vortex mixer
4. Magnetic stirrer
5. Tube shaker
6. 5 ml automatic syringe (Cornwall type) for washing
7. Aspiration system (optional).
8. Any gamma counter capable of measuring ^{125}I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- A. **Calibrators** : Reconstitute the zero calibrator with 2 ml distilled water and the other calibrators with 1 ml distilled water..
- B. **Controls** : Reconstitute the controls with 0.5 ml distilled water.
- C. **Working Wash solution** : Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, calibrators and controls are stable for 3 days at 2-8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 3 months. Avoid subsequent freeze-thaw cycles.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum and plasma must be kept at 2 – 8°C.
 - If the test is not run within 24 hours, storage in aliquots at -20°C is recommended.
 - Avoid subsequent freeze-thaw cycles.
 - Serum, heparinized or EDTA plasma provide similar results.
- $Y \text{ (serum)} = 1.00x \text{ (hep. plasma)} + 0.3 \quad r = 0.99 \quad n = 39$
 $Y \text{ (serum)} = 1.02x \text{ (EDTA plasma)} + 0.3 \quad r = 0.99 \quad n = 39$

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature prior to use. Thoroughly mix all reagents and samples by gentle agitation or swirling. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample. High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times. Prepare a calibration curve for each run, do not use data from previous runs.

B. Procedure

1. Label coated tubes in duplicate for each calibrator, control and sample. For determination of total counts, label 2 normal tubes.
2. Briefly vortex calibrators, samples and controls and dispense 100 µl of each into the respective tubes.
3. Dispense 50 µl of anti-FSH- ^{125}I tracer into each tube, including the uncoated tubes for total counts.
4. Shake the rack containing the tubes gently by hand to liberate any trapped air bubbles.
5. Incubate for 1 hour at room temperature on a tube shaker (700 rpm).
6. Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
7. Wash the tubes with 2 ml Wash Solution (except total counts). Avoid foaming during the addition of the Working Wash Solution.
8. Aspirate (or decant) the content of each tube (except total counts).
9. Wash again the tubes with 2 ml Wash Solution (except total counts) and aspirate (or decant).
10. After the last washing, let the tubes standing upright for two minutes and aspirate the remaining drop of liquid.
11. Count the tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. On semi-logarithmic or linear graph paper plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of FSH (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
3. Read the concentration for each control and sample by interpolation on the calibration curve.
4. Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

FSH-IRMA		cpm	B/T (%)
Total count		362200	100
Calibrator	0.0 mIU/ml	119	0.03
	0.7 mIU/ml	1265	0.35
	3.5 mIU/ml	2783	0.77
	10.4 mIU/ml	8400	2.32
	30.0 mIU/ml	28008	7.73
	100.0 mIU/ml	102364	28.26
	152.0 mIU/ml	160536	44.32

XIII. PERFORMANCE AND LIMITATIONS

A. Detection Limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations above the average counts at zero binding, was 0.1 mIU/ml.

B. Specificity

Cross-reactive hormones were added to a low and to a high FSH value calibrator. The apparent FSH response was measured.

added Hormone	FSH CAL 1		FSH CAL 5	
	mIU/ml	mIU/ml	mIU/ml	mIU/ml
-	0.65	100		
LH 300 mIU/ml	1.31	99.16		
hCG 300000 mIU/ml	1.46	100.31		
TSH 300 µIU/ml	0.98	98.31		

C. Precision

INTRA ASSAY				INTER ASSAY			
Serum	N	$\bar{X} \pm S.D.$ mIU/ml	CV (%)	Serum	N	$\bar{X} \pm S.D.$ (mIU/ml)	CV (%)
A	10	4.0 ± 0.1	1.8	A	20	15.0 ± 3.0	4.4
B	10	9.0 ± 0.2	2.0	B	20	41.9 ± 1.0	2.4
C	10	50.7 ± 0.5	1.1				

D. Accuracy

RECOVERY TEST			
Sample	Added FSH (mIU/ml)	Recovered FSH (mIU/ml)	Recovery (%)
Serum	6.4	6.2	97
	25.1	23.6	94
	88.8	86.1	97
Heparin plasma	6.4	7	109
	25.1	25.1	100
	88.8	96	108
EDTA plasma	6.4	7.3	114
	25.1	23.6	94
	88.8	84.1	95

DILUTION TEST			
Sample	Dilution	Theoretical Concent. (mIU/ml)	Measured Concent. (mIU/ml)
Serum 1	1/1	-	8.6
	1/2	4.3	4.1
	1/4	2.2	2
	1/8	1.1	1.1
	1/16	0.5	0.5
Serum 2	1/1	-	59.6
	1/2	29.8	31.3
	1/4	14.9	16.5
	1/8	7.5	8.2
	1/16	3.7	4.0
	1/32	1.9	2.0
	1/64	0.9	1.0

Samples were diluted with zero calibrator.

E. Time Delay

As shown hereafter, assay results remain accurate even when a sample is dispensed 45 minutes after the calibrator has been added to coated tubes.

TIME DELAY

	0'	15'	30'	45'
Serum 1 (mIU/ml)	1.9	1.7	1.8	1.7
Serum 2 (mIU/ml)	4.2	4.3	4.2	4.2
Serum 3 (mIU/ml)	9.6	9.5	9.3	9.3
Serum 4 (mIU/ml)	52.0	52.0	52.0	53.0

F. Hook effect

A sample spiked with FSH up to 1000 mIU/ml gives higher counts than the last calibrator point.

XIV. LIMITATIONS

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays. Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed in case of the presence of heterophilic antibodies. Carefully evaluate the results of patients suspected of having these antibodies. If results are not consistent with other clinical observations, additional information should be required before diagnosis.

XV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises

XVI. REFERENCE INTERVALS

The values provided below are given only for guidance; each laboratory should establish its own normal range of values.

The range is expressed as 2.5% to 97.5% percentiles.

Identification	Number of subjects	Mean (mIU/ml)	Range (mIU/ml)
Children (0 to 12 years)			
▪ Boys	20	0.9	0.1 – 2.3
▪ Girls	20	2.7	0.4 – 12.7
Pubers (12 to 18 years)	18	4.5	0.3 – 9.0
Adult males	69	4.2	1.3 – 8.1
Women			
▪ Ovulatory cycles			
- Follicular phase (day -12 to -6)	34	6.1	1.8 – 9.4
- Ovulatory peak (day 0)	48	11.9	3.4 – 33.1
- Luteal phase (day +6 to +12)	63	4.3	1.2 – 13.4
▪ Postmenopausal	53	64.7	27.7 – 93.3

XVII. PRECAUTIONS AND WARNINGS

Safety

For in vitro diagnostic use only.

This kit contains ^{125}I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVIII. BIBLIOGRAPHY

- CLARKE J.J., CUMMINS J.T. (1987)
Pulsatility of reproductive hormones: physiological basis and clinical implications.
Baillière's Clin. Endocrinol. Metab., 1:1
- FILICORI M., SANTORO N., MERRIAN G.R., CROWLEY W.F. Jr., (1986)
Characterization of the physiological pattern of episodic gonadotropin secretion throughout the human menstrual cycle.
J. Clin. Endocrinol. Metab., 62:1136
- MAIS V., CETEL N.S., MUSE K.N., QUIGLEY M.E., REID R.L., YEN S.S.C. (1987)
Hormonal dynamics during luteal-follicular transition.
J. Clin. Endocrinol. Metab., 64:1109
- SAWYER-STEFFAN J.E., LASLEY B.L., HOFF J.D., YEN S.S.C. (1982)
Comparison of in-vitro bioactivity and immunoreactivity of serum LH in normal cyclic and hypogonadal women treated with low doses of LH-RH.

J. Reprod. Fertil., 65:45

- THOMAS C.M.G., SEGERS M.F.G. (1988)
Measurement of Serum hLH: hCG interference evaluated for two hLH-Specific IRMA kits.
Clin. Chem., 34:768
- LOUMAYE E., VANKRIEKEN L., DEPREESTER S., PSALTI I., de COOMAN S., THOMAS K. (1989)
Hormonal changes induced by short-term administration of a gonadotropin-releasing hormone agonist during ovarian hyperstimulation for in vitro fertilization and their consequences for embryo development.
Fertil. and Steril., 51:105
- VERMES I., BONTE H.A., SLUIS VEER G., SCHOEMAKER J. (1991)
Interpretations of five monoclonal immunoassays of Lutropin and Follitropin: effects of normalization with WHO standard.
Clin. Chem., 37:415
- DEMOULIN A., DUBOIS M., GERDAY C., GILLAIN D., LAMBOTTE R., FRANCHIMONT P. (1991)
Variation of luteinizing hormone serum concentration after exogenous human chorionic gonadotropin administration during ovarian stimulation.
Fertil. and Steril., 55:796
- MASSA G., de ZECHER F., VANDERSCHUREN-LODEWKY K. (1992)
Serum levels of immunoreactive Inhibine, FSH and LH in human infants at Preterm and Term Birth.
Biol. of the Neonat., 61:150
- DE HERTOGH R., VANKRIEKEN L., THOMAS K. de GASPARO M. (1992)
Circhoral fluctuations of serum total renin, inhibin and related hormones around the mid-cycle in normal human females.
Hum. Reprod., 7:337

XIX. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS ml	CALIBRATORS ml	SAMPLE(S) CONTROLS ml
Calibrators (0-6) Samples, Controls Tracer	- 0.05	0.1 0.05	- 0.1 0.05
Incubation	1 hour at room temperature with shaking at 700 rpm		
Separation Washing solution Separation Washing solution Separation	- - - - -	aspirate (or decant) 2.0 aspirate (or decant) 2.0 aspirate (or decant)	
Counting	Count tubes for 60 seconds		

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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XVIII. BIBLIOGRAPHIE

- CLARKE J.J., CUMMINS J.T. (1987)
Pulsatility of reproductive hormones: physiological basis and clinical implications.
Baillière's Clin. Endocrinol. Metab., 1:1
- FILICORI M., SANTORO N., MERRIAN G.R., CROWLEY W.F. Jr., (1986)
Characterization of the physiological pattern of episodic gonadotropin secretion throughout the human menstrual cycle.
J. Clin. Endocrinol. Metab., 62:1136
- MAIS V., CETEL N.S., MUSE K.N., QUIGLEY M.E., REID R.L., YEN S.S.C. (1987)
Hormonal dynamics during luteal-follicular transition.
J. Clin. Endocrinol. Metab., 64:1109
- SAWYER-STEFFAN J.E., LASLEY B.L., HOFF J.D., YEN S.S.C. (1982)
- THOMAS C.M.G., SEGERS M.F.G. (1988)
Measurement of Serum hLH: hCG interference evaluated for two hLH-Specific IRMA kits.
Clin. Chem., 34:768
- LOUMAYE E., VANKRIEKEN L., DEPREESTER S., PSALTI I., de COOMAN S., THOMAS K. (1989)
Hormonal changes induced by short-term administration of a gonadotropin-releasing hormone agonist during ovarian hyperstimulation for in vitro fertilization and their consequences for embryo development.
Fertil. and Steril., 51:105
- VERMES I., BONTE H.A., SLUIS VEER G., SCHOEMAKER J. (1991)
Interpretations of five monoclonal immunoassays of Lutropin and Follitropin: effects of normalization with WHO standard.
Clin. Chem., 37:415
- DEMOULIN A., DUBOIS M., GERDAY C., GILLAIN D., LAMBOTTE R., FRANCHIMONT P. (1991)
Variation of luteinizing hormone serum concentration after exogenous human chorionic gonadotropin administration during ovarian stimulation.
Fertil. and Steril., 55:796
- MASSA G., de ZECHER F., VANDERSCHUREN-LODEWK K. (1992)
Serum levels of immunoreactive Inhibine, FSH and LH in human infants at Preterm and Term Birth.
Biol. of the Neonat., 61:150
- DE HERTOGH R., VANKRIEKEN L., THOMAS K. de GASPARO M. (1992)
Circhoral fluctuations of serum total renin, inhibin and related hormones around the mid-cycle in normal human females.
Hum. Reprod., 7:337

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1. CLARKE J.J., CUMMINS J.T. (1987)
Pulsatility of reproductive hormones: physiological basis and clinical implications.
Bailliére's Clin. Endocrinol. Metab., 1:1
2. FILICORI M., SANTORO N., MERRIAN G.R., CROWLEY W.F. Jr., (1986)
Characterization of the physiological pattern of episodic gonadotropin secretion throughout the human menstrual cycle.
J. Clin. Endocrinol. Metab., 62:1136
3. MAIS V., CETEL N.S., MUSE K.N., QUIGLEY M.E., REID R.L., YEN S.S.C. (1987)
Hormonal dynamics during luteal-follicular transition.
J. Clin. Endocrinol. Metab., 64:1109
4. SAWYER-STEFFAN J.E., LASLEY B.L., HOFF J.D., YEN S.S.C. (1982)
Comparison of in-vitro bioactivity and immunoreactivity of serum LH in normal cyclic and hypogonadal women treated with low doses of LH-RH.
J. Reprod. Fertil., 65:45
5. THOMAS C.M.G., SEGERS M.F.G. (1988)
Measurement of Serum hLH: hCG interference evaluated for two hLH-Specific IRMA kits.
Clin. Chem., 34:768
6. LOUMAYE E., VANKRIEKEN L., DEPREESTER S., PSALTI I., de COOMAN S., THOMAS K. (1989)
Hormonal changes induced by short-term administration of a gonadotropin-releasing hormone agonist during ovarian hyperstimulation for in vitro fertilization and their consequences for embryo development.
Fertil. and Steril., 51:105
7. VERMES I., BONTE H.A., SLUIS VEER G., SCHOE MAKER J. (1991)
Interpretations of five monoclonal immunoassays of Lutropin and Follitropin: effects of normalization with WHO standard.
Clin. Chem., 37:415
8. DEMOULIN A., DUBOIS M., GERDAY C., GILLAIN D., LAMBOTTE R., FRANCHIMONT P. (1991)
Variation of luteinizing hormone serum concentration after exogenous human chorionic gonadotropin administration during ovarian stimulation.
Fertil. and Steril., 55:796
9. MASSA G., de ZECHER F., VANDERSCHUREN-LODEWYK K. (1992)
Serum levels of immunoreactive Inhibine, FSH and LH in human infants at Preterm and Term Birth.
Biol. of the Neonat., 61:150
10. DE HERTOGH R., VANKRIEKEN L., THOMAS K. de GASPARO M. (1992)
Circhoral fluctuations of serum total renin, inhibin and related hormones around the mid-cycle in normal human females.
Hum. Reprod., 7:337

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- MAIS V., CETEL N.S., MUSE K.N., QUIGLEY M.E., REID R.L., YEN S.S.C. (1987)
Hormonal dynamics during luteal-follicular transition.
J. Clin. Endocrinol. Metab., 64:1109
- SAWYER-STEFFAN J.E., LASLEY B.L., HOFF J.D., YEN S.S.C. (1982)
Comparison of in-vitro bioactivity and immunoreactivity of serum LH in normal cyclic and hypogonadal women treated with low doses of LH-RH.
J. Reprod. Fertil., 65:45
- THOMAS C.M.G., SEGERS M.F.G. (1988)
Measurement of Serum hLH: hCG interference evaluated for two hLH-Specific IRMA kits.
Clin. Chem., 34:768
- LOUMAYE E., VANKRIEKEN L., DEPREESTER S., PSALTI I., de COOMAN S., THOMAS K. (1989)
Hormonal changes induced by short-term administration of a gonadotropin-releasing hormone agonist during ovarian hyperstimulation for in vitro fertilization and their consequences for embryo development.
Fertil. and Steril., 51:105
- VERMES I., BONTE H.A., SLUIS VEER G., SCHOEMAKER J. (1991)
Interpretations of five monoclonal immunoassays of Lutropin and Follitropin: effects of normalization with WHO standard.
Clin. Chem., 37:415
- DEMOULIN A., DUBOIS M., GERDAY C., GILLAIN D., LAMBOTTE R., FRANCHIMONT P. (1991)
Variation of luteinizing hormone serum concentration after exogenous human chorionic gonadotropin administration during ovarian stimulation.
Fertil. and Steril., 55:796
- MASSA G., de ZECHER F., VANDERSCHUREN-LODEWKY K. (1992)
Serum levels of immunoreactive Inhibine, FSH and LH in human infants at Preterm and Term Birth.
Biol. of the Neonat., 61:150
- DE HERTOGH R., VANKRIEKEN L., THOMAS K. de GASPARO M. (1992)
Circhoral fluctuations of serum total renin, inhibin and related hormones around the mid-cycle in normal human females.
Hum. Reprod., 7:337

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- THOMAS C.M.G., SEGERS M.F.G. (1988)
Measurement of Serum hLH: hCG interference evaluated for two hLH-Specific IRMA kits.
Clin. Chem., 34:768
- LOUMAYE E., VANKRIEKEN L., DEPREESTER S., PSALTI I., de COOMAN S., THOMAS K. (1989)
Hormonal changes induced by short-term administration of a gonadotropin-releasing hormone agonist during ovarian hyperstimulation for in vitro fertilization and their consequences for embryo development.
Fertil. and Steril., 51:105
- VERMES I., BONTE H.A., SLUIS VEER G., SCHOOEMAKER J. (1991)
Interpretations of five monoclonal immunoassays of Lutropin and Follitropin: effects of normalization with WHO standard.
Clin. Chem., 37:415
- DEMOULIN A., DUBOIS M., GERDAY C., GILLAIN D., LAMBOTTE R., FRANCHIMONT P. (1991)
Variation of luteinizing hormone serum concentration after exogenous human chorionic gonadotropin administration during ovarian stimulation.
Fertil. and Steril., 55:796
- MASSA G., de ZECHER F., VANDERSCHUREN-LODEWYK K. (1992)
Serum levels of immunoreactive Inhibine, FSH and LH in human infants at Preterm and Term Birth.
Biol. of the Neonat., 61:150
- DE HERTOGH R., VANKRIEKEN L., THOMAS K. de GASPARO M. (1992)
Circhoral fluctuations of serum total renin, inhibin and related hormones around the mid-cycle in normal human females.
Hum. Reprod., 7:337

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1. CLARKE J.J., CUMMINS J.T. (1987) **Pulsatility of reproductive hormones: physiological basis and clinical implications.** Baillière's Clin. Endocrinol. Metab., 1:1
2. FILICORI M., SANTORO N., MERRIAN G.R., CROWLEY W.F. Jr., (1986) **Characterization of the physiological pattern of episodic gonadotropin secretion throughout the human menstrual cycle.** J. Clin. Endocrinol. Metab., 62:1136
3. MAIS V., CETEL N.S., MUSE K.N., QUIGLEY M.E., REID R.L., YEN S.S.C. (1987) **Hormonal dynamics during luteal-follicular transition.** J. Clin. Endocrinol. Metab., 64:1109
4. SAWYER-STEFFAN J.E., LASLEY B.L., HOFF J.D., YEN S.S.C. (1982) **Comparison of in-vitro bioactivity and immunoreactivity of serum LH in normal cyclic and hypogonadal women treated with low doses of LH-RH.** J. Reprod. Fertil., 65:45
5. THOMAS C.M.G., SEGERS M.F.G. (1988) **Measurement of Serum hLH: hCG interference evaluated for two hLH-Specific IRMA kits.** Clin. Chem., 34:768
6. LOUMAYE E., VANKRIEKEN L., DEPREESTER S., PSALTI I., de COOMAN S., THOMAS K. (1989) **Hormonal changes induced by short-term administration of a gonadotropin-releasing hormone agonist during ovarian hyperstimulation for in vitro fertilization and their consequences for embryo development.** Fertil. and Steril., 51:105
7. VERMES I., BONTE H.A., SLUIS VEER G., SCHOE MAKER J. (1991) **Interpretations of five monoclonal immunoassays of Lutropin and Follitropin: effects of normalization with WHO standard.** Clin. Chem., 37:415
8. DEMOULIN A., DUBOIS M., GERDAY C., GILLAIN D., LAMBOTTE R., FRANCHIMONT P. (1991) **Variation of luteinizing hormone serum concentration after exogenous human chorionic gonadotropin administration during ovarian stimulation.** Fertil. and Steril., 55:796
9. MASSA G., de ZECHER F., VANDERSCHUREN-LODEWYK K. (1992) **Serum levels of immunoreactive Inhibine, FSH and LH in human infants at Preterm and Term Birth.** Biol. of the Neonat., 61:150
10. DE HERTOGH R., VANKRIEKEN L., THOMAS K. de GASPARO M. (1992) **Circhoral fluctuations of serum total renin, inhibin and related hormones around the mid-cycle in normal human females.** Hum. Reprod., 7:337

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- CLARKE J.J., CUMMINS J.T. (1987) *Pulsatility of reproductive hormones: physiological basis and clinical implications.*

Baillière's Clin. Endocrinol. Metab., 1:1

- FILICORI M., SANTORO N., MERRIAN G.R., CROWLEY W.F. Jr., (1986) *Characterization of the physiological pattern of episodic gonadotropin secretion throughout the human menstrual cycle.* J. Clin. Endocrinol. Metab., 62:1136
- MAIS V., CETEL N.S., MUSE K.N., QUIGLEY M.E., REID R.L., YEN S.S.C. (1987) *Hormonal dynamics during luteal-follicular transition.* J. Clin. Endocrinol. Metab., 64:1109
- SAWYER-STEFFAN J.E., LASLEY B.L., HOFF J.D., YEN S.S.C. (1982) *Comparison of in-vitro bioactivity and immunoreactivity of serum LH in normal cyclic and hypogonadal women treated with low doses of LH-RH.* J. Reprod. Fertil., 65:45
- THOMAS C.M.G., SEGERS M.F.G. (1988) *Measurement of Serum hLH: hCG interference evaluated for two hLH-Specific IRMA kits.* Clin. Chem., 34:768
- LOUMAYE E., VANKRIEKEN L., DEPREESTER S., PSALTI I., de COOMAN S., THOMAS K. (1989) *Hormonal changes induced by short-term administration of a gonadotropin-releasing hormone agonist during ovarian hyperstimulation for in vitro fertilization and their consequences for embryo development.* Fertil. and Steril., 51:105
- VERMES I., BONTE H.A., SLUIS VEER G., SCHOE MAKER J. (1991) *Interpretations of five monoclonal immunoassays of Lutropin and Follitropin: effects of normalization with WHO standard.* Clin. Chem., 37:415
- DEMOULIN A., DUBOIS M., GERDAY C., GILLAIN D., LAMBOTTE R., FRANCHIMONT P. (1991) *Variation of luteinizing hormone serum concentration after exogenous human chorionic gonadotropin administration during ovarian stimulation.* Fertil. and Steril., 55:796
- MASSA G., de ZECHER F., VANDERSCHUREN-LODEWYK K. (1992) *Serum levels of immunoreactive Inhibine, FSH and LH in human infants at Preterm and Term Birth.* Biol. of the Neonat., 61:150
- DE HERTOGH R., VANKRIEKEN L., THOMAS K. de GASPARO M. (1992) *Circhoral fluctuations of serum total renin, inhibin and related hormones around the mid-cycle in normal human females.* Hum. Reprod., 7:337

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