RENIN-IRMA

KIP1531

Version : 230123

Date of issue : 23/01/2023

Revision date: 23/01/2023

History

Summary of change:

Current Version:
230123
New logo

en

Read entire protocol before use.

RENIN-IRMA

I. INTENDED USE

Immunoradiometric assay kit for the in vitro quantitative determination of active Renin in human EDTA plasma.

II. GENERAL INFORMATION

- A. Proprietary name : DIAsource Renin-IRMA Kit
- B. Catalog number : KIP1531 : 96 tests
- C. Manufactured by : DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

Renin , a polypeptidic enzyme (MW~ 40000) (1) also known as angiotensinogenase, is a circulating protease secreted by juxtaglomerular cells in the juxtaglomerular apparatus of the kidneys in response to low blood volume or low body NaCl content.

Renin activates the renin-angiotensin system by cleaving angiotensinogen produced in the liver into angiotensin I (inactive) which is further converted into angiotensin II (active) in the vascular epithelium of the lung. Angiotensin II can cause vasoconstriction by stimulating the central nervous system , in addition it stimulates ADH (antidiuretic hormone) secretion and aldosterone secretion from the adrenal gland .(6)

Regulation of blood pressure and renal glomerular filtration control (2) are the most important functions of renin -angiotensin system .

Plasmatic concentration of renin is influenced by concentration of circulating angiotensinogen and subsequently the concentration of angiotensin II. High plasmatic levels of angiotensin II reduce renin secretion.(negative feed back)

Determination of renin plasma levels is useful in the diagnosis of hypertension and in the therapeutic follow up of hypertensive patients (3).

Plasmatic concentration of renin decreases in patients with hypertension due to a primary hyperaldosteronism (4), contrary to renovascular hypertension (5) where concentrations of renin and aldosterone are both elevated.

IV. PRINCIPLES OF THE METHOD

The DIAsource Renin-IRMA is an immunoradiometric assay based on coated-tube. Mab1, the capture antibody, is attached to the lower and inner surface of the plastic tube. Calibrators or samples added to the tubes will at first show low affinity for Mab1. Addition of Mab2, the signal antibody labelled with ^{125}I , will complete the system and trigger the immunological reaction. After washing, the remaining radioactivity bound to the tube reflects the antigen concentration.

V. REAGENTS PROVIDED

Reagents	96 Test Kit	Colour Code	Reconstitution
Tubes coated with anti Renin (monoclonal antibody)	2 x 48	black	Ready for use
Anti-Renin ^{125}I (monoclonal antibody) in Phosphate buffer with bovine serum, azide (<0.1%)	1 vial 10.5 ml 760 kBq	red	Ready for use
Calibrators 0-6 in human plasma and thymol. See exact value on vial labels.	7 vials lyophilised	yellow	Add 2 ml distilled water
Wash solution (Tween 20-NaCl)	1 vial 40 ml	brown	Dilute 20 x with distilled water (use a magnetic stirrer).
Controls - N = 1 or 2 in human plasma and thymol	2 vials lyophilised	silver	Add 2 ml distilled water

Note: 1. Use the zero calibrator for sample dilution.

2. 1 pg of the calibrator preparation is equivalent to 2.2 +/- 0.2 μIU of NIBSC 68/356.

Values obtained in pg/ml must be multiplied by 2.2 to obtain results in $\mu\text{IU}/\text{ml}$ or mIU/l .

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Pipettes for delivery of: 100 μl , 300 μl , and 2 ml (the use of accurate pipettes with disposable plastic tips is recommended)
3. Plastic tubes for total counts
4. Vortex mixer
5. Tube shaker (400 rpm)
6. Magnetic stirrer
7. 5 ml automatic syringe (Cornwall type) for washing
8. Aspiration system (optional)
9. Any gamma counter capable of measuring ^{125}I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- A. **Calibrators:** Reconstitute the calibrators with 2 ml distilled water
! For a complete solubilisation : after reconstitution, let the vials 30 min on a shaker, then vortex them.
- B. **Controls:** Reconstitute the controls with 2 ml distilled water.
- C. **Working Wash solution:** Prepare an adequate volume of Working Wash solution by adding 19 volumes of distilled water to 1 volume of Wash Solution (20x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, the calibrators and controls are unstable, use them immediately after reconstitution, freeze them immediately in aliquots and keep them at -20°C for maximum 6 weeks. They are stable after 1 freeze-thaw cycle.
- Avoid subsequent freeze-thaw cycles.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Samples must be EDTA plasma.
- If the test is not run within 4 h., plasma should be aliquoted and stored at -20°C.
- Avoid subsequent freeze-thaw cycles.

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date.
Do not mix materials from different kit lots.
Bring all the reagents to room temperature prior to use.
Thoroughly mix all reagents and samples by gentle agitation or swirling.
Use a clean disposable pipette tip for addition of each different reagent and sample in order to avoid cross-contamination. High precision pipettes or automated pipetting equipment will improve the precision.
Respect the incubation times.
Prepare a calibration curve for each run, do not use data from previous runs.

B. Procedure

1. Label coated tubes in duplicate for each calibrator, sample and control. For the determination of total counts, label 2 normal tubes.
2. Briefly vortex calibrators, controls and samples and dispense 300 μl of each into the respective tubes.
3. Dispense 100 μl of ^{125}I odine labelled anti Renin into each tube, including the uncoated tubes for total counts.
4. Incubate for 180 minutes at room temperature on a tube shaker (400 rpm).
5. Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
6. Wash tubes with 2 ml Working Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
7. Wash tubes again with 2 ml Wash solution (except total counts) and aspirate (or decant).
8. Let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
9. Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. On log-log, semi logarithmic or linear graph paper plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of Renin (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
3. Read the concentration for each control and sample by interpolation on the calibration curve.
4. Computer assisted data reduction will simplify these calculations. If automatic result processing is to be used, a 4-parameter logistic function curve fitting is recommended

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

Renin-IRMA		cpm	B/T (%)	
Total count			100	
Calibrator	0 pg/ml 4 pg/ml 9 pg/ml 47 pg/ml 95 pg/ml 250 pg/ml 520 pg/ml	0 µIU/ml 8.8 µIU/ml 19.8 µIU/ml 103.4 µIU/ml 209 µIU/ml 550 µIU/ml 1144 µIU/ml	152 579 984 3826 8161 20851 55190	0.14 0.27 1.18 2.58 6.67 17.73

1 pg of the calibrator preparation is equivalent to 2.2 +/- 0.2 µIU of NIBSC 68/356

XIII. PERFORMANCE AND LIMITATIONS

A. Detection limit

Twelve zero calibrators were assayed along with a set of the other calibrators. The detection limit, defined as the apparent concentration of the average count at zero binding plus two standard deviations, was 0.78 pg/ml

B. Specificity

The Cross-reactivity of Prorenin in this Renin IRMA assay was determined by adding various concentrations of Prorenin to a plasma matrix and by measuring the apparent Renin response. Prorenin Cross-reactivity was found to be 0.3 %.

The potentially interfering effects of hemoglobin at 7.5 mg/ml and of bilirubin at 0.2 mg/ml have been evaluated. The results of this test (see the table below) show a decrease of approximatively 10% of plasma values. The recommendation is to avoid hemolized samples and bilirubin containing samples.

Sample tested	Renin value (pg/ml)	+ Human Hb at 7.5 mg/ml (pg/ml)
1.	21.2	18.5
2.	57.1	50
Sample tested	Renin value (pg/ml)	+ bilirubine at 0.2 mg/ml (pg/ml)
3..	22.6	20
4.	58.1	53.4

C. Precision

INTRA-ASSAY PRECISION

Sample	N	$\bar{X} \pm SD$ (pg/ml)	CV (%)
A	10	20.2 ± 1.7	8.5
B	10	67.7 ± 2.0	3.0

INTER-ASSAY PRECISION

Sample	N	$\bar{X} \pm SD$ (pg/ml)	CV (%)
A	6	15.5 ± 1.7	11
B	6	59.3 ± 2.4	4

SD: Standard Deviation; CV: Coefficient of variation

D. Accuracy

DILUTION TEST

Sample	Dilution	Theoretical Concent. (pg/ml)	Measured Concent. (pg/ml)
	1/4	258.6	
	1/8	129.3	125.5
	1/16	64.6	62.3
	1/32	32.3	29.4
	1/64	16.1	14.2
	1/128	8	6.7

Sample was diluted with zero calibrator

Value of undiluted sample : 1034 pg/ml

RECOVERY TEST		
Added Renin (pg/ml)	Recovered Renin (pg/ml)	Recovery (%)
11.8	10.5	89
24.2	21.4	88
57	53	93
117	96	82

E. Hook effect

A sample spiked with human Renin up to 90 000 pg/ml gives a signal above the highest calibrator concentration.

F. Reference Intervals

Normal range has been determined on EDTA plasma of 60 apparently healthy subjects, fasting, aged from 20 to 60:40 males and 20 females, without oestrogen-progesterone treatment and without treatment for hypertension. For each subject, two blood samples were taken: one after 1 hour of activity in an upright position and another after remaining in a supine position for one hour.

Range limits are fixed from 5 to 95 percentile

Upright 1.3 - 13.8 pg/ml

Supine 1.0 - 8.2 pg/ml

Each laboratory should establish its own normal range of values.

Be careful many factors can influence renin levels (age, posture, estrogen treatment, antihypertensive medication,...)

XIV. LIMITATIONS

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays.
- Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed in case of the presence of heterophilic antibodies. Carefully evaluate the results of patients suspected of having these antibodies.
- If results are not consistent with other clinical observations, additional information should be required before diagnosis.

XV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises.

XVI. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

This kit contains ^{125}I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by Europe and approved and/or FDA approved methods and found negative for HbsAg, anti HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with the local safety procedures. All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVII. BIBLIOGRAPHY

- (1) Primary structure of the human Renin gene : Hardman J.A; and al ; DNA (1984): 3(6):457-468.
- (2) Control of glomerular filtration rate by rennin-angiotensin system : Hall J. E. and al; Am.J.Physiol. (1977) : 233(5) :F366-F372
- (3) Raised aldosterone to renine ratio predicts antihypertensive efficacy of spironolactone: a prospective cohort follow-up study : Lim P.O. and al ; Br. J. Clin. Pharmacol. (1999) : 48(5) :756-760 .
- (4) Screening of primary aldosteronism :Schirpenbach C. and al ; Best Pract. Res. Clin. Endocrinol. Metab.(2006) : 20(3) : 369-384 .
- (5) Diagnostic procedure in renovascular hypertension :Distler A. and al. Clinical nephrology (1991) :36(4):174-180
- (6) Circulating and tissue angiotensin systems :Campbell D.J. ; J. Clin . Invest. (1987) :79(1) :1-6

XVIII SUMMARY OF THE PROTOCOL

	TOTAL COUNTS µl	CALIBRATORS CONTROLS µl	SAMPLE(S) µl
INCUBATION Calibrators (0 to 6), controls Samples	- -	300 -	- 300
Tracer	100	100	100
Incubation	180 minutes at room temperature with shaking at 400 rpm		
Separation Working Wash solution Separation Working Wash solution Separation	-	Aspirate (or decant) 2 ml Aspirate (or decant) 2 ml Aspirate (or decant)	
Counting	Count tubes for 60 seconds		

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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XVII. BIBLIOGRAPHIE

- (1) Primary structure of the human Renin gene : Hardman J.A; and al. ;
DNA (1984): 3(6):457-468.
- (2) Control of glomerular filtration rate by rennin-angiotensin system : Hall J. E.
and al; Am.J.Physiol. (1977) : 233(5) :F366-F372
- (3) Raised aldosterone to renine ratio predicts antihypertensive efficacy of
spironolactone: a prospective cohort follow-up study : Lim P.O. and al ;
Br. J. Clin. Pharmacol. (1999) : 48(5) :756-760 .
- (4) Screening of primary aldosteronism :Schirpenbach C. and al ;
Best Pract. Res. Clin. Endocrinol. Metab.(2006) : 20(3) : 369-384 .
- (5) Diagnostic procedure in renovascular hypertension :Distler A. and al.
Clinical nephrology (1991) :36(4):174-180
- (6) Circulating and tissue angiotensin systems :Campbell D.J. ;
J. Clin . Invest. (1987) :79(1) :1-6

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XVII. LITERATUR

- (1) Primary structure of the human Renin gene : Hardman J.A; and al ; DNA (1984): 3(6):457-468.
- (2) Control of glomerular filtration rate by rennin-angiotensin system : Hall J. E. and al; Am.J.Physiol. (1977) : 233(5) :F366-F372
- (3) Raised aldosterone to renin ratio predicts antihypertensive efficacy of spironolactone: a prospective cohort follow-up study : Lim P.O. and al ; Br. J. Clin. Pharmacol. (1999) : 48(5) :756-760 .
- (4) Screening of primary aldosteronism :Schirpenbach C. and al ; Best Pract. Res. Clin. Endocrinol. Metab.(2006) : 20(3) : 369-384 .
- (5) Diagnostic procedure in renovascular hypertension :Distler A. and al. Clinical nephrology (1991) :36(4):174-180
- (6) Circulating and tissue angiotensin systems :Campbell D.J. ; J. Clin . Invest. (1987) :79(1) :1-6

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XVII. BIBLIOGRAFIA

- (1) Primary structure of the human Renin gene : Hardman J.A; and al. ; DNA (1984): 3(6):457-468.
- (2) Control of glomerular filtration rate by rennin-angiotensin system : Hall J. E. and al; Am.J.Physiol. (1977) : 233(5) :F366-F372
- (3) Raised aldosterone to renine ratio predicts antihypertensive efficacy of spironolactone: a prospective cohort follow-up study : Lim P.O. and al ; Br. J. Clin. Pharmacol. (1999) : 48(5) :756-760 .
- (4) Screening of primary aldosteronism :Schirpenbach C. and al ; Best Pract. Res. Clin. Endocrinol. Metab.(2006) : 20(3) : 369-384 .
- (5) Diagnostic procedure in renovascular hypertension :Distler A. and al. Clinical nephrology (1991) :36(4):174-180
- (6) Circulating and tissue angiotensin systems :Campbell D.J. ; J. Clin . Invest. (1987) :79(1) :1-6

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XVII. BIBLIOGRAFIA

- (1) Primary structure of the human Renin gene : Hardman J.A; and al ; DNA (1984): 3(6):457-468.
- (2) Control of glomerular filtration rate by rennin-angiotensin system : Hall J. E. and al; Am.J.Physiol. (1977) : 233(5):F366-F372
- (3) Raised aldosterone to renine ratio predicts antihypertensive efficacy of spironolactone: a prospective cohort follow-up study : Lim P.O. and al ; Br. J. Clin. Pharmacol. (1999) : 48(5):756-760 .
- (4) Screening of primary aldosteronism :Schirpenbach C. and al ; Best Pract. Res. Clin. Endocrinol. Metab.(2006) : 20(3) : 369-384 .
- (5) Diagnostic procedure in renovascular hypertension :Distler A. and al. Clinical nephrology (1991) :36(4):174-180
- (6) Circulating and tissue angiotensin systems :Campbell D.J. ; J. Clin . Invest. (1987) :79(1) :1-6

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- (1) Primary structure of the human Renin gene : Hardman J.A; and al ; DNA (1984): 3(6):457-468.
- (2) Control of glomerular filtration rate by rennin-angiotensin system : Hall J. E. and al; Am.J.Physiol. (1977) : 233(5):F366-F372
- (3) Raised aldosterone to renine ratio predicts antihypertensive efficacy of spironolactone: a prospective cohort follow-up study : Lim P.O. and al ; Br. J. Clin. Pharmacol. (1999) : 48(5) :756-760 .
- (4) Screening of primary aldosteronism :Schirpenbach C. and al ; Best Pract. Res. Clin. Endocrinol. Metab.(2006) : 20(3) : 369-384 .
- (5) Diagnostic procedure in renovascular hypertension :Distler A. and al. Clinical nephrology (1991) :36(4):174-180
- (6) Circulating and tissue angiotensin systems :Campbell D.J. ; J. Clin . Invest. (1987) :79(1) :1-6

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- (1) Primary structure of the human Renin gene : Hardman J.A; and al ; DNA (1984): 3(6):457-468.
- (2) Control of glomerular filtration rate by rennin-angiotensin system : Hall J. E. and al; Am.J.Physiol. (1977) : 233(5) :F366-F372
- (3) Raised aldosterone to renin ratio predicts antihypertensive efficacy of spironolactone: a prospective cohort follow-up study : Lim P.O. and al ; Br. J. Clin. Pharmacol. (1999) : 48(5) :756-760 .
- (4) Screening of primary aldosteronism :Schirpenbach C. and al ; Best Pract. Res. Clin. Endocrinol. Metab.(2006) : 20(3) : 369-384 .
- (5) Diagnostic procedure in renovascular hypertension :Distler A. and al. Clinical nephrology (1991) :36(4):174-180
- (6) Circulating and tissue angiotensin systems :Campbell D.J. ; J. Clin . Invest. (1987) :79(1):1-6

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