IVD

CE

11-DESOXYCORTISOL-RIA-CT

KIPI20000

Version : 230123

Date of issue : 23/01/2023

Revision date: 23/01/2023

History

Summary of change:

| |
|-------------------------|
| Current Version: |
| 230123 |
| New logo |

11-DESOXYCORTISOL-RIA-CT

Radioimmunoassay for the Quantitative Determination of

11- Desoxycortisol (11-DOC) in Human serum

KIPI20000

FOR IN VITRO DIAGNOSTIC USE ONLY

en

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 91

1. INTENDED USE :

The 11-Desoxycortisol (Compound S) is an intermediate steroid in the glucocorticoids biosynthesis. Precursor of the cortisol, it comes from 17-hydroxyprogesterone after action of the 21-hydroxylase. This parameter is interesting for the diagnosis and follow up of treatment in case of surrenal enzymatic deficiency in 11 beta hydroxylase which is responsible for congenital surrenal hyperplasia in children and hyperandrogenics in adult women

Under hypothalamic-pituitary control via the Adrenocorticotrophic hormone (ACTH), the secretion of 11-Desoxycortisol follows a nycthemeral variation: it reaches a maximum in the morning (around 8 a.m.) and a minimal during the night (between 0 to 4 a.m.).


The measure of the 11-Desoxycortisol can be performed in serum using immunological competition methods (RIA).

The metyrapone inhibits the 11 beta hydroxylase and the conversion of the 11-Desoxycortisol in cortisol.

The metyrapone (Metopirone®) test is an indicator of the ACTH reserve.

2. PRINCIPLE OF THE METHOD :

The 11-Desoxycortisol RIA obeys the law of mass action according to the following equation :

Since the concentrations of ^{125}I - 11-DOC and coated antibodies are constant, the advancing state of the equation depends on the concentration of 11-DOC. The amount of ^{125}I - 11-DOC bound to the coated tube is inversely proportional to the concentration of 11 DOC in the sample.

Following the incubation, the tube is aspirated to remove excess unbound ^{125}I labelled 11-DOC.

Patient sample concentrations are read from a calibration curve.

3. MATERIAL PROVIDED AND STORAGE :

Stored at 2 - 8°C, the material can be used up to the expiration date printed on each label.

3.1. **2 x 48 Polystyrene tubes** (12 x 75 mm)
coated with anti-11 Doc polyclonal antibodies.

Systematically, allow the coated tubes to reach room temperature (18-25°C) before use. Store at 2-8°C.

3.2. **^{125}I - 11 DOC E analog tracer** (yellow, 52 mL) :
1 bottle, containing < 5 μCi (105 kBq) of ^{125}I -labelled 11 DOC analog in protein-based buffer containing < 0.1% NaN_3 as a preservative. Store at 2-8°C.

3.3. **11-DOC Calibrators**

5 vials, 0.5 ml each, containing 11-DOC calibrators in human serum (**see exact values on vial labels**) with < 0.1% sodium azide as a preservative. Store at 2-8°C. Can be used up to the expiration date printed on each label.

3.4. **11-DOC Zero Calibrator**

One vial, 1 ml, containing 0 ng/ml 11-DOC zero calibrator in human serum with < 0.1% sodium azide as a preservative. Store at 2-8°C. Can be used up to the expiration date printed on the label.

3.5.

Control 1 - 2

2 vials, lyophilised, 0.5 ml each, level I and level II containing low and high 11-DOC levels in human plasma with < 0.1% sodium azide as a preservative. (**see exact values on vial labels**)

Store at 2-8°C. Can be used up to the expiration date printed on each label.

Before use, reconstitute the contents of the controls with 0.5 ml of distilled water.

Concentrated Wash Buffer

1 bottle concentrated buffered solution containing sodium azide (NaN_3 < 0.1%). Pour the solution in 700 ml of distilled water.

The reconstituted washing solution is stable for 2 weeks at 2-8°C if covered with adhesive film to avoid contamination.

Note: Conversion factor: 1 ng/ml = 2.886 nmol/L

1 nmol/L = 0.3465 ng/ml

4. MATERIAL REQUIRED BUT NOT PROVIDED :

- bench surfaces protected by absorbent paper to reduce the effects of radioactive spillage.
- waste disposal containers appropriately labelled and designed as suitable for solid or liquid radioactive materials.
- manual or automated precision micropipettes for dispensing samples or reagents without cross-contamination.
- absorbent paper.
- vacuum pump connected through a trap for aspiration.
- water bath
- a gamma scintillation counter.

5. METHODOLOGY

5.1. Collection and handling of blood samples :

The blood sample may be collected into a dry tube.

After separation from the red blood cells, serum samples may be assayed immediately, within 24 hours if stored at 2 - 8°C, or later, after period up to several months if stored at -20°C. Repeating freezing and thawing must be avoided.

For samples from children of less than 2 years old, the 11-Desoxycortisol must be assayed immediately after an extraction step.

5.2. Assay procedure :

Reagents stored at 2°- 8° C. must be brought at room temperature (18-25°C) prior to use. Do not mix reagents of different lots. Label the tubes for T (« Total Counts » do not use coated tubes) calibrators, samples and control sera.

Perform the assay in duplicate. Calibrators, controls and samples must be assayed at the same time.

Each tube can only be used once.

1. Calibration curve:

Pipette 25 μl of each calibrator into the corresponding tubes.

2. Unknowns and control sera:

Pipette 25 μl of each sample or control sera into the corresponding tubes.

3. Add 500 μl of ^{125}I - 11-DOC tracer to each tube.

4. Vortex, cover and incubate 2.5 hours at $37 \pm 2^\circ\text{C}$.

5. Carefully aspirate or decant the solution of all tubes (Except total counts tubes).

6. Add 2 ml of washing solution to each tube. Aspirate or decant carefully (Except total counts tubes).

7. Add 2 ml of washing solution to each tube. Aspirate or decant carefully (Except total counts tubes).

8. Count the radioactivity fixed in each tube for at least 60 seconds.

Data processing :

Determine the mean count rate for each set of duplicate tubes.

Calculate the ratio B/B₀ as follows :

$$B/B_0 \% = [\text{CAL or Sample cpm} / B_0 (\text{CAL } 0) \text{ cpm}] \times 100$$

Draw the calibration curve by plotting the ratio B/B₀ % (linear scale) obtained for each calibrator versus its respective concentration expressed in ng/ml. 11-DOC concentrations in samples may be read directly from the calibration curve.

If a computer is used to calculate the results, the data can be fitted to the appropriate equation, 4-parameter logistic function curve fitting is the preferred method.

5.3. Example of a typical assay:

| | Content (ng/ml) | cpm 1st duplicate | cpm 2nd duplicate | Mean cpm | B/B ₀ (%) | 11- DOC (ng/ml) |
|---------------------|--------------------|----------------------|----------------------|-------------|-------------------------|-----------------------|
| Total counts | - | 63015 | 62510 | 62762 | - | - |
| CAL 0 | 0 | 31865 | 31445 | 31655 | 100 | |
| CAL 1 | 0.3 | 27429 | 28004 | 27711 | 87.6 | |
| CAL 2 | 1.5 | 20079 | 19803 | 19941 | 62.9 | |
| CAL 3 | 5 | 13246 | 13152 | 13199 | 41.7 | |
| CAL 4 | 20 | 6369 | 6509 | 6439 | 20.3 | |
| CAL 5 | 65 | 3464 | 3480 | 3472 | 11 | |
| C ₁ Low | - | 22977 | 22569 | 22773 | 72 | 1 |
| C ₂ High | - | 13711 | 13699 | 13705 | 43.3 | 4.5 |
| Sample | - | 4172 | 4242 | 4214 | 13.3 | 45.2 |

Example of a typical assay, do not use for calculations

6. PERFORMANCE CHARACTERISTICS:

6.1. Specificity

| Steroid | % Cross-reactivity |
|--|--------------------|
| 11-desoxycortisol | 100 |
| 17 α Hydroxyprogesterone | 5.6 |
| Desoxycorticosterone | 0.46 |
| Progesterone | 0.59 |
| Cortisol | 0.09 |
| Estradiol-17 β | 0.03 |
| Androstenedione, Corticosterone, Dexamethasone, Cortisone, Testosterone, DHEA-S, Aldosterone | N.D. |

6.2. Minimum detectable concentration of 11-DOC

The minimum detectable concentration has been assayed at 0.11 ng/ml and corresponds to the concentration given by two standard deviations below the mean cpm of 20 replicate determinations of the zero calibrators.

6.3. Recovery test:

When sera of known 11-DOC contents have their 11-DOC supplemented by addition of 11-DOC in equal volumes (1/1), a satisfactory correlation between theoretical and assayed 11-DOC is obtained.

| | | | | |
|---|------|------|------|------|
| Added 11-DOC (ng/ml) (1:1 in serum sample) | 1.5 | 5 | 20 | 65 |
| Theoretical (ng/ml) | 0.85 | 2.6 | 10.1 | 32.6 |
| Assayed 11-DOC (ng/ml) | 0.70 | 2.5 | 10.2 | 30.2 |
| % recovery | 82.4 | 96.2 | 101 | 92.6 |

6.4. Dilution test:

The dilution test indicates that there is immunological identity between the 11-DOC present in the sample and the 11-DOC used to calibrate the test.

| Dilution Factor | 1 | 1/2 | 1/4 | 1/8 | 1/16 | 1/32 | 1/64 |
|----------------------------|------|-------|-----|-----|------|------|------|
| Assayed 11-DOC (ng/ml) | 40.1 | 18.4 | 9.8 | 4.7 | 2.5 | 1.1 | 0.6 |
| Expected 11-DOC (ng/ml) | - | 20.05 | 10 | 5 | 2.5 | 1.25 | 0.63 |
| % Recovery | - | 92 | 98 | 94 | 100 | 88 | 96 |

6.5. Reproducibility:

| | Within assay variation 5 replicates | | Between assay variation 8 replicates | |
|--------|--|------|---|------|
| | Mean (ng/mL) | % CV | Mean (ng/mL) | % CV |
| Pool 1 | 1.1 | 6.2 | 1.95 | 8.7 |
| Pool 2 | 3.7 | 5.2 | 5.48 | 11.5 |
| Pool 3 | 28.3 | 7.7 | 36.85 | 15.1 |

7. LIMITATION OF THE PROCEDURE

- The results obtained from this or any other diagnostic kit should be used and interpreted only in the context of an overall clinical picture.
- Do not use lipemic, haemolyzed, icteric or turbid specimens.

8. EXPECTED VALUES

These values are given only for guidance; each laboratory should establish its own normal range of values.

The concentration range (2.5 to 97.5% percentile) obtained on 714 patients without stimulation : 0.43 to 7.56 ng/ml (1.24 to 21.82 nmol/L). According to literature, ICE (UK) and J.endocrinol.Invest.2003; 26:25-30, after stimulation with Metopirone® : > 200 nmol/L (> 69.3 ng/ml).

9. WARNING AND PRECAUTION

For in vitro diagnostic use only

Safety

This kit contains ¹²⁵I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

For more information, see Material Safety Data Sheet (MSDS).

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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