IVD

CE

Tg-S IRMA

R-CM-100

Version : 230123

Date of issue : 23/01/2023

Revision date: 23/01/2023

History

Summary of change:

Current Version: 230123
New logo

ENGLISH

IMMUNORADIOMETRIC ASSAY FOR THE QUANTITATIVE DETERMINATION OF THE HUMAN THYROGLOBULIN IN HUMAN SERUM AND PLASMA SENSITIVE METHOD.

R-CM-100 - 96 Determinations

FOR IN VITRO DIAGNOSTIC USE ONLY

1. CLINICAL APPLICATIONS

Thyroglobulin (hTg) is the major component of colloid in thyroid follicles. hTg is a glycoprotein of 660 kda. hTg being specific of the thyroid cell, the determination of its blood concentration is thus of great significance for the diagnosis of thyroid diseases.

In areas of endemic goiter, a majority of affected persons had substantially higher thyroglobulin levels than normal controls. However elevated thyroglobulin levels are also observed among patients with sporadic goiter in area where iodine intake is substantially higher. In both situations, these elevated values are correlated with the development of the goiter itself.

During the initial phase of subacute thyroiditis, the thyroglobulin levels are markedly elevated. Persistent high hTg levels at the end of a cure by antithyroid drug are considered predictive for a relapse of hyperthyroidism at the withdrawal of therapy. hTg determination is thus part of the monitoring of Graves' disease and its management.

2. PRINCIPLE OF THE ASSAY

This test kit is an immunoradiometric assay (IRMA) based on coated-tubes with monoclonal antibodies directed against distinct epitopes of the molecule of Tg.

Three capture antibodies are coated on the inner wall of tubes. Tg of the calibrators or of the samples is captured by these antibodies. Addition of the fourth antibody labeled with ¹²⁵Iodine completes the system, allowing the formation of a bridge between the coated antibodies and the labeled antibody.

After washing, the remaining radioactivity bound to the tubes is directly related to the concentration of Tg in the calibrators or samples.

The careful choice of the four monoclonal antibodies allows high specificity and high sensitivity and avoids excess of specificity which is sometimes reproached to immunometric assays with only two monoclonals.

3. REAGENTS PROVIDED WITH THE KIT

- The reagents are sufficient for 96 determinations.
- Store the kit and reagents at 2-8°C.
- The expiration date of each reagent is shown on the label.

Reagents	96 Tests Kit	Reconstitution		
 Tubes coated with three anti-human thyroglobulin mouse monoclonal antibodies.	2 x 48	Ready to Use		
<table border="1" data-bbox="76 1675 252 1729"> <tr> <td>Ab</td> <td>¹²⁵I</td> </tr> </table> TRACER : anti-human thyroglobulin mouse monoclonal antibody. Preservative: NaN ₃ (<0.1%),	Ab	¹²⁵ I	1 Vial 27.5 ml ±480 kBq	Ready to Use
Ab	¹²⁵ I			
[CAL]0] Calibrator Zero Preservative: NaN ₃ (<0.1%)	1 Vial 3 ml	Ready to Use		
[CAL]N] Calibrator - N = 7 Human thyroglobulin in phosphate buffer. Preservative: NaN ₃ (<0.1%)	7 Vials 1 ml	Ready to Use		

[CONTROL] Control human thyroglobulin in phosphate buffer containing proteins and protein stabilizers. Preservative: NaN ₃ (<0.1%).	1 Vial 1 ml	Ready to Use			
<table border="1" data-bbox="842 235 992 264"> <tr> <td>DIL</td> <td>BUF</td> </tr> </table> Diluent Buffer phosphate buffer solution with BSA and preservative NaN ₃ (< 0.1 %).	DIL	BUF	1 Vial 25 ml	Ready to Use	
DIL	BUF				
<table border="1" data-bbox="842 347 1088 376"> <tr> <td>RECOVERY</td> <td>SOLN</td> </tr> </table> Recovery Solution human thyroglobulin in phosphate buffer containing proteins and protein stabilizers. The concentration of thyroglobulin is variable from lot to lot. Please refer to the label of the bottle for the exact concentration. Preservatives : NaN ₃ < 0.1 %.	RECOVERY	SOLN	1 Vial 4.2 ml	Ready to Use	
RECOVERY	SOLN				
<table border="1" data-bbox="842 548 1120 577"> <tr> <td>WASH</td> <td>SOLN</td> <td>CONC</td> </tr> </table> Washing Solution (50 x concentrated) Tris-HCl buffer with detergent and preservative NaN ₃ (<0.1%).	WASH	SOLN	CONC	1 Vial 20 ml	Bring to 1000 ml with distilled water
WASH	SOLN	CONC			

- Control & Calibrators :** For the exact value, refer to the value written on the vial label. Calibrator 1 may be used only with the sensitive method.
The calibrators are calibrated against the hTg european reference preparation MRC457.
- Coated Tubes:** Unused tubes must be stored at 2-8°C , protected from moisture.
- Washing Solution (50 x concentrated):** The diluted washing solution is stable for 2 months at 2-8°C.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Plastic test tubes
- Test tube racks.
- Adjustable, automatic micropipettes with disposable tips.
- Vortex mixer.
- Graduated cylinder.
- Aspiration pump or automated washing device.
- Scintillation gamma counter.
- Distilled water.
- Orbital shaker adjustable at 150-450 rpm.

5. WARNINGS AND PRECAUTIONS FOR USERS

For *in vitro* diagnostic use.

Only experienced laboratory personnel should use this test and handling should be in agreement with GLP.

Radioactive Material - Not for Internal or External Use in Humans or Animals.

In order to obtain reproducible results, the following rules must be observed :

- Do not mix reagents of different lots.
- Do not use reagents beyond their expiry date.
- Use thoroughly clean glassware.
- Use distilled water, stored in clean containers.
- Avoid any contamination among samples; to this purpose disposable tips should be used for each sample and reagent.

Safety :

This kit contains ¹²⁵I (half-life: 60 days) ,emitting ionizing X (28 keV) and γ (35.5 keV) radiations. This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area. away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative

for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up. Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

For more information, see Material Safety Data Sheet (MSDS).

6. SPECIMEN COLLECTION

Lipemic sera may give aberrant results. Samples may be stored at 2°-8°C for up to 24 hours. For long term storage, samples should be divided in aliquots and frozen at -20°C. Avoid repeated freezing and thawing. After thawing and before use, mix samples by inversion.

Although interferences from anti-thyroglobulin autoantibodies (TgAb) have been greatly reduced by a careful choice of the monoclonal antibodies involved in the assay, it is advisable to verify the validity of the thyroglobulin result by performing a recovery test as described below

Serum or EDTA plasma provide similar results.

$$Y (\text{plasma}) = 0.99x (\text{serum}) + 0.44 \quad r = 0.99 \quad n = 20$$

7. ASSAY PROCEDURE

- Bring all reagents and samples to warm up at room temperature (18-25°C).
- Mix samples by gentle agitation before use.
- For all calibrators, a duplicate measure is recommended before use.

A. SENSITIVE ASSAY PROCEDURE

1. Prepare tubes to accommodate Calibrators, test Sera/Plasmas, Recoveries and Control in duplicate. Use non sensitized polystyrene tubes for the measurement of the Total activity.
2. Pipette 100 µl of the Calibrators, Samples, Spiked samples and Control into the appropriate coated tubes. Pipette directly to the bottom of the tubes.
3. Add 100 µl of Diluent to all these tubes, except the tubes for Total activity.
4. Mix gently the tubes on vortex and incubate for 2 hours on an orbital shaker set at 150-450 rpm.
5. At the end of the incubation, aspirate completely the contents of the tubes. Wash once the tubes with 2 ml washing solution. Aspirate completely and remove any residual moisture.
6. Add 250 µl of radioactive tracer (red colored) to all the tubes.
7. Incubate the tubes for 18-24 hours at room temperature (18-25°C).
8. Aspirate carefully the incubation mixture from all tubes except those of Total activity.
9. Wash the tubes twice with 2 ml of washing solution to all tubes except total tubes. Aspirate completely the contents of the tubes and remove any residual moisture.
10. Count the radioactivity bound to the tubes for 1 minute in a gamma counter. We suggest to control the background of the instrument before counting the assay. In order to avoid variations in the sensitivity of the system, the background must be reduced to a minimum or adjusted properly.

SCHEME OF THE ASSAY WITH SENSITIVE PROCEDURE

Tubes Reagent	Total Activity	Calibrators	Control	Samples	Recovery
Calibrators	-	100 µl	-	-	-
Control	-	-	100 µl	-	-
Samples	-	-	-	100 µl	-
Spiked samples*	-	-	-	-	100 µl
Diluent	-	100 µl	100 µl	100 µl	100 µl
- Mix on vortex and incubate for 2 h at room temperature (18-25°C) shaking (150 – 450 rpm). - Aspirate - wash 1 x 2 ml - aspirate					
Tracer	250 µl	250 µl	250 µl	250 µl	250 µl
- Incubate: 18 – 24 h at room temperature (18-25°C) - Aspirate - wash 2 x 2 ml - aspirate - Count					

B. NON-SENSITIVE ASSAY PROCEDURE

1. Prepare tubes to accommodate Calibrators, test Sera/Plasmas, Recoveries and Control in duplicate. Use non sensitized polystyrene tubes for the measurement of the Total activity.
2. Pipette 50 µl of the Calibrators, Samples, Spiked samples and Control into the appropriate coated tubes. Pipette directly to the bottom of the tubes.
Do not use the calibrator 1.
3. Add 200 µl of Diluent to all these tubes, except the tubes for Total activity.
4. Mix gently the tubes on vortex and incubate for 2 hours on an orbital shaker set at 150-450 rpm.
5. At the end of the incubation, aspirate completely the contents of the tubes. Wash once the tubes with 2 ml of washing solution. **Aspirate completely and remove any residual moisture.**
6. Add 250 µl of radioactive tracer (red colored) to all the tubes.
7. Incubate the tubes for 18-24 hours at room temperature (18-25°C).
8. Aspirate carefully the incubation mixture from all tubes except those of Total activity.
9. Wash the tubes twice with 2 ml of washing solution to all tubes except total tubes. **Aspirate completely the contents of the tubes and remove any residual moisture.**
10. Count the radioactivity bound to the tubes for 1 minute in a gamma counter. We suggest to control the background of the instrument before counting the assay. In order to avoid variations in the sensitivity of the system, the background must be reduced to a minimum or adjusted properly.

SCHEME OF THE ASSAY WITH NON-SENSITIVE PROCEDURE

Tubes Reagent	Total Activity	Calibrators	Control	Samples	Recovery
Calibrators	-	50 µl	-	-	-
Control	-	-	50 µl	-	-
Samples	-	-	-	50 µl	-
Spiked samples*	-	-	-	-	50 µl
Diluent	-	200 µl	200 µl	200 µl	200 µl
- Mix on vortex and incubate for 2 h at room temperature shaking (18-25°C) (150-450 rpm). - Aspirate - wash 1 x 2 ml - aspirate					
Tracer	250 µl	250 µl	250 µl	250 µl	250 µl
- Incubate: 18 – 24 h at room temperature (18-25°C) - Aspirate - wash: 2 x 2 ml - aspirate - Count					

C. RECOVERY TEST

1. **Preparation of the recovery samples (spiked samples*).**
Dilute 1/2 the samples with the contents of the recovery vial, for example 125 µl sample (75 µl with the non-sensitive method) + 125 µl (75 µl with the non-sensitive method) recovery solution. The thyroglobulin concentration in the recovery solution is variable from lot to lot. Please refer to the label of the bottle for the exact concentration. The so prepared samples are the "spiked samples".

2. Assay of the unspiked samples and of the spiked samples.

In the same serie, measure thyroglobulin in **both** the unspiked and spiked samples.

3. Calculation.

The recovery percentage is calculated as follow :

$$\text{Recovery (\%)} = \frac{\text{Value of the spiked sample}}{\text{(Value of the unspiked sample + concentration of the recovery solution)/2}} \times 100$$

If there is no interference, the recovery will be near 100 %. If the recovery is lower than 80 %, an interference of TgAb may be supposed and should be verified by measuring TgAb in the sample. If the recovery is higher than 120 %, the observed interference is not due to TgAb.

8. CALCULATION OF RESULTS

Draw the calibration curve on log/lin graph by plotting the B/T (%) obtained for each calibrator (y-axis) against the relative concentration (x-axis). Calculate the B/T (%) of each sample and read the concentration by interpolating on the calibration curve.

EXAMPLE OF CALCULATION

The values reported below must be considered as an example and may not be used in place of experimental data.

Description	Average cpm.	B/T (%)	hTg conc. (ng/ml)
Total Activity (T)	247160	-	-
CAL 0	110	0.00	0
CAL 1	297	0.08	0.75
CAL 2	525	0.17	1.5
CAL 3	1351	0.50	5
CAL 4	4022	1.58	15
CAL 5	12023	4.82	50
CAL 6	45802	18.49	200
CAL 7	110743	44.76	600
CONTROL	6890	3.1	27.0

INTERNAL QUALITY CONTROL

- If the results obtained for Control are not within the range specified, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises

9. REFERENCE VALUES

It is recommended that each laboratory determines its own reference interval. Values reported below are only indicative.

The normal values for Thyroglobulin have been determined on 99 samples and Thyroglobulin levels have been found **less than 50 ng/ml**.

10. PERFORMANCE OF THE ASSAY

SPECIFICITY

Samples were processed in assays which included one to one dilution with the Thyroglobulin Recovery. The percentages of recovery were calculated. This evaluation was made with 100 samples containing no anti-thyroglobulin auto-antibodies (TgAb negative) and 100 samples containing various levels of anti-thyroglobulin auto-antibodies (TgAb positive) ranging from 50 to more than 5,000 IU/ml.

The mean percentages of recovery of "TgAb negative" and "TgAb positive" are 91 and 85%, respectively. Therefore, no cross-reaction is observed with TgAb present in human serum.

However, a recovery test is recommended for each studied sample to validate the results since TgAb vary from patients in terms of quality and quantity (cfr. Section 7.C).

SENSITIVITY

Analytical sensitivity

The LoB (limit of blank) was calculated by measuring the blank several times and was calculated as the mean + 2 standard deviations of the distribution of the test values ; the LoB was calculated to be 0.18 ng/ml. The LoD (limit of detection) was calculated as the LoB + 1.645 standard deviations of a low concentration sample tested in 11 different run. The LoD was calculated to be 0.42 ng/ml.

The LoQ (limit of quantification) was calculated by testing 5 low values samples, 11 times. The LoQ was calculated to be 0.69 ng/ml.

PRECISION

Precision was evaluated upon intra- and inter-assay variability at different analyte concentrations.

Intra-assay

Serum	Mean	±	S.D.	C.V. (%)	N
1	9.4	(ng/ml)	0.2	2.4	20
2	51.9	±	0.8	1.6	20
3	112.8	±	1.7	1.5	20

Inter-assay

Serum	Mean	±	S.D.	C.V. (%)	N
1	9.2	(ng/ml)	0.2	2.2	9
2	51.6	±	1.7	3.3	9
3	215.5	±	4.8	2.2	9

ACCURACY

The accuracy of the method has been evaluated by analytical recovery and parallelism tests.

Analytical recovery

Sample	Expected (ng/ml)	Measured (ng/ml)	Recovery (%)
1	27.8	28.1	103.2
	28.2	29.6	105.0
	28.6	29.3	102.4
	30.3	31.5	104.0
	35.3	36.5	103.4
	52.8	51.4	97.3
	127.8	113.9	89.1
	327.8	327.6	99.9
2	27.0	30.9	114.4
	27.4	29.5	107.7
	27.8	29.3	105.4
	29.5	31.5	106.8
	34.5	35.5	102.9
	52.0	51.8	99.6
	127	113.8	89.6
	327	312.5	95.6

Parallelism

Sera with high analyte concentration were tested at different dilutions with the Zero Calibrator.

Dilution	Expected (ng/ml)	Measured (ng/ml)	Recovery (%)
S1 undiluted	-	161	-
1/2	80.5	88.1	109.4
1/4	40.2	44.7	111.0
1/8	20.1	24.5	121.6
1/16	10.1	12.0	119.3
1/32	5.0	6.4	127.2
S2 undiluted	-	241.5	-
1/2	120.8	121.2	100.4
1/4	60.4	64.0	106.0
1/8	30.2	33.3	110.3
1/16	15.1	15.6	103.4
1/32	7.5	8.7	115.1

HOOK EFFECT

Samples spiked with purified human thyroglobulin up to 100,000 ng/ml gave values higher than the last calibrator in both procedures.

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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