# TPO Ab RIA

***R-CO-100***

---

Version : 230123

Date of issue : 23/01/2023

Revision date: 23/01/2023

# **History**

---

## **Summary of change:**

<b>Current Version:</b>
<b>230123</b>
New logo

# ENGLISH

## RADIOIMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF ANTI-THYROPEROXIDASE AUTOANTIBODIES IN HUMAN SERUM R-CO-100 - 96 Determinations

FOR IN VITRO DIAGNOSTIC USE ONLY

### 1. CLINICAL APPLICATIONS

Thyroid autoimmunity is more frequently registered in women. Antibody prevalence in women increases with age, rising from approximately 10% at the age of 18-24 up to 30% at the age of 55-65 for TgAb and from 15% at the age of 18-24 up to 24% at the age of 55-65 for TPOAb.

### 2. PRINCIPLE OF THE ASSAY

The present method is based on a competitive radioimmunoassay (RIA). During the incubation, the monoclonal anti-TPO antibodies in solid phase competes with the sample/calibrator anti-TPO autoantibodies for the specific sites of the  $^{125}\text{I}$  labelled recombinant TPO (tracer). After aspiration and washing, the radioactivity in the tubes is measured in a gamma counter. The degree of binding will be inversely proportional to the sample anti-TPO autoantibodies concentration. The calibrators of anti-TPO auto-antibodies are calibrated on the TMAb international reference preparation MRC66/387.

### 3. REAGENTS PROVIDED WITH THE KIT

- The reagents are sufficient for 96 determinations.
- Store the kit and reagents at 2-8°C.
- The expiration date of each reagent is shown on the label.

Reagents	96 Tests Kit	Colour Code	Reconstitution
 Tubes coated with mouse monoclonal anti-TPO antibodies	2 x 48	Green	<b>Ready to Use</b>
Ag $^{125}\text{I}$ TRACER : Recombinant TPO tracer. Preservative: NaN3 (<0.1%)	1 Vial 22 ml $\pm 114 \text{ kBq}$	Red	<b>Ready to Use</b>
[CAL]0 Calibrator zero Preservative: NaN3 (<0.1%)	1 Vial 4.2 ml	Yellow	<b>Ready to Use</b>
[CAL]N Calibrator – N = 5 Anti-TPO (human serum). Preservative: NaN3 (<0.1%) (see exact values on vial labels)	5 Vials 0.7 ml	Yellow	<b>Ready to Use</b>
[CONTROL N] Control - N = 2 Anti-TPO (human serum). Preservative: NaN3 (<0.1%). (see exact values on vial labels)	2 Vials 0.7 ml	Silver	<b>Ready to Use</b>
WASH SOLN CONC Washing Solution (50 x concentrated) Tris-HCl buffer with detergent and preservative NaN3 (<0.1%).	1 Vial 20 ml	Brown	Bring to 1000 ml with distilled water

- **Controls & Calibrators :** For the exact value, refer to the value written on the vial label.
- **Coated Tubes:** Unused tubes must be stored at 2-8°C , protected from moisture.
- **Washing Solution (50 x concentrated):** The diluted washing solution is stable for 2 months at 2-8°C.

### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Plastic test tubes
- Test tube racks.
- Adjustable, automatic micropipettes with disposable tips.
- Vortex mixer.
- Graduated cylinder.
- Aspiration pump or automated washing device.
- Scintillation gamma counter.
- Distilled water.
- Orbital shaker adjustable at 150 rpm.

### 5. WARNINGS AND PRECAUTIONS FOR USERS

#### In vitro diagnostic use.

Only experienced laboratory personnel should use this test and handling should be in agreement with GLP.

Radioactive Material - Not for Internal or External Use in Humans or Animals.

In order to obtain reproducible results, the following rules must be observed :

- Do not mix reagents of different lots.
- Do not use reagents beyond their expiry date.
- Use thoroughly clean glassware.
- Use distilled water, stored in clean containers.
- Avoid any contamination among samples; to this purpose disposable tips should be used for each sample and reagent.

#### Safety :

This kit contains  $^{125}\text{I}$  (half-life: 60 days) ,emitting ionizing X (28 keV) and  $\gamma$  (35.5 keV) radiations. This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up. Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

For more information, see Material Safety Data Sheet (MSDS)

## 6. SPECIMEN COLLECTION

It is recommended to use serum. Highly lipemic or hemolyzed samples must be discarded. Keep samples at 2-8°C for 1 day; for longer periods it is advisable to freeze samples in aliquots at -20°C. Repeated freezing and thawing of samples should be avoided.

## 7. ASSAY PROCEDURE

- Bring all reagents and samples to warm up at room temperature (18-25°C).
- Mix samples by gentle agitation before use.
- For all calibrators, a duplicate measurement is recommended
- Each well can only be used once.

1. Prepare plain tubes for Total Counts, and coated tubes for Calibrators, Samples and Controls.
2. Pipette **20 µl** of each Calibrator, Control Serum and Sample into the corresponding tubes.
3. Add **200 µl** of Radioactive Tracer into all tubes.
4. Incubate for **90 minutes** at room temperature (18-25°C) on an orbital shaker set at 150 rpm.
5. Aspirate the contents of each tube, except the tubes for Total Counts.
6. Add **2 ml** of diluted washing solution to every tube, except tubes for Total Counts and aspirate thoroughly or decant the contents of all tubes on absorbent paper.
7. Count the radioactivity bound to the tubes for 1 minute in a gamma counter. We suggest to control the background of the instrument before counting the assay. In order to avoid variations in the sensitivity of the system, the background must be reduced to a minimum or adjusted properly.

### Note

If high anti-TPO values (> 1000 IU/ml) are expected for any patient sample, the original sample should be diluted. This dilution should be done using the zero calibrator.

**Do not use any other buffers for this purpose.**

## ASSAY SCHEME

Tubes Reagent	Total Activity	Calibrators	Controls	Samples
Calibrators	----	20 µl	----	----
Controls	----	----	20 µl	----
Samples	----	----	----	20 µl
Tracer	200 µl	200 µl	200 µl	200 µl
- Incubate: 90 min R.T. (18-25°C) shaking (150 rpm)				
- Aspirate - wash: 1 x 2 ml - aspirate				
- Count				

## 8. CALCULATION OF RESULTS

$$\text{Binding capacity} = \frac{B_0}{T} \times 100 = \frac{B_0 \text{ Cpm}}{T \text{ Cpm}} \times 100$$

Percent binding for Calibrators, Controls and Samples =

$$\text{Binding (\%)} = \frac{B}{B_0} \times 100 = \frac{\text{Calibrator or Sample Cpm}}{B_0 \text{ Cpm}} \times 100$$

Draw the calibration curve on log/logit graph by plotting the  $B/B_0 (\%)$  obtained for each calibrator (y-axis) against the relative concentration (x-axis). Calculate the  $B/B_0 (\%)$  of each sample and read the TPO Ab concentration in IU/ml, by interpolating on the calibration curve.

### EXAMPLE OF CALCULATION

The values reported below must be considered as an example and may not be used in place of experimental data.

Description	Average cpm.	B/B0 (%)	TPO Ab conc. (IU/ml)
Total Activity (T)	46591	-	-
CAL 0	15003	100	0
CAL 1	10978	73.2	40
CAL 2	6798	45.3	80
CAL 3	4408	29.4	150
CAL 4	1693	11.3	400
CAL 5	606	4.0	1000
CONTROL 1	5481	36.5	108
CONTROL 2	1926	12.8	353
P1	10736	71.6	40.7
P2	4085	27.2	157
P3	1037	6.9	629

### INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises

## 9. REFERENCE VALUES

It is recommended that each laboratory determines its own reference interval. Values reported below are only indicative.

Lower than 60 IU/ml	Negative for anti-TPO auto-antibodies.
Between 60-80 IU/ml	Borderline
Higher than 80 IU/ml	Positive for anti-TPO auto-antibodies.

## 10. PERFORMANCE OF THE ASSAY

### SPECIFICITY

No cross-reactions have been observed with anti-thyroglobulin autoantibodies present in human serum.

### SENSITIVITY

#### Analytical sensitivity

The sensitivity was calculated based upon the calibration curve and expressed as the minimal dose showing a significant difference from the Zero Calibrator (mean value - 2 S.D.). This dose is 7.4 IU/ml.

#### Functional sensitivity

The functional assay sensitivity is the lowest value which is measured with a precision of max. 20% inter-assay variance. For the anti-TPO autoantibodies, this value is lower than 20 IU/ml.

### PRECISION

Precision was evaluated upon intra- and inter-assay variability, at different analyte concentrations.

#### Intra-assay

Serum	Mean	± (IU/ml)	S.D.	C.V. (%)	N
1	40.1	± 4.4	4.4	11.0	95
2	105	± 5.5	5.5	5.2	19
3	553	± 43.8	43.8	7.9	19

#### Inter-assay

Serum	Mean	± (IU/ml)	S.D.	C.V. (%)	N
1	36.7	± 4.2	4.2	11.4	8
2	103	± 8.1	8.1	7.9	9
3	720	± 79.1	79.1	11.0	9

### ACCURACY

Accuracy of the method has been checked by the recovery and parallelism tests:

**Recovery Test**

Samples, mixed with equal volumes of each calibrator, were tested.

	<b>Expected (IU/ml)</b>	<b>Measured (IU/ml)</b>	<b>Recovery (%)</b>
S1	-	23.7	-
S1 + CAL 1	31.9	28.1	88.1
S1 + CAL 2	51.9	49.1	94.6
S1 + CAL 3	86.9	81.2	93.4
S1 + CAL 4	211.9	243	114.7
S2	-	32.9	-
S2 + CAL 1	36.5	28.9	79.2
S2 + CAL 2	56.5	45.9	81.2
S2 + CAL 3	91.5	76.2	83.3
S2 + CAL 4	216.5	240	110.9

**Parallelism Test**

Serums with high analyte concentration were tested at different dilutions with the Zero Calibrator.

<b>Dilution</b>	<b>Expected (IU/ml)</b>	<b>Measured (IU/ml)</b>	<b>Recovery (%)</b>
S1 undiluted	-	857	-
1/2	428	418	97.7
1/4	214	248	115.9
1/8	107	126	117.8
1/16	53.5	64	119.6
S2 undiluted	-	273	-
1/2	136	136	100.0
1/4	68.1	60.0	88.1
1/8	34.1	39.8	116.7

NB. Due to the heterogeneity of autoantibodies, for some patient samples a nonlinear dilution is possible

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com>

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1. Prentice L.M. et al. Geographical distribution of subclinical autoimmune thyroid disease in Britain : A study using highly sensitive direct assays for autoantibodies to thyroglobulin and thyroid peroxidase. *Acta Endocrinologica*, 1990, 123, 493-8.
2. Jaume J.C. et al. Thyroid peroxidase autoantibody epitopic « fingerprints » in juvenile Hashimoto's thyroiditis : evidence for conservation over time and in families. *Clin. Exp. Immunol.*, 1996, 104, 115-123.
3. Czarnocka B. et al. Immunochemical properties of hTPO. *Thyroperoxidase and Thyroid Autoimmunity*, Ed. P. Carayon, J. Ruf. 1990, 207, 59-67.
4. Czarnocka B. et al. Majority of thyroid peroxidase autoantibodies in patients with autoimmune thyroid disease are directed to a single TPO domain. *Autoimmunity*, 1996, 23, 145-154.
5. Takamatsu J. et al. Correlation of antithyroglobulin and antithyroid-peroxidase antibody profiles with clinical and ultrasonic characteristics of chronic thyroiditis. *Thyroid*, 1998, 8, 1101-6.
6. Smyth P. et al. Serum thyroid peroxidase autoantibodies, thyroid volume, and outcome in breast carcinoma. *J. Endocrinol. Metab.*, 1998, 83, 2711-6.
7. Gauna A. et al. Immunological aspects of Graves' disease patients in different clinical stages. *J. Endocrinol. Invest.*, 1989, 12, 671-7.
8. Stagnaro-Green A. et al. Detection of at-risk pregnancy by means of highly sensitive assays for thyroid autoantibodies. *JAMA*, 1990, 264, 1422-5.

Physikalisches Daten of  $^{125}\text{I}$

Physical characteristics of  $^{125}\text{I}$

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