# **Cortisol Saliva Elisa**

***KAPDB290***

# History

---

Summary of change :

<b>Previous Version :</b> 220421	<b>Current Version :</b> 230629
Old Diasource logo	New DiaSource logo on the front page

# Cortisol Saliva Elisa

en

For the quantitative determination of Cortisol by enzyme immunoassay in human saliva.

**KAPDB290**

**IN VITRO DIAGNOSTIC**

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 90

## INTENDED USE

For the quantitative determination of Cortisol by enzyme immunoassay in human saliva.

For *in vitro* diagnostic use only.

## PRINCIPLE OF THE TEST

The cortisol saliva ELISA is a competitive immunoassay. Competition occurs between an unlabeled antigen present in calibrators, controls and patient samples and an enzyme-labelled antigen (conjugate) for a limiting number of anti-cortisol antibody binding sites on the microplate wells. After a washing step that removes unbound materials, the enzyme substrate is added and approximately 20 minutes later the enzymatic reaction is terminated by addition of stopping solution. The resulting optical density (OD), measured with a microplate reader, is inversely proportional to the concentration of cortisol in the sample. A calibrator curve is plotted with a provided set of calibrators to calculate directly the concentration of cortisol in patient samples and controls.

## CLINICAL APPLICATIONS

Cortisol is the most abundant circulating steroid and the major glucocorticoid secreted by the adrenal cortex. Cortisol is physiologically effective in blood pressure maintenance and anti-inflammatory activity. It is also involved in calcium absorption, gluconeogenesis as well as the secretion of gastric acid and pepsin. The levels of cortisol are increased under stress situations, physical exercise and external administration of ACTH.

Most circulating cortisol is bound to corticosteroid-binding globulin (transporting) and albumin. The amount of unbound or free cortisol (which is considered the active fraction in blood), represents approximately 1–2% of the total amount in blood.

The measurement of salivary cortisol is considered to represent the free fraction, due to the absence of appreciable amounts of cortisol-binding proteins in saliva. The level of salivary cortisol shows a diurnal rhythm with the highest levels in the morning and the lowest levels at night. Studies consistently report high correlations between serum and salivary cortisol, indicating that salivary cortisol levels reliably estimate serum cortisol levels.

The measurement of salivary cortisol levels can be used as an indicator of adrenal function and the differential diagnosis of Addison's and Cushing's diseases.

## PROCEDURAL CAUTIONS AND WARNINGS

1. This kit is intended for *in vitro* use only.
2. Practice good laboratory practices when handling kit reagents and specimens. This includes:
  - Do not pipette by mouth.
  - Do not smoke, drink, or eat in areas where specimens or kit reagents are handled.
  - Wear protective clothing and disposable gloves.
  - Wash hands thoroughly after performing the test.
  - Avoid contact with eyes; use safety glasses; in case of contact with eyes, flush eyes with water immediately and contact a doctor.
3. Users should have a thorough understanding of this protocol for the successful use of this kit. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
4. Avoid microbial contamination of reagents.
5. A calibrator curve must be established for every run.
6. It is recommended to all customers to prepare their own control materials or saliva pools which should be included in every run at a high and low level for assessing the reliability of results.
7. The controls (included in kit) must be included in every run and their results must fall within the ranges stated in the quality control certificate; a failed control result might indicate improper procedural techniques or pipetting, incomplete washing or improper reagent storage.
8. When the use of water is specified for dilution or reconstitution, use deionized or distilled water.
9. All kit reagents and specimens should be brought to room temperature and mixed gently but thoroughly before use.
10. When reading the microplate, the presence of bubbles in the wells will affect the optical densities (ODs). Carefully remove any bubbles before performing the reading step.

11. The substrate solution (TMB) is sensitive to light and should remain colourless if properly stored. Instability or contamination may be indicated by the development of a blue colour, in which case it should not be used.
12. When dispensing the substrate and stopping solutions, do not use pipettes in which these liquids will come into contact with any metal parts.
13. To prevent contamination of reagents, use a new disposable pipette tip for dispensing each reagent, sample, calibrator and control.
14. Do not use kit components from different kit lots within a test and do not use any component beyond the expiration date printed on the label.
15. Kit reagents must be regarded as hazardous waste and disposed of according to local and/or national regulations.

## LIMITATIONS

1. The kit is calibrated for the determination of cortisol in human saliva. The kit is not calibrated for the determination of cortisol in other specimens of human or animal origin.
2. Do not use blood contaminated saliva samples.
3. Patients being treated with prednisone or prednisolone should not be tested with this cortisol saliva test until at least 24 hours after the last treatment.
4. Samples or controls containing azide or thimerosal are not compatible with this kit, they may lead to false results.
5. Only calibrator 0 may be used to dilute high saliva samples. The use of any other reagent will lead to false results.
6. The results obtained with this kit shall never be used as the sole basis for a clinical diagnosis. For example, some drugs and the occurrence of heterophilic antibodies in patients regularly exposed to animals or animal products have the potential of causing interferences in immunological tests. Consequently, the clinical diagnosis should comprise all aspects of a patient's background including the frequency of exposure to animals/products.

## SAFETY CAUTIONS AND WARNINGS POTENTIAL BIOHAZARDOUS MATERIAL

The reagents should be considered a potential biohazard and handled with the same precautions applied to saliva specimens. All human specimens should be considered a potential biohazard and handled as if capable of transmitting infections and in accordance with good laboratory practices.

## CHEMICAL HAZARDS

Avoid contact with reagents containing TMB, hydrogen peroxide and sulfuric acid. If contacted with any of these reagents, wash with plenty of water. TMB is a suspected carcinogen.

## SPECIMEN COLLECTION, PRE-TREATMENT AND STORAGE

Avoid sample collection within 1 hour after eating a major meal or within 12 hours after consuming alcohol. Acidic or high sugar foods can compromise assay performance by lowering sample pH and influencing bacterial growth. To minimize these factors, rinse mouth thoroughly with water 10 minutes before the sample is collected. Do not use blood-contaminated specimens.

### Specimen Collection

Approximately 0.1 mL of saliva is required per duplicate determination. Rinse mouth thoroughly with water 10 minutes before the sample is collected. Collect 1–2 mL of saliva into a clean polypropylene tube without force or inducement.

### Specimen Pre-Treatment

Following collection, the sample must be pretreated according to the following procedure:

1. Freeze the sample for a minimum of 2 hours.
2. Thaw the sample.
3. Vortex to mix and centrifuge the sample at 2000x g for 10 minutes.
4. Carefully remove the supernatant and transfer to a new labeled tube. The supernatant will be used in the assay procedure of the test.

## Specimen Storage

Store pretreated samples at 4°C for up to 24 hours or freeze at or below -20°C for up to 6 months. Samples that have been stored should be inspected to ensure they are free from precipitates before being used in the assay. If there are precipitates present, follow steps 3–4 in the specimen pretreatment section. Consider all human specimens as possible biohazardous materials.

## REAGENTS AND EQUIPMENT NEEDED BUT NOT PROVIDED

1. Precision pipette to dispense 25, 50, 100, 150, 300 and 1000 µL
2. Disposable pipette tips
3. Distilled or deionized water
4. Benchtop centrifuge
5. Vortex
6. Microplate reader with a filter set at 450 nm and an upper OD limit of 3.0 or greater
7. Microplate washer (recommended)

## REAGENTS PROVIDED

**Anti-Cortisol Antibody-Coated Break-Apart Well Microplate** - Ready To Use.

Contents: One 96 well (12x8) polyclonal antibody-coated microwell plate in a resealable pouch with desiccant.

Storage: Refrigerate at 2–8°C

**Cortisol Horseradish Peroxidase (HRP) Conjugate Concentrate** – Requires Preparation **X50**

Contents: Cortisol-HRP conjugate in a protein-based buffer with a non-mercury preservative.

Volume: 450 µL/vial

Storage: Refrigerate at 2–8°C

Preparation of working conjugate: Dilute 1:50 in assay buffer before use (eg. 40 µL of conjugate concentrate in 2 mL of assay buffer). If the whole plate is to be used dilute 340 µL of conjugate concentrate in 17 mL of assay buffer. Discard any that is left over.

**Cortisol Saliva Calibrators** - Ready To Use. N = 0 to 5

Contents: Six vials containing cortisol in a protein-based buffer with a non-mercury preservative. Prepared by spiking buffer with a defined quantity of Cortisol.

\*Listed below are approximate concentrations, please refer to vial labels for exact concentrations.

Calibrator	Concentration ng/mL	Volume/Vial
Calibrator 0	0	2.0
Calibrator 1	0.1	1.0
Calibrator 2	0.5	1.0
Calibrator 3	2.5	1.0
Calibrator 4	10	1.0
Calibrator 5	50	1.0

Storage: Refrigerate at 2–8°C

Stability: Once opened, calibrators should be used within 14 days or aliquoted and stored frozen. Avoid multiple freezing and thawing cycles.

**Controls** - Ready To Use.

Contents: Two vials containing cortisol in a protein-based buffer with a non-mercury preservative. Prepared by spiking buffer with a defined quantity of cortisol. Refer to vial labels for the acceptable ranges.

Volume: 1.0 mL/vial

Storage: Refrigerate at 2–8°C

Stability: Once opened, the controls should be used within 14 days or aliquoted and stored frozen. Avoid multiple freezing and thawing cycles.

**Wash Buffer Concentrate** – **X10**

Contents: One bottle containing buffer with a non-ionic detergent and a non-mercury preservative.

Volume: 50 mL/bottle

Storage: Refrigerate at 2–8°C

Preparation of working wash buffer: Dilute 1:10 in distilled or deionized water before use. If the whole plate is to be used dilute 50 mL of the wash buffer concentrate in 450 mL of water.

**Assay Buffer** - Ready To Use.

Contents: One bottle containing a protein-based buffer with a non-mercury preservative.

Volume: 20 mL/bottle

Storage: Refrigerate at 2–8°C

**TMB Substrate** - Ready To Use.

Contents: One bottle containing tetramethylbenzidine and hydrogen peroxide in a non-DMF or DMSO containing buffer.

Volume: 16 mL/bottle

Storage: Refrigerate at 2–8°C

**Stopping Solution** - Ready To Use.

Contents: One bottle containing 1M sulfuric acid.

Volume: 6 mL/bottle

Storage: Refrigerate at 2–8°C

## ASSAY PROCEDURE

Specimen Pretreatment: See Specimen Collection, Pre-Treatment and Storage Section.

All kit components, controls and specimen samples must reach room temperature prior to use. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.

1. After all kit components have reached room temperature, mix gently by inversion. Prepare the working conjugate and working wash buffer (see Cortisol-Horseradish Peroxidase (HRP) Conjugate Concentrate and Wash Buffer Concentrate reagents in the REAGENTS PROVIDED section).
2. Remove the required number of strips from the microplate and assemble into a plate frame. Reseal the bag and return any unused strips to the refrigerator.
3. Pipette 25 µL of each calibrator, control and specimen sample into correspondingly labelled microplate wells.
4. Pipette 150 µL of the Cortisol-HRP working conjugate solution into each microplate well (the use of a multi-channel pipette is recommended).
5. Gently tap the microplate frame for 10 seconds to mix the contents of the wells and incubate the microplate at room temperature (no shaking) for 45 minutes.
6. Wash the microplate wells 3 times with working wash buffer (350 µL/well for each wash) and tap the plate firmly against absorbent paper to ensure that it is dry. The use of a microplate washer is highly recommended. If a microplate washer is not available, ensure that the wash buffer reaches the top edge of the wells and that no liquid remains in the microplate after the final washing, avoid splashing.
7. Pipette 150 µL of TMB substrate into each microplate well at timed intervals.
8. Gently tap the microplate frame for 10 seconds to mix the contents of the wells and incubate the microplate at room temperature (no shaking) for 15–20 minutes.
9. Pipette 50 µL of stopping solution into each microplate well at the same timed intervals as in step 7 and gently tap the microplate frame to mix the contents of the wells.
10. Read the microplate in a microplate reader at 450 nm, within 20 minutes after addition of the stopping solution.

## CALCULATIONS

1. Calculate the mean optical density of each calibrator, control and sample duplicate.
2. Use a 4-parameter or 5-parameter curve fit with immunoassay software to generate a calibrator curve.
3. Read the values of the unknowns directly off the calibrator curve.


4. If a sample reads more than 50 ng/mL dilute it with calibrator 0 not more than 10-fold. The result obtained must be multiplied by the dilution factor.

#### TYPICAL TABULATED DATA

Calibrator	Mean OD 450 nm	% Binding	Value (ng/mL)
0	2.555	100	0
1	2.219	87	0.1
2	1.695	66	0.5
3	0.963	38	2.5
4	0.498	20	10
5	0.204	8	50
Unknown	0.988	—	2.4

#### TYPICAL CALIBRATION CURVE

Sample curve only. **Do not** use to calculate results.

#### PERFORMANCE CHARACTERISTICS

##### SENSITIVITY

The lower detection limit is calculated from the standard curve by determining the resulting concentration of the mean OD of Calibrator 0 (based on 10 replicate analyses) minus 2 SD. Therefore, the sensitivity of the Cortisol Saliva ELISA kit is 0.033 ng/mL.

##### SPECIFICITY (CROSS REACTIVITY)

The following compounds were tested for cross-reactivity with cortisol cross-reacting at 100%.

Compound	% Cross Reactivity
Cortisol	100
Progesterone	1.9
Corticosterone	1.4
11-Deoxycorticosterone	< 0.04
Cortisone	6.3
Prednisone	3.6
DHEAS	< LoD
Prednisolone	17.8

Prednisolone does cross react with most assays (~30% with most commercially available assays). Prednisolone has a half-life of ~3.5 hours and needs 16–18 hours to be adequately cleared from the circulation.

Since prednisone is converted to prednisolone in vivo, caution must be exercised when assaying the cortisol levels of patients undergoing either therapy. Patients being treated with prednisone or prednisolone should not be tested with this cortisol saliva test until at least 24 hours after the last treatment.


##### PRECISION

The experimental protocol used a nested components-of-variance design with 10 testing days, two lots and two scientists per day. Each scientist ran two tests with two lots per day and two replicate measurements per run (a 10 x 2 x 2 x 2 design) for each sample. The results were analyzed with a two-way nested ANOVA and summarized in the table below.

Sample	Mean	Within Run SD	Within Run CV	Total SD	Total CV
1	1.870	0.213	11.4%	0.218	11.7%
2	6.637	0.462	7.0%	0.589	8.9%
3	10.196	0.834	8.2%	1.061	10.4%
4	18.756	1.383	7.4%	1.763	9.4%
5	3.035	0.287	9.4%	0.291	9.6%
6	0.679	0.114	16.9%	0.118	17.4%
7	24.800	2.103	8.5%	2.567	10.4%
8	0.356	0.056	15.8%	0.058	16.4%
9	2.309	0.214	9.2%	0.228	9.9%
10	12.939	0.991	7.7%	1.090	8.4%

#### LINEARITY

The linearity study was performed with five human saliva samples with concentrations covering the range of the assay and following CLSI guideline EP06-A. The samples were diluted in calibrator 0 up to ten-fold (1:10), tested in duplicate, and the results (y) compared to the predicted concentration. The statistical analysis shows that the assay is sufficiently linear throughout the dynamic range of the kit up to a 1:10 dilution, when using calibrator 0 as the diluent.

#### COMPARATIVE STUDIES

The DIIAsource Cortisol Saliva ELISA kit assay (y) was compared to a commercial Cortisol Saliva ELISA assay (x). The comparison of 47 saliva samples yielded the following linear regression results:  $y = 0.844x + 0.634$ ,  $R^2 = 0.929$

#### REFERENCE RANGES

Reference ranges (95%) were established using samples obtained from individuals of diverse races. Each laboratory shall establish their own range of reference values.

Group	N	Median (ng/mL)	95% Range (ng/mL)	95% Range (ng/mL)	Total Range (ng/mL)	Total Range (ng/mL)
Males AM	82	1.638	0.21–6.33	AM range (ng/mL) 0.21–7.27	0.07–7.63	AM range (ng/mL) 0.07–7.90
Females AM	78	2.265	0.23–7.27		0.113–7.90	
Males PM	76	0.897	0.16–3.42	PM range (ng/mL) 0.074–3.42	0.045–4.34	PM range (ng/mL) 0.02–4.34
Females PM	77	0.396	0.074–2.38		0.02–3.78	

## LITERATURE

7. Rose JQ, Yurchak AM, Jusko JW. Dose dependent pharmacokinetics of prednisone and prednisolone in man. *J Pharmacokinet Biopharm.* 1981; 9:389–417.
8. Vining RF, McGinley RA. The measurement of hormones in saliva: possibilities and pitfalls. *J Steroid Biochem.* 1987; 27:81–94.
9. Aardal E, Holm, AC. Cortisol in saliva – Reference ranges in relation to cortisol in serum. *Eur J Clin Chem Clin Biochem.* 1995; 33(12):927–32.
10. Nicholson N, Storms C, Ponds R, Sulon J. Salivary Cortisol Levels and Stress Reactivity in Human Aging . *J Gerontol A Biol Sci Med Sci.* 1997; 52:68–75
11. Yaneva M, et al. Midnight Salivary Cortisol for the Initial Diagnosis of Cushing's Syndrome of Various Cause. *J Clin Endocrinol Metab.* 2004; 89(7):3345–51.
12. Dorn LD, Lucke JF, Loucks TL, Berga SL. Salivary cortisol reflects serum cortisol: Analysis of circadian profiles. *Ann Clin Biochem.* 2007; 44(pt 3):281–4.
13. Dickmeis T. Glucocorticoids and the circadian clock. *J Endocrinol.* 2009; 200(1):3–22.
14. VanBruggen MD, et al. The Relationship Between Serum and Salivary Cortisol Levels in Response to Different Intensities of Exercise. *Int J Sports Physiol Perform.* 2011; 6(3):396–407.
15. Prednisolone, Dexamethasone and Cortisol assays. Biochemical Investigations in Laboratory Medicine. [http://www.pathology.leedsth.nhs.uk/dnn\\_bilm/Misc/Syntheticglucocorticoidsandcortisolassays.aspx](http://www.pathology.leedsth.nhs.uk/dnn_bilm/Misc/Syntheticglucocorticoidsandcortisolassays.aspx) Published 17 Nov 2009. Accessed 21 Aug 2019.

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