IVD

CE

17 α -HYDROXYPROGESTERONE (17-OHP)-RIA-CT

KIP1409

Version : 230911

Date of issue : 11/09/2023

Revision date: 11/09/2023

History

Summary of change:

Current Version: 230123	Current Version: 230911
	Correction in the Korean translation only : Addition of the Extraction protocol for newborn samples Other languages remain unchanged

Read entire protocol before use.

17 α -HYDROXYPROGESTERONE (17-OHP)-RIA-CT

I. INTENDED USE

Radioimmunoassay for the *in vitro* quantitative measurement of human 17 α -hydroxyprogesterone (17-OHP) in serum and plasma.

II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource 17 α -HYDROXYPROGESTERONE (17-OHP)-RIA-CT Kit
- B. **Catalog number :** KIP1409 : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel: +32 (0) 10 84.99.11 Fax: +32 (0) 10 84.99.91

III. CLINICAL BACKGROUND

A. *Biological activities of 17-hydroxyprogesterone*

17 α -hydroxyprogesterone (17-OHP) is a C-21 steroid hormone (molecular weight 330.3) which is produced from 17 α -hydroxypregnenolone in the adrenals and also in the ovaries, testes and placenta. 17-OHP is hydroxylated at positions 11 and 21 to produce cortisol via 11-deoxycortisol.

B. *Clinical applications of 17 α -hydroxyprogesterone determination*



As a rule, serum or amniotic fluid 17-OHP dosages are relevant to diagnose congenital adrenal hyperplasia (CAH). This CAH is due to a specific enzyme defect (six distinct enzyme deficiencies have been described). As a result of these deficiencies, ACTH increases and produces adrenal hyperplasia and the raise of many steroid precursors. But it is also very interesting to know the value of 17-OHP in patients with varicocele. (17-OHP and testosterone represent markers of Leydig cell function) and in ageing male patients to detect Benign Prostatic Hypertrophy (BPH) and carcinoma of the prostate (PCA) (plasma 17-OHP is significantly lower in PCA and BPH groups than in normal men).

There are other domains for 17-OHP investigations as : male infertility, girls with peripubertal virilization, children with premature adrenache (in these cases, the values of 17-OHP are increased without or after ACTH stimulation).

IV. PRINCIPLES OF THE METHOD

A fixed amount of ¹²⁵I labelled steroid competes with the steroid to be measured present in the sample or in the calibrator for a fixed amount of antibody sites immobilized on the wall of a polystyrene tube. Neither extraction nor chromatography is required because of the high specificity of the coated antibodies. After 3 hours incubation at 37°C, an aspiration step terminates the competition reaction. The tubes are then washed with 3 ml of washing solution and aspirated. A calibration curve is plotted and the 17-OHP concentrations of the samples are determined by dose interpolation from the calibration curve.

V. REAGENTS PROVIDED

Reagents	96 Test Kit	Reconstitution
 Tubes coated with anti 17-OHP	2 x 48	Ready for use
 TRACER: ¹²⁵ Iodine labelled 17-OHP (HPLC grade) in phosphate/citric acid buffer with bovine casein and azide (<0.1%)	1 vial 55 ml 190 kBq	Ready for use
 Zero calibrator in human serum and azide (0.5%)	1 vial 3 ml	Ready for use
 Calibrators 17-OHP N = 1 to 5 (see exact values on vial labels) in human serum and azide (0.5%)	5 vials 0.5 ml	Ready for use
 Wash solution (TRIS-HCl)	1 vial 10 ml	Dilute 70 x with distilled water (use a magnetic stirrer).
 Controls - N = 1 or 2 in human plasma and thymol	2 vials lyophilized	Add 0.5 ml distilled water

Note : Use the zero calibrator for sample dilutions.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- Distilled water
- Pipettes for delivery of: 25 µl, 500 µl and 3 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Vortex mixer
- Magnetic stirrer
- Water bath at 37°C
- 5 ml automatic syringe (Cornwall type) for washing
- Aspiration system (optional)
- Any gamma counter capable of measuring ¹²⁵I may be used (minimal yield 70%).
- Incubator at 37°C
- Optional (for extraction)**
- Diethyl ether (analytical grade; purity > 98%)
- Disposable glass tubes (12 x 75 mm), with stoppers.
- DIAsource Reconstitution Solution for (17-OHP)-RIA-CT. ref: 32.140.74.

VII. REAGENT PREPARATION

- Controls** : Reconstitute the controls with 0.5 ml distilled water.
- Working Wash solution** : Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, controls are stable for one week at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for 3 months.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum or plasma samples must be kept at 2-8°C.
- If the test is not run within 48 hrs., storage at -20°C is recommended.
- Avoid successive freezing and thawing.
- Heparinized plasma provide 10% lower results than serum :
Y (hep. plasma) = 0.9 x (Serum) + 0.20 r = 0.92 n = 20
- EDTA plasma provides 10% lower results than serum :
Y (EDTA plasma) = 0.9 x (Serum) + 0.18 r = 0.87 n = 20

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date.
Do not mix materials from different kit lots.
Bring all the reagents to room temperature prior to use.
Thoroughly mix all reagents and samples by gentle agitation or swirling.
Use a clean disposable pipette tip for addition of each different reagent and sample in order to avoid cross-contamination. High precision pipettes or automated pipetting equipment will improve the precision.
Respect the incubation times.
Prepare a calibration curve for each run, do not use data from previous runs.

B. Extraction of newborn samples (optional) – FOR SERUM SAMPLE

- Label one glass tube for each sample (**do not extract calibrators or controls**).
- Pipette 75 µl of serum, followed by 2.25 ml of diethyl ether.
- Vortex all the tubes vigorously (2 x 1 minute)
- Let stand on the table for 5 minutes
- Place the tubes at -20°C in order to freeze the aqueous phase.
- Prepare a second series of glass tubes and, for each sample, transfer the organic phase into these new tubes. Avoid contamination by the aqueous phase.
- Evaporate the organic (ether) phase completely under a stream of air or by placing the tubes at 37°C (water bath). Manipulate under a hood.
- Dissolve the dry ether extract with 150 µl of Reconstitution solution (**not provided, see VI**). Vortex vigorously for 1 minute.
- Let stand for 10-15 minutes and vortex again for 1 minute. These volumes allow to perform the determination in duplicate.

C. Procedure

- Label coated tubes in duplicate for each calibrator, sample and control. For the determination of total counts, label 2 normal tubes
- Briefly vortex calibrators, samples and controls and dispense 25µl of each into respective tubes. **For extracted newborn samples, pipette 50 µl instead of 25 µl.**
- Dispense 0.5 ml of ¹²⁵Iodine labelled 17-OHP into each tube, including the uncoated tubes for total counts.
- Shake the tube rack gently by hand, to liberate any trapped bubbles.
- Incubate for 3 hours at 37°C.
- Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
- Wash tubes with 3 ml Working Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
- Let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
- Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

- Calculate the mean of duplicate determinations.
- Calculate the bound radioactivity as a percentage of the binding determined at the zero calibrator point (0) according to the following formula :

$$B/B_0(\%) = \frac{\text{Counts (Calibrator or sample)}}{\text{Counts (Zero Calibrator)}} \times 100$$

- Using a 3 cycle semi-logarithmic or logit-log graph paper, plot the (B/B₀(%)) values for each calibrator point as a function of the 17-OHP concentration of each calibrator point, reject obvious outliers.
- Computer assisted methods can also be used to construct the calibration curve. If automatic result processing is used, a 4 parameter logistic function curve fitting is recommended.
- By interpolation of the sample (B/B₀(%)) values, determine the 17-OHP concentrations of the samples from the reference curve.
- The obtained concentration for extracted newborn samples (see X.B.) can be read directly on the curve (no correction factor needed)
- For each assay, the percentage of total tracer bound in the absence of unlabelled 17-OHP (B₀/T) must be checked.

XII. TYPICAL DATA

The following data, obtained with a 6 weeks old tracer, are for illustration only and should never be used instead of the real time calibration curve.

(17-OHP)-RIA-CT	cpm	B/B ₀ (%)
Total count	32363	
Calibrator		
0.00 ng/ml	12596	100.0
0.17 ng/ml	9769	77.6
0.59 ng/ml	6529	51.8
1.4 ng/ml	3752	29.8
4.3 ng/ml	1396	11.1
14.0 ng/ml	479	3.8

XIII. PERFORMANCE AND LIMITATIONS

A. Sensitivity

The LoB (Limit of Blank) was tested by measuring 20 times the blank and was calculated as the mean ± 2 Standard Deviation of the distribution of test values.

LoB was calculated at 0.03 ng/ml.

The LoD (Limit of Detection) was calculated as the LoB ± 1.645 Standard Deviation of a low concentration sample measured in 10 different tests.

The LoD was calculated at 0.09 ng/ml.

The LoQ (Limit of Quantification) was calculated by testing 4 low concentration samples, in 10 different assays.

The LoQ was calculated at 0.14 ng/ml.

B. Specificity

The specificity was estimated by spiking a pool of 17-OHP samples (± 0.5 ng/ml) with steroids which might be present in patient samples.

Compound	Added amount (ng/ml)	Cross-Reactivity (%)
17α-hydroxyprogesterone	-	100
Progesterone	500	0.7
17α-hydroxypregnenolone	1000	0.3
21-deoxycortisol	1000	1.07
Pregnenolone	500	0.02
11-deoxycortisol	500	0.21
Corticosterone	1000	0.008
11-deoxycorticosterone	500	0.24
Cortisol	5000	0.005
Testosterone	5000	0.004
Androstenedione	5000	0.01
Estradiol	1000	0.008
17α-hydroxypregnenolone sulfate	5000	0.34

The effect of potential interfering substances on samples using the DIASource 17OHP-RIA-CT was evaluated. Different levels of Hemoglobin, Bilirubin, Triglyceride and Bilirubin Conjugate in serum samples were tested on samples with different 17α-hydroxyprogesterone concentrations. Our acceptance criteria was to have interference of less than 10%. The tested substances did not affect the performance of the DIASource 17OHP-RIA-CT test.

Substance	17OHP (ng/ml)	Concentration of interferent (mg/dl)	Mean Variation
Hemoglobin	0.57	250	98.6%
		500	
	1.8	250	
Bilirubin Conjugated	0.57	50	95.1%
	1.8	50	
Bilirubin unconjugated	0.57	50	99.2%
		100	
	1.8	50	
Triglyceride	0.57	50	97.0%
		100	
		250	
		500	
		1000	
	1.54	50	
		100	
		250	
		500	
		1000	

C. Precision

INTRA-ASSAY

INTRA-ASSAY

AFTER EXTRACTION

Serum	N	<X> ± SD (ng/ml)	CV (%)	Serum	N	<X> ± SD (ng/ml)	CV (%)
A	20	0.73 ± 0.05	6.8	C	24	0.56 ± 0.06	10.7
B	20	1.94 ± 0.09	4.6	D	24	1.81 ± 0.17	9.4

SD : Standard Deviation; CV: Coefficient of variation

INTER-ASSAY

INTER-ASSAY

AFTER EXTRACTION

Serum	N	<X> ± SD (ng/ml)	CV (%)	Serum	N	<X> ± SD (ng/ml)	CV (%)
A	12	0.69 ± 0.06	8.7	A	8	0.72 ± 0.13	18
B	12	2.08 ± 0.16	7.6	B	8	2.03 ± 0.39	19.2

SD : Standard Deviation; CV: Coefficient of variation

D. Accuracy

DILUTION TEST

Serum	Dilution	Theoretical Concent. (ng/ml)	Measured Concent. (ng/ml)
	1/1	-	3.3
	1/2	1.65	1.79
	1/4	0.83	1.0
	1/8	0.42	0.5
	1/16	0.21	0.27

Samples were diluted with the zero calibrator.

RECOVERY TEST

Sample	added 17-OHP (ng/ml)	Recovered 17-OHP (ng/ml)	Recovered (%)
1	6	5.7	95
	3	2.99	99
	1.5	1.64	109
	0.75	0.81	108
2	6	5.27	87
	3	2.86	95
	1.5	1.49	99
	0.75	0.74	98

Conversion factor :

From ng/ml to nmol/L : x 3.03

From nmol/L to ng/ml : x 0.33

To the best of our knowledge, no international reference material exists for this parameter.

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 40 minutes after the calibrator has been added to coated tubes.

TIME DELAY

Sample (ng/ml)	0'	10'	30'	40'
Sample 1	0.63	0.65	0.73	0.74
Sample 2	2.03	2.4	2.37	2.2

F. Limitations for newborn samples

Very young children (< 3 months old) present high levels of 17 α -hydroxypregnenolone sulphate. Despite the fact that there is almost no cross reaction with this molecule, we recommend to apply the extraction procedure (X. B.)

XIV. INTERNAL QUALITY CONTROL

- The percentage of total tracer bound in the absence of unlabelled 17-OHP (B0/T) must be > 25%.
- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Do not freeze-thaw more than twice.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises.

XV. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

	Concentration range (2.5 to 97.5% percentiles) (ng/ml)	Number of subjects
Normal males	0.67 – 3.32	79
Normal Females		
. Follicular phase	0.41 – 2.27	50
. Luteal phase	0.23 – 3.87	28
Pregnancy		
. First trimester	2.60 – 7.95	30
. Second trimester	1.83 – 9.74	30
. Third trimester	3.54 – 18.97	15
Newborns (after extraction)	0.34 – 1.96	23 (0 to 3 months)

XVI. PRECAUTIONS AND WARNINGS

Safety

For in vitro diagnostic use only.

This kit contains ¹²⁵I (half-life: 60 days) ,emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons. Purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiosafety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVII. BIBLIOGRAPHY

1. ANDO S. et al (1983)
Plasma levels of 17-OH-progesterone and testosterone in patients with varicoceles.
Acta Endocrinologica, 102,463-469.
2. BONACCORSI A. et al (1987)
Male infertility due to congenital adrenal hyperplasia: testicular biopsy findings, hormonal evaluation, and therapeutic results in three patients.
Fertility and Sterility, 47, 4, 664-670.
3. CACCIARE et al (1987)
Neonatal screening program for congenital adrenal hyperplasia in a homogeneous Caucasian population.
Annals New York Academy of sciences 85-88.
4. DRAFTA D. et al (1982)
Plasma steroids in benign prostatic hypertrophy and carcinoma of the prostate.
J. Steroid Biochem. 17, 689-693.
5. GRANOFF A. et al (1985)
17-hydroxyprogesterone responses to adrenocorticotropin in children with premature adrenarche.
Journal of Clinical Endocrinology and Metabolism, 3, 409-415.
6. NEW M. et al (1983)
Genotyping steroid 21-hydroxylase deficiency hormonal reference data.
Journal of Clinical Endocrinology and Metabolism 57, 320-326.
7. SOLYOM J. et al (1987)
Detection of late-onset adrenal hyperplasia in girls with peripubertal virilization.
Acta Endocrinologica (Copenh.) 115, 413-418.

XVIII. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS (μl)	CALIBR ATORS (μl)	SAMPLE(S) CONTROLS (μl)	EXTRACTED SAMPLES (μl)
Calibrators (0-5)	-	25	-	-
Samples, controls	-	-	25	-
Extracted samples	-	-	-	50
Tracer	500	500	500	500
Incubation	3 hours at 37°C			
Separation Working Wash solution Separation	-	aspirate 3.0 ml aspirate carefully		
Counting	Count tubes for 60 seconds			

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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XVII. BIBLIOGRAPHIE

1. ANDO S. et al (1983)
Plasma levels of 17-OH-progesterone and testosterone in patients with varicoceles.
Acta Endocrinologica, 102,463-469.
2. BONACCORSI A. et al (1987)
Male infertility due to congenital adrenal hyperplasia: testicular biopsy findings, hormonal evaluation, and therapeutic results in three patients.
Fertility and Sterility, 47, 4, 664-670.
3. CACCIARE et al (1987)
Neonatal screening program for congenital adrenal hyperplasia in a homogeneous Caucasian population.
Annals New York Academy of sciences 85-88.
4. DRAFTA D. et al (1982)
Plasma steroids in benign prostatic hypertrophy and carcinoma of the prostate.
J. Steroid Biochem. 17, 689-693.
5. GRANOFF A. et al (1985)
17-hydroxyprogesterone responses to adrenocorticotropin in children with premature adrenceance.
Journal of Clinical Endrocrinology and Metabolism, 3, 409-415.
6. NEW M. et al (1983)
Genotyping steroid 21-hydroxylase deficiency hormonal reference data.
Journal of Clinical Endocrinology and Metabolism 57, 320-326.
7. SOLYOM J. et al (1987)
Detection of late-onset adrenal hyperplasia in girls with peripubertal virilization.
Acta Endocrinologica (Copenh.) 115, 413-418.

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XVII. BIBLIOGRAFIE

1. ANDO S. et al (1983)
Plasma levels of 17-OH-progesterone and testosterone in patients with varicoceles.
Acta Endocrinologica, 102,463-469.
2. BONACCORSI A. et al (1987)
Male infertility due to congenital adrenal hyperplasia: testicular biopsy findings, hormonal evaluation, and therapeutic results in three patients.
Fertility and Sterility, 47, 4, 664-670.

3. CACCIARE et al (1987)
Neonatal screening program for congenital adrenal hyperplasia in a homogeneous Caucasian population.
Annals New York Academy of sciences 85-88.
4. DRAFTA D. et al (1982)
Plasma steroids in benign prostatic hypertrophy and carcinoma of the prostate.
J. Steroid Biochem. 17, 689-693.
5. GRANOFF A. et al (1985)
17-hydroxyprogesterone responses to adrenocorticotropin in children with premature adrenance.
Journal of Clinical Endocrinology and Metabolism, 3, 409-415.
6. NEW M. et al (1983)
Genotyping steroid 21-hydroxylase deficiency hormonal reference data.
Journal of Clinical Endocrinology and Metabolism 57, 320-326.
7. SOLYOM J. et al (1987)
Detection of late-onset adrenal hyperplasia in girls with peripubertal virilization.
Acta Endocrinologica (Copenh.) 115, 413-418.

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XVII. LITERATUR

1. ANDO S. et al (1983)
Plasma levels of 17-OH-progesterone and testosterone in patients with varicoceles.
 Acta Endocrinologica, 102,463-469.
2. BONACCORSI A. et al (1987) **Male infertility due to congenital adrenal hyperplasia: testicular biopsy findings, hormonal evaluation, and therapeutic results in three patients.**
 Fertility and Sterility, 47, 4, 664-670.
3. CACCIARE et al (1987)
Neonatal screening program for congenital adrenal hyperplasia in a homogeneous Caucasian population.
 Annals New York Academy of sciences 85-88.
4. DRAFTA D. et al (1982)
Plasma steroids in benign prostatic hypertrophy and carcinoma of the prostate.
 J. Steroid Biochem. 17, 689-693.
5. GRANOFF A. et al (1985)
17-hydroxyprogesterone responses to adrenocorticotropin in children with premature adreanence.
 Journal of Clinical Endocrinology and Metabolism, 3, 409-415.
6. NEW M. et al (1983)
Genotyping steroid 21-hydroxylase deficiency hormonal reference data.
 Journal of Clinical Endocrinology and Metabolism 57, 320-326.
7. SOLYOM J. et al (1987)
Detection of late-onset adrenal hyperplasia in girls with peripubertal virilization.
 Acta Endocrinologica (Copenh.) 115, 413-418.

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1. ANDO S. et al (1983)
Plasma levels of 17-OH-progesterone and testosterone in patients with varicoceles.
Acta Endocrinologica, 102,463-469.
2. BONACCORSI A. et al (1987)
Male infertility due to congenital adrenal hyperplasia: testicular biopsy findings, hormonal evaluation, and therapeutic results in three patients.
Fertility and Sterility, 47, 4, 664-670.
3. CACCIARE et al (1987)
Neonatal screening program for congenital adrenal hyperplasia in a homogeneous Caucasian population.
Annals New York Academy of sciences 85-88.
4. DRAFTA D. et al (1982)
Plasma steroids in benign prostatic hypertrophy and carcinoma of the prostate.
J. Steroid Biochem. 17, 689-693.
5. GRANOFF A. et al (1985)
17-hydroxyprogesterone responses to adrenocorticotropin in children with premature adrenance.
Journal of Clinical Endocrinology and Metabolism, 3, 409-415.
6. NEW M. et al (1983)
Genotyping steroid 21-hydroxylase deficiency hormonal reference data.
Journal of Clinical Endocrinology and Metabolism 57, 320-326.
7. SOLYOM J. et al (1987)
Detection of late-onset adrenal hyperplasia in girls with peripubertal virilization.
Acta Endocrinologica (Copenh.) 115, 413-418.

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XVII. BIBLIOGRAFIA

1. ANDO S. et al (1983)
Plasma levels of 17-OH-progesterone and testosterone in patients with varicoceles.
Acta Endocrinologica, 102,463-469.
2. BONACCORSI A. et al (1987)
Male infertility due to congenital adrenal hyperplasia: testicular biopsy findings, hormonal evaluation, and therapeutic results in three patients.
Fertility and Sterility, 47, 4, 664-670.
3. CACCIARE et al (1987)
Neonatal screening program for congenital adrenal hyperplasia in a homogeneous Caucasian population.
Annals New York Academy of sciences 85-88.

4. DRAFTA D. et al (1982)
Plasma steroids in benign prostatic hypertrophy and carcinoma of the prostate.
J. Steroid Biochem. 17, 689-693.
5. GRANOFF A. et al (1985)
17-hydroxyprogesterone responses to adrenocorticotropin in children with premature adrenance.
Journal of Clinical Endocrinology and Metabolism, 3, 409-415.
6. NEW M. et al (1983)
Genotyping steroid 21-hydroxylase deficiency hormonal reference data.
Journal of Clinical Endocrinology and Metabolism 57, 320-326.
7. SOLYOM J. et al (1987)
Detection of late-onset adrenal hyperplasia in girls with peripubertal virilization.
Acta Endocrinologica (Copenh.) 115, 413-418.

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1. ANDO S. et al (1983)
Plasma levels of 17OHprogesterone and testosterone in patients with varicoceles.
Acta Endocrinologica, 102,463469.
2. BONACCORSI A. et al (1987)
Male infertility due to congenital adrenal hyperplasia: testicular biopsy findings, hormonal evaluation, and therapeutic results in three patients.
Fertility and Sterility, 47, 4, 664670.
3. CACCIARE et al (1987)
Neonatal screening program for congenital adrenal hyperplasia in a homogeneous Caucasian population.
Annals New York Academy of sciences 8588.
4. DRAFTA D. et al (1982)
Plasma steroids in benign prostatic hypertrophy and carcinoma of the prostate.
J. Steroid Biochem. 17, 689693.
5. GRANOFF A. et al (1985)
17hydroxyprogesterone responses to adrenocorticotropin in children with premature adrenance.
Journal of Clinical Endocrinology and Metabolism, 3, 409415.
6. NEW M. et al (1983)
Genotyping steroid 21hydroxylase deficiency hormonal reference data.
Journal of Clinical Endocrinology and Metabolism 57, 320326.
7. SOLYOM J. et al (1987)
Detection of lateonset adrenal hyperplasia in girls with peripubertal virilization.
Acta Endocrinologica (Copenh.) 115, 413418.

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