# **Fecal Adenovirus Antigen Test Kit**

***RAPEPKT918***

---

Version: 230313

Date of issue : 13/03/2023

Revision date: 13/03/2023

# History

**Summary of change :**

Previous version: 200224-1	Current Version: 230313
	New logo
<b>REAGENTS: Preparation and Storage</b>	<b>REAGENTS: Preparation and Storage</b> Information reviewed
<b>MATERIALS REQUIRED BUT NOT PROVIDED</b> Positive Control	<b>MATERIALS REQUIRED BUT NOT PROVIDED</b> Timer or clock
<b>SAFTY PRECAUTIONS</b> The reagents are for professional use only. Source material from which reagents of bovine serum was derived in the contiguous 48 United States. It was obtained only from donor health animals maintained under veterinary supervision and found free of contagious diseases. Wear gloves while performing this test and handle these reagents and patient samples as if they are infectious. Do not get in eyes, on skin, or on clothing. Do not ingest. On contact, flush with copious amounts of water for at lease 15 minutes. When the assay procedure is completed, dispose of specimens (biohazard materials) carefully after autoclaving for at least one hour. Alternatively, treat with a 0.5 or 1% solution of sodium hypochlorite for one hour before disposal. Use Good Laboratory Practices	<b>SAFTY PRECAUTIONS</b> 1. For in-vitro diagnostic use only. Not to be taken internally. 2. Do not use product beyond the expiration date. 3. Handle all specimens as potentially infectious. 4. Do not reuse the test.
<b>No existing section "PATIENT PREPARATION"</b>	<b>PATIENT PREPARATION</b> Dietary restrictions are not necessary.
<b>SPECIMEN COLLECTION</b>	<b>SPECIMEN COLLECTION</b> Information reviewed
<b>No existing section "PROCEDURAL NOTES"</b>	<b>PROCEDURAL NOTES</b> 1. After the test strip tube is screwed completely into the sampling tube, you should see a minimum 5 mm extraction buffer liquid in the bottom of the strip tube. 2. You should see liquid migrating across the membrane area right after the screw in process. If not, take the tube and tap against the table several times, and the migration of the liquid should be observed.

<p><b>INTERPRETATION OF RESULTS</b></p> <ul style="list-style-type: none"> <li>• <b>Positive:</b> If two red/pink colored bands are visible within 10 minutes, the test result is positive and valid.</li> <li>• <b>Negative:</b> If the test area has no red/pink colored band and the control area displays a red/pink colored band, the test result is negative.</li> </ul>	<p><b>INTERPRETATION OF RESULTS</b></p> <ul style="list-style-type: none"> <li>• <b>Positive:</b> If two red/pink colored bands are visible within 5 minutes, the test result is positive and valid.</li> <li>• <b>Negative:</b> If the test area has no red/pink colored band and the control area displays a red/pink colored band, the test result is negative.</li> </ul>
<p><b>1. Internal procedural control:</b> Each fecal adenovirus antigen test contains a build in procedural control. A coloured control line will appear if the test has been performed correctly occurred and the reagents are reactive. This control does not ensure that the test line antibody is accurately detecting the presence or absence of viral antigen in the test fecal sample.</p>	<p><b>1. Internal procedural control:</b> Each fecal adenovirus antigen test contains a build in procedural control. A coloured control line will appear if the test has been performed correctly, sample wicking has occurred and the reagents are reactive. This control does not ensure that the test line antibody is accurately detecting the presence or absence of viral antigen in the test fecal sample.</p>
<p><b>LIMITATION OF THE PROCEDURE</b></p> <p>4. No cross-reactions with other viruses or substances were observed during the evaluation of the test. A negative result does not totally exclude a possible adenovirus infection. The significance of the results must be evaluated in relation to the patient's clinical symptoms.</p>	<p><b>LIMITATION OF THE PROCEDURE</b></p> <p>Information removed</p>

# Fecal Adenovirus Antigen Test Kit

en

Rapid Immunochromatographic Test Device for the detection of  
Adenovirus Antigen In Feces

**RAPEPKT918**

**IN VITRO DIAGNOSTIC**

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 91

## INTENDED USE

This Adenovirus antigen test kit is intended for the direct qualitative detection of the presence of adenovirus antigen in patient fecal samples. The test might be used as an aid for the diagnosis of an adenovirus infection resulting in acute gastroenteritis. For professional use only.

## SUMMARY OF PHYSIOLOGY

Adenoviruses are one of the main causes of acute gastroenteritis and diarrhea, especially in children under the age of two years. Adenoviruses have been identified in almost 12% of the feces of children with gastroenteritis. It was reported that adenovirus is the second leading cause of the hospitalized cases of diarrhea in infant and young children. If not treated, the infection may result in severe dehydration and disorders of body electrolyte balance. Therefore, it can be mortal in risk populations such as children, the elderly or immunosuppressed individuals. Adenovirus is transmitted by oral-fecal contact, but can result from the inhalation of aerosols as well. Its incubation period lasts 5 to 8 days. Characteristic symptoms include vomiting, hydrodiarrhoea, high temperature and stomach pains. There are 41 known human adenoviruses primarily differentiated by serology and DNA analysis. Morphologically, the viruses are non-enveloped icosahedral structured with a diameter of about 80 nm.

## ASSAY PRINCIPLE

The Adenovirus Rapid Test Strip employs dye-conjugated monoclonal antibodies specifically against human adenovirus antigens and a solid-phase coated monoclonal antibody specifically against adenovirus hexon antigens. In this test the specimen is first treated with an extraction solution to extract adenovirus antigens from the feces. Following extraction, the only step required is to screw the adenovirus test strip into the sample collection device. As the sample extraction flows through the chamber and reach the test strip, the colored particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the colored particles. Different colored lines will be visible, depending upon the virus content of the sample. These lines, after 5 minutes of incubation at room temperature, are used to interpret the result.

## REAGENTS: Preparation and Storage

- **TUBE** 30 Fecal specimen collection device: containing sampling tube, sampling lid and pre-added extraction solution (1.1 ml Tris Buffer containing 0.1% bovine serum albumin) in the sampling tube. This device should be stored at 2 to 8°C. Do not freeze.
- **STRIP** 30 Test strip: one dipstick for the Fecal Adenovirus Antigen test is assembled in a transparent housing and sealed in a foil pouch with desiccant. It should remain in its original sealed pouch until ready for use. The test strip should be stored at 2 to 30°C. Do not freeze.
- Single use green pipette.
- Patient label for sample collection device.
- Instruction for use.

## SAFETY PRECAUTIONS

1. For in-vitro diagnostic use only. Not to be taken internally.
2. Do not use product beyond the expiration date.
3. Handle all specimens as potentially infectious.
4. Do not reuse the test.

## MATERIALS REQUIRED BUT NOT PROVIDED

Timer or clock

## PATIENT PREPARATION

Dietary restrictions are not necessary.

## SPECIMEN COLLECTION

Stool samples must be taken as soon as the symptoms appear. Viral particles decrease in number after one week, making the diagnosis more difficult. The samples can be stored at 2-8 °C for 1 to 2 days. For longer storage they must be kept frozen at -20°C. In this case, the sample should be totally thawed, and brought to room temperature and homogenized before testing.

For each solid specimen: Unscrew the sampling lid and keep the sampling tube in a vertical position to prevent the loss of any extraction solution. Insert and twist the tip of the sampling lid into the stool specimen at two or more different sites. Collect fecal sample that is stuck to the surface of the sampling lid. Do not intentionally collect any separate and large pieces of fecal sample into the tube. Replace the sampling lid into the tube and secure tightly.

For each liquid specimen: Unscrew the sampling lid and keep the sampling tube in a vertical position to prevent the loss of any extraction solution. Using the green pipette to collect three drops of the liquid stool sample and add to the sampling tube. Replace the sampling lid into the tube and secure tightly.

*Note: Two specimens from three consecutive bowel movements are recommended from American Cancer Society. Specimen should not be collected during digital rectal examination.*

## TEST PROCEDURE

1. Bring the sealed foil pouch test strips and collected specimens to room temperature.
2. Shake the sampling tube vigorously to ensure a good liquid suspension.
3. Position the sampling tube upside down vertically and let it settle for about 1 minute.
4. Remove the test strip from the sealed foil pouch (figure A).
5. Screw the test strip in a **vertical position** into the sampling tube by **breaking** the bottom seal of the sampling tube. Secure tightly! (Figure B)
6. Allow the solution to flow into the bottom space of the test strip and keeping the device in a **vertical position**.
7. Read test result at 5 minutes. Do not interpret test result after 10 minutes.

## PROCEDURAL NOTES

1. After the test strip tube is screwed completely into the sampling tube, you should see a minimum 5 mm extraction buffer liquid in the bottom of the strip tube.
2. You should see liquid migrating across the membrane area right after the screw in process. If not, take the tube and tap against the table several times, and the migration of the liquid should be observed.

## INTERPRETATION OF RESULTS

- **Positive:**  
If two red/pink colored bands are visible within 5 minutes, the test result is positive and valid.
- **Negative:**  
If the test area has no red/pink colored band and the control area displays a red/pink colored band, the test result is negative.

• **Invalid:**

If a colored band does not appear in the control area regardless whether there is any band in the test area, the test result is invalid.

Note: The control line is located at the upper half of the test membrane area and the test line is located at the lower half of it.

**QUALITY CONTROL**

Good laboratory practices recommend the use of appropriate controls. There are two types of controls for the Adenovirus antigen test, the internal procedural control and external controls.

- 1. Internal procedural control:** Each fecal adenovirus antigen test contains a build in procedural control. A coloured control line will appear if the test has been performed correctly, sample wicking has occurred and the reagents are reactive. This control does not ensure that the test line antibody is accurately detecting the presence or absence of viral antigen in these test fecal sample.
- 2. External controls:** It is recommended to use external positive controls. The external positive controls are not provided with this kit. External controls are used to assure that the test line antibody is reactive. However, external controls will not detect an error in the performance of a particular test with a patient sample. It is recommended that the external control be tested once per kit.

You should always follow local, state, federal guidelines for running quality control.

**LIMITATION OF THE PROCEDURE**

1. The test should be used only for the detection of adenovirus antigen in fecal samples.
2. The test is qualitative and no quantitative interpretation should be made with respect to the intensity of the positive line, when reporting the result
3. Two hundred samples were evaluated to assure the correct performance of the test. The correlation of the results with other techniques (ELISA) was satisfactory. However, interferences in the performance of the tests should not be excluded.
4. As with all diagnostic tests, the definitive clinical diagnosis must not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated. Fecal adenovirus antigen test is designed for the aid of clinical diagnosis and should not replace other diagnostic procedures.

**PERFORMANCE**

The sensitivity and specificity of this adenovirus antigen test device are studied with 212 clinical samples and compared with an adenovirus antigen ELISA test. The results are shown below.

Strip \ ELISA	Positive	Negative	Total
Positive	61	1	62
Negative	1	149	150
Total	62	150	212

**Sensitivity:** 98 % (61/62 = 98.4%)

**Specificity:** 99 % (149/150 = 99.3%)

**Accuracy:** 99% (210/212 = 99.1%)

**Inter-series and intra-series accuracy:** 100 %

**Interference:** Cross reactivity has been evaluated and found to be negative compared to positive specimens of *Cryptosporidium parvum* and rotavirus.

**LITERATURES**

1. Nishio O, Ooseto M, Takagi K, Yamasita Y, Ishihara Y, Isomura S. Enzyme-linked immunosorbent assay employing monoclonal antibodies for direct identification of enteric adenoviruses (Ad40,41) in feces. Microbiol Immunol. 1990;34(10):871-7.
2. Vizzi E, Ferraro D, Cascio A, Di Stefano R, Arista S. Detection of enteric adenoviruses 40 and 41 in stool specimens by monoclonal antibody-based enzyme immunoassays. Res Virol. 1996 Nov-Dec;147(6):333-9.
3. Singh-Naz N, Rodriguez WJ, Kidd AH, Brandt CD. Monoclonal antibody enzyme-linked immunosorbent assay for specific identification and typing of subgroup F adenoviruses. J Clin Microbiol. 1988 Feb;26(2):297-300.
4. Uhnoo I, Wadell G, Svensson L, Johansson ME. Importance of enteric adenoviruses 40 and 41 in acute gastroenteritis in infants and young children. J Clin Microbiol. 1984 Sep;20(3):365-72.
5. Uhnoo I, Svensson L, Wadell G. Enteric adenoviruses. Baillieres Clin Gastroenterol. 1990 Sep;4(3):627-42.
6. Shinozaki T, Araki K, Fujita Y, Kobayashi M, Tajima T, Abe T. Epidemiology of enteric adenoviruses 40 and 41 in acute gastroenteritis in infants and young children in the Tokyo area. Scand J Infect Dis. 1991;23(5):543-7.

Consult instructions for use

Manufacturer

Storage temperature

Contains sufficient for n tests

Use by

In vitro diagnostic medical device

Batch code

Test Strip

Catalogue number

Collection device

Other translations of this Instructions for Use can be downloaded from our website:

<https://www.diasource-diagnostics.com/>

# Fecal Adenovirus Antigen Test Kit

Rapid Immunochromatographic Test Device for the detection of Adenovirus Antigen In Feces – Instructions for Fecal Sample Collection


**RAPEPKT918**

*IN VITRO DIAGNOSTIC*

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 91


**1**

1. Take the sampling tube and unscrew the sampling lid, keeping the sampling tube in a **vertical position** to prevent loss of solution.


**2**

1. Hold the sampling lid by the **Thumb Grip**.
2. Use the **tip** of the sampling lid to collect a small amount of fecal sample at two or more sites. Only take the fecal sample that sticks to the sampling lid tip (never intentionally place any separate piece of fecal sample into the tube).

**3**

1. Insert and screw the sampling lid back into the sampling tube in a **vertical position**. Do not spill any solution from the tube.
2. Tightly seal the lid with the tube.

Rapid Immunochromatographic Test Device for the detection of Adenovirus Antigen In Feces – Instructions for Test Procedures


**RAPEPKT918**

***IN VITRO DIAGNOSTIC***

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 91


**A**

1. Shake to dissolve the stool into solution.
2. Turn the sampling tube **upside down vertically**.
3. Remove the test strip from foil pouch.

**B**

1. Insert and screw the test strip **in a vertical position** into the sampling tube by breaking the bottom seal of the sampling tube.
2. Allow the solution to flow into the bottom space of test strip, keeping the device **in a vertical position**.
3. You may soon see a red fluid moving across the white area of the test strip. Read test result after 5 minutes.

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