# Trypsin RIA

***KIRCE07***

---

Version : 230123

Date of issue : 23/01/2023

Revision date: 23/01/2023

# **History**

---

## **Summary of change:**

<b>Current Version:</b>
<b>230123</b>
New logo

en

Read entire protocol before use.

# Trypsin RIA

## I. INTENDED USE

Radioimmunoassay for the quantitative determination of human trypsin (h trypsin) in serum.

**For Research use only. Not for use in diagnostic procedures.**

## II. GENERAL INFORMATION

- A. Proprietary name : DIAsource Trypsin RIA Kit
- B. Catalog number : KIRCE07 : 100 tests
- C. Manufactured by : DIAsource ImmunoAssays S.A.

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## III. BACKGROUND

Trypsin is one of the enzymes involved in the digestion of dietary protein in the small bowel. It functions as an endopeptidase with a marked preferential specificity for peptide bonds arising from the carbonyl groups of L-arginyl and L-lysyl residues. Chemically, trypsin is a protein consisting of 201 amino acids with a molecular weight of 22,900 Da.

In common with other proteolytic enzymes of pancreatic origin, trypsin is synthesized and secreted from the acinar cells of the exocrine pancreas as an inactive precursor or proenzyme, trypsinogen. This occurs in two forms, both of which are converted to enzymatically active trypsin by the action of enterokinase produced by cells in the duodenum. This kinase, in the presence of Ca++ ions, removes the hexapeptide Val-(Asp)4-Lys from trypsinogen, thereby unmasking the active center of the trypsin molecule. In addition trypsin itself can activate the trypsinogen by autocatalysis.

Some trypsin appears to be excreted intact from the gut, as it can be detected and measured in stool. In addition, the trypsinogen, trypsin and probably trypsin bound to a pancreatic inhibitor may find their way from the pancreatic cells into the blood via the interstitial fluid. In blood all trypic activity is blocked by three specific inhibitors ( $\alpha$ 1-antitrypsin,  $\alpha$ 2-macroglobulin and inter- $\alpha$ -trypsin inhibitor). Nevertheless, normal serum still retains some proteolytic activity as shown by the hydrolysis of synthetic substrates, but this cannot be due solely to trypsin.

#### **IV. PRINCIPLES OF THE METHOD**

Trypsin can be measured in duodenal or pancreatic juice by spectrophotometric methods, which depend upon its ability to hydrolyze synthetic substrates. However, as the trypsins are inactive and the enzymatic activity of trypsin in serum is blocked by the inhibitors mentioned above, this approach cannot be used to determine blood levels of the enzyme.

If the immunological rather than the enzymatic properties of trypsin could be used in a method for its determination in serum, the foregoing problem would not arise.

The Trypsin kit permits radioimmunochemical determination of the enzyme trypsin in human serum and other biological fluids. The kit utilizes the principle of competitive protein binding analysis, using a double antibody radioimmunoassay method, free trypsin being separated from antibody-bound trypsin by an anti-rabbit precipitating serum.

During the first incubation, serum trypsin and  $^{125}\text{I}$ -labelled trypsin compete for highly specific antibodies against trypsin, with the result that  $^{125}\text{I}$ -labelled trypsin is bound to the antibodies in inverse proportion to the amount of serum trypsin present.

During the second incubation the trypsin or 125I-trypsin antibody complex reacts with an anti-rabbit-gamma globulin, causing precipitation. The precipitate can be centrifuged down, the supernatant decanted and the iodine-125 activity in the precipitate measured in a gamma scintillation counter. The activity measured is compared with a calibrator curve prepared under the same conditions and provides an indication of the trypsin concentration in the serum.

## V. REAGENTS PROVIDED

Reagents	100 Tests Kit	Reconstitution		
<table border="1" style="display: inline-table; vertical-align: middle;"><tr><td>Ag</td><td><math>^{125}\text{I}</math></td></tr></table> $^{125}\text{I}$ odine labelled human trypsin	Ag	$^{125}\text{I}$	1 vial <95 kBq 22 ml	<b>Ready</b> for use
Ag	$^{125}\text{I}$			
<table border="1" style="display: inline-table; vertical-align: middle;"><tr><td>CAL</td><td>0</td></tr></table> Calibrator 0: in horse serum and sodium azide.	CAL	0	1 vial Lyophilized	<b>Add</b> 0.5 ml distilled water
CAL	0			
<table border="1" style="display: inline-table; vertical-align: middle;"><tr><td>CAL</td><td>N</td></tr></table> Calibrators 1-6: in horse serum and sodium azide (see exact values on vial labels)	CAL	N	6 vials Lyophilized	<b>Add</b> 0.5 ml distilled water
CAL	N			
<table border="1" style="display: inline-table; vertical-align: middle;"><tr><td>CONTROL</td><td>N</td></tr></table> Control in human serum and sodium azide. (see exact value on vial label)	CONTROL	N	1 vial 0.5 ml	<b>Ready</b> for use
CONTROL	N			
<table border="1" style="display: inline-table; vertical-align: middle;"><tr><td>ANTISERUM</td><td></td></tr></table> Rabbit anti human trypsin in buffer with bovine albumin and sodium azide.	ANTISERUM		1 vial 10.5 ml	<b>Ready</b> for use
ANTISERUM				
<table border="1" style="display: inline-table; vertical-align: middle;"><tr><td>PREC</td><td>AGENT</td></tr></table> Precipitation reagent: buffer with goat anti-rabbit IgG , bovine albumin and sodium azide.	PREC	AGENT	1 vial 55 ml	<b>Ready</b> for use
PREC	AGENT			

#### **VI. SUPPLIES NOT PROVIDED**

The following material is required but not provided in the kit:

- The following material is required but not provided in the kit:

  1. Distilled water: for reagents reconstitution and for washing step.
  2. Pipettes for delivery of 100 $\mu$ l, 200 $\mu$ l and 500 $\mu$ l (the use of accurate pipettes with disposable plastic tips is recommended).
  3. Incubation tubes 3 to 5ml ; e.g.: polystyrene tubes 12 x 75 mm.
  4. Vortex mixer.
  5. 1 ml automatic syringe (Cornwall type) for washing.
  6. Aspiration system.
  7. Centrifuge  $\geq$  1500 g.
  8. Any gamma counter capable of measuring  $^{125}\text{I}$  may be used (minimal yield 70%).

## VII. REAGENT PREPARATION

**Calibrators :** Carefully dissolve the calibrators in distilled water, proceeding as follows: Gently tap the vials of calibrator to dislodge any particles which may be adhering to the stopper. Then carefully remove the stoppers, placing them upside down on the work surface. Using a suitable micropipette, dispense precisely 500  $\mu$ l distilled water into each vial, replace the stoppers and allow approximately 10 minutes for the freeze-dried material to dissolve completely. Reconstitution and mixing may be accelerated by swirling the contents and/or rolling the vials between the hands.

#### **VIII. STORAGE AND EXPIRATION DATING OF REAGENTS**

- Before opening, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
  - After its first use, tracer is stable, if kept in the original well-closed vial at 2-8°C until the tracer expiry date.
  - Alterations in physical appearance of kit reagents may indicate instability or deterioration.

## *IX. SPECIMEN COLLECTION AND PREPARATION*

After blood sampling, serum is obtained by the usual methods. It is pipetted off, assayed directly or stored up to 3 days at 2-8°C for later use, or alternatively frozen at -20°C.

The Trypsin radioimmunoassay kit permits measurement of trypsin as a protein in the serum, and is not affected by the presence of serum inhibitors. The radioimmunoassay is specific for trypsin, trypsinogen and the enzyme's inhibited forms.

## X. PROCEDURE

#### A. Handling notes

**Handling note:**  
Do not use the kit or components beyond expiry date.  
Do not mix materials from different kit lots.

Do not mix materials from different kit lots.  
Bring all other reagents to room temperature (18-25 °C) prior to use.  
Thoroughly mix all reagents and samples by gentle agitation or swirling.  
Use a clean disposable pipette tip for addition of each different reagent and sample in order to avoid cross-contamination. High precision pipettes or automated pipetting equipment will improve the precision.  
Prepare a calibration curve for each run, do not use data from previous runs.

## B. Procedure

- 2. Procedure**

  - 1- It is recommended to perform the assays in duplicate for the calibrators, the control and the samples.
  - 2- Label sufficient incubation tubes (3-5 mL), as given in *VI. Supplies not provided*. (calibrators, control serum, serum samples and "total activity" tubes).
  - 3- Pipette 100 µL calibrators, control serum and serum samples into test tubes that have been prepared.
  - 4- Dispense 200 µL  $^{125}\text{I}$ -trypsin into each test tube (including "total activity" tubes).
  - 5- Dispense 100 µL Antiserum into each test tube (excluding "total activity" tubes), mix the contents of the tubes on the vortex mixer, and incubate for 3h (3-5h) at room temperature (18-25 °C) away from direct light.
  - 6- Dispense 500 µL Precipitation reagent into each test tube (excluding "total activity" tubes), mix the contents of the tubes on the vortex mixer, and incubate for 30 min (30-60 min) at room temperature (18-25 °C) away from direct light.
  - 7- Dispense 1 mL distilled water into each test tube (excluding "total activity" tubes), centrifuge the tubes at  $\geq 1500$  g for 15 min at 20 °C and decant the supernatants.
  - 8- Measure the test tubes for 1 minute in a gamma scintillation counter.

## XI CALCULATION OF RESULTS

There is no internationally accepted calibrator for trypsin at present. The Trypsin values cannot be readily compared with the results of other trypsin assays.

The calibrators and control values are assigned against normalized reference controls and the value of previous batch of calibrators and control.

reference controls and the value of previous batch of calibrators and control. The trypsin control serum included in the kit provides an analytical means checking the accuracy of the assay procedure performed in the laboratory. In order to establish the calibration curve, the following procedure must be followed:

1. Calculate the mean of duplicate determinations.
  2. Calculate the bound radioactivity as a percentage of the binding determined at the zero calibrator point (0) according to the following formula :

$$B/B_0 (\%) = \frac{\text{Counts (Calibrator or sample)}}{\text{Counts (Zero Calibrator)}} \times 100$$

3. Plot the (B/B<sub>0</sub>(%)) values for each calibrator point as a function of trypsin concentration of each calibrator point. Reject obvious outliers.
4. Computer assisted methods can also be used to construct the calibration curve. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.
5. By interpolation of the sample (B/B<sub>0</sub> (%)) values, determine trypsin concentrations of the samples from the calibration curve.
6. For each assay, the percentage of total tracer bound in the absence of unlabelled trypsin (B<sub>0</sub>/T) must be checked.

## XII. TYPICAL CALIBRATION CURVE

The following data are for illustration only and should never be used instead of the real time calibration curve.

Trypsin RIA		cpm	B/B <sub>0</sub> (%)
Total count		29722	
Calibrator	0.0 ng/ml	15091	100
	36 ng/ml	13067	86.6
	73 ng/ml	11503	76.2
	160 ng/ml	9242	61.2
	300 ng/ml	6904	45.7
	630 ng/ml	4900	32.5
	1560 ng/ml	3842	26.5

## XIII. PERFORMANCE AND LIMITATIONS

### B. Detection limit

The LoB (Limit of blank) was calculated by measuring the blank several times and was calculated as the mean + 1.65 standard deviation of the distribution of the test values . The LoB was calculated to be 4 ng/ml.

### B. Interferences

No interference has been observed :  
with haemoglobin up to 5 g/l ,  
with triglycerides up to 2.5 g/l ,  
with bilirubin non conjugated up to 1 g/l ,  
with conjugated bilirubin up to 0.5 g/l.

### C. Precision

INTRA ASSAY			
Sample	Replicate	$\bar{X} \pm SD$ (ng/ml)	CV (%)
A	20	$382.9 \pm 25.6$	6.7
B	20	$230.3 \pm 16.2$	7.0

INTER ASSAY			
Sample	Replicate	$\bar{X} \pm SD$ (ng/ml)	CV (%)
1	10	$151.6 \pm 7.5$	4.9
2	10	$406 \pm 30.9$	7.6

SD : Standard Deviation; CV: Coefficient of variation

## XIV. INTERNAL QUALITY CONTROL

If the results obtained for the kit control are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.

If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Do not freeze-thaw more than twice.

Acceptance criteria for the difference between the duplicates results of the samples should rely on Good Laboratory Practises.

## XV. PRECAUTIONS AND WARNINGS

### Safety

For *in vitro* diagnostic use only.

This kit contains <sup>125</sup>I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area. away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

## XVI. BIBLIOGRAPHY

- (1) Adrian TE, Bestermann HS, Mallinson CN, Galarotis C, Bloom SR. Clin Sci Mol Med. 1978;54:24.
- (2) Bray PT, Clark GCF, Moody GJ, Thomas JDR. Clin Chim. 1977;80:333.
- (3) Crossley JR, Smith PA, Elliott RB. Lancet. March 3. 1979;472.
- (4) Elias E, Redshaw M, Wood T. Lancet II. 1977;66.
- (5) Fedail SS, Salmon PR, Harvey RF, Brown PB, Read AE. Abstract. 1978;19:A445.
- (6) Goldberg JM. Clin Chem. 1976;22:638.
- (7) Lesser PB, Marshaw AL. Gastroenterology. 1975;68:630.
- (8) Nagasawa S, Han BH, Sugihari H, Suzuki T. J Biochem. 1970;67:821.
- (9) Temler RS, Felber JP. Biochim Biophys. 1976;445:720.

**XVII. SUMMARY OF THE PROTOCOL**

	TOTAL COUNTS	CALIBRATORS CONTROL	SAMPLES
	µl	µl	µl
Calibrators/control Samples	-	100	-
Tracer : $^{125}\text{I}$ human Trypsin:	200	200	200
Antiserum	-	100	100
Incubation	Vortex tubes. Incubate 3 hours at 18-25°C.		
Precipitation reagent	-	500	
Incubation	Vortex tubes. Incubate for 30 minutes at 18-25°C.		
Washing step	-	Add 1.0 ml distilled water	
Centrifugation	-	Centrifuge for 15 minutes at $\geq 1500\text{g}$ , at 20°C. Decant the supernatant.	
Counting	Count tubes for 1 minute		

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