# CanAg CA242 EIA

REF

101-10

IVD

Instructions for use. 2022-06

Read highlighted changes

|    |                             |
|----|-----------------------------|
| EN | EXPLANATION OF SYMBOLS      |
| BG | ОБЯСНЕНИЕ НА СИМВОЛИТЕ      |
| CS | VÝZNAM SYMBOLŮ              |
| DA | SYMBOLFORKLARING            |
| DE | ERKLÄRUNG DER SYMBOLE       |
| EL | ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ      |
| ES | SIGNIFICADO DE LOS SÍMBOLOS |
| ET | SÜMBOLITE SELGITUS          |
| FR | EXPLICATION DES SYMBOLES    |
| HR | OBJAŠNJENJE SIMBOLA         |
| HU | JELMAGYARÁZAT               |
| IT | SPIEGAZIONE DEI SIMBOLI     |
| LT | SIMBOLIŲ PAAIŠKINIMAI       |
| LV | SIMBOLU SKAIDROJUMS         |
| NL | VERKLARING DER SYMBOLEN     |
| NO | SYMBOLFORKLARING            |
| PL | OBJAŚNIENIE SYMBOLI         |
| PT | EXPLICAÇÃO DOS SÍMBOLOS     |
| RO | SEMNIȚAȚIA SIMBOLURILOR     |
| RU | ОБОЗНАЧЕНИЯ                 |
| SV | SYMBOLFÖRKLARING            |
| SK | VÝZNAM SYMBOLOV             |
| SL | RAZLAGA SIMBOLOV            |
| SR | OBJAŠNJENJE SIMBOLA         |
| TR | SEMBOLLERİN AÇIKLAMALARI    |

Use By/Годно до/Použitelné do/  
Holdbar til/Verwendbar bis/  
Ημερομηνία λήξης/Fecha  
de caducidad/Kölblik kuni/  
Utiliser jusque/Rok valjanosti/  
Felhasználható/Utilizzare entro/  
Sunautoti iki/Izlietot līdz/Houdbaar  
tot/Brukes innen/Użyç przed/  
Prazo de validade/Expirã la/  
Использовать до/Använd före/  
Použite né do/Uporabno do/  
Upotrebljivo do/Son Kullanna Tarihi

LOT

Batch code/Номер на партида/  
Číslo šarže/Lotnummer/  
Chargenbezeichnung/Αριθμός  
Παρτίδας/Código de lote/Partii  
kood/Code du lot/Kod serije/  
Sarzsám/Codice del lotto/  
Partijos kodas/Partijas kods/Lot  
nummer/Partikode/Kod partii/  
Código do lote/Număr de lot/  
Номер лота/Lotnummer/Číslo  
šarže/Številka serije/Kod partije/  
Parti Kodu

Date of manufacture/Дата на производство/Datum výroby/  
Produktionsdato/Herstellungsdatum/  
Ημερομηνία παραγωγής/Fecha de fabricación/Valmistamise kuupäev/  
Date de fabrication/Datum proizvodnje/  
Gyártási idő/Data di produzione/  
Pagaminimo data/Ražošanas datums/  
Productiedatum/Fremstillingsdato/  
Data produkcji/Data de fabrico/Data fabricației/Дата производства/  
Tillverkningsdatum/Dátum výroby/Datum izdelave/Datum proizvodnje/Úretim tarihi

Temperature limitation/  
Температурни граници/  
Теплотни омеzeи/  
Temperaturbegrænsning/  
Temperaturbegrenzung/  
Περιορισμοί θερμοκρασίας/  
Limites de temperatura/  
Temperatuuri piirang/  
Limite de température/  
Temperaturno ograničenje/  
Hőmérsékletre vonatkozó korlátozás/  
Limiti di temperatura/  
Temperatūriniai apribojimai/  
Temperatūras ierobežojums/  
Temperaturbepæring/  
Temperaturbegrensninger/  
Temperaturey granične/  
Limite de temperatura/  
Limite de temperatură/  
Температурный режим/  
Temperaturbegrænsning/  
Теплотне обмеzenie  
Omejitve temperature/  
Temperaturno ograničenje/  
Sıcaklık sınırlaması/

**IVD**

In Vitro Diagnostic Medical Device/  
Медицински уред за диагностика  
ин vitro/Diagnostický zdravotnícký  
prostředek in vitro/Medicinsk udstyr til  
in vitro-diagnostik/In-vitro-Diagnostikum/  
Ιατροτεχνολογικό προϊόν για διάγνωση  
In Vitro/Dispositivo médico para  
diagnóstico in vitro/In vitro diagnostiline  
meditsiiniseade/Dispositif médical de  
diagnostic in vitro/Diagnostički medicinski  
uređaj In Vitro/In vitro orvosdiagnostikai  
eszköz/Dispositivo medico per test  
diagnostici in vitro/In Vitro Diagnostinė  
Medicinos Priemonė/Medicínska ierice  
in vitro diagnostikai/In vitro-diagnostisch  
medisch instrument/In vitro diagnostisk  
medisinsk utstyr/Wyrób medyczny do  
diagnostyki in vitro/Dispositivo Médico  
de Diagnóstico In Vitro/Dispozitiv medical  
pentru diagnostic in vitro/Только для  
диагностики In Vitro/Endast för in  
vitro-diagnostik/ Zdravotnícka pomôcka na  
diagnostiku in vitro/In vitro diagnostični  
pripomoček/Diagnostički medicinski  
uređaj In Vitro/<96> testleri için yeterlilik  
içerir

Contains sufficient for <96> tests/Съдържа  
достатъчно количество за тестове  
<96>/Lze použít pro <96> testů/Ineholder  
tilstrækkeligt/Inhalt ausreichend für <96>  
Prüfungen/Περεχόμενο επαρκές για  
«96» εξετάσεις/Contenido suficiente para  
<96> ensayos/Kogusest piisab <96> testi  
lääbiviimiseks/Contenu suffisant pour «96»  
tests/Sadržaj dovoljno za <96> testova/A  
doboz tartalma <96> vizsgálat elvégzéséhez  
elegendő/Contenuto sufficiente per «96»  
saggi/Turiny's skirtas atlikti <96> tyrimus/  
Saturis pietiekams <96> testiem/Inhoud  
voldoende voor «96» testen/til «96» test/  
Tilstrækkelig innhold for <96> prøver/  
Wystarczy na wykonanie <96> testów/  
Conteúdo suficiente para «96» ensaios/  
Conținut suficient pentru 96 de teste/  
Содержит достаточные количества для  
«96» определений/Innehåller tillräckligt  
till «96» antal tester/Obsah postačuje na  
tento počet testov: <96>/Vsebinsa zadostuje  
za <96> testov/Sadržina dovoljna za <96>  
testova/<96> testleri için yeterlilik içerir

**REF**

Catalogue number/Каталожен номер/  
Katalogové číslo/Katalognummer/  
Bestellnummer/Αριθμός καταλόγου/  
Número de catálogo/Katalogoi number/  
Numéro de catalogue/Kataloški broj/  
Katalógusszám/Numero di catalogo/  
Katalogo numeris/Numurs katalogā/  
Catalogusnummer/Katalognummer/  
Numer katalogowy/Número do catálogo/  
Număr de catalog/Номер по каталогу/  
Produktnummer/Katalógové číslo/  
Kataloška številka/Kataloški broj/  
Katalog numarası

Consult Instructions for Use/  
Прочетете инструкцията за  
употреба/Konzultujte s návodem  
k použití/Se brugsanvisning/Siehe  
Gebrauchsanweisung/Συμβουλευτείτε  
της Οδηγίας σχετικά με τη χρήση/  
Consulte las instrucciones de uso/  
Vt kasutusjuhendit/Consulter le mode  
d'emploi/Pročítajte upute za uporabu/  
Olvassa el a használati utasítást/  
Consultare le istruzioni per l'uso/Dél  
naudojimo žiūrėkite instrukcijas/Izlasiet  
lietošanas instrukciju/Aadpleeg de  
instructies voor gebruik/Les instruksene  
for bruk/Sprawdzić w instrukcji użycia/  
Consulte as Instruções de Utilização/  
Consultați instrucțiunile de utilizare/  
Обратитесь к инструкции по  
применению/Se bruksanvisning/  
Prečítajte si návod na používanie/  
Pročítajte uputstvo za upotrebu/  
Kullanım Talimatlarını Bakınız

Contents of kit/Съдържание на набора/  
Obsah soupravy/Kittets indhold/Inhalt  
des Kits/Περιεχόμενα του κιτ/Contenido  
del kit/Komplekt sisaldab/Contenu du  
kit/Sadržaj opreme/A készlet tartalma/  
Contenuto del kit/Rinkinio turinys/  
Komplekta saturs/Inhoud van de set/  
Settets innhold/Zawartość zestawu/  
Conteúdo do kit/Conținutul setului/  
Компоненты набора/Kit innehåll/  
Obsah súpravy/Vsebina kompleta/Sadržaj  
opreme/Kitin içindekiler

Biological risks/Биологическа  
опасност/Biológická rizika/Biologisk  
fare/Biologische Gefahren/Βιολογικοί  
κίνδυνοι/Riesgos biológicos/  
Biolooilised ohud/Risques biologiques/  
Biolóskli rizici/Biológiai kockázatok/Rischi  
biologici/Biologinis pavojus/Biológiskais  
risks/Biologische risico's/Biologische  
risikoer/Zagroženie biologiczne/Riscos  
biológicos/ Biologisk risk/Pericole  
biologice/Биологическая опасность/  
Biologicky rizikové/Biológické riziká/  
Biolóskli rizici/Biyolojik riskler

Human/C човешки производ/Ľidské/  
Human/Human/ἄνθρωπος αναφοράς/  
Humano/Inimāritolu/Humaine/Ljudskog  
porjekla/Humán/Origine Umana/  
Žmogaus kilmės/Cilvēku izcelsmes/  
Human/Menneske/Ludzka/Humano/  
Origine umână/Человеческого  
происхождения/Human/Ludské/  
Humanega izvora/Ljudskog porekla/Ľnsan

From mouse/C миши производ/Myši/  
Fra mus/Maus/από ποντίκι/de ratón/  
Hiirtelt/De souris/Mišijeg porjekla/  
Egérböli/Murino/Pelés kilmés/No peles/  
Van muizen/Fra mus/Mysia/Do rato/De  
la șoareci/Мышиного происхождения/  
Från mus/Myšie/Mišjega izvora/Mišijeg  
porekla/Fareden

Bovine/C говежди производ/  
Hovēži/Bovin/Rind/από βοοειδή/  
Bovino/Veistelt/Bovine/Rogate stoke/  
Szarvasmarha/Bovino/Jaučio/No  
liellopa/Bovien/Bovin/Wolowy/Bovino/  
Origine bovină/крупного рогатого  
скота/Från ko/Hovädzie/Rogvega  
izvora/Rogate krupne stoke/Bovin

Reconstitute with/Пазтваряне с/  
Rozfeďte pomoci/Rekonstitueres med/  
Rekonstituieren mit/Ανασύσταση με/  
Reconstituir con/Lahjendamine/  
Reconstituer avec/Rekonstituiraite s/  
Feloldáshoz/Ricostituire con/Atkurti,  
ištirpdant su/Atšķaidīt ar/Reconstituite  
met/Rekonstitueres med/Odtworzyć  
za pomocą/Reconstituir com/A  
se reconstitui cu/Пастворить в/  
Rekonstituera med/Rozriedte pomocou/  
Rekonstituiraite z/s/Ponovno formiranje  
sa/Yeniden oluşturalur

Manufacturer/Производитель/Výrobce/  
Producent/Hersteller/Κατασκευαστής/  
Fabricante/Tootja/Fabricant/Proizvođač/  
Gyártó/Fabbricante/Gamintojas/  
Ražotājs/Fabrikant/Produsent/  
Producent/Fabricante/Producător/  
Производитель/Tilverkare/ Výrobca/  
Izdelovalec/Proizvođač/Üretici

# INSTRUCTIONS FOR USE

EN

## INSTRUCTIONS FOR USE

Please visit our website [www.fdi.com/ifu](http://www.fdi.com/ifu) to obtain the Instructions For Use (IFU) in additional languages.

To ensure that you download the correct IFU for your kit lot, please select the revision corresponding to the issue date printed on the front page of the IFU provided with this kit.

Please follow the IFU carefully. Instructions for safe handling are found in the WARNINGS AND PRECAUTIONS section. Material Safety Data Sheets (MSDS) are available on our website [www.fdi.com](http://www.fdi.com). If you do not have access to the internet, please contact your local distributor, or Fujirebio Diagnostics AB for assistance.

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# CanAg CA242 EIA

Instructions for use

Enzyme immunometric assay kit  
For 96 determinations

## INTENDED USE

The CanAg CA242 EIA kit is intended for the quantitative determination of CA242 cancer antigen in serum.

## SUMMARY AND EXPLANATION OF THE ASSAY

The tumor marker CA242 is defined by the monoclonal antibody C242. The chemical structure of the antigenic determinant is not exactly known, but the determinant have been shown to be a sialylated carbohydrate structure. In serum, CA242 is found on the same mucin-complex as CA50 and sialylated Lewis<sup>a</sup> (CA19-9). Thus, CA242 is related, but not identical to the epitope of CA19-9 (1, 2). Serum levels of CA242 are low in healthy subjects and subjects with benign diseases, while elevated levels are commonly found in serum from patients with gastro-intestinal cancer (3).

The CA242 marker may be used as an aid in the diagnosis and management of patients with known or suspected gastro-intestinal carcinomas (4-9). The CanAg CA242 EIA should not be used as a substitute for any established clinical examination of malignancy, but may be used as a complement to existing clinical and laboratory methods.

## PRINCIPLE OF THE TEST

The CanAg CA242 EIA is a solid-phase, non-competitive immunoassay based on the direct sandwich technique. Calibrators, controls and patient samples are incubated together with biotinylated anti-CA242 monoclonal antibody (MAb) C241 in Streptavidin coated microstrips. CA242 present in calibrators, controls or samples is adsorbed to the Streptavidin coated microstrips by the biotinylated anti-CA242 MAb during the incubation. The strips are then washed and incubated with horseradish peroxidase (HRP) labelled Anti-CA242 MAb C242. After washing, buffered Substrate/Chromogen reagent (hydrogen peroxide and 3, 3', 5, 5' tetra-methylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction a blue colour will develop if antigen is present. The intensity of the colour is proportional to the amount of CA242 antigen present in the samples.

The colour intensity is determined in a microplate spectrophotometer at 620 nm (or optionally at 405 nm after addition of Stop Solution). Calibration curves are

constructed for each assay by plotting absorbance value versus the concentration for each calibrator. The CA242 concentrations of patient samples are then read from the calibration curve.

## REAGENTS

- Each CanAg CA242 EIA kit contains reagents for 96 tests.
- The expiry date of the kit is stated on the label on the outside of the kit box.
- Do not use the kit beyond the expiry date.
- Do not mix reagents from different kit lots.
- Store the kit at 2–8° C. Do not freeze.
- Opened reagents are stable according to the table below provided they are not contaminated, stored in resealed original containers and handled as prescribed. Return to 2–8° C immediately after use.

| Component | Quantity | Storage and stability after first opening |
|-----------|----------|---|
|-----------|----------|---|

### MICROPLA

|                   |         |  |
|-------------------|---------|--|
| <b>Microplate</b> | 1 Plate | 2–8° C until expiry date stated on the plate |
|-------------------|---------|--|

12 x 8 wells coated with Streptavidin. After opening, immediately return unused strips to the aluminium pouch, containing desiccant. Reseal carefully to keep dry.

|                          |         |  |
|--------------------------|---------|--|
| <b>CA242 Calibrators</b> | 5 vials | 2–8° C until expiry date stated on the vials |
|--------------------------|---------|--|

|            |              |          |        |             |
|------------|--------------|----------|--------|-------------|
| <b>CAL</b> | <b>CA242</b> | <b>0</b> | 0 U/mL | 1 x 0.75 mL |
|------------|--------------|----------|--------|-------------|

|            |              |           |         |             |
|------------|--------------|-----------|---------|-------------|
| <b>CAL</b> | <b>CA242</b> | <b>15</b> | 15 U/mL | 1 x 0.75 mL |
|------------|--------------|-----------|---------|-------------|

|            |              |           |         |             |
|------------|--------------|-----------|---------|-------------|
| <b>CAL</b> | <b>CA242</b> | <b>50</b> | 50 U/mL | 1 x 0.75 mL |
|------------|--------------|-----------|---------|-------------|

|            |              |            |          |             |
|------------|--------------|------------|----------|-------------|
| <b>CAL</b> | <b>CA242</b> | <b>100</b> | 100 U/mL | 1 x 0.75 mL |
|------------|--------------|------------|----------|-------------|

|            |              |            |          |             |
|------------|--------------|------------|----------|-------------|
| <b>CAL</b> | <b>CA242</b> | <b>150</b> | 150 U/mL | 1 x 0.75 mL |
|------------|--------------|------------|----------|-------------|

Human CA242 antigen in a Tris-HCl buffered salt solution containing bovine serum albumin, an inert yellow dye, and 0.05% NaN<sub>3</sub> as preservative. Ready for use.

| Component | Quantity | Storage and stability after first opening |
|-----------|----------|---|
|-----------|----------|---|

**CA242 Controls** 2 vials 2–8° C until expiry date stated on the vials

|                |              |          |
|----------------|--------------|----------|
| <b>CONTROL</b> | <b>CA242</b> | <b>1</b> |
|----------------|--------------|----------|

1 x 0.75 mL

|                |              |          |
|----------------|--------------|----------|
| <b>CONTROL</b> | <b>CA242</b> | <b>2</b> |
|----------------|--------------|----------|

1 x 0.75 mL

Human CA242 antigen in a Tris-HCl buffered salt solution containing bovine serum albumin, and 0.05% NaN<sub>3</sub> as preservative. Ready for use.

|               |                   |
|---------------|-------------------|
| <b>BIOTIN</b> | <b>Anti-CA242</b> |
|---------------|-------------------|

**Biotin Anti-CA242** 1 x 15 mL 2–8° C until expiry date stated on the vial

Biotin Anti-CA242 monoclonal antibody from mouse, approximately 1.5 µg/mL. Contains Tris-HCl buffered saline (pH 7.75) with bovine serum albumin, bovine immunoglobulin, blocking agents, detergent, an inert red dye and 0.05% NaN<sub>3</sub> as preservative. Ready for use.

|             |                   |
|-------------|-------------------|
| <b>CONJ</b> | <b>Anti-CA242</b> |
|-------------|-------------------|

**Tracer, HRP Anti-CA242** 1 x 0.75 mL 2–8° C until expiry date stated on the vial

Stock solution of HRP Anti-CA242 monoclonal antibody from mouse, approximately 40 µg/mL. Contains preservatives. To be diluted with Tracer Diluent prior to use.

|            |             |
|------------|-------------|
| <b>DIL</b> | <b>CONJ</b> |
|------------|-------------|

**Tracer Diluent** 1 x 15 mL 2–8° C until expiry date stated on the vial

Phosphate buffered saline (pH 7.2) with bovine serum albumin, bovine immunoglobulin, blocking agents, detergent, an inert blue dye, and 0.01 % methyl-isothiazolone (MIT) as preservative. Ready for use.

| Component | Quantity | Storage and stability after first opening |
|-----------|----------|---|
|-----------|----------|---|

|      |     |
|------|-----|
| SUBS | TMB |
|------|-----|

|                          |           |   |
|--------------------------|-----------|---|
| <b>TMB HRP-Substrate</b> | 1 x 12 mL | 2–8° C until expiry date stated on the vial |
|--------------------------|-----------|---|

Contains buffered hydrogen peroxide and 3, 3', 5, 5' tetramethyl-benzidine (TMB). Ready for use.

|      |
|------|
| STOP |
|------|

|                      |           |   |
|----------------------|-----------|---|
| <b>STOP Solution</b> | 1 x 15 mL | 2–8° C until expiry date stated on the vial |
|----------------------|-----------|---|

Contains 0.12 M hydrochloric acid. Ready for use.

|         |     |
|---------|-----|
| WASHBUF | 25X |
|---------|-----|

|                         |           |   |
|-------------------------|-----------|---|
| <b>Wash Concentrate</b> | 1 x 50 mL | 2–8° C until expiry date stated on the bottle |
|-------------------------|-----------|---|

A Tris-HCl buffered salt solution with Tween 20. Contains Germall II as preservative. To be diluted with water 25 times before use.

### Indications of instability

The TMB HRP-Substrate should be colourless or slightly bluish. A blue colour indicates that the reagent has been contaminated and should be discarded.

### WARNINGS AND PRECAUTIONS

#### For in vitro diagnostic use.

- For professional use only.
- Please refer to the US Department of Health and Human Services (Bethesda, Md., US) publication No. (CDC) 88-8395 on laboratory safety or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Reagents contain sodium azide ( $\text{NaN}_3$ ) as a preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.
- Follow local guidelines for disposal of all waste material.

# Protocol Sheet

**CanAg CA242 EIA** REF **101-10**

Mix the components directly before use. Use shaking conditions according to the Instructions.

| Step  | Vial/Plate  | Procedure  |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
|---|---|--|---------------|-----------------------------|---------------------|---|----|---|---|-----|---|---|-----|---|---|-----|---|---|-----|---|---|-----|---|---|-----|---|---|-----|---|---|-----|---|----|-----|----|----|-----|----|----|-----|----|
| 1. Prepare Wash Solution<br><br>Prepare Tracer working solution | <b>WASHBUF</b> <b>25X</b>                               | Dilute 50 mL of Wash Concentrate with 1200 mL of distilled or deionized water.   |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
|   | <b>CONJ</b> <b>Anti-CA242</b><br><b>DIL</b> <b>CONJ</b> | Mix 50 µL of Tracer, HRP Anti-CA242 with 1 mL of Tracer Diluent per strip:   |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
|   |   | <table border="1"><thead><tr><th>No. of Strips</th><th>Tracer, HRP Anti-CA242 (µL)</th><th>Tracer Diluent (mL)</th></tr></thead><tbody><tr><td>1</td><td>50</td><td>1</td></tr><tr><td>2</td><td>100</td><td>2</td></tr><tr><td>3</td><td>150</td><td>3</td></tr><tr><td>4</td><td>200</td><td>4</td></tr><tr><td>5</td><td>250</td><td>5</td></tr><tr><td>6</td><td>300</td><td>6</td></tr><tr><td>7</td><td>350</td><td>7</td></tr><tr><td>8</td><td>400</td><td>8</td></tr><tr><td>9</td><td>450</td><td>9</td></tr><tr><td>10</td><td>500</td><td>10</td></tr><tr><td>11</td><td>550</td><td>11</td></tr><tr><td>12</td><td>600</td><td>12</td></tr></tbody></table> | No. of Strips | Tracer, HRP Anti-CA242 (µL) | Tracer Diluent (mL) | 1 | 50 | 1 | 2 | 100 | 2 | 3 | 150 | 3 | 4 | 200 | 4 | 5 | 250 | 5 | 6 | 300 | 6 | 7 | 350 | 7 | 8 | 400 | 8 | 9 | 450 | 9 | 10 | 500 | 10 | 11 | 550 | 11 | 12 | 600 | 12 |
| No. of Strips   | Tracer, HRP Anti-CA242 (µL)                             | Tracer Diluent (mL)  |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
| 1   | 50  | 1  |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
| 2   | 100   | 2  |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
| 3   | 150   | 3  |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
| 4   | 200   | 4  |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
| 5   | 250   | 5  |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
| 6   | 300   | 6  |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
| 7   | 350   | 7  |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
| 8   | 400   | 8  |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
| 9   | 450   | 9  |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
| 10  | 500   | 10   |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
| 11  | 550   | 11   |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
| 12  | 600   | 12   |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |
| 2. Wash   | <b>MICROPLA</b>   | Wash each well once with Wash Solution. Use manual or automatic washer.  |               |                             |                     |   |    |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |   |     |   |    |     |    |    |     |    |    |     |    |

|  |  |                          |  |                          |
|--|--|--------------------------|--|--------------------------|
| 3. Add calibrators, controls and samples | <table border="1"> <tr> <td>CAL</td> <td>CA242</td> </tr> </table>         | CAL                      | CA242  | 25 $\mu$ L in each well  |
| CAL                                      | CA242  |                          |  |                          |
|  | 0, 15, 50, 100, 150  |                          |  |                          |
|  | <table border="1"> <tr> <td>CONTROL</td> <td>CA242</td> </tr> </table>     | CONTROL                  | CA242  |                          |
| CONTROL                                  | CA242  |                          |  |                          |
|  | 1, 2   |                          |  |                          |
| 4. Add Biotin Anti-CA242                 | <table border="1"> <tr> <td>BIOTIN</td> <td>Anti-CA242</td> </tr> </table> | BIOTIN                   | Anti-CA242   | 100 $\mu$ L in each well |
| BIOTIN                                   | Anti-CA242   |                          |  |                          |
| 5. Incubate                              | <table border="1"> <tr> <td>MICROPLA</td> </tr> </table>                   | MICROPLA                 | 2 hours shaking at room temperature  |                          |
| MICROPLA                                 |  |                          |  |                          |
| 6. Wash                                  | <table border="1"> <tr> <td>MICROPLA</td> </tr> </table>                   | MICROPLA                 | Wash each well three times with Wash Solution. Use manual or automatic washer. |                          |
| MICROPLA                                 |  |                          |  |                          |
| 7. Add Tracer working solution           | TRACER WORKING SOLUTION  | 100 $\mu$ L in each well |  |                          |
| 8. Incubate                              | <table border="1"> <tr> <td>MICROPLA</td> </tr> </table>                   | MICROPLA                 | 1 hour shaking at room temperature   |                          |
| MICROPLA                                 |  |                          |  |                          |
| 9. Wash                                  | <table border="1"> <tr> <td>MICROPLA</td> </tr> </table>                   | MICROPLA                 | Wash each well six times with Wash Solution. Use manual or automatic washer.   |                          |
| MICROPLA                                 |  |                          |  |                          |
| 10. Add TMB HRP-Substrate                | <table border="1"> <tr> <td>SUBS</td> <td>TMB</td> </tr> </table>          | SUBS                     | TMB  | 100 $\mu$ L in each well |
| SUBS                                     | TMB  |                          |  |                          |
| 11. Incubate                             | <table border="1"> <tr> <td>MICROPLA</td> </tr> </table>                   | MICROPLA                 | 30 min shaking at room temperature   |                          |
| MICROPLA                                 |  |                          |  |                          |
| 12. Read absorbance                      | <table border="1"> <tr> <td>MICROPLA</td> </tr> </table>                   | MICROPLA                 | 620 nm   |                          |
| MICROPLA                                 |  |                          |  |                          |
| Alt.12 Add Stop Solution                 | <table border="1"> <tr> <td>STOP</td> </tr> </table>                       | STOP                     | 100 $\mu$ L in each well   |                          |
| STOP                                     |  |                          |  |                          |
| Alt.13 Incubate                          | <table border="1"> <tr> <td>MICROPLA</td> </tr> </table>                   | MICROPLA                 | 1 min shaking at room temperature  |                          |
| MICROPLA                                 |  |                          |  |                          |
| Alt.14 Read absorbance                   | <table border="1"> <tr> <td>MICROPLA</td> </tr> </table>                   | MICROPLA                 | Read at 405 nm within 5 min  |                          |
| MICROPLA                                 |  |                          |  |                          |

## Caution

Material used in the preparation of human source reagent has been tested and found to be Non Reactive for HIV-1/2 Antibody, HCV Antibody and Hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they were potentially infectious.

## CLP (1272/2008) HAZARD CLASSIFICATION

Information about CLP (1272/2008) HAZARD CLASSIFICATION can be found at the end of this document.

## SPECIMEN COLLECTION AND HANDLING

The CanAg CA242 EIA is intended for use with serum. Collect blood by venipuncture and separate the serum according to common procedures. Samples can be stored at 2–8° C for 24 hours. For longer periods it is recommended to store the samples at –20° C or below. Avoid repeated freezing and thawing of the samples. Allow frozen samples to thaw slowly, preferably at 2–8° C over night and then bring the samples to room temperature before analysis.

## PROCEDURE

### Materials required but not supplied with the kit

- 1. Microplate shaker**  
Shaking should be medium to vigorous, approximately 700-1100 oscillations/min.
- 2. Microplate wash device**  
Automatic plate washer capable of performing 1, 3 and 6 washing cycles with a minimal fill volume of 350 µL/well/washcycle.  
An 8-channel pipette with disposable plastic tips for delivery of 350 µL is recommended if an automatic microplate washer is not used.
- 3. Microplate spectrophotometer**  
With a wavelength of 620 nm and/or 405 nm and an absorbance range of 0 to 3.0.
- 4. Precision pipettes**  
With disposable plastic tips to deliver microlitre and millilitre volumes.  
An 8-channel pipette or respenser pipette with disposable plastic tips for delivery of 100 µL is useful but not essential.
- 5. Distilled or deionized water**  
For preparation of Wash Solution.

## Procedural notes

1. A thorough understanding of this package insert is necessary to ensure proper use of the CanAg CA242 EIA kit. The reagents supplied with the kit are intended for use as an integral unit. Do not mix identical reagents from kits having different lot numbers. Do not use the kit reagents after the expiry date printed on the outside of the kit box.
2. Reagents should be allowed to reach room temperature (20–28° C) prior to use. The assay should only be performed at temperatures between 20–28° C to obtain accurate results. Frozen specimens should be brought to room temperature slowly and must be gently but thoroughly mixed after thawing.
3. Before starting to pipette calibrators, controls and patient specimens it is advisable to mark the strips to be able to clearly identify the samples during and after the assay.
4. The requirement for efficient and thorough washing for separation of bound and unbound antigen and reagents from the solid-phase bound antibody-antigen complexes is one of the most important steps in an EIA. In order to ensure efficient washing make sure that all wells are completely filled to the top edge with wash solution during each wash cycle, that wash solution is dispensed at a good flow rate, that the aspiration of the wells between and after the wash cycles is complete and that the wells are empty. If there is liquid left, invert the plate and tap it carefully against absorbent paper.
  - Automatic strip washer: Follow the manufacturer's instructions for cleaning and maintenance diligently and wash the required number of wash cycles prior to and after each incubation step. It's highly recommended to use *strip* process mode and *overflow* wash mode with a dispensing volume of 800 µL. The aspiration/wash device should not be left standing with the Wash Solution for long periods, as the needles may get clogged resulting in poor liquid delivery and aspiration.
5. The TMB HRP-Substrate is very sensitive for contamination. For optimal stability of the TMB HRP-Substrate, pour the required amount from the vial to a carefully cleaned reservoir or preferably a disposable plastic tray to avoid contamination of the reagent. Be sure to use clean disposable plastic pipette tips (or respenser pipette tip).
6. Be sure to use clean disposable plastic pipette tips and a proper pipetting technique when handling samples and reagents. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid. A proper pipetting technique is of particular importance when handling the TMB HRP-Substrate.

| Preparation of reagents | Stability of prepared reagent            |
|-------------------------|--|
| <b>Wash Solution</b>    | 2 weeks at 2–25° C in a sealed container |

Pour the 50 mL Wash Concentrate into a clean container and dilute 25- fold by adding 1200 mL of distilled or deionized water to give a buffered Wash Solution.

|                                |   |
|--------------------------------|---|
| <b>Tracer Working Solution</b> | 3 weeks at 2–8° C in a sealed container |
|--------------------------------|---|

Prepare the required quantity of Tracer working solution by mixing 50 µL of Tracer, HRP Anti-CA242 with 1 mL of Tracer Diluent per strip (see table below):

| No. of Strips | Tracer, HRP Anti-CA242 (µL) | Tracer Diluent (mL) |
|---------------|-----------------------------|---------------------|
| 1             | 50                          | 1                   |
| 2             | 100                         | 2                   |
| 3             | 150                         | 3                   |
| 4             | 200                         | 4                   |
| 5             | 250                         | 5                   |
| 6             | 300                         | 6                   |
| 7             | 350                         | 7                   |
| 8             | 400                         | 8                   |
| 9             | 450                         | 9                   |
| 10            | 500                         | 10                  |
| 11            | 550                         | 11                  |
| 12            | 600                         | 12                  |

Be sure to use a clean plastic or glass bottle for preparation of the Tracer Working Solution.

**Alternative:** Pour the content of the Tracer, HRP Anti-CA242 into the vial of Tracer Diluent and mix gently. Make sure that all of the Tracer, HRP Anti-CA242 is transferred to the vial of Tracer Diluent.

**NOTE:** The Tracer Working Solution is stable for 3 weeks at 2–8° C. Do not prepare more Tracer Working Solution than will be used within this period and make sure that it is stored properly.

## Assay procedure

Perform each determination in duplicate for calibrators, controls and patient samples. A calibration curve should be run with each assay. All reagents and samples must be brought to room temperature (20–28°C) before use.

1. Start to prepare Wash Solution and Tracer Working Solution. It is important to use clean containers. Follow the instructions carefully.
2. Transfer the required number of microplate strips to a strip frame. (Immediately return the remaining strips to the aluminium pouch containing a desiccant and reseal carefully). Wash each strip once with the Wash Solution. Do not wash more strips than can be handled within 30 min.
3. Pipette 25 µL of the CA242 Calibrators (CAL 0, 15, 50, 100, 150), Controls (C1, C2) and patient samples (unknowns Unk) into the strip wells according to the following scheme:

|   | 1          | 2          | 3     | 4 | 5 | 6 | 7 etc |
|---|------------|------------|-------|---|---|---|-------|
| A | Cal<br>0   | Cal<br>150 | Unk 2 |   |   |   |       |
| B | Cal<br>0   | Cal<br>150 | Unk 2 |   |   |   |       |
| C | Cal<br>15  | C1         | etc.  |   |   |   |       |
| D | Cal<br>15  | C1         |       |   |   |   |       |
| E | Cal<br>50  | C2         |       |   |   |   |       |
| F | Cal<br>50  | C2         |       |   |   |   |       |
| G | Cal<br>100 | Unk 1      |       |   |   |   |       |
| H | Cal<br>100 | Unk 1      |       |   |   |   |       |

4. Add 100 µL of Biotin Anti-CA242 to each well using a 100 µL precision pipette (or an 8-channel 100 µL precision pipette). Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or the surface of the liquid.
5. Incubate the frame containing the strips for 2 hours (± 10 min) at room temperature (20–28°C) with constant shaking of the plate using a microplate shaker.
6. After the first incubation aspirate and wash each strip 3 times using the wash procedure described in Procedural notes, item 4.

7. Add 100  $\mu$ L of Tracer working solution to each well. Use the same pipetting procedure as in item 4 above.
8. Incubate the frame for 1 hour ( $\pm$  5 min) at room temperature (20–28°C) with constant shaking.
9. After the second incubation aspirate and wash each strip 6 times, using the wash procedure described in Procedural notes, item 4.
10. Add 100  $\mu$ L of TMB HRP-Substrate to each well using the same pipetting procedure as in item 4. The TMB HRP-Substrate should be added to the wells as quickly as possible and the time between the addition to the first and last well should not exceed 5 min.
11. Incubate for 30 min ( $\pm$  5 min) at room temperature with constant shaking. Avoid direct sunlight.
12. Immediately read the absorbance at 620 nm in a microplate spectrophotometer.

### Option

If the laboratory does not have access to a microplate spectrophotometer capable of reading at 620 nm, the absorbance can be determined as follows:

Alt. 12. Add 100  $\mu$ L of Stop Solution. Mix and read the absorbance at 405 nm in a microplate spectrophotometer within 5 minutes after addition of Stop Solution.

### Measurement range

The CanAg CA242 EIA measures concentrations between 1 and 150 U/mL. If CA242 concentrations above the measuring range are to be expected, it is recommended to dilute samples with normal human serum prior to analysis. **NOTE:** The serum used for dilution must also be measured in order to determine the endogenous CA242 concentration (see “Calculation of results”).

### Quality control

CA242 Control 1 and 2 may be used for validation of the assay series. Ranges of expected results are indicated on the vial labels. If values outside of the specified range are obtained, a complete check of reagents and reader performance should be made and the analysis repeated. Each laboratory may in addition prepare its own serum pools at different levels, which can be used as internal controls in order to assure the precision of the assay.

## Reference material

Since no common reference material is available for CA242 antigen, CanAg CA242 EIA Calibrator values are assigned against a set of in-house reference standards.

## CALCULATION OF RESULTS

If a microplate spectrophotometer reader with built-in data calculation program is used, refer to the manual for the plate reader and create a program using the concentration stated on the labels of each of the CA242 Calibrators.

For automatic calculation of CA242 results it is recommended to use either of the following methods:

- Cubic spline curve fit method. Calibrator 0 should be included in the curve with the value 0 U/mL.
- Spline smoothed curve fit method. Calibrator 0 should be used as plate blank.
- Interpolation with point-to-point evaluation. Calibrator 0 should be included in the curve with the value 0 U/mL.
- Quadratic curve fit method. Calibrator 0 should be included in the curve with the value 0 U/mL.

**NOTE:** 4-parametric or linear regression should not be used.

For manual evaluation, a calibration curve is constructed by plotting the absorbance (A) values obtained for each CA242 Calibrator against the corresponding CA242 concentration (in U/mL), see figure below. The unknown CA242 concentrations can then be read from the calibration curve using the mean absorbance value of each patient specimen.


If samples in an initial analysis give CA242 levels higher than 150 U/mL, the samples should be diluted 1/10 with normal human serum and reanalysed to obtain the accurate CA242 concentration. **NOTE:** The sample used for dilution must also be measured in order to determine the endogenous CA242 concentration.

The CA242 concentration of the undiluted sample is calculated as:

$$\text{Dilution 1/10: } 10 \times ([\text{CA242}]_{\text{Diluted sample}}^{-0.9 \times [\text{CA242}]_{\text{Normal serum}}})$$

## Example of results

| Specimen   |       |     | Calibrator values | Mean abs value (A) | CA242 (U/mL) |
|------------|-------|-----|-------------------|--------------------|--------------|
| CAL        | CA242 | 0   | 0 U/mL            | 0.050              |              |
| CAL        | CA242 | 15  | 15 U/mL           | 0.390              |              |
| CAL        | CA242 | 50  | 50 U/mL           | 1.107              |              |
| CAL        | CA242 | 100 | 100 U/mL          | 1.922              |              |
| CAL        | CA242 | 150 | 150 U/mL          | 2.617              |              |
| Specimen A |       |     |                   | 0.410              | 16.1         |
| Specimen B |       |     |                   | 1.636              | 80.9         |

*Example (do not use this curve or table above to determine actual assay results).*

## LIMITATIONS OF THE PROCEDURE

The level of CA242 cannot be used as absolute evidence for the presence or absence of malignant disease, and the CA242 test should not be used in cancer screening. The results of the test should be interpreted only in conjunction with other investigations and procedures in the diagnosis of disease and the management of patients, and the CA242 test should not replace any established clinical examination.

Anti-reagent antibodies (human anti-mouse antibody (HAMA) or heterophilic antibodies) in the patient sample may occasionally interfere with the assay, even though specific blocking agents are included in the buffer.

Biotin may interfere with the assay giving false low results. This should be taken into consideration for patients taking dietary supplements or receiving therapy containing high (> 5mg/day) or extremely high (300 mg/day) biotin doses. Peak serum levels have been reported to occur 1-3h post ingestion and a physiological half-life of 8-16h depending on renal function (12,13). In a study by Grimsey et al. (12) a specimen concentration of 30 ng/mL was reached 8h following intake of 10 mg of biotin. For patients on very high doses of biotin it is recommended to stop taking biotin for at least 2 days before blood draw (14).

## EXPECTED VALUES

Serum samples from 184 apparently healthy blood donors, 97 women and 87 men were analysed resulting in a mean value of  $8.5 \pm 7.6$  U/mL. The lower and upper extremes of the normal range were examined using IFCC recommended non-parametric statistical treatment. The reference interval contains the central 95% fraction of the reference distribution. Reference limits may accordingly be estimated as the 2.5% (lower) and 97.5% (upper) fractiles. These limits cut off a fraction of 2.5% of the values in each tail of the reference distribution. Non-parametric estimates:

| Fraction                   | Reference limit (U/mL) | 95% confidence interval |
|----------------------------|------------------------|-------------------------|
| 2.5 <sup>th</sup> (lower)  | 0                      | 0-0                     |
| 97.5 <sup>th</sup> (upper) | 29                     | 25-44                   |

93% of the healthy subjects had assay values  $\leq 20$  U/mL.

It is recommended that each laboratory establish their own normal range to account for such local environmental factors as diet, climate, living conditions, patient selection, etc. It should also be borne in mind that the individual patient's own baseline result provides the most important reference point for interpretation of marker results (9).

## PERFORMANCE CHARACTERISTICS

### Precision

Total precision was determined according to NCCLS guideline EP5-A (10) using four levels of frozen pooled human serum containing added CA242 from patients with gastro-intestinal cancer. Each sample was randomly pipetted (n=2/analysis) and analysed twice each day over 20 days. The analyses were undertaken during a period of 42 months, by  $\geq$  two different technicians and using 20 different CanAg CA242 EIA kit batches.

| Sample  | Replicates | Mean (U/mL) | Within-run SD (U/mL) | Within-run CV % | Between-day SD (U/mL) | Between-day CV % |
|---------|------------|-------------|----------------------|-----------------|-----------------------|------------------|
| CA242 1 | 80         | 16.2        | 0.67                 | 4.1             | 0.39                  | 2.4              |
| CA242 2 | 80         | 48.4        | 1.93                 | 4.0             | 1.82                  | 3.8              |
| CA242 3 | 80         | 79.5        | 2.99                 | 3.8             | 2.46                  | 3.1              |
| CA242 4 | 80         | 125         | 5.81                 | 4.7             | 2.74                  | 2.2              |

### Detection limit

The detection limit of the CanAg CA242 EIA is  $< 1$  U/mL defined as the concentration corresponding to the mean of the absorbance values of the CA242 calibrator 0 plus 2 standard deviations according to formula:

$$\frac{2 \times \text{SD CAL } 0}{\text{OD CAL } 15 - \text{OD CAL } 0} \times 15 \text{ U/mL}$$

### Recovery

Spiked serum samples were prepared by adding human CA242 antigen to normal serum samples. The recovery of the added antigen was in the range 87 –107 %.

**NOTE:** recovery studies should **not** be performed using the kit calibrators.

### Hook effect

No hook effect has been noticed with samples up to 150 000 U/mL. **NOTE:** In very high samples the colour of the substrate will change from blue to greenish (and eventually yellow in extremely high samples). This will lead to a falsely low absorbance at 620 nm, and in extreme cases the absorbance may fall within the calibration curve range and noticed as a hook.

### Linearity

Patient samples were serially diluted with normal lipid stripped human serum and analysed. The obtained values were in the range 97–108%.

## Specificity

The CanAg CA242 EIA is based on two mouse monoclonal antibodies, the catching MAb C241 targeting sialylated Lewis<sup>a</sup> and the detecting MAb C242 specific for the CA242 epitope. Thus the assay determines S-Le<sup>a</sup> containing mucin antigen(s) expressing the CA242 epitope. The NCCLS guideline EP7-P (11) was followed to determine possible sources of interference. The following substances and concentrations were tested and found not to interfere with the test.

|                         | Concentration with no significant<br>(± 10%) interference |
|-------------------------|---|
| Lipemia (Intralipid®)   | 8 mg/mL   |
| Bilirubin, unconjugated | 0.6 mg/mL   |
| Hemoglobin              | 5 mg/mL   |

## Biotin interference

A study was conducted to evaluate biotin interference. Low and high serum control samples were spiked to final biotin concentrations of 15, 30, 60, and 200 ng/mL. The mean CA242 concentration was determined for each sample and the percent recovery for each biotin concentration was calculated using the formula: Recovery (%) = 100x (Mean CA242 concentration w. biotin added/Mean CA242 concentration w. diluent only added).

| CA242 analyte level | Biotin test conc. (ng/mL) | Expected CA242 conc. (U/mL) | Observed CA242 conc. (U/mL) | Recovery (%) |
|---------------------|---------------------------|-----------------------------|-----------------------------|--------------|
| Low                 | 15                        | 16,5                        | 16,3                        | 99           |
| High                | 15                        | 122                         | 121                         | 99           |
| Low                 | 30                        | 15,3                        | 14,9                        | 97           |
| High                | 30                        | 118                         | 115                         | 98           |
| Low                 | 60                        | 15,2                        | 13,8                        | 91           |
| High                | 60                        | 117                         | 109                         | 94           |
| Low                 | 200                       | 16,8                        | 11,1                        | 66           |
| High                | 200                       | 125                         | 86,7                        | 70           |

Based on linear regression analysis, the lowest concentration of biotin found to influence test results (≥10%) was 72 ng/mL.

## Method comparison

The CanAg CA242 EIA was compared to the CA242 Delfia. 145 serum samples from healthy blood donors and from patients with malignant and non-malignant diseases, ranging in values from 0-250 U/mL were analysed and linear regression analyses of the results yielded (4):

$$\text{CanAg CA242 EIA} = 1.02 \times \text{CA242 Delfia} - 1.1 \quad r=0.99$$

## CLP (1272/2008) HAZARD CLASSIFICATION

The following warnings and precautions apply to

SUBS

TMB

### Hazard pictograms:

### Signal word:

Danger

### Hazard Statement:

Repr. 1B: H360D May damage the unborn child.

### Prevention statement:

P202 Do not handle until all safety precautions have been read and understood.

### Prevention:

P280 Wear protective gloves / protective clothing / eye protection / face protection.

### Precautionary statement response:

P308+P313 IF exposed or concerned get medical advice/attention.

### Precautionary statement disposal:

P501 Dispose of contents / container to an approved hazardous / special waste disposal facility in accordance with local and national regulations.

### Restricted to professional users.

**Hazardous substances:** 2- Pyrrolidone

### Other hazards

None of the mixtures in the kit contains any substances considered to meet the criteria classifying them as PBT and/or vPvB.

## WARRANTY

The performance data presented here were obtained using the assay procedure indicated. Any change or modification of the procedure not recommended by Fujirebio Diagnostics may affect the results, in which event Fujirebio Diagnostics disclaims all warranties expressed, implied or statutory including the implied warranty of merchantability and fitness for use.

## REFERENCES

1. Johansson, C., Nilsson, O., Baeckström, D., Jansson, E.-L. and Lindholm, L. (1991) Novel Epitopes on the CA50-Carrying Antigen: Chemical and Immunochemical Studies, *Tumor Biol.*, 12, 159-179.
2. Johansson, C., Nilsson, O. and Lindholm, L., (1991) Comparison of Serological Expression of Different Epitopes on the CA50-Carrying Antigen CanAg, *Int. J. Cancer*, 48, 757-763.
3. Nilsson, O., Johansson, C., Glimelius, B., Persson B., Nørgaard-Pedersen, B., Andrén-Sandberg, Å., and Lindholm L. (1992): Sensitivity and specificity of CA242 in gastro-intestinal cancer. A comparison with CEA, CA50 and CA19-9. *Br J Cancer* 65, 215-221.
4. Dahlén U., Karlsson B., Lindholm L., Nilsson O., (1993) Development of an enzyme immuno-assay for determination of the tumour associated antigen CA242, *J. Tumor Marker Oncology* 8, 3, p 111.
5. Kawa, S., Tokoo, M., Hasebe, O., Hayashi, K., Imai, H., Oguchi, H., Kiyosawa, K., Furuta, S., and Homma, T.,(1994) Comparative study of CA242 and CA19-9 for the diagnosis of pancreatic cancer, *Br. J. Cancer*, 70, 481-486.
6. Von Kleist, Hesse, Kananeeh, (1996) Comparative Evaluation of Four Tumor Markers, CA 242, CA 19-9, TPA, and CEA in Carcinomas of the Colon, *Anticancer Research* 16: 2325-2332.
7. Spila A, Ferroni P, Cosimelli M, D'Alessandro R, Abbolito MR, Mariotti S, Aloe S, Carone MD, Graziano F, Tedesco M, Martini F, Mancini R, Stigliano V, Roselli M, Guadagni F. (2001) Comparative analysis of CA 242 and CA 19-9 serum tumor markers in colorectal cancer patients. A longitudinal evaluation. *Anticancer Res. Mar-Apr;21(2B):1263-70.*
8. Carpelan-Holmström, M., Haglund, K., Lundin, J., Järvinen, H and Roberts, P., (1996) Pre-operative serum levels of CA 242 and CEA predict outcome in colorectal cancer. *Eur. J. Cancer* 32(7), 1156-1161.
9. Engarås B. (2003 ) Individual cutoff levels of carcinoembryonic antigen and CA 242 indicate recurrence of colorectal cancer with high sensitivity. *Dis Colon Rectum.* Mar;46(3):313-21.
10. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
11. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).
12. Grimsp P et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. *J Pharmacokinet.* 2017; 2(4); 247-256.

13. Jenkins Colon P, Greene D.N. Biotin Interference in Clinical Immunoassays. *JALM* 2018; 2(6); 941-951.
14. Chun KY. Biotin Interference in Diagnostic Tests. *Clinical Chemistry* 2017; 63; 619-620.

---

**CanAg<sup>®</sup>** is a registered trademark of Fujirebio Diagnostics AB

**Fujirebio Diagnostics AB**

**Elof Lindälvs gata 13**

**SE-414 58 Göteborg**

**Sweden**

**Phone + 46 31 85 70 30**

**Fax + 46 31 85 70 40**

**info@fdab.com**

**www.fdab.com**