# CanAg PSA EIA

REF

340-10

IVD

CE<sub>0197</sub>

Instructions for use. 2022-06

Read highlighted changes

EN	EXPLANATION OF SYMBOLS
BG	ОБЯСНЕНИЕ НА СИМВОЛИТЕ
CS	VÝZNAM SYMBOLŮ
DA	SYMBOLFORKLARING
DE	ERKLÄRUNG DER SYMBOLE
EL	ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
ES	SIGNIFICADO DE LOS SÍMBOLOS
ET	SÜMBOLITE SELGITUS
FR	EXPLICATION DES SYMBOLES
HR	OBJAŠNJENJE SIMBOLA
HU	JELMAGYARÁZAT
IT	SPIEGAZIONE DEI SIMBOLI
LT	SIMBOLIŲ PAAIŠKINIMAI
LV	SIMBOLU SKAIDROJUMS
NL	VERKLARING DER SYMBOLEN
NO	SYMBOLFORKLARING
PL	OBJAŚNIENIE SYMBOLI
PT	EXPLICAÇÃO DOS SÍMBOLOS
RO	SEMNIȚAȚIA SIMBOLURILOR
RU	ОБОЗНАЧЕНИЯ
SV	SYMBOLFÖRKLARING
SK	VÝZNAM SYMBOLOV
SL	RAZLAGA SIMBOLOV
SR	OBJAŠNJENJE SIMBOLA
TR	SEMBOLLERİN AÇIKLAMALARI

Use By/Годно до/Použitelné do/  
Holdbar til/Verwendbar bis/  
Ημερομηνία λήξης/Fecha  
de caducidad/Kölblik kuni/  
Utiliser jusque/Rok valjanosti/  
Felhasználható/Utilizzare entro/  
Sunautoti iki/Izlietot līdz/Houdbaar  
tot/Brukes innen/Użyç przed/  
Prazo de validade/Expířã la/  
Использовать до/Använd före/  
Použite'né do/ Uporabno do/  
Upotrebljivo do/Son Kullanna Tarihi

LOT

Batch code/Номер на партида/  
Číslo šarže/Lotnummer/  
Chargenbezeichnung/Αριθμός  
Παρτίδας/Código de lote/Partii  
kood/Code du lot/Kod serije/  
Sarzsám/Codice del lotto/  
Partijos kodus/Partijas kods/Lot  
nummer/Partikode/Kod partii/  
Código do lote/Număr de lot/  
Номер лота/Lotnummer/Číslo  
šarže/Številka serije/Kod partije/  
Parti Kodu

Date of manufacture/Дата на производство/Datum výroby/  
Produktionsdato/Herstellungsdatum/  
Ημερομηνία παραγωγής/Fecha de fabricación/Valmistamise kuupäev/  
Date de fabrication/Datum proizvodnje/  
Gyártási idő/Data di produzione/  
Pagaminimo data/Ražošanas datums/  
Productiedatum/Fremstillingsdato/  
Data produkcyj/Data de fabrico/Data fabricației/Дата производства/  
Tillverkningsdatum/Dátum výroby/Datum izdelave/Datum proizvodnje/Üretim tarihi

Temperature limitation/  
Температурни граници/  
Teplotní omezení/  
Temperaturbegrensning/  
Temperaturbegrenzung/  
Περιορισμοί θερμοκρασίας/  
Limites de temperatura/  
Temperatuuri piirang/  
Limite de température/  
Temperaturno ograničenje/  
Hőmérsékletre vonatkozó korlátozás/  
Limiti di temperatura/  
Temperatūriņa ierobežojumi/  
Temperatūras ierobežojums/  
Temperaturbepërking/  
Temperaturbegrensninger/  
Temperatury graniczne/  
Limite de temperatura/  
Limite de temperatură/  
Температурный режим/  
Temperaturbegrensning/  
Teplotné obmedzenie  
Omejitve temperature/  
Temperaturno ograničenje/  
Sıcaklık sınırlaması/

## IVD

In Vitro Diagnostic Medical Device/  
Медицински уред за диагностика  
ин витро/Diagnostický zdravotnický  
prostředek in vitro/Medicinsk udstyr til  
in vitro-diagnostik/In-vitro-Diagnostikum/  
Ιατροτεχνολογικό προϊόν για διάγνωση  
In Vitro/Dispositivo médico para  
diagnóstico in vitro/In vitro diagnostiline  
meditsiiniseade/Dispositif médical de  
diagnostic in vitro/Diagnostički medicinski  
uređaj In Vitro/In vitro orvosdiagnostikai  
eszköz/Dispositivo medico per test  
diagnostici in vitro/In Vitro Diagnostiné  
Medicinos Priemonė/Medicīniska ierīce  
in vitro diagnostikai/In vitro-diagnostisch  
medisch instrument/In vitro diagnostisk  
medisinsk utstyr/Wyrób medyczny do  
diagnostyki in vitro/Dispositivo Médico  
de Diagnóstico In Vitro/Dispozitiv medical  
pentru diagnostic in vitro/Только для  
диагностики In Vitro/Endast för in  
vitro-diagnostik/ Zdravotnička pomôcka na  
diagnostiku in vitro/In vitro diagnostični  
pripomoček/Diagnostički medicinski  
uređaj In Vitro/<96> testleri için yeterlilik  
içerir

Contains sufficient for <96> tests/Съдържа  
достатъчно количество за тестове  
<96>/Lze použít pro <96> testů/Ineholder  
tilstrækkeligt/Inhalt ausreichend für <96>  
Prüfungen/Περιεχόμενο επαρκές για  
«96» εξετάσεις/Contenido suficiente para  
<96> ensayos/Kogusest pisab <96> testi  
lābiviimeseks/Contenu suffisant pour «96»  
tests/Sadržaj dovoljno za <96> testova/A  
doboz tartalma <96> vizsgálat elvégzéséhez  
elegendő/Contenuto sufficiente per «96»  
saggi/Turinys skirtas atliki <96> tyrimus/  
Saturis pietiekams <96> testiem/Inhoud  
voldoende voor «96» testen/til «96» test/  
Tilstrækkelig innhold for <96> prøver/  
Wystarczy na wykonanie <96> testów/  
Conteúdo suficiente para «96» ensaios/  
Conținut suficient pentru 96 de teste/  
Содержит достаточные количества для  
«96» определений/Innehåller tillräckligt  
till «96» antal tester/Obsah postačuje na  
tento počet testov: <96>/Vsebina zadostuje  
za <96> testov/Sadržina dovoljna za <96>  
testova/<96> testleri için yeterlilik içerir

## REF

Catalogue number/Каталожен номер/  
Katalogové číslo/Katalognummer/  
Bestellnummer/Αριθμός καταλόγου/  
Número de catálogo/Kataloogi number/  
Numéro de catalogue/Kataloški broj/  
Katalógusszám/Numero di catalogo/  
Katalogo numeris/Numurs katalogā/  
Catalogusnummer/Katalognummer/  
Numer katalogowy/Número do catálogo/  
Număr de catalog/Номер по каталогу/  
Produktnummer/Katalógové číslo/  
Kataloška številka/Kataloški broj/  
Katalog numarası

Consult Instructions for Use/  
 Прочетете инструкцията за  
 употреба/Konzultujte s návodem  
 k použití/Se brugsanvisning/Siehe  
 Gebrauchsanweisung/Συμβουλευτείτε  
 τις Οδηγίες σχετικά με τη χρήση/  
 Consulte las instrucciones de uso/  
 Vt kasutusjuhendit/Consulter le mode  
 d'emploi/Pročitajte upute za uporabu/  
 Olvassa el a használati utasítást/  
 Consultare le istruzioni per l'uso/Dél  
 naudojimo žiūrėkite instrukcijas/Izlasiet  
 lietošanas instrukciju/Raadpleeg de  
 instructies voor gebruik/Les instruksene  
 før bruk/Sprawdź w instrukcji użycia/  
 Consulte as Instruções de Utilização/  
 Consultați instrucțiunile de utilizare/  
 Обратитесь к инструкции по  
 применению/Se brugsanvisning/  
 Prečítajte si návod na používanie/  
 Pročitajte uputstvo za upotrebu/  
 Kullanim Talimatlarına Bakınız

**CONT**

Contents of kit/Съдържание на набора/  
 Obsah soupravy/Kittets indhold/Inhalt  
 des Kits/Περιεχόμενα του κιτ/Contenido  
 del kit/Komplekt sisaldab/Contenu du  
 kit/Sadržaj opreme/A készlet tartalma/  
 Contenido del kit/Rinkinio turinys/  
 Komplekta saturs/Inhoud van de set/  
 Settets indhold/Zawartość zestawu/  
 Conteúdo do kit/Conținutul setului/  
 Компоненты набора/Kit innehåll/  
 Obsah súpravy/Vsebina kompleta/Sadržaj  
 opreme/Kitin içindekiler

Biological risks/Биологическа  
 опасност/Biológická rizika/Biologisk  
 fare/Biologische Gefahren/Bioloγικό  
 κίνδυνος/Riesgos biológicos/  
 Bioloogilised ohud/Risques biologiques/  
 Biološkii rizici/Biológiai kockázatok/Rischi  
 biologicci/Biologinis pavojus/Biolóģiskais  
 risks/Biologische risico's/Biologiske  
 risikoer/Zagrozenie biologiczne/Riscos  
 biológicos/ Biologisk risk/Pericole  
 biologice/Биологическая опасность/  
 Biologicky rizikové/Biológické riziká/  
 Biološkii rizici/Bijolojik riskler

**Rx ONLY**

US: For use by or on the order of a  
 physician only

**ORIG HUM**

Human/C човешки произход/Lidské/  
 Human/Human/δείγματα αναφοράς/  
 Humano/Inimpăritolu/Humaine/Ljudskog  
 porjekla/Human/Origine Umana/  
 Žmogaus kilmės/Cilvēku izcelsmes/  
 Human/Menneske/Ludzka/Humano/  
 Origine umaná/Человеческого  
 происхождения/Human/L'udské/  
 Humanega izvora/Ljudskog porekla/İnsan

**ORIG MOU**

From mouse/C миши произход/Мыši/  
 Fra mus/Maus/από ποντίκι/de ratón/  
 Hiirtelt/De souris/Mišijeg porjekla/  
 Egérböi/Murino/Pelės kilmės/No peles/  
 Van muizen/Fra mus/Mysia/Do rato/De  
 la șoareci/Мышиного происхождения/  
 Frán mus/Myšie/Mišjega izvora/Mišijeg  
 porekla/Fareden

**ORIG BOV**

Bovine/C говежди произход/  
 Novézi/Bovin/Rind/από βοοειδή/  
 Bovino/Veistelt/Bovine/Rogate stoke/  
 Szarvasmarha/Bovino/Jaučio/No  
 liellopa/Bovien/Bovin/Wolowy/Bovino/  
 Origine bovină/крупного рогатого  
 скота/Frán ko/Hovádzie/Govejega  
 izvora/Rogate krupne stoke/Bovin

Reconstitute with/Пазтваряне с/  
 Rozfedte pomocí/Rekonstitueres med/  
 Rekonstituieren mit/Ανασύσταση με/  
 Reconstituir con/Lahjendamine/  
 Reconstituer avec/Rekonstituiraite s/  
 Feloldáshoz/Ricostituire con/Atkurti,  
 ištirpdant su/Atšķaidīt ar/Reconstitutie  
 met/Rekonstitueres med/Odtworzyć  
 za pomocą/Reconstituir com/A  
 se reconstitui cu/Растворить в/  
 Rekonstituera med/Rozriedte pomocou/  
 Rekonstituiraite z/s/Ponovno formiranje  
 sa/Yeniden oluşturur

Manufacturer/Производител/Výrobce/  
 Producent/Hersteller/Κατασκευαστής/  
 Fabricante/Tootja/Fabricant/Proizvođač/  
 Gyártó/Fabbricante/Gamintojas/  
 Ražotājs/Fabrikant/Produsent/  
 Producent/Fabricante/Producător/  
 Производитель/Тилverkare/ Výrobca/  
 Izdelovalec/Proizvođač/Üretici

# INSTRUCTIONS FOR USE

EN

## INSTRUCTIONS FOR USE

Please visit our website [www.fdi.com/ifu](http://www.fdi.com/ifu) to obtain the Instructions For Use (IFU) in additional languages.

To ensure that you download the correct IFU for your kit lot, please select the revision corresponding to the issue date printed on the front page of the IFU provided with this kit.

Please follow the IFU carefully. Instructions for safe handling are found in the WARNINGS AND PRECAUTIONS section. Material Safety Data Sheets (MSDS) are available on our website [www.fdi.com](http://www.fdi.com). If you do not have access to the internet, please contact your local distributor, or Fujirebio Diagnostics AB for assistance.

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# CanAg PSA EIA

Instructions for use

Enzyme immunometric assay kit  
For 96 determinations

## INTENDED USE

The CanAg PSA EIA kit is intended for the quantitative determination of Total PSA (Prostate Specific Antigen) in human serum.

## SUMMARY AND EXPLANATION OF THE ASSAY

PSA is a 32 kDa single chain glycoprotein serine protease with a chymotrypsin like specificity produced by the secretory epithelium of the prostate gland (1). PSA is normally secreted into the seminal fluid and plays a functional role in the cleavage of the seminal vesicle proteins and the liquefaction of the seminal coagulum (2). Only low levels of PSA are normally present in the blood stream, and increasing serum concentrations indicate prostatic pathology, including benign prostatic hyperplasia and cancer of the prostate. Determination of PSA is now widely used for detection and management of patients with prostatic cancer and considered as the superior serological marker for cancer of the prostate (3, 8).

PSA has been shown to form stable complexes with different antiproteases and the dominating portion of PSA in patient serum occurs in complex with  $\alpha_1$ -antichymotrypsin (PSA-ACT) (4). However there are large variations in the relation between free PSA and PSA-ACT complex between different individuals. Studies have also indicated that the proportion of free PSA is higher in benign prostatic disease as compared to prostatic cancer (5). The antibodies in the CanAg PSA EIA have been carefully evaluated and selected to give the same molar response for free PSA as for the PSA-ACT complex (7). The CanAg PSA EIA thus gives a true "total" PSA value independent of the individual variations of free and ACT-complexed PSA.

## PRINCIPLE OF THE TEST

The CanAg PSA EIA is a solid-phase, non-competitive immunoassay based upon the direct sandwich technique. Calibrators, controls and patient samples are incubated together with biotinylated anti-PSA monoclonal antibody and horseradish peroxidase (HRP) labelled anti-PSA monoclonal antibody in Streptavidin coated microstrips. After washing, buffered Substrate/Chromogen reagent (hydrogen peroxide and 3, 3', 5, 5' tetra-methylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction a blue colour will develop if antigen is present. The intensity of the colour is proportional to the amount of PSA present in the samples.

The colour intensity is determined in a microplate spectrophotometer at 620 nm (or optionally at 405 nm after addition of Stop Solution).

Calibration curves are constructed for each assay by plotting absorbance value versus the concentration for each calibrator. The PSA concentrations of patients samples are then read from the calibration curve.

## REAGENTS

- Each CanAg PSA EIA kit contains reagents for 96 tests.
- The expiry date of the kit is stated on the label on the outside of the kit box.
- Do not use the kit beyond the expiry date.
- Do not mix reagents from different kit lots.
- Store the kit at 2–8° C. Do not freeze.
- Opened reagents are stable according to the table below provided they are not contaminated, stored in resealed original containers and handled as prescribed. Return to 2–8° C immediately after use.

Component	Quantity	Storage and stability after first opening
-----------	----------	---

### MICROPLA

<b>Microplate</b>	1 Plate	2–8° C until expiry date stated on the plate
-------------------	---------	--

12 x 8 wells coated with Streptavidin. After opening, immediately return unused strips to the aluminium pouch, containing desiccant. Reseal carefully to keep dry.

<b>PSA Calibrators</b>	6 vials	2–8° C until expiry date stated on the vials
------------------------	---------	--

<b>CAL</b>	<b>PSA</b>	<b>0</b>	0 µg/L	1 x 0.75 mL
------------	------------	----------	--------	-------------

<b>CAL</b>	<b>PSA</b>	<b>1</b>	1 µg/L	1 x 0.75 mL
------------	------------	----------	--------	-------------

<b>CAL</b>	<b>PSA</b>	<b>2</b>	2 µg/L	1 x 0.75 mL
------------	------------	----------	--------	-------------

<b>CAL</b>	<b>PSA</b>	<b>10</b>	10 µg/L	1 x 0.75 mL
------------	------------	-----------	---------	-------------

<b>CAL</b>	<b>PSA</b>	<b>30</b>	30 µg/L	1 x 0.75 mL
------------	------------	-----------	---------	-------------

<b>CAL</b>	<b>PSA</b>	<b>60</b>	60 µg/L	1 x 0.75 mL
------------	------------	-----------	---------	-------------

Human PSA in a Tris-HCl buffered salt solution containing bovine serum albumin, an inert yellow dye and 0.01 % methyl-isothiazolone (MIT) as preservative. Ready for use.

Component	Quantity	Storage and stability after first opening
<b>PSA Controls</b>	2 vials	2–8° C until expiry date stated on the vials
<b>CONTROL</b>   <b>PSA</b>   <b>1</b>	1 x 0.75 mL	
<b>CONTROL</b>   <b>PSA</b>   <b>2</b>	1 x 0.75 mL	

Human PSA in a Tris-HCl buffered salt solution containing bovine serum albumin, and 0.01 % methyl-isothiazolone (MIT) as preservative. Ready for use.

<b>BIOTIN</b>   <b>Anti-PSA</b>		
<b>Biotin Anti-PSA</b>	1 x 15 mL	2–8° C until expiry date stated on the vial

Biotin Anti-PSA monoclonal antibody from mouse, approximately 1 µg/mL. Contains phosphate buffered saline (pH 7.2), bovine serum albumin, bovine immunoglobulin, blocking agents, Tween 20, an inert blue dye and 0.01 % methyl-isothiazolone (MIT) as preservative. Ready for use.

<b>CONJ</b>   <b>Anti-PSA</b>		
<b>Tracer, HRP Anti-PSA</b>	1 x 0.75 mL	2–8° C until expiry date stated on the vial

Stock solution of HRP Anti-PSA monoclonal antibody from mouse, approximately 20 µg/mL. Contains preservatives. To be mixed with Biotin Anti-PSA prior to use.

<b>SUBS</b>   <b>TMB</b>		
<b>TMB HRP-Substrate</b>	1 x 12 mL	2–8° C until expiry date stated on the vial

Contains buffered hydrogen peroxide and 3, 3', 5, 5' tetramethyl-benzidine (TMB). Ready for use.

Component	Quantity	Storage and stability after first opening
-----------	----------	---

**STOP**

<b>STOP Solution</b>	1 x 15 mL	2–8° C until expiry date stated on the vial
----------------------	-----------	---

Contains 0.12 M hydrochloric acid. Ready for use.

**WASHBUF 25X**

<b>Wash Concentrate</b>	1 x 50 mL	2–8° C until expiry date stated on the bottle
-------------------------	-----------	---

A Tris-HCl buffered salt solution with Tween 20. Contains Germall II as preservative. To be diluted with water 25 times before use.

### Indications of instability

The TMB HRP-Substrate should be colourless or slightly bluish. A blue colour indicates that the reagent has been contaminated and should be discarded.

### WARNINGS AND PRECAUTIONS

#### For in vitro diagnostic use.

- For professional use only.
- Please refer to the US Department of Health and Human Services (Bethesda, Md., US) publication No. (CDC) 88-8395 on laboratory safety or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Follow local guidelines for disposal of all waste material.

#### Caution

Material used in the preparation of human source reagent has been tested and found to be Non Reactive for HIV-1/2 Antibody, HCV Antibody and Hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they were potentially infectious.

### CLP (1272/2008) HAZARD CLASSIFICATION

Information about CLP (1272/2008) HAZARD CLASSIFICATION can be found at the end of this document.

# Protocol Sheet

## CanAg PSA EIA REF 340-10

Mix the components directly before use. Use shaking conditions according to the Instructions.

Step	Vial/Plate	Procedure																																	
1. Prepare Wash Solution	WASHBUF 25X	Dilute 50 mL of Wash Concentrate with 1200 mL of distilled water or deionized water.																																	
Prepare Antibody Solution	CONJ Anti-PSA	Mix 50 $\mu$ L of Tracer, HRP Anti-PSA, with 1 mL of Biotin Anti-PSA per strip:																																	
	BIOTIN Anti-PSA																																		
		<table border="1"> <thead> <tr> <th>No. of Strips</th> <th>Tracer, HRP Anti-PSA (<math>\mu</math>L)</th> <th>Biotin Anti-PSA (mL)</th> </tr> </thead> <tbody> <tr><td>1</td><td>50</td><td>1</td></tr> <tr><td>2</td><td>100</td><td>2</td></tr> <tr><td>3</td><td>150</td><td>3</td></tr> <tr><td>4</td><td>200</td><td>4</td></tr> <tr><td>5</td><td>250</td><td>5</td></tr> <tr><td>6</td><td>300</td><td>6</td></tr> <tr><td>7</td><td>350</td><td>7</td></tr> <tr><td>8</td><td>400</td><td>8</td></tr> <tr><td>9</td><td>450</td><td>9</td></tr> <tr><td>10</td><td>500</td><td>10</td></tr> </tbody> </table>	No. of Strips	Tracer, HRP Anti-PSA ( $\mu$ L)	Biotin Anti-PSA (mL)	1	50	1	2	100	2	3	150	3	4	200	4	5	250	5	6	300	6	7	350	7	8	400	8	9	450	9	10	500	10
No. of Strips	Tracer, HRP Anti-PSA ( $\mu$ L)	Biotin Anti-PSA (mL)																																	
1	50	1																																	
2	100	2																																	
3	150	3																																	
4	200	4																																	
5	250	5																																	
6	300	6																																	
7	350	7																																	
8	400	8																																	
9	450	9																																	
10	500	10																																	

## SPECIMEN COLLECTION AND HANDLING

The CanAg PSA EIA is intended for use with serum. Collect blood by venipuncture and separate the serum according to common procedures. Samples can be stored at 2–8°C for 48 hours. For longer periods store samples at -20°C or below. Samples should not be stored in a self-defrosting freezer. Allow frozen samples to thaw slowly, preferably at 2–8°C over night and then bring the samples to room temperature before analysis.

Elevated levels of PSA may be expected following manipulation of the prostate. It is therefore recommended that blood be drawn before digital rectal examination. Following surgical manipulation of the prostate, such as needle biopsy or transurethral resection it is recommended to wait  $\geq$  than 6 weeks before drawing blood for PSA testing (8). It should be taken into account that Finasteride treatment of BPH have been shown to decrease PSA levels (8).

## PROCEDURE

### Materials required but not supplied with the kit

**1. Microplate shaker**

Shaking should be medium to vigorous, approximately 700-1100 oscillations/min .

**2. Microplate wash device**

Automatic plate washer capable of performing 1 and 6 washing cycles with a minimal fill volume of 350  $\mu$ L/well/washcycle.

An 8-channel pipette with disposable plastic tips for delivery of 350  $\mu$ L is recommended if an automatic microplate washer is not used.

**3. Microplate spectrophotometer**

With a wavelength of 620 nm and/or 405 nm and an absorbance range of 0 to 3.0.

**4. Precision pipettes**

With disposable plastic tips to deliver microlitre and millilitre volumes. An 8-channel pipette or respenser pipette with disposable plastic tips for delivery of 100  $\mu$ L is useful but not essential.

**5. Distilled or deionized water**

For preparation of Wash Solution.

## Procedural notes

1. A thorough understanding of this package insert is necessary to ensure proper use of the CanAg PSA EIA kit. The reagents supplied with the kit are intended for use as an integral unit. Do not mix identical reagents from kits having different lot numbers. Do not use the kit reagents after the expiry date printed on the outside of the kit box.
2. Reagents should be allowed to reach room temperature (20–25° C) prior to use. The assay should only be performed at temperatures between 20–25° C to obtain accurate results. Frozen specimens should be brought to room temperature slowly and must be gently but thoroughly mixed after thawing.
3. Before starting to pipette calibrators, controls and patient specimens it is advisable to mark the strips to be able to clearly identify the samples during and after the assay.
4. The requirement for efficient and thorough washing for separation of bound and unbound antigen and reagents from the solid-phase bound antibody-antigen complexes is one of the most important steps in an EIA. In order to ensure efficient washing make sure that all wells are completely filled to the top edge with wash solution during each wash cycle, that wash solution is dispensed at a good flow rate, that the aspiration of the wells between and after the wash cycles is complete and that the wells are empty. If there is liquid left, invert the plate and tap it carefully against absorbent paper.
  - Automatic strip washer: Follow the manufacturer's instructions for cleaning and maintenance diligently and wash the required number of wash cycles prior to and after each incubation step. It's highly recommended to use *strip* process mode and *overflow* wash mode with a dispensing volume of 800 µL. The aspiration/wash device should not be left standing with the Wash Solution for long periods, as the needles may get clogged resulting in poor liquid delivery and aspiration.
5. 

CAL	PSA	1
-----	-----	---

 shall be used when ultrasensitive measurement is desired. Please refer to Option 2 at page 22 for full instructions.
6. The TMB HRP-Substrate is very sensitive for contamination. For optimal stability of the TMB HRP-Substrate, pour the required amount from the vial to a carefully cleaned reservoir or preferably a disposable plastic tray to avoid contamination of the reagent. Be sure to use clean disposable plastic pipette tips (or respenser pipette tip).
7. Be sure to use clean disposable plastic pipette tips and a proper pipetting technique when handling samples and reagents. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid. A proper pipetting technique is of particular importance when handling the TMB HRP-Substrate solution.

Preparation of reagents	Stability of prepared reagent
<b>Wash Solution</b>	2 weeks at 2–25° C in a sealed container

Pour the 50 mL Wash Concentrate into a clean container and dilute 25- fold by adding 1200 mL of distilled or deionized water to give a buffered Wash Solution.

<b>Antibody Solution</b>	3 weeks at 2–8° C
Prepare the required quantity of Antibody Solution by mixing 50 µL of Tracer, HRP-Anti PSA with 1 mL of Biotin Anti-PSA per strip (see table below and the Protocol Sheet).	

No. of Strips	Tracer, HRP Anti-PSA (µL)	Biotin Anti-PSA (mL)
1	50	1
2	100	2
3	150	3
4	200	4
5	250	5
6	300	6
7	350	7
8	400	8
9	450	9
10	500	10
11	550	11
12	600	12

Be sure to use a clean plastic or glass bottle for preparation of the Antibody Solution.

**Alternative:** Pour the content of the Tracer, HRP Anti-PSA into the vial of Biotin Anti-PSA and mix gently. Make sure that all of the Tracer is transferred to the vial of Biotin Anti-PSA.

**Note:** The Antibody Solution is stable for 3 weeks at 2–8° C. Do not prepare more Antibody Solution than will be used within this period and make sure that it is stored properly.

## ASSAY PROCEDURE

Perform each determination in duplicate for calibrators, controls and patient samples. A calibration curve should be run with each assay. All reagents and samples must be brought to room temperature (20–25°C) before use.

1. Start to prepare Wash Solution and Antibody Solution. It is important to use clean containers. Follow the instructions carefully.
2. Transfer the required number of microplate strips to a strip frame. (Immediately return the remaining strips to the aluminium pouch containing a desiccant and reseal carefully). Wash each strip once with the Wash Solution. Do not wash more strips than can be handled within 30 min.
3. Pipette 25 µL of the PSA Calibrators (CAL 0, 2, 10, 30, 60), controls (c) and patient samples (unknowns-Unk) into the strip wells according to the following scheme:

	1	2	3	4	5	6	7 etc
A	CAL 0	CAL 60	Unk 2				
B	CAL 0	CAL 60	etc.				
C	CAL 2	C 1					
D	CAL 2	C 1					
E	CAL 10	C 2					
F	CAL 10	C 2					
G	CAL 30	Unk 1					
H	CAL 30	Unk 1					

4. Add 100 µL of Antibody Solution to each well using a 100 µL precision pipette (or an 8-channel 100 µL precision pipette). Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or the surface of the liquid.
5. Incubate the frame containing the strips for 1 hour ( $\pm$  10 min) at room temperature (20–25°C) with constant shaking of the plate using a microplate shaker.
6. Wash each strip 6 times, using the wash procedure described in Procedural notes item 4.

7. Add 100  $\mu\text{L}$  of TMB HRP-Substrate to each well using the same pipetting procedure as in item 4. The TMB HRP-Substrate should be added to the wells as quickly as possible and the time between the addition to the first and last well should not exceed 5 min.
8. Incubate for 30 min ( $\pm$  5 min) at room temperature with constant shaking. Avoid direct sunlight.
9. Immediately read the absorbance at 620 nm in a microplate spectrophotometer.

### Option 1

If the laboratory does not have access to a microplate spectrophotometer capable of reading at 620 nm, the absorbance can be determined as follows:

Add 100  $\mu\text{L}$  of Stop Solution. Mix and read the absorbance at 405 nm in a micro-plate spectrophotometer within 5 minutes after addition of Stop Solution.

### Option 2

For ultrasensitive determination of PSA in the low range (0-10  $\mu\text{g/L}$ ),

CAL	PSA	1
-----	-----	---

 (1  $\mu\text{g/L}$ ) should be included in the calibration curve. It is then necessary to exclude the two highest calibrators (30 and 60  $\mu\text{g/L}$ ). Measurements of absorbance should be done according to Option 1 but reading at 450 nm.

### Measurement range

The CanAg PSA EIA measures concentrations between 0.1 and 60  $\mu\text{g/L}$ . If PSA concentrations above the measuring range are to be expected, it is recommended to dilute samples with normal human male serum prior to analysis. **NOTE:** The serum used for dilution should also be measured in order to determine the endogenous PSA concentration (see "Calculation of results").

### Quality control

PSA Control 1 and 2 may be used for validation of the assay series. Ranges of expected results are indicated on the vial labels. If values outside of the specified range are obtained, a complete check of reagents and reader performance should be made and the analysis repeated. Each laboratory may in addition prepare its own serum pools at different levels, which can be used as internal controls in order to assure the precision of the assay.

### Reference material

The 1<sup>st</sup> International Standard 96/670 may be used as a reference standard. Values for PSA Calibrators and Controls were assigned against a set of in-house reference standards whose values are traceable to the 1<sup>st</sup> International Standard.

## CALCULATION OF RESULTS

If a microplate spectrophotometer reader with built-in data calculation program is used refer to the manual for the plate reader and create a program using the concentration stated on the labels of each of the PSA Calibrators.

For automatic calculation of PSA results it is recommended to use either of the following methods:


- Cubic spline curve fit method. Calibrator 0 should be included in the curve with the value 0  $\mu\text{g/L}$ .
- Spline smoothed curve fit method. Calibrator 0 should be used as plate blank.
- Interpolation with point-to-point evaluation. Calibrator 0 should be included in the curve with the value 0  $\mu\text{g/L}$ .
- Quadratic curve fit method. Calibrator 0 should be included in the curve with the value 0  $\mu\text{g/L}$ .

**Note:** 4-parametric or linear regression should not be used.

For manual evaluation, a calibration curve is constructed by plotting the absorbance (A) values obtained for each PSA calibrator against the corresponding PSA concentration (in  $\mu\text{g/L}$ ), see figure below. The unknown PSA concentrations can then be read from the calibration curve using the mean absorbance value of each patient specimen.

### Example of results

Specimen	Calibrator values	Mean abs value (A)	PSA ( $\mu\text{g/L}$ )
CAL PSA 0	0 $\mu\text{g/L}$	0.036	
CAL PSA 2	2 $\mu\text{g/L}$	0.174	
CAL PSA 10	10 $\mu\text{g/L}$	0.705	
CAL PSA 30	30 $\mu\text{g/L}$	1.807	
CAL PSA 60	60 $\mu\text{g/L}$	2.864	
Specimen A		0.453	6.1
Specimen B		1.739	28.6

**Example (do not use this curve or table above to determine actual assay results).**

### Calculation of results with diluted samples

If samples in an initial analysis give PSA levels higher than 60 µg/L the samples should be diluted 1/10 with normal male human serum and reanalysed to obtain the accurate PSA concentration. **NOTE:** The sample used for dilution should also be measured in order to determine the endogenous PSA concentration.

The PSA concentration of the undiluted sample is calculated as:

$$\text{Dilution 1/10: } 10 \times ([\text{PSA}]_{\text{Diluted sample}} - (0.9 \times [\text{PSA}]_{\text{Normal male serum}}))$$

### LIMITATIONS OF THE PROCEDURE

The level of PSA cannot be used as absolute evidence for the presence or absence of malignant disease. The results of the test should be interpreted only in conjunction with other investigations and procedures in the diagnosis of disease and the management of patients. The PSA test should not replace any established clinical examination.

The calibrators of the CanAg PSA EIA kit should not be used for recovery studies of PSA. For recovery studies it is recommended to use a highly elevated patient sample.

Anti-reagent antibodies (human anti-mouse antibody (HAMA) or heterophilic antibodies) in the patient sample may occasionally interfere with the assay, even though specific blocking agents are included in the buffer.

Biotin may interfere with the assay giving false low results. This should be taken into consideration for patients taking dietary supplements or receiving therapy containing high (> 5mg/day) or extremely high (300 mg/day) biotin doses. Peak

serum levels have been reported to occur 1-3h post ingestion and a physiological half-life of 8-16h depending on renal function (11,12). In a study by Grimsey et al. (11) a specimen concentration of 30 ng/mL was reached 8h following intake of 10 mg of biotin. For patients on very high doses of biotin it is recommended to stop taking biotin for at least 2 days before blood draw (13).

## EXPECTED VALUES

Healthy males are expected to have PSA values below 4 µg/L. However, as PSA levels increase with age, the use of age-specific reference ranges have been suggested in order to increase the sensitivity in younger men and increase the specificity in older men (6, 8). It is recommended that each laboratory establish its own normal range to account for such local environmental factors as diet, climate, living conditions, patient selection, etc.

Results from the CanAg PSA EIA assay should not be interpreted as being definitive for the presence or absence of prostate cancer. Individuals with PSA levels within the reference interval found in apparently healthy subjects may have prostate cancer and individuals with levels above the reference interval found in apparently healthy subjects, may be free of prostate cancer.

### Reference range

The distribution of PSA values was determined in a study of 43 apparently healthy male subjects in the age interval 62- 65 years. Result is shown in the table below:

Age	No of subjects	Mean (µg/L)	Range	PSA (µg/L)	
				0-4	> 4
62 – 65	43	1,61	0,22-8,16	93%	7%

### Diagnostic sensitivity and specificity

The diagnostic accuracy of CanAg PSA EIA was evaluated in a study using samples from seventy-seven (77) men with confirmed prostate cancer and fifty-two (52) men with benign prostate hypertrophy (BPH). The sensitivity and specificity for discrimination between prostate cancer and BPH are shown in the table below:

	CanAg PSA EIA < 4,0 µg/L	CanAg PSA EIA > 4,0 µg/L
BPH n=52	28	24
Prostate cancer, n=77	11	66

The use of a cut-point of 4 µg/L resulted in 85,7% sensitivity and 53,8% specificity for detection of prostate cancer.

## PERFORMANCE CHARACTERISTICS

### Precision

Total precision was determined according to NCCLS guideline EP5-A (9) using four levels of frozen pooled human serum containing added PSA and six different CanAg PSA EIA reagent combinations. Each sample was randomly pipetted (n=2/analysis) and analysed twice each day over 20 days.

Sample	Replicates	Mean (µg/L)	Within-run SD (µg/L)	Within-run CV %	Between-day SD (µg/L)	Between-day CV %
PSA 1	80	1.42	0.04	2.7	0.03	2.2
PSA 2	80	5.92	0.13	2.2	0.06	1.0
PSA 3	80	14.2	0.35	2.5	0.12	0.8
PSA 4	80	39.2	0.60	1.5	0.60	1.5

### Detection limit

The detection limit of the CanAg PSA EIA assay is < 0.1 µg/L defined as the concentration corresponding to the mean of the absorbance values for the PSA Calibrator 0 plus 2 standard deviations according to the formula:

$$\frac{2 \times \text{SD CAL } 0}{\text{OD CAL } 2 - \text{OD CAL } 0} \times 2 \mu\text{g/L}$$

### Recovery

Spiked serum samples were prepared by adding aliquots of samples with highly elevated PSA to normal male serum samples. The recovery of the antigen was within  $\pm 10\%$  of the expected values. **Note:** Recovery studies should **not** be performed using the kit calibrators.

### Hook effect

No hook effect has been noticed with samples up to 23000 µg/L. **Note:** In very high samples the colour of the substrate will change from blue to greenish (and eventually yellow in extremely high samples). This will lead to a falsely low absorbance at 620 nm, and in extreme cases the absorbance may fall within the calibration curve range and noticed as a hook.

### Linearity

Patient samples were diluted with normal male human serum and analysed. The obtained values were within  $\pm 10\%$  of the expected values.

### Specificity

The CanAg PSA EIA is based on two mouse monoclonal antibodies, PSA10 and PSA66, directed against two distinct epitopes exposed both in the PSA-ACT complex and in free PSA. This antibody combination provides an equimolar response for free PSA and the PSA-ACT complex (7). The NCCLS guideline EP7-P (10) was followed to determine possible sources of interference. The following substances and concentrations were tested and found not to interfere with the test.

	Concentration with no significant ( $\pm 10\%$ ) interference
Lipemia (Intralipid®)	10 mg/mL
Bilirubin, unconjugated	0.6 mg/mL
Hemoglobin	5 mg/mL

### Biotin interference

A study was conducted to evaluate biotin interference. Low and high serum control samples were spiked to final biotin concentrations of 15, 30, 60, and 200 ng/mL. The mean PSA concentration was determined for each sample and the percent recovery for each biotin concentration was calculated using the formula: Recovery (%) = 100x (Mean PSA concentration w. biotin added/Mean PSA concentration w. diluent only added).

PSA analyte level	Biotin test conc. (ng/mL)	Expected PSA conc. (µg/L)	Observed PSA conc. (µg/L)	Recovery (%)
Low	15	2,98	2,99	100
High	15	27,9	27,1	97
Low	30	2,98	2,83	95
High	30	27,1	25,8	95
Low	60	3,00	2,68	89
High	60	27,4	23,9	87
Low	200	3,01	0,75	25
High	200	28,0	7,32	26

Based on linear regression analysis, the lowest concentration of biotin found to influence test results ( $\geq 10\%$ ) was 44 ng/mL.

#### Method comparison

The CanAg PSA EIA (Prod. No. 340-10) was compared to the two-step CanAg PSA EIA (300-10). Three hundred eight male human serum samples ranging in values from 0-53 µg/L were measured and linear regression analyses of the results yielded:

$$[\text{PSA}]_{\text{Prod. No. 340-10}} = 0.92 \times [\text{PSA}]_{\text{Prod. No. 300-10}} - 0.036 \quad r=1.00$$

## CLP (1272/2008) HAZARD CLASSIFICATION

The following warnings and precautions apply to

SUBS	TMB
------	-----

### Hazard pictograms:

<b>Signal word:</b>	Danger
<b>Hazard Statement:</b>	Repr. 1B: H360D May damage the unborn child.
<b>Prevention statement:</b>	P202 Do not handle until all safety precautions have been read and understood.
<b>Prevention:</b>	P280 Wear protective gloves / protective clothing / eye protection / face protection.
<b>Precautionary statement response:</b>	P308+P313 IF exposed or concerned get medical advice/attention.
<b>Precautionary statement disposal:</b>	P501 Dispose of contents / container to an approved hazardous / special waste disposal facility in accordance with local and national regulations.

### Restricted to professional users.

**Hazardous substances:** 2- Pyrrolidone

### Other hazards

None of the mixtures in the kit contains any substances considered to meet the criteria classifying them as PBT and/or vPvB.

## WARRANTY

The performance data presented here were obtained using the assay procedure indicated. Any change or modification of the procedure not recommended by Fujirebio Diagnostics may affect the results, in which Fujirebio Diagnostics disclaims all warranties expressed, implied or statutory including the implied warranty of merchantability and fitness for use.

## REFERENCES

1. Wang MC, Valenzuela LA, Murphy GP, Chu TM (1979). Purification of a human prostate specific antigen. *Invest Urol* 17: 159–163.
2. Lilja H. (1985). A kallikrein-like serine protease in prostatic fluid cleaves the predominant seminal vesicle protein. *J Clin Invest* 76: 1899–1903.
3. Oesterling JE (1991). Prostate specific antigen: A critical assessment of the most useful tumor marker for adenocarcinoma of the prostate. *J Urology* 145: 907–923.
4. Lilja H., Christensson A., Dahlén U., Matikainen M-T, Nilsson O., Pettersson K., Lövgren T. (1991). Prostate-specific antigen in serum occurs predominantly in complex with  $\alpha_1$ -antichymotrypsin. *Clin Chem* 37: 1618–1625.
5. Christensson A., Björk T., Nilsson O., Dahlén U., Matikainen M-T., Cockett ATK, Abrahamsson PA, Lilja H. (1993). Serum prostate specific antigen complexed to  $\alpha_1$ -antichymotrypsin as an indicator of prostate cancer. *J Urology* 150: 100–105.
6. Oesterling JE., Cooner WH., Jacobsen SJ., Guess HA., Lieber MM. (1993). The influence of patient age on the serum prostate specific antigen concentration: An important clinical observation. *Urol Clin North Am* 20: 671–80, 1993a.
7. Nilsson O., Peter A. Andersson I., Nilsson K., Grundström B., and Karlsson B. (1997) Antigenic determinants of prostatespecific antigen (PSA) and development of assays specific for different forms of PSA. *Br J Cancer* 75(6): 789–797.
8. P Price C., Allard J., Davies G., Dawney A., J Duffy M., France M., Mandarino G., Milford Ward A., Patel B., Sibley P. and Sturgeon C. (2001) Pre-and post-analytical factors that may influence use of serum prostate specific antigen and its isoforms in a screening programme for prostate cancer. *Ann Clin Biochem* 38: 188–216.
9. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
10. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).
11. Grimsey P et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. *J Pharmacokinet.* 2017; 2(4); 247-256.
12. Jenkins Colon P, Greene D.N. Biotin Interference in Clinical Immunoassays. *JALM* 2018; 2(6); 941-951.
13. Chun KY. Biotin Interference in Diagnostic Tests. *Clinical Chemistry* 2017; 63; 619-620.

---

**CanAg<sup>®</sup>** is a registered trademark of Fujirebio Diagnostics AB

**Fujirebio Diagnostics AB**

**Elof Lindälv's gata 13**

**SE-414 58 Göteborg**

**Sweden**

**Phone + 46 31 85 70 30**

**Fax + 46 31 85 70 40**

**info@fdab.com**

**www.fdab.com**