CanAg CA125 EIA

REF

400-10

IVD

Instructions for use. 2022-06

Read highlighted changes

EN	EXPLANATION OF SYMBOLS
BG	ОБЯСНЕНИЕ НА СИМВОЛИТЕ
CS	VÝZNAM SYMBOLŮ
DA	SYMBOLFORKLARING
DE	ERKLÄRUNG DER SYMBOLE
EL	ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
ES	SIGNIFICADO DE LOS SÍMBOLOS
ET	SÜMBOLITE SELGITUS
FR	EXPLICATION DES SYMBOLES
HR	OBJAŠNJENJE SIMBOLA
HU	JELMAGYARÁZAT
IT	SPIEGAZIONE DEI SIMBOLI
LT	SIMBOLIŲ PAAIŠKINIMAI
LV	SIMBOLU SKAIDROJUMS
NL	VERKLARING DER SYMBOLEN
NO	SYMBOLFORKLARING
PL	OBJAŚNIENIE SYMBOLI
PT	EXPLICAÇÃO DOS SÍMBOLOS
RO	SEMNIȚAȚIA SIMBOLURILOR
RU	ОБОЗНАЧЕНИЯ
SV	SYMBOLFÖRKLARING
SK	VÝZNAM SYMBOLOV
SL	RAZLAGA SIMBOLOV
SR	OBJAŠNJENJE SIMBOLA
TR	SEMBOLLERİN AÇIKLAMALARI

Use By/Годно до/Použitelné do/
Holdbar til/Verwendbar bis/
Ημερομηνία λήξης/Fecha
de caducidad/Kölblik kuni/
Utiliser jusque/Rok valjanosti/
Felhasználható/Utilizzare entro/
Sunautoti iki/Izlietot līdz/Houdbaar
tot/Brukes innen/Użyć przed/
Prazo de validade/Expirã ła/
Использовать до/Använd före/
Použite né do/ Uporabno do/
Upotrebljivo do/Son Kullanna Tarihi

LOT

Batch code/Номер на партида/
Číslo šarže/Lotnummer/
Chargenbezeichnung/Αριθμός
Παρτίδας/Código de lote/Partii
kood/Code du lot/Kod serije/
Sarzsám/Codice del lotto/
Partijas kods/Partijas kods/Lot
nummer/Partikode/Kod partii/
Código do lote/Număr de lot/
Номер лота/Lotnummer/Číslo
šarže/Številka serije/Kod partije/
Parti Kodu

Date of manufacture/Дата на производство/Datum výroby/
Produktionsdato/Herstellungsdatum/
Ημερομηνία παραγωγής/Fecha de fabricación/Valmistamise kuupäev/
Date de fabrication/Datum proizvodnje/
Gyártási idő/Data di produzione/
Pagaminimo data/Ražošanas datums/
Productiedatum/Fremstillingsdato/
Data produkcji/Data de fabrico/Data fabricației/Дата производства/
Tillverkningsdatum/Dátum výroby/Datum izdelave/Datum proizvodnje/Úretim tarihi

Temperature limitation/
Температурни граници/
Теплотни омеzeи/
Temperaturbegrænsning/
Temperaturbegrenzung/
Περιορισμοί θερμοκρασίας/
Limites de temperatura/
Temperatuuri piirang/
Limite de température/
Temperaturno ograničenje/
Hőmérsékletre vonatkozó korlátozás/
Limiti di temperatura/

IVD

In Vitro Diagnostic Medical Device/
Медицински уред за диагностика
ин vitro/Diagnostický zdravotnícký
prostředek in vitro/Medicinsk udstyr til
in vitro-diagnostik/In-vitro-Diagnostikum/
Ιατροτεχνολογικό προϊόν για διάγνωση
In Vitro/Dispositivo médico para
diagnóstico in vitro/In vitro diagnostiline
meditsiiniseade/Dispositif médical de
diagnostic in vitro/Diagnostički medicinski
uređaj In Vitro/In vitro orvosdiagnostikai
eszköz/Dispositivo medico per test
diagnostici in vitro/In Vitro Diagnostinė
Medicinos Priemonė/Medicínska ierice
in vitro diagnostikai/In vitro-diagnostisch
medisch instrument/In vitro diagnostisk
medisinsk utstyr/Wyrób medyczny do
diagnostyki in vitro/Dispositivo Médico
de Diagnóstico In Vitro/Dispozitiv medical
pentru diagnostic in vitro/Только для
диагностики In Vitro/Endast för in
vitro-diagnostik/ Zdravotnícka pomôcka na
diagnostiku in vitro/In vitro diagnostični
pripomoček/Diagnostički medicinski
uređaj In Vitro/<96> testleri için yeterlilik
içerir

Temperatúriniai apribojimai/
Temperatūras ierobežojums/
Temperatuurbepierking/
Temperaturbegrensninger/
Temperatury graniczne/
Limite de temperatura/
Limite de temperatură/
Температурный режим/
Temperaturbegrænsning/
Teplotné obmedzenie
Omejitve temperature/
Temperaturno ograničenje/
Sıcaklık sınırlaması/

Contains sufficient for <96> tests/Съдържа
достатъчно количество за тестове
<96>/Lze použít pro <96> testů/Ineholder
tilstrækkeligt/Inhalt ausreichend für <96>
Prüfungen/Περεχόμενο επαρκές για
«96» εξετάσεις/Contenido suficiente para
<96> ensayos/Kogusest piisab <96> testi
läbiviimiseks/Contenu suffisant pour «96»
tests/Sadržaj dovoljno za <96> testova/A
doboz tartalma <96> vizsgálat elvégzéséhez
elegendő/Contenuto sufficiente per «96»
saggi/Turiny's skirtas atlikti <96> tyrimus/
Saturis pietiekams <96> testiem/Inhoud
voldoende voor «96» testen/til «96» test/
Tilstrækkelig innhold for <96> prøver/
Wystarczy na wykonanie <96> testów/
Conteúdo suficiente para «96» ensaios/
Conținut suficient pentru 96 de teste/
Содержит достаточные количества для
«96» определений/Innehåller tillräckligt
till «96» antal tester/Obsah postačuje na
tento počet testov: <96>/Vsebinsa zadostuje
za <96> testov/Sadržina dovoljna za <96>
testova/<96> testleri için yeterlilik içerir

REF

Catalogue number/Каталожен номер/
Katalogové číslo/Katalognummer/
Bestellnummer/Αριθμός καταλόγου/
Número de catálogo/Katalogoi number/
Numéro de catalogue/Kataloški broj/
Katalógusszám/Numero di catalogo/
Katalogo numeris/Numurs katalogā/
Catalogusnummer/Katalognummer/
Numer katalogowy/Número do catálogo/
Număr de catalog/Номер по каталогу/
Produktnummer/Katalógové číslo/
Kataloška številka/Kataloški broj/
Katalog numarası

Consult Instructions for Use/
Прочетете инструкцията за
употреба/Konzultujte s návodem
k použití/Se brugsanvisning/Siehe
Gebrauchsweisung/Συμβουλευτείτε
της Οδηγίας σχετικά με τη χρήση/
Consulte las instrucciones de uso/
Vt kasutusjuhendit/Consulter le mode
d'emploi/Pročítajte upute za uporabu/
Olvassa el a használati utasítást/
Consultare le istruzioni per l'uso/Dél
naudojimo žiūrėkite instrukcijas/Izlasiet
lietošanas instrukciju/Raadpleeg de
instructies voor gebruik/Les instruksene
for bruk/Sprawdzić w instrukcji użycia/
Consulte as Instruções de Utilização/
Consultați instrucțiunile de utilizare/
Обратитесь к инструкции по
применению/Se bruksanvisning/
Prečítajte si návod na používanie/
Pročítajte uputstvo za upotrebu/
Kullanım Talimatlarını Bakınız

CONT

Contents of kit/Съдържание на набора/
Obsah soupravy/Kittets indhold/Inhalt
des Kits/Περιεχόμενα του κιτ/Contenido
del kit/Komplekt sisaldab/Contenu du
kit/Sadržaj opreme/A készlet tartalma/
Contenuto del kit/Rinkinio turinys/
Komplekta saturs/Inhoud van de set/
Settets innhold/Zawartość zestawu/
Conteúdo do kit/Conținutul setului/
Компоненты набора/Kit innehåll/
Obsah súpravy/Vsebina kompleta/Sadržaj
opreme/Kitin içindekiler

Biological risks/Биологическа
опасност/Biológická rizika/Biologisk
fare/Biologische Gefahren/Βιολογικοί
κίνδυνοι/Riesgos biológicos/
Biolooigilised ohud/Risques biologiques/
Biolóskli rizici/Biológiai kockázatok/Rischi
biologici/Biologinis pavojus/Biológiskais
risks/Biologische risico's/Biologiske
risikoer/Zagroženie biologiczne/Riscos
biológicos/ Biologisk risk/Pericole
biologice/Биологическая опасность/
Biologický rizikové/Biológické riziká/
Biolóskli rizici/Biyolojik riskler

ORIG HUM

Human/C човешки производ/Lidské/
Human/Human/δείγματα αναφοράς/
Humano/Inimpãritolu/Humaine/Ljudskog
porjekla/Humán/Origine Umana/
Žmogaus kilmės/Cilvēku izcelsmes/
Human/Menneske/Ludzka/Humano/
Origine umaná/Человеческого
происхождения/Human/Ludské/
Humanega izvora/Ljudskog porekla/İnsan

ORIG MOU

From mouse/C миши производ/Мыši/
Fra mus/Maus/από πογτίκι/de ratón/
Hiirtelt/De souris/Mišijeg porjekla/
Egérbdli/Murino/Pelés kilmés/No peles/
Van muizen/Fra mus/Mysia/Do rato/De
la șoareci/Мышиного происхождения/
Från mus/Мыšie/Mišjega izvora/Mišijeg
porekla/Fareden

ORIG BOV

Bovine/C говежди производ/Hovézi/
Bovin/Rind/από βοοειδή/Bovino/
Veistelt/Bovine/Rogate stoke/Szarvas-
marha/Bovino/Jaučio/No liellopa/
Bovien/Bovin/Wolowy/Bovino/Origine
bovină/крупного рогатого скота/Från
ko/Hovädzie/Govejega izvora/Rogate
krupne stoke/Bovin

Reconstitute with/Разтваряне с/
Rozfeđe pomoci/Rekonstitueres med/
Rekonstituieren mit/Ανασύσταση με/
Reconstituir con/Lahjendamine/
Reconstituer avec/Rekonstituiraite s/
Feloldáshoz/Ricostituire con/Atkurti,
ištirpdant su/Atšķaidīt ar/Reconstituite
met/Rekonstitueres med/Odtworzyć
za pomocą/Reconstituir com/A
se reconstitui cu/Разтворит в/
Rekonstituera med/Roziedte pomocou/
Rekonstituiraite z/s/Ponovno formiranje
sa/Yeniden oluşturulur

Manufacturer/Производител/Výrobce/
Producent/Hersteller/Κατασκευαστής/
Fabricante/Tootja/Fabricant/Proizvođač/
Gyártó/Fabbricante/Gamintojas/
Ražotājs/Fabrikant/Producent/
Producent/Fabricante/Producător/
Производитель/Tilverkare/ Výrobca/
Izdelovalec/Proizvođač/Üretici

INSTRUCTIONS FOR USE

EN

INSTRUCTIONS FOR USE

Please visit our website www.fdi.com/ifu to obtain the Instructions For Use (IFU) in additional languages.

To ensure that you download the correct IFU for your kit lot, please select the revision corresponding to the issue date printed on the front page of the IFU provided with this kit.

Please follow the IFU carefully. Instructions for safe handling are found in the WARNINGS AND PRECAUTIONS section. Material Safety Data Sheets (MSDS) are available on our website www.fdi.com. If you do not have access to the internet, please contact your local distributor, or Fujirebio Diagnostics AB for assistance.

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CanAg CA125 EIA

Instructions for use

Enzyme immunometric assay kit
For 96 determinations

INTENDED USE

The CanAg CA125 EIA kit is intended for the quantitative determination of the cancer associated antigen CA125 in serum.

SUMMARY AND EXPLANATION OF THE ASSAY

CA125 is a high molecular weight mucin type glycoprotein, originally defined by the Oc125 monoclonal antibody (MAb) established by Bast et al. (1). Different epitopes, co-expressed with the Oc125 epitope on the CA125 antigen, have been used for the development of heterologous assays for determination of the CA125 antigen (2). The CanAg CA125 EIA is based on two mouse monoclonal antibodies, Ov197 and Ov185, directed against two independent epitopes of the protein core of the CA125 antigen (3, 4).

Assays for CA125 are frequently used to monitor patients with gynecological malignancies such as epithelial ovarian cancer (5).

PRINCIPLE OF THE TEST

The CanAg CA125 EIA is a solid-phase, non-competitive immunoassay based upon the direct sandwich technique. Calibrators, controls and patient samples are incubated together with a biotinylated Anti-CA125 monoclonal antibody (MAb) Ov197 (derived from mice) in streptavidin coated microstrips. CA125 present in calibrators or samples is adsorbed to the streptavidin coated microstrips by the biotinylated Anti-CA125 MAb during the incubation. The strips are then washed and incubated with HRP labeled Anti-CA125 MAb Ov185. After washing, buffered Substrate/Chromogen reagent (hydrogen peroxide and 3, 3', 5, 5' tetra-methyl-benzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction a blue color will develop if antigen is present. The intensity of the color is proportional to the amount of CA125 present in the samples.

The color intensity is determined in a microplate spectrophotometer at 620 nm (or optionally at 405 nm after addition of Stop Solution).

Calibration curves are constructed for each assay by plotting absorbance value versus the concentration for each calibrator. The CA125 concentrations of patient samples are then read from the calibration curve.

REAGENTS

- Each CanAg CA125 EIA kit contains reagents for 96 tests.
- The expiry date of the kit is stated on the label on the outside of the kit box.
- Do not use the kit beyond the expiry date.
- Do not mix reagents from different kit lots.
- Store the kit at 2–8 °C . Do not freeze.
- Opened reagents are stable according to the table below provided they are not contaminated, stored in resealed original containers and handled as prescribed. Return to 2–8 °C immediately after use.

Component	Quantity	Storage and stability after first opening
-----------	----------	---

MICROPLA

Microplate	1 Plate	2–8 °C until expiry date stated on the plate
-------------------	---------	--

12 x 8 breakable wells coated with streptavidin. After opening, immediately return unused strips to the aluminium pouch, containing desiccant. Reseal carefully to keep dry.

CA125 Calibrators	5 vials	2–8 °C until expiry date stated on the vials
--------------------------	---------	--

CAL	CA125	0	0 U/mL	1 x 8 mL
CAL	CA125	10	10 U/mL	1 x 0.75 mL
CAL	CA125	40	40 U/mL	1 x 0.75 mL
CAL	CA125	200	200 U/mL	1 x 0.75 mL
CAL	CA125	500	500 U/mL	1 x 0.75 mL

CA125 antigen in a Tris-HCl buffered salt solution containing bovine serum albumin, detergent, an inert yellow dye, and 0.05 % sodium azide as preservative. Ready for use.

CAL	CA125	0
-----	-------	---

 should also be used for dilution of samples.

Component	Quantity	Storage and stability after first opening
CA125 Controls	2 vials	2–8 °C until expiry date stated on the vials
CONTROL CA125 1	1 x 0.75 mL	
CONTROL CA125 2	1 x 0.75 mL	

CA125 antigen in a Tris-HCl buffered salt solution containing bovine serum albumin, detergent, and 0.05 % sodium azide as preservative. Ready for use.

BIOTIN Anti-CA125		
Biotin Anti-CA125	1 x 15 mL	2–8 °C until expiry date stated on the vial

Biotin Anti-CA125 monoclonal antibody from mouse, approximately 2 µg/mL. Contains Tris-HCl buffered saline (pH 7.75), bovine serum albumin, blocking agents, detergents, an inert red dye, and 0.05 % sodium azide as preservative. Ready for use.

CONJ Anti-CA125		
Tracer, HRP Anti-CA125	1 x 0.75 mL	2–8 °C until expiry date stated on the vial

Stock Solution of HRP Anti-CA125 monoclonal antibody from mouse, approximately 30 µg/mL. Contains preservatives. To be diluted with Tracer Diluent prior to use.

DIL CONJ		
Tracer Diluent	1 x 15 mL	2–8 °C until expiry date stated on the vial

Phosphate buffered saline (pH 7.2) with bovine serum albumin, blocking agents, detergents, an inert blue dye, and 0.01 % methyl-isothiazolone (MIT) as preservative. Ready for use.

Component	Quantity	Storage and stability after first opening
-----------	----------	---

SUBS	TMB
------	-----

TMB HRP-Substrate	1 x 12 mL	2–8 °C until expiry date stated on the vial
--------------------------	-----------	---

Contains buffered hydrogen peroxide and 3, 3', 5, 5' tetra-methylbenzidine (TMB). Ready for use.

STOP

Stop Solution	1 x 15 mL	2–8 °C until expiry date stated on the vial
----------------------	-----------	---

Contains 0.12 M hydrochloric acid. Ready for use.

WASHBUF	25X
---------	-----

Wash Concentrate	1 x 50 mL	2–8 °C until expiry date stated on the bottle
-------------------------	-----------	---

A Tris-HCl buffered salt solution with Tween 20. Contains Germall II as preservative. To be diluted with water 25 times before use.

Indications of instability

The TMB HRP-Substrate should be colorless or slightly bluish. A blue color indicates that the reagent has been contaminated and should be discarded.

WARNINGS AND PRECAUTIONS

For in vitro diagnostic use

- Please refer to the U.S. Department of Health and Human Services (Bethesda, Md., USA) publication No. (CDC) 88–8395 on laboratory safety procedures or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Reagents contain sodium azide (NaN_3) as a preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.
- Follow local guidelines for disposal of all waste material.

Protocol Sheet

CanAg CA125 EIA REF **400-10**

Mix the components directly before use. Use shaking conditions according to the Instructions.

Step	Vial/Plate	Procedure																																							
1. Prepare Wash Solution	WASHBUF 25X	Dilute 50 mL of Wash Concentrate with 1 200 mL of distilled or deionized water.																																							
Prepare Tracer working solution	CONJ Anti-CA125	Mix 50 µL of Tracer, HRP Anti-CA125 with 1mL of Tracer Diluent per strip:																																							
	DIL CONJ																																								
		<table border="1"><thead><tr><th>No. of Strips</th><th>Tracer, HRP Anti-CA125 (µL)</th><th>Tracer Diluent (mL)</th></tr></thead><tbody><tr><td>1</td><td>50</td><td>1</td></tr><tr><td>2</td><td>100</td><td>2</td></tr><tr><td>3</td><td>150</td><td>3</td></tr><tr><td>4</td><td>200</td><td>4</td></tr><tr><td>5</td><td>250</td><td>5</td></tr><tr><td>6</td><td>300</td><td>6</td></tr><tr><td>7</td><td>350</td><td>7</td></tr><tr><td>8</td><td>400</td><td>8</td></tr><tr><td>9</td><td>450</td><td>9</td></tr><tr><td>10</td><td>500</td><td>10</td></tr><tr><td>11</td><td>550</td><td>11</td></tr><tr><td>12</td><td>600</td><td>12</td></tr></tbody></table>	No. of Strips	Tracer, HRP Anti-CA125 (µL)	Tracer Diluent (mL)	1	50	1	2	100	2	3	150	3	4	200	4	5	250	5	6	300	6	7	350	7	8	400	8	9	450	9	10	500	10	11	550	11	12	600	12
No. of Strips	Tracer, HRP Anti-CA125 (µL)	Tracer Diluent (mL)																																							
1	50	1																																							
2	100	2																																							
3	150	3																																							
4	200	4																																							
5	250	5																																							
6	300	6																																							
7	350	7																																							
8	400	8																																							
9	450	9																																							
10	500	10																																							
11	550	11																																							
12	600	12																																							
2. Wash	MICROPLA	Wash each well once with Wash Solution. Use manual or automatic washer.																																							

3. Add calibrators, controls and samples	<table border="1"> <tr> <td data-bbox="10 713 41 888">CAL</td> <td data-bbox="41 713 165 888">CA125</td> </tr> <tr> <td colspan="2" data-bbox="10 655 165 713">0, 10, 40, 200, 500</td> </tr> <tr> <td data-bbox="10 728 41 888">CONTROL</td> <td data-bbox="41 728 165 888">CA125</td> </tr> </table>	CAL	CA125	0, 10, 40, 200, 500		CONTROL	CA125	25 µL in each well
CAL	CA125							
0, 10, 40, 200, 500								
CONTROL	CA125							
1, 2								
4. Add Biotin Anti-CA125	<table border="1"> <tr> <td data-bbox="186 713 227 888">BIOTIN</td> <td data-bbox="227 713 383 888">Anti-CA125</td> </tr> </table>	BIOTIN	Anti-CA125	100 µL in each well				
BIOTIN	Anti-CA125							
5. Incubate	<table border="1"> <tr> <td data-bbox="248 713 290 888">MICROPLA</td> </tr> </table>	MICROPLA	2 hour shaking at room temperature					
MICROPLA								
6. Wash	<table border="1"> <tr> <td data-bbox="310 713 352 888">MICROPLA</td> </tr> </table>	MICROPLA	Wash each well three times with Wash Solution. Use manual or automatic washer.					
MICROPLA								
7. Add Tracer working solution	<table border="1"> <tr> <td data-bbox="372 713 414 888">TRACER WORKING SOLUTION</td> </tr> </table>	TRACER WORKING SOLUTION	100 µL in each well					
TRACER WORKING SOLUTION								
8. Incubate	<table border="1"> <tr> <td data-bbox="435 713 476 888">MICROPLA</td> </tr> </table>	MICROPLA	1 hour shaking at room temperature					
MICROPLA								
9. Wash	<table border="1"> <tr> <td data-bbox="497 713 538 888">MICROPLA</td> </tr> </table>	MICROPLA	Wash each well six times with Wash Solution. Use manual or automatic washer.					
MICROPLA								
10. Add TMB HRP-Substrate	<table border="1"> <tr> <td data-bbox="559 713 600 888">SUBS</td> <td data-bbox="600 713 652 888">TMB</td> </tr> </table>	SUBS	TMB	100 µL in each well				
SUBS	TMB							
11. Incubate	<table border="1"> <tr> <td data-bbox="621 713 663 888">MICROPLA</td> </tr> </table>	MICROPLA	30 min shaking at room temperature					
MICROPLA								
12. Read absorbance	<table border="1"> <tr> <td data-bbox="683 713 725 888">MICROPLA</td> </tr> </table>	MICROPLA	620 nm					
MICROPLA								
Alt.12 Add Stop Solution	<table border="1"> <tr> <td data-bbox="745 713 787 888">STOP</td> </tr> </table>	STOP	100 µL in each well					
STOP								
Alt.13 Incubate	<table border="1"> <tr> <td data-bbox="808 713 849 888">MICROPLA</td> </tr> </table>	MICROPLA	1 min shaking at room temperature					
MICROPLA								
Alt.14 Read absorbance	<table border="1"> <tr> <td data-bbox="870 713 911 888">MICROPLA</td> </tr> </table>	MICROPLA	Read at 405 nm within 5 min					
MICROPLA								

Caution

Material used in the preparation of human source reagent has been tested and found to be Non Reactive for HIV 1 and 2 Antibody, HCV Antibody and Hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they were potentially infectious.

CLP (1272/2008) HAZARD CLASSIFICATION

Information about CLP (1272/2008) HAZARD CLASSIFICATION can be found at the end of this document.

SPECIMEN COLLECTION AND HANDLING

The CanAg CA125 EIA is intended for use with serum. Collect blood by venipuncture and separate the serum according to common procedures. Samples can be stored at 2–8 °C for 24 hours. For longer periods store samples at -70 °C or below. Samples should not be stored in a self-defrosting freezer and not be thawed and refrozen before analysis. Allow frozen samples to thaw slowly at 2–8 °C over night and then bring the samples to room temperature before analysis.

PROCEDURE

Materials required but not supplied with the kit

- 1. Microplate shaker**
Shaking should be medium to vigorous, approximately 700-1100 oscillations/min.
- 2. Microplate wash device**
Automatic plate wash capable of performing 1, 3 and 6 washing cycles, with a minimal fill volume of 350 µL/well/washcycle.
An 8-channel pipette with disposable plastic tips for delivery of 350 µL is recommended if an automatic microplate washer is not used.
- 3. Microplate spectrophotometer**
With a wavelength of 620 nm and/or 405 nm, and an absorbance range of 0 to 3.0.
- 4. Precision pipettes**
With disposable plastic tips for dispensing microlitre volumes. An 8-channel pipette or resenser pipette with disposable plastic tips for delivery of 100 µL is useful but not essential. Pipettes for dispensing millilitre volumes.
- 5. Distilled or deionized water**
For preparation of Wash Solution.

Procedural notes

1. A thorough understanding of this package insert is necessary to ensure proper use of the CanAg CA125 EIA kit. The reagents supplied with the kit are intended for use as an integral unit. Do not mix identical reagents from kits having different lot numbers. Do not use the kit reagents after the expiry date printed on the outside of the kit box.
2. Reagents should be allowed to reach room temperature (20–25 °C) prior to use. The assay should only be performed at temperatures between 20–25 °C to obtain accurate results. Frozen specimens must be gently but thoroughly mixed after thawing.
3. Before starting to pipette calibrators and unknown specimens it is advisable to mark the strips to be able to clearly identify the samples during and after the assay.
4. The requirement for efficient and thorough washing for separation of bound and unbound antigen and reagents from the solid-phase bound antibody-antigen complexes is one of the most important steps in an EIA. In order to ensure efficient washing make sure that all wells are completely filled to the top edge with wash solution during each wash cycle, that wash solution is dispensed at a good flow rate, that the aspiration of the wells between and after the wash cycles is complete and that the wells are empty. If there is liquid left, invert the plate and tap it carefully against absorbent paper.
 - Automatic strip washer: Follow the manufacturer's instructions for cleaning and maintenance diligently and wash the required number of wash cycles prior to and after each incubation step. It's highly recommended to use *plate* process mode and *overflow* wash mode with a dispensing volume of 800 µL. The aspiration/wash device should not be left standing with the Wash Solution for long periods, as the needles may get clogged resulting in poor liquid delivery and aspiration.
5. The TMB HRP-Substrate is very sensitive to contamination. For optimal stability of the TMB HRP-Substrate, pour the required amount from the vial into a carefully cleaned reservoir or preferably a disposable plastic tray to avoid contamination of the reagent. Be sure to use clean disposable plastic pipette tips (or dispenser pipette tip).
6. Be sure to use clean disposable plastic pipette tips and a proper precision pipetting technique when handling samples and reagents. Do not allow the pipette tip to touch the surface of the liquid in order to avoid carry-over. A diligent pipetting technique is of particular importance when handling the samples and the TMB HRP-Substrate solution.

Preparation of reagents	Stability of prepared reagent
Wash Solution	2 weeks at 2–25 °C in a sealed container
Tracer working solution	3 weeks at 2–8 °C in a sealed container

Pour the 50 mL Wash Concentrate into a clean container and dilute 25-fold by adding 1200 mL of distilled or deionised water to give a buffered Wash Solution.

Prepare the required quantity of Tracer working solution by mixing 50 µL of Tracer, HRP Anti-CA125 with 1 mL of Tracer Diluent per strip (see table below):

No. of Strips	Tracer, HRP Anti-CA125 (µL)	Tracer Diluent (mL)
1	50	1
2	100	2
3	150	3
4	200	4
5	250	5
6	300	6
7	350	7
8	400	8
9	450	9
10	500	10
11	550	11
12	600	12

Be sure to use a clean plastic or glass bottle for preparation of Tracer working solution.

Alternative: Pour the content of the Tracer, HRP Anti-CA125 into the vial of Tracer Diluent and mix gently. Make sure that the entire content of the Tracer, HRP Anti-CA125 is transferred to the vial of Tracer Diluent.

NOTE: The Tracer working solution is stable for 3 weeks at 2–8 °C. Do not prepare more Tracer working solution than will be used within this period and make sure that it is stored properly.

Assay procedure

Perform each determination in duplicate for both calibrators and patient samples. A calibration curve should be run with each assay. All reagents and samples must be brought to room temperature (20–25 °C) before use.

1. Start to prepare Wash Solution and Tracer working solution. It is important to use clean containers. Follow the instructions carefully.
2. Transfer the required number of microplate strips to a strip frame. (Immediately return the remaining strips to the aluminium pouch containing a desiccant and reseal carefully). Wash each strip once with the Wash Solution. Do not wash more strips than can be handled within 30 min.
3. Pipette 25 μL of the CA125 Calibrators (CAL 0, 10, 40, 200, 500), CA125 Controls (C 1, C 2) and patient specimens (unknowns-Unk) into the strip wells according to the following scheme:

	1	2	3	4	5	6	7 etc
A	Cal 0	Cal 500	2nd Unk				
B	Cal 0	Cal 500	2nd Unk				
C	Cal 10	C1	etc.				
D	Cal 10	C1					
E	Cal 40	C2					
F	Cal 40	C2					
G	Cal 200	1st Unk					
H	Cal 200	1st Unk					

4. Add 100 μL of Biotin Anti-CA125 to each well using a 100 μL precision pipette (or an 8-channel 100 μL precision pipette). Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid.
5. Incubate the plate for 2 hours (\pm 10 min) at room temperature (20-25 °C) with constant shaking of the plate using a microplate shaker.

6. After the first incubation aspirate and wash each strip 3 times using the wash procedure described in Procedural notes, item 4.
7. Add 100 μ L of Tracer working solution to each well. Use the same pipetting procedure as in item 4 above.
8. Incubate the frame for 1 hour (\pm 5 min) at room temperature (20–25 °C) with constant shaking.
9. After the second incubation aspirate and wash each strip 6 times, using the wash procedure described in Procedural notes, item 4.
10. Add 100 μ L of TMB HRP-Substrate to each well using the same pipetting technique as in item 4. The TMB HRP-Substrate should be added to the wells as quickly as possible and the time between addition to the first and last well should not exceed 5 min.
11. Incubate for 30 min (\pm 5 min) at room temperature with constant shaking. Avoid exposure to direct sunlight.
12. Immediately read the absorbance at 620 nm in a microplate spectrophotometer.

Option

If the laboratory does not have access to a microplate reader capable of reading at 620 nm, the absorbance can be determined as in item alternative 12:

Alt. 12. Add 100 μ L of Stop Solution, mix and read the absorbance at 405 nm in a microplate spectrophotometer within 5 min after addition of Stop Solution.

Measurement range

The CanAg CA125 EIA measures concentrations between 1.5 and 500 U/mL. If CA125 concentrations above the measuring range are to be expected, it is recommended to dilute samples with CA125 Calibrator 0 prior to analysis.

Quality control

CA125 Control 1 and 2 should be used for validation of the assay series. Ranges of expected results are indicated on the vial labels. If values outside of the specified range are obtained, a complete check of reagents and reader performance should be made and the analysis repeated. Each laboratory may also prepare its own serum pools at different levels, which can be used as internal controls in order to assure the precision of the assay.

Reference materials

Since no common reference material is available for CA125 antigen, CanAg CA125 EIA Calibrator values are assigned against a set of in-house reference standards.

CALCULATION OF RESULTS

If a microplate spectrophotometer with built-in data calculation program is used, refer to the manual for the spectrophotometer and create a program using the concentration stated on the label of each of the CA125 Calibrators.

For automatic calculation of CA125 results it is recommended to use either of the following methods:

- Cubic spline curve fit method. Calibrator 0 should be included in the curve with the value 0 U/mL.
- Spline smoothed curve fit method. Calibrator 0 should be used as plate blank.
- Interpolation with point-to-point evaluation. Calibrator 0 should be included in the curve with the value 0 U/mL.
- Quadratic curve fit method. Calibrator 0 should be included in the curve with the value 0 U/mL.

NOTE: 4-Parametric or Linear regression evaluation methods should not be used.

For manual evaluation, a calibration curve is constructed by plotting the absorbance (A) values obtained for each CA125 Calibrator against the corresponding CA125 concentration (in U/mL), see figure below. The unknown CA125 concentrations can then be read from the calibration curve using the mean absorbance value of each patient specimen.

If samples in an initial analysis give CA125 levels higher than 500 U/mL the samples should be diluted 1/10 and 1/100 with CA125 Calibrator 0 to obtain the accurate CA125 concentration of the samples.

1/10 dilution = 50 μ L of specimen + 450 μ L of CA125 Calibrator 0

1/100 dilution = 50 μ L of 1/10 dilution + 450 μ L of CA125 Calibrator 0

The CA125 concentration of the undiluted sample is then calculated as:

Dilution 1/10: 10 x measured value

Dilution 1/100: 100 x measured value

Example of results

Specimen			Calibrator values	Mean abs value (A)	CA125 U/mL
CAL	CA125	0	0 U/mL	0.047	
CAL	CA125	10	10 U/mL	0.116	
CAL	CA125	40	40 U/mL	0.298	
CAL	CA125	200	200 U/mL	1.269	
CAL	CA125	500	500 U/mL	2.218	

Specimen A


0.490

69.8

Specimen B

1.650

325

Example, do not use this curve to determine assay results.

LIMITATIONS OF THE PROCEDURE

The level of CA125 cannot be used as absolute evidence for the presence or absence of malignant disease and the CA125 test should not be used in cancer screening. The results of the test should be interpreted only in conjunction with other investigations and procedures in the diagnosis of disease and the management of patients, and the CA125 test should not replace any established clinical examination.

Anti-reagent antibodies (human anti-mouse antibody (HAMA) or heterophilic antibodies) in the patient sample may occasionally interfere with the assay, even though specific blocking agents are included in the buffers.

Biotin may interfere with the assay giving false low results. This should be taken into consideration for patients taking dietary supplements or receiving therapy containing high (> 5mg/day) or extremely high (300 mg/day) biotin doses. Peak serum levels have been reported to occur 1-3h post ingestion and a physiological half-life of 8-16h depending on renal function (9,10). In a study by Grimsey et al. (9) a specimen concentration of 30 ng/mL was reached 8h following intake of 10 mg of biotin. For patients on very high doses of biotin it is recommended to stop taking biotin for at least 2 days before blood draw (11).

EXPECTED VALUES

CanAg CA125 was measured in 100 healthy female blood donors. The mean value obtained was 14.7 U/mL with a standard deviation of 7.7. The median value was 13.1 U/mL, range 5.06– 47.9 U/mL. The lower and upper extremes of the normal range were examined using IFCC recommended non-parametric statistical treatment. The reference interval contains the central 95% fraction of the reference distribution. Reference limits may accordingly be estimated as the 2.5% (lower) and 97.5% (upper) fractiles. These limits cut off a fraction of 2.5% of the values in each tail of the reference distribution. Non-parametric estimates:

N=100

Fraction	Reference limit (U/mL)
2.5 th (lower)	5
97.5 th (upper)	39

96% of the healthy women had assay values below 35 U/mL.

It is recommended that each laboratory establish its own normal range to account for such local environmental factors as diet, climate, living conditions, patient selection, etc. It should also be borne in mind that the individual patient's own baseline result provides the most important reference point for interpretation of marker results (6).

PERFORMANCE CHARACTERISTICS

Precision

Total precision was determined according to NCCLS guideline EP5-A (7) using four levels of frozen pooled human serum containing added ascitespool. Each sample was randomly pipetted (n=2/analysis) and analysed twice each day over 10 days. The analyses were undertaken during a period of 8 months, by ≥ two different technicians and using 10 different CanAg CA125 EIA kit batches.

Sample	Replicates	Mean U/mL	Within-run SD (U/mL)	Within-run CV %	Between-day SD (U/mL)	Between-day CV %
CA125 1	40	16.8	0.74	4.4	0.53	3.1
CA125 2	40	75.7	3.26	4.3	2.42	3.2
CA125 3	40	201	8.55	4.3	7.58	3.8
CA125 4	40	392	11.4	2.9	15.5	4.0

Detection limit

The detection limit of the CanAg CA125 EIA assay is < 1.5 U/mL defined as the concentration corresponding to the mean of the absorbance values for the CA125 Calibrator 0 plus 2 standard deviations according to the formula:

$$\frac{2 \times \text{SD CAL 0}}{\text{OD CAL 10} - \text{OD CAL 0}} \times 10 \text{ U/mL}$$

Recovery

Spiked serum samples were prepared by adding human CA125 antigen to normal serum samples. The recovery of the antigen was 100 ±15%.

Hook effect

No hook effect has been noticed with samples up to 50 000 U/mL. **NOTE:** In very high samples the color of the substrate will change from blue to greenish (and eventually yellow in extremely high samples). This will lead to a falsely low absorbance at 620 nm, and in extreme cases the absorbance may fall within the calibration curve range and noticed as a hook.

Linearity

Patient samples were serially diluted with CA125 Calibrator 0 and analyzed. The obtained values were 100 ±15% of the expected values.

Specificity

The CanAg CA125 EIA is based on two mouse monoclonal antibodies, Ov197 and Ov185, directed against two independent epitopes of the protein core of the CA125 antigen (4). The NCCLS guideline EP7-P (8) was followed to determine possible sources of interference. The following substances and concentrations were tested and found not to interfere with the test.

	Concentration with no significant (\pm 10 %) interference
Lipemia (Intralipid®)	4 mg/mL
Bilirubin, unconjugated	0.6 mg/mL
Hemoglobin	5 mg/mL

Biotin interference

A study was conducted to evaluate biotin interference. Low, medium and high serum control samples were spiked to final biotin concentrations of 15, 30, 60, and 200 ng/mL. The mean CA125 concentration was determined for each sample and the percent recovery for each biotin concentration was calculated using the formula: Recovery (%) = 100x (Mean CA125 concentration w. biotin added/Mean CA125 concentration w. diluent only added).

CA125 analyte level	Biotin test conc. (ng/mL)	Expected CA125 conc. (U/mL)	Observed CA125 conc. (U/mL)	Recovery (%)
Low	15	27,6	27,2	98
Medium	15	112,1	110,8	99
High	15	408,0	410,6	101
Low	30	27,6	26,6	96
Medium	30	112,1	110,1	98
High	30	408,0	414,0	101
Low	60	27,6	26,3	95
Medium	60	112,1	106,3	95
High	60	408,0	419,3	103
Low	200	27,8	17,0	61
Medium	200	112,4	65,9	59
High	200	411,5	301,5	73

Based on linear regression analysis, the lowest concentration of biotin found to influence test results (\geq 10%) was 76 ng/mL.

CLP (1272/2008) HAZARD CLASSIFICATION

The following warnings and precautions apply to

SUBS TMB

Hazard pictograms:

Signal word:	Danger
Hazard Statement:	Repr. 1B: H360D May damage the unborn child.
Prevention statement:	P202 Do not handle until all safety precautions have been read and understood.
Prevention:	P280 Wear protective gloves / protective clothing / eye protection / face protection.
Precautionary statement response:	P308+P313 IF exposed or concerned get medical advice/attention.
Precautionary statement disposal:	P501 Dispose of contents / container to an approved hazardous / special waste disposal facility in accordance with local and national regulations.

Restricted to professional users.

Hazardous substances: 2- Pyrrolidone

Other hazards

None of the mixtures in the kit contains any substances considered to meet the criteria classifying them as PBT and/or vPvB.

WARRANTY

The performance data presented here were obtained using the assay procedure indicated. Any change or modification of the procedure not recommended by Fujirebio Diagnostics may affect the results, in which Fujirebio Diagnostics disclaims all warranties expressed, implied or statutory including the implied warranty of merchantability and fitness for use.

REFERENCES

1. Bast R., Feeney M., Lazarus H., Nadler L., Colvin R., Knapp R., (1981). Reactivity of a monoclonal antibody with human ovarian carcinoma. *J Clin Invest* 68: 1331 - 1337.
2. O'Brien T, Raymond L, Bannon G, Ford D, Hardardottir H, Miller F, Quirk J (1991). New monoclonal antibodies identify the glycoprotein carrying the CA125 epitope. *Am J Obstet Gynecol* 165: 1857 - 1864.
3. Nilsson O, Jansson E-L, Dahlen U, Nilsson K, Nustad K, Högberg T, Lindholm L (1994). In: Current Tumor Diagnosis: Applications, Clinical Relevance, Research, Trends. Ed R. Klapdor. Zucksverdt Verlag München, pp 401-405.
4. Nustad K. et al., (1996) Specificity and affinity of 26 monoclonal antibodies against the CA125 antigen: First report from the ISOBM TD-1 Workshop. *Tumor Biol* 17:196-219.
5. Bonfrer J., Duffy M., Radtke, M., Segurado O., Torre G., van Dalen A., Zwirner M. (1999) Tumor markers in gynaecological cancers-EGTM recommendations. *Anticancer Research* 19:2785-2820.
6. Tuxen M., Sölétormos, G., Rustin, G., Nelstrop A., Dombrowsky, P. (2000) Biological variation and analytical imprecision of CA125 in patients with ovarian cancer. *Scand J Clin Invest* 60:713-722.
7. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
8. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).
9. Grimsey P et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. *J Pharmacokinet.* 2017; 2(4); 247-256.
10. Jenkins Colon P, Greene D.N. Biotin Interference in Clinical Immunoassays. *JALM* 2018; 2(6); 941-951.
11. Chun KY. Biotin Interference in Diagnostic Tests. *Clinical Chemistry* 2017; 63; 619-620.

CanAg[®] is a registered trademark of Fujirebio Diagnostics AB

Fujirebio Diagnostics AB

Elof Lindälvs gata 13

SE-414 58 Göteborg

Sweden

Phone + 46 31 85 70 30

Fax + 46 31 85 70 40

info@fdab.com

www.fdab.com