CanAg NSE EIA

REF

420-10

IVD

Instructions for use. 2022-06

Read highlighted changes

EN	EXPLANATION OF SYMBOLS
BG	ОБЯСНЕНИЕ НА СИМВОЛИТЕ
CS	VÝZNAM SYMBOLŮ
DA	SYMBOLFORKLARING
DE	ERKLÄRUNG DER SYMBOLE
EL	ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
ES	SIGNIFICADO DE LOS SÍMBOLOS
ET	SÜMBOLITE SELGITUS
FR	EXPLICATION DES SYMBOLES
HR	OBJAŠNJENJE SIMBOLA
HU	JELMAGYARÁZAT
IT	SPIEGAZIONE DEI SIMBOLI
LT	SIMBOLIŲ PAAIŠKINIMAI
LV	SIMBOLU SKAIDROJUMS
NL	VERKLARING DER SYMBOLEN
NO	SYMBOLFORKLARING
PL	OBJAŚNIENIE SYMBOLI
PT	EXPLICAÇÃO DOS SÍMBOLOS
RO	SEMNIȚAȚIA SIMBOLURILOR
RU	ОБОНАЧЕНИЯ
SV	SYMBOLFÖRKLARING
SK	VÝZNAM SYMBOLOV
SL	RAZLAGA SIMBOLOV
SR	OBJAŠNJENJE SIMBOLA
TR	SEMBOLLERİN AÇIKLAMALARI

Use By/Годно до/Použitelné do/
Holdbar til/Verwendbar bis/
Ημερομηνία λήξης/Fecha
de caducidad/Kölblik kuni/
Utiliser jusque/Rok valjanosti/
Felhasználható/Utilizzare entro/
Sunautoti iki/Izlietot līdz/Houdbaar
tot/Brukes innen/Użyç przed/
Prazo de validade/Expirã la/
Использовать до/Använd före/
Použite né do/ Uporabno do/
Upotrebljivo do/Son Kullanna Tarihi

LOT

Batch code/Номер на партида/
Číslo šarže/Lotnummer/
Chargenbezeichnung/Αριθμός
Παρτίδας/Código de lote/Partii
kood/Code du lot/Kod serije/
Sarzsám/Codice del lotto/
Partijos kudas/Partijas kods/Lot
nummer/Partikode/Kod partii/
Código do lote/Număr de lot/
Номер лота/Lotnummer/Číslo
šarže/Številka serije/Kod partije/
Parti Kodu

Date of manufacture/Дата на производство/Datum výroby/
Produktionsdato/Herstellungsdatum/
Ημερομηνία παραγωγής/Fecha de fabricación/Valmistamise kuupäev/
Date de fabrication/Datum proizvodnje/
Gyártási idő/Data di produzione/
Pagaminimo data/Ražošanas datums/
Productiedatum/Fremstillingsdato/
Data produkcji/Data de fabrico/Data fabricației/Дата производства/
Tillverkningsdatum/Dátum výroby/Datum izdelave/Datum proizvodnje/Úretim tarihi

Temperature limitation/
Температурни граници/
Теплотни омеzeи/
Temperaturbegrænsning/
Temperaturbegrenzung/
Περιορισμοί θερμοκρασίας/
Limites de temperatura/
Temperatuuri piirang/
Limite de température/
Temperaturno ograničenje/
Hőmérsékletre vonatkozó korlátozás/
Limiti di temperatura/
Temperatūriniai apribojimai/
Temperatūras ierobežojums/
Temperaturbepæring/
Temperaturbegrensninger/
Temperaturey granicne/
Limite de temperatura/
Limite de temperatură/
Температурный режим/
Temperaturbegrænsning/
Теплотне обмеzenie
Omejitve temperature/
Temperaturno ograničenje/
Sıcaklık sınırlaması/

IVD

In Vitro Diagnostic Medical Device/
Медицински уред за диагностика
ин vitro/Diagnostický zdravotnícký
prostředek in vitro/Medicinsk udstyr til
in vitro-diagnostik/In-vitro-Diagnostikum/
Ιατροτεχνολογικό προϊόν για διάγνωση
In Vitro/Dispositivo médico para
diagnóstico in vitro/In vitro diagnostiline
meditsiiniseade/Dispositif médical de
diagnostic in vitro/Diagnostički medicinski
uređaj In Vitro/In vitro orvosdiagnostikai
eszköz/Dispositivo medico per test
diagnostici in vitro/In Vitro Diagnostinė
Medicinos Priemonė/Medicínska ierice
in vitro diagnostikai/In vitro-diagnostisch
medisch instrument/In vitro diagnostisk
medisinsk utstyr/Wyrób medyczny do
diagnostyki in vitro/Dispositivo Médico
de Diagnóstico In Vitro/Dispozitiv medical
pentru diagnostic in vitro/Только для
диагностики In Vitro/Endast för in
vitro-diagnostik/ Zdravotnicka pomôcka na
diagnostiku in vitro/In vitro diagnostični
pripomoček/Diagnostički medicinski
uređaj In Vitro/<96> testleri için yeterlilik
içerir

Contains sufficient for <96> tests/Съдържа
достатъчно количество за тестове
<96>/Lze použít pro <96> testů/Ineholder
tilstrækkeligt/Inhalt ausreichend für <96>
Prüfungen/Περιοχόμενο επαρκές για
«96» εξετάσεις/Contenido suficiente para
<96> ensayos/Kogusest piisab <96> testi
lääbiviimiseks/Contenu suffisant pour «96»
tests/Sadržaj dovoljno za <96> testova/A
doboz tartalma <96> vizsgálat elvégzéséhez
elegendő/Contenuto sufficiente per «96»
saggi/Turiny's skirtas atlikti <96> tyrimus/
Saturis pietiekams <96> testiem/Inhoud
voldoende voor «96» testen/til «96» test/
Tilstrækkelig innhold for <96> prøver/
Wystarczy na wykonanie <96> testów/
Conteúdo suficiente para «96» ensaios/
Conținut suficient pentru 96 de teste/
Содержит достаточные количества для
«96» определений/Innehåller tillräckligt
till «96» antal tester/Obsah postačuje na
tento počet testov: <96>/Vsebinsa zadostuje
za <96> testov/Sadržina dovoljna za <96>
testova/<96> testleri için yeterlilik içerir

REF

Catalogue number/Каталожен номер/
Katalogové číslo/Katalognummer/
Bestellnummer/Αριθμός καταλόγου/
Número de catálogo/Katalogoi number/
Numéro de catalogue/Kataloški broj/
Katalógusszám/Numero di catalogo/
Katalogo numeris/Numurs katalogā/
Catalogusnummer/Katalognummer/
Numer katalogowy/Número do catálogo/
Număr de catalog/Номер по каталогу/
Produktnummer/Katalógové číslo/
Kataloška številka/Kataloški broj/
Katalog numarası

Consult Instructions for Use/
Прочетете инструкцията за
употреба/Konzultujte s návodem
k použití/Se brugsanvisning/Siehe
Gebrauchsanweisung/Συμβουλευτείτε
της Οδηγίας σχετικά με τη χρήση/
Consulte las instrucciones de uso/
Vt kasutusjuhendit/Consulter le mode
d'emploi/Pročítajte upute za uporabu/
Olvassa el a használati utasítást/
Consultare le istruzioni per l'uso/Dél
naudojimo žiūrėkite instrukcijas/Izlasiet
lietošanas instrukciju/Raadpleeg de
instructies voor gebruik/Les instruksene
for bruk/Sprawdzić w instrukcji użycia/
Consulte as Instruções de Utilização/
Consultați instrucțiunile de utilizare/
Обратитесь к инструкции по
применению/Se bruksanvisning/
Prečítajte si návod na používanie/
Pročítajte uputstvo za upotrebu/
Kullanım Talimatlarını Bakınız

Contents of kit/Съдържание на набора/
Obsah soupravy/Kittets indhold/Inhalt
des Kits/Περιεχόμενα του κιτ/Contenido
del kit/Komplekt sisaldab/Contenu du
kit/Sadržaj opreme/A készlet tartalma/
Contenuto del kit/Rinkinio turinys/
Komplekta saturs/Inhoud van de set/
Settets innhold/Zawartość zestawu/
Conteúdo do kit/Conținutul setului/
Компоненты набора/Kit innehåll/
Obsah súpravy/Vsebina kompleta/Sadržaj
opreme/Kitin içindekiler

Biological risks/Биологическа
опасност/Biológická rizika/Biologisk
fare/Biologische Gefahren/Βιολογικοί
κίνδυνοι/Riesgos biológicos/
Biolooigilised ohud/Risques biologiques/
Biolóškli rizici/Biológiai kockázatok/Rischi
biologici/Biologinis pavojus/Biolóškiskais
risks/Biologische risico's/Biologiske
risikoer/Zagroženie biologiczne/Riscos
biológicos/ Biologisk risk/Pericole
biologice/Биологическая опасность/
Biologický rizikové/Biológické riziká/
Biolóškli rizici/Biyolojik riskler

Human/C човешки производ/Ľidské/
Human/Human/ἄνθρωπος αναφοράς/
Humano/Inimāritolu/Humaine/Ljudskog
porjekla/Humán/Origine Umana/
Žmogaus kilmės/Cilvēku izcelsmes/
Human/Menneske/Ludzka/Humano/
Origine umână/Человеческого
происхождения/Human/Ludské/
Humanega izvora/Ljudskog porekla/Ľnsan

From mouse/C миши производ/Myši/
Fra mus/Maus/από ποντίκι/de ratón/
Hiirtelt/De souris/Mišijeg porjekla/
Egérböli/Murino/Pelés kilmés/No peles/
Van muizen/Fra mus/Mysia/Do rato/De
la șoareci/Мышиного происхождения/
Från mus/Myšije/Mišjega izvora/Mišijeg
porekla/Fareden

Bovine/C говежди производ/
Hovēži/Bovin/Rind/από βοοειδή/
Bovino/Veistelt/Bovine/Rogate stoke/
Szarvasmarha/Bovino/Jaučio/No
liellopa/Bovien/Bovin/Wolowy/Bovino/
Origine bovină/крупного рогатого
скота/Från ko/Hovädzie/Rogvega
izvora/Rogate krupne stoke/Bovin

Reconstitute with/Пазтваряне с/
Rozfeďte pomoci/Rekonstitueres med/
Rekonstituieren mit/Ανασύσταση με/
Reconstituir con/Lahjendamine/
Reconstituer avec/Rekonstituiraite s/
Feloldáshoz/Ricostituire con/Atkurti,
ištirpdant su/Atšķaidīt ar/Reconstituire
met/Rekonstitueres med/Odtworzyć
za pomocą/Reconstituir com/A
se reconstitui cu/Пастворить в/
Rekonstituera med/Rozriedte pomocou/
Rekonstituiraite z/s/Ponovno formiranje
sa/Yeniden oluşturalur

Manufacturer/Производитель/Výrobce/
Producent/Hersteller/Κατασκευαστής/
Fabricante/Tootja/Fabricant/Proizvođač/
Gyártó/Fabbricante/Gamintojas/
Ražotājs/Fabrikant/Produsent/
Producent/Fabricante/Producător/
Производитель/Tilverkare/ Výrobca/
Izdelovalec/Proizvođač/Üretici

INSTRUCTIONS FOR USE

EN

INSTRUCTIONS FOR USE

Please visit our website www.fdi.com/ifu to obtain the Instructions For Use (IFU) in additional languages.

To ensure that you download the correct IFU for your kit lot, please select the revision corresponding to the issue date printed on the front page of the IFU provided with this kit.

Please follow the IFU carefully. Instructions for safe handling are found in the WARNINGS AND PRECAUTIONS section. Material Safety Data Sheets (MSDS) are available on our website www.fdi.com. If you do not have access to the internet, please contact your local distributor, or Fujirebio Diagnostics AB for assistance.

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CanAg NSE EIA

Instructions for use

Enzyme immunometric assay kit
For 96 determinations

INTENDED USE

The CanAg NSE EIA kit is intended for the quantitative determination of NSE in human serum.

SUMMARY AND EXPLANATION OF THE ASSAY

The glycolytic enzyme enolase (2-phospho-D-glycerate hydrolase, EC 4.2.1.11) exists as several dimeric isoenzymes ($\alpha\alpha$, $\alpha\beta$, $\alpha\gamma$, $\beta\beta$ and $\gamma\gamma$) composed of three distinct subunits α , β and γ . The γ unit is found either in a homologous $\gamma\gamma$ - or in a heterologous $\alpha\gamma$ -isoenzyme and is known as neuron-specific enolase (NSE). The monoclonal antibodies used in the CanAg NSE EIA bind to the γ -subunit of the enzyme and thereby detects both the $\gamma\gamma$ and the $\alpha\gamma$ forms (1, 2). The NSE levels are low in healthy subjects and subjects with benign diseases. Elevated levels are commonly found in patients with malignant tumours with neuroendocrine differentiation, especially small cell lung cancer (SCLC) (3) and neuroblastoma (4). Quantitative determination of NSE in serum may be valuable in the management of patients with suspected or diagnosed SCLC or neuroblastoma, to aid in the differential diagnosis and to monitor the effect of treatment (5, 6).

PRINCIPLE OF THE TEST

The CanAg NSE EIA is a solid phase, non-competitive immunoassay based on two monoclonal antibodies (derived from mice) directed against two separate antigenic determinants of the NSE molecule. The monoclonal antibodies (MAb) used bind to the γ -subunit of the enzyme and thereby detects both the $\gamma\gamma$ and the $\alpha\gamma$ form. Calibrators and patient samples are incubated together with biotinylated Anti-NSE MAb E21 and horseradish peroxidase (HRP) labelled Anti-NSE MAb E17 in streptavidin coated micro strips. After washing, buffered Substrate/Chromogen reagent (hydrogen peroxide and 3, 3', 5, 5' tetramethylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction a blue colour will develop if antigen is present. The intensity of the colour development is proportional to the amount of NSE present in the samples.

The colour intensity is determined in a microplate spectrophotometer at 620 nm (or optionally at 405 nm after addition of Stop Solution).

Calibration curves are constructed for each assay by plotting absorbance value versus the concentration for each calibrator. The NSE concentrations of patient samples are then read from the calibration curve.

REAGENTS

- Each CanAg NSE EIA kit contains reagents for 96 tests.
- The expiry date of the kit is stated on the label on the outside of the kit box.
- Do not use the kit beyond the expiry date.
- Do not mix reagents from different kit lots.
- Store the kit at 2–8 °C. Do not freeze.
- Opened reagents are stable according to the table below provided they are not contaminated, stored in resealed original containers and handled as prescribed. Return to 2–8 °C immediately after use.

Component	Quantity	Storage and stability after first opening
-----------	----------	---

MICROPLA

Microplate	1 Plate	2–8 °C until expiry date stated on the plate
-------------------	---------	--

12 x 8 breakable wells coated with streptavidin. After opening, immediately return unused strips to the aluminium pouch containing desiccant and reseal carefully to keep dry.

NSE Calibrators	5 vials, lyophilised	4 weeks at 2–8 °C 3 months at –20 °C
------------------------	----------------------	---

CAL	NSE	A	1 x 0.75 mL
-----	-----	---	-------------

CAL	NSE	B	1 x 0.75 mL
-----	-----	---	-------------

CAL	NSE	C	1 x 0.75 mL
-----	-----	---	-------------

CAL	NSE	D	1 x 0.75 mL
-----	-----	---	-------------

CAL	NSE	E	1 x 0.75 mL
-----	-----	---	-------------

The lyophilised calibrators contain human NSE in a protein matrix with 0.01 % of a non-azide preservative. To be reconstituted with 0.75 mL distilled water before use.

NOTE: The exact NSE concentration is lot specific and is indicated on the label of each vial.

Component	Quantity	Storage and stability after first opening
-----------	----------	---

BIOTIN	Anti-NSE
---------------	-----------------

Biotin Anti-NSE	1 x 15 mL	2–8 °C until expiry date stated on the vial
-----------------	-----------	---

Biotin Anti-NSE monoclonal antibody from mouse, approximately 2 µg/mL. Contains phosphate buffer (pH 7.1), bovine serum albumin, blocking agents, an inert blue dye and 0.01 % methyl-isothiazolone (MIT) as preservative. To be mixed with Tracer, HRP Anti-NSE before use.

CONJ	Anti-NSE
-------------	-----------------

Tracer, HRP Anti-NSE	1 x 0.75 mL	2–8 °C until expiry date stated on the vial
----------------------	-------------	---

Stock solution of HRP Anti-NSE monoclonal antibody from mouse, approximately 40 µg/mL. To be mixed with Biotin Anti-NSE prior to use. Contains 0.02 % methyl-isothiazolone (MIT), 0.02 % bromonitrodioxane and 20 ppm Proclin™ 300 as preservatives.

SUBS	TMB
-------------	------------

TMB HRP-Substrate	1 x 12 mL	2–8 °C until expiry date stated on the vial
-------------------	-----------	---

Ready for use. Contains buffered hydrogen peroxide and 3, 3', 5, 5' tetramethylbenzidine (TMB).

Component	Quantity	Storage and stability after first opening
-----------	----------	---

STOP

Stop Solution	1 x 15 mL	2–8 °C until expiry date stated on the vial
---------------	-----------	---

Ready for use. Contains 0.12 M hydrochloric acid.

WASHBUF	25X
----------------	------------

Wash Concentrate	1 x 50 mL	2–8 °C until expiry date stated on the bottle
------------------	-----------	---

To be diluted with water 25 times before use. A Tris-HCl buffered salt solution with Tween 20. Contains Germall II as preservative.

Indications of instability

The TMB HRP-Substrate should be colourless or slightly bluish. A blue colour indicates that the reagent has been contaminated and should be discarded.

WARNINGS AND PRECAUTIONS

For in vitro diagnostic use

- For Professional Use Only
- Please refer to the U.S. Department of Health and Human Services (Bethesda, Md., US) publication No. (CDC) 88-8395 on laboratory safety or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Follow local guidelines for disposal of all waste material.

Caution

Each donor unit used in the preparation of human source reagent has been tested and found to be Non Reactive for HIV-1/2 Antibody, HCV Antibody and Hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they were potentially infectious.

Protocol Sheet

CanAg NSE EIA REF 420-10

Mix the components directly before use. Use shaking conditions according to the Instructions.

Step	Vial/Plate	Procedure																						
1. Prepare NSE Calibrators	<table border="1"><tr><td>CAL</td><td>NSE</td></tr></table> A, B, C, D, E	CAL	NSE	Add 0.75 mL of distilled water to each vial and mix gently. Allow to stand for at least 15 minutes. NOTE: The exact concentration of each calibrator is stated on the label. This value of the calibrators should be used for calculations.																				
CAL	NSE																							
2. Prepare Wash Solution	<table border="1"><tr><td>WASHBUF</td><td>25X</td></tr></table>	WASHBUF	25X	Dilute 50 mL of Wash Concentrate with 1200 mL of distilled water or deionized water.																				
WASHBUF	25X																							
3. Prepare Antibody Solution	<table border="1"><tr><td>CONJ</td><td>Anti-NSE</td></tr><tr><td>BIOTIN</td><td>Anti-NSE</td></tr></table>	CONJ	Anti-NSE	BIOTIN	Anti-NSE	Mix 50 μ L of Tracer, HRP Anti-NSE with 1 mL of Biotin Anti-NSE per strip: <table border="1"><thead><tr><th>No. of Strips</th><th>HRP Anti-NSE (μL)</th><th>Biotin Anti-NSE (mL)</th></tr></thead><tbody><tr><td>1</td><td>50</td><td>1</td></tr><tr><td>2</td><td>100</td><td>2</td></tr><tr><td>3</td><td>150</td><td>3</td></tr><tr><td>4</td><td>200</td><td>4</td></tr><tr><td>5</td><td>250</td><td>5</td></tr></tbody></table>	No. of Strips	HRP Anti-NSE (μ L)	Biotin Anti-NSE (mL)	1	50	1	2	100	2	3	150	3	4	200	4	5	250	5
CONJ	Anti-NSE																							
BIOTIN	Anti-NSE																							
No. of Strips	HRP Anti-NSE (μ L)	Biotin Anti-NSE (mL)																						
1	50	1																						
2	100	2																						
3	150	3																						
4	200	4																						
5	250	5																						

250
0

6
7
8
9
10
11
12

300
350
400
450
500
550
600

5
6
7
8
9
10
11
12

4. Wash	MICROPLA	Wash each well once with Wash Solution
5. Add calibrators and samples	CAL NSE A, B, C, D, E	25 μ L in each well
6. Add Antibody Solution	ANTIBODY SOLUTION	100 μ L in each well
7. Incubate	MICROPLA	1 hour shaking at room temperature
8. Wash	MICROPLA	Wash each well six times with Wash Solution
9. Add TMB HRP-Substrate	SUBS TMB	100 μ L in each well
10. Incubate	MICROPLA	30 min shaking at room temperature
11. Read absorbance	MICROPLA	620 nm
Alt.11 Add Stop Solution	STOP	100 μ L in each well
Alt.12 Incubate	MICROPLA	1 min shaking at room temperature
Alt.13 Read absorbance	MICROPLA	Read at 405 nm within 5 min

CLP (1272/2008) HAZARD CLASSIFICATION

Information about CLP (1272/2008) HAZARD CLASSIFICATION can be found at the end of this document.

SPECIMEN COLLECTION AND HANDLING

The CanAg NSE EIA is intended for use with serum. Collect blood by venipuncture and separate the serum according to common procedures. Serum should be separated from the clot within 60 minutes of collection to avoid leaking of NSE from blood cells. Do not use haemolysed samples. Plasma is not recommended since significant amounts of NSE can be released from platelets. Samples can be stored at 2–8 °C for 24 hours. For longer periods store samples at -70 °C or below. Samples should not be stored in a self-defrosting freezer and not be thawed and refrozen before analysis. Bring frozen samples to room temperature and mix THOROUGHLY by gently inverting multiple times before analysis. Samples that contain gross particulates should be centrifuged at 10.000 x g for 10 minutes, prior to use to eliminate any particulate matter that may have developed from the thawing process. Analyze thawed samples within one hour.

PROCEDURE

Materials required but not supplied with the kit

1. Microplate shaker

Shaking should be medium to vigorous, approximately 700-1100 oscillations/min.

2. Microplate wash device

Automatic plate washer capable of performing 1 and 6 washing cycles with a minimal fill volume of 350 µL/well/washcycle.

An 8-channel pipette with disposable plastic tips for delivery of 350 µL is recommended if an automatic microplate washer is not used.

3. Microplate spectrophotometer

With a wavelength of 620 nm and/or 405 nm, and an absorbance range of 0 to 3.0.

4. Precision pipettes

With disposable plastic tips for dispensing microlitre volumes. An 8-channel pipette or respenser pipette with disposable plastic tips for delivery of 100 µL is useful but not essential. Pipettes for dispensing millilitre volumes.

5. Distilled or deionized water

For reconstitution of NSE Calibrators and for preparation of diluted wash solution.

Procedural notes

1. A thorough understanding of this package insert is necessary to ensure proper use of the CanAg NSE EIA kit. The reagents supplied with the kit are intended for use as an integral unit. Do not mix identical reagents from kits having different lot numbers. Do not use the kit reagents after the expiry date printed on the outside of the kit box.
2. Reagents should be allowed to reach room temperature (20–25 °C) prior to use. The assay should only be performed at temperatures between 20–25 °C to obtain accurate results. Frozen sera must be gently but thoroughly mixed after thawing.
3. Before starting to pipette calibrators and patient specimens it is advisable to mark the strips to be able to clearly identify the samples during and after the assay.
4. The requirement for efficient and thorough washing for separation of bound and unbound antigen and reagents from the solid-phase bound antibody-antigen complexes is one of the most important steps in an EIA. In order to ensure efficient washing make sure that all wells are completely filled to the top edge with wash solution during each wash cycle, that wash solution is dispensed at a good flow rate, that the aspiration of the wells between and after the wash cycles is complete and that the wells are empty. If there is liquid left, invert the plate and tap it carefully against absorbent paper.
 - Automatic strip washer: Follow the manufacturer's instructions for cleaning and maintenance diligently and wash the required number of wash cycles prior to and after each incubation step. It's highly recommended to use *strip* process mode and *overflow* wash mode with a dispensing volume of 800 µL. The aspiration/wash device should not be left standing with the Wash Solution for long periods, as the needles may get clogged resulting in poor liquid delivery and aspiration.
5. The TMB HRP-Substrate is very sensitive for contamination. For optimal stability of the TMB HRP-Substrate, pour the required amount from the vial to a carefully cleaned reservoir or preferably a disposable plastic tray to avoid contamination of the reagent. Be sure to use clean disposable plastic pipette tips (or respenser pipette tip).
6. Be sure to use clean disposable plastic pipette tips and a proper pipetting technique when handling samples and reagents. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid. A proper pipetting technique is of particular importance when handling the TMB HRP-Substrate solution.

Preparation of reagents	Stability of prepared reagent
-------------------------	-------------------------------

NSE Calibrators

4 weeks at 2–8 °C
3 months at -20 °C

Add exactly 0.75 mL of distilled water to each vial and mix gently. Allow standing for at least 15 minutes to reconstitute. **NOTE:** The concentration of the calibrators is stated on the labels and should be used for calculation of the results.

Wash Solution

2 weeks at 2–25 °C
in a sealed container

Pour the 50 mL Wash Concentrate into a clean container and dilute 25-fold by adding 1200 mL of distilled or deionised water to give a buffered Wash Solution.

Antibody Solution

3 weeks at 2–8 °C

Prepare the required quantity of Antibody Solution by mixing 50 µL of Tracer, HRP Anti-NSE with 1 mL of Biotin Anti-NSE per strip (see table below):

No. of Strips	Tracer, HRP Anti-NSE (µL)	Biotin Anti-NSE (mL)
1	50	1
2	100	2
3	150	3
4	200	4
5	250	5
6	300	6
7	350	7
8	400	8
9	450	9
10	500	10
11	550	11
12	600	12

Be sure to use a clean plastic or glass bottle for preparation of Antibody Solution.

Alternative: Pour the content of the Tracer, HRP Anti-NSE into the vial of Biotin Anti-NSE and mix gently. Be sure that all content of the Tracer is transferred to the vial of Biotin Anti-NSE.

NOTE: The Antibody Solution is stable for 3 weeks at 2–8 °C. Do not prepare more Antibody Solution than will be used within this period and make sure that it is stored properly.

Assay procedure

Perform each determination in duplicate for both calibrators and patient samples. A calibration curve should be run with each assay. All reagents and samples must be brought to room temperature (20–25 °C) before use.

1. Start to prepare NSE Calibrators, Wash Solution and Antibody Solution. It is important to use clean containers. Follow the instructions carefully.
2. Transfer the required number of microplate strips to a strip frame. (Immediately return the remaining strips to the aluminium pouch containing a desiccant and reseal carefully). Wash each strip once with the Wash Solution. Do not wash more strips than can be handled within 30 min.
3. Pipette 25 μ L of the NSE Calibrators (CAL A, B, C, D, E) and patient specimens (unknowns-Unk) into the strip wells according to the following scheme:

	1	2	3	4	5	6	7 etc
A	Cal A	Cal E	4th Unk				
B	Cal A	Cal E	etc				
C	Cal B	1st Unk					
D	Cal B	1st Unk					
E	Cal C	2nd Unk					
F	Cal C	2nd Unk					
G	Cal D	3rd Unk					
H	Cal D	3rd Unk					

4. Add 100 μ L of Antibody Solution to each well using a 100 μ L precision pipette (or an 8-channel 100 μ L precision pipette). Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid.
5. Incubate the plate for 1 hour (\pm 10 min) at room temperature (20-25 °C) with constant shaking of the plate using a microplate shaker.
6. After the incubation aspirate and wash each strip 6 times.

7. Add 100 μL of TMB HRP-Substrate to each well using the same procedure as in item 4. The TMB HRP-Substrate should be added to the wells as quickly as possible and the time between addition to the first and last well should not exceed 5 min.
8. Incubate for 30 min (\pm 5 min) at room temperature with constant shaking. Avoid exposure to direct sunlight.
9. Immediately read the absorbance at 620 nm in a microplate spectrophotometer.

Option

If the laboratory does not have access to a microplate spectrophotometer capable of reading at 620 nm the absorbance can be determined as in item 10.

10. Add 100 μL of Stop Solution, mix and read the absorbance at 405 nm in a microplate spectrophotometer within 5 min after addition of Stop Solution.

Measurement range

The CanAg NSE EIA measures concentrations between 1 and approximately 150 $\mu\text{g/L}$. If NSE concentrations above the measuring range are to be expected, it is recommended to dilute samples with normal human serum prior to analysis.

NOTE: The serum used for dilution should also be measured in order to determine the endogenous NSE concentration (see "Calculation of results").

Quality control

CanChek Tumor Marker Control Sera Levels 1 and 2 (available separately, REF 107-20) are recommended for validation of the assay series. If values outside of the specified range are obtained, a complete check of reagents and reader performance should be made and the analysis repeated.

Reference materials

Since no common reference material is available for NSE antigen, CanAg NSE EIA Calibrator values are assigned against a set of in-house reference standards.

CALCULATION OF RESULTS

If a microplate spectrophotometer with built-in data calculation program is used refer to the manual for the spectrophotometer and create a program using the concentration stated on the label of each of the NSE calibrators.

For automatic calculation of NSE results it is recommended to use either of the following methods:

- Cubic spline curve fit method. Calibrator A should be included in the curve

with the value 0 µg/L.

- Spline smoothed curve fit method. Calibrator A should be used as plate blank.
- Interpolation with point-to-point evaluation. Calibrator A should be included in the curve with the value 0 µg/L.
- Quadratic curve fit method. Calibrator A should be included in the curve with the value 0 µg/L.


NOTE: 4-Parametric or Linear regression evaluation methods should not be used.

For manual evaluation, a calibration curve is constructed by plotting the absorbance (A) values obtained for each NSE Calibrator against the corresponding NSE concentration (in µg/L), see figure. The unknown NSE concentrations can then be read from the calibration curve using the mean absorbance value of each patient specimen. If samples in an initial analysis give NSE levels above the concentration of calibrator E, it is necessary to dilute the sample 1/10 with normal human serum in order to obtain accurate results. The result should then be calculated according to the following procedure:

$$\text{Dilution 1/10: } 10 \times ([\text{NSE}]_{\text{Diluted sample}} - (0.9 \times [\text{NSE}]_{\text{Normal human serum}}))$$

Example of results

Specimen	Calibrator values	Mean abs value (A)	NSE µg/L
CAL NSE A	0 µg/L	0.037	
CAL NSE B	7.5 µg/L	0.238	
CAL NSE C	22.9 µg/L	0.663	
CAL NSE D	68.4 µg/L	1.688	
CAL NSE E	138.0 µg/L	2.720	
Specimen 1		0.518	17.5
Specimen 2		1.474	57.8

Example, do not use this curve to determine assay results.

The exact NSE concentration is indicated on the label of each calibrator vial.

LIMITATIONS OF THE PROCEDURE

The level of NSE cannot be used as absolute evidence for the presence or absence of malignant disease and the NSE test should not be used in cancer screening. The results of the test should be interpreted only in conjunction with other investigations and procedures in the diagnosis of disease and the NSE test should not replace any established clinical examination.

Elevated NSE values not due to tumours may occur in dialysis patients and patients with leukaemic diseases.

Serum should not contain visible haemolysis (the absorbance at 500 nm for non-turbid sample should not exceed 0.3) since erythrocytes contain significant amounts of NSE (7). Prolonged storage of whole blood can cause release of NSE from the blood cells.

Anti-reagent antibodies (human anti-mouse antibody (HAMA) or heterophilic antibodies) in the patient sample may occasionally interfere with the assay, even though specific blocking agents are included in the buffer.

Biotin may interfere with the assay giving false low results. This should be taken into consideration for patients taking dietary supplements or receiving therapy containing high (> 5mg/day) or extremely high (300 mg/day) biotin doses. Peak serum levels have been reported to occur 1-3h post ingestion and a physiological half-life of 8-16h depending on renal function (10,11). In a study by Grimsey et al. (10) a specimen concentration of 30 ng/mL was reached 8h following intake of 10 mg of biotin. For patients on very high doses of biotin it is recommended to stop taking biotin for at least 2 days before blood draw (12).

EXPECTED VALUES

CanAg NSE EIA was used to measure NSE antigen levels in specimens from 495 apparently healthy blood donors. In this study 97.5 % of the subjects had a NSE concentration at or below 10.5 µg/L and 95 % of the subjects had a NSE concentration at or below 9.9 µg /L. Median serum level was 6.5 µg/L.

It is recommended that each laboratory establish its own reference range with the population of interest and sample collection procedures used within the laboratory.

PERFORMANCE CHARACTERISTICS

Precision

Total precision was determined according to NCCLS guideline EP5-A (8) using four levels of frozen pooled human serum containing added NSE. Each sample was randomly pipetted in duplicates and analysed twice each day over 20 days. The analyses were undertaken during a period of 40 months, by \geq three different technicians and using 20 different CanAg NSE EIA kit batches.

Sample	Replicates	Mean µg/L	Within-run SD (µg/L)	Within-run CV %	Between-day SD (µg/L)	Between-day CV %
NSE 1	80	10.3	0.24	2.3	0.57	5.5
NSE 2	80	23.7	0.82	3.5	0.97	4.1
NSE 3	80	48.2	1.02	2.1	1.93	4.0
NSE 4	80	92.7	1.60	1.7	3.44	3.7

Detection limit

The detection limit of the CanAg NSE EIA assay is < 1 µg/L defined as the concentration corresponding to the mean of the absorbance values for the NSE Calibrator A plus 2 standard deviations according to the formula:

$$\frac{2 \times \text{SD CAL A}}{\text{OD CAL B} - \text{OD CAL A}} \times [\text{CAL B}] \mu\text{g/L}$$

Hook effect

No hook effect has been noticed for NSE concentrations up to 200 000 µg/L.

Linearity

Patient samples were diluted with normal serum and analysed. The obtained values were in the range 93–101 % of the expected values.

Specificity

The monoclonal antibodies used are specific for the γ -subunit of enolase. No measurable cross-reactions with other enolase have been observed.

The NCCLS guideline EP7-P (9) was followed to determine possible sources of interference. The following substances and concentrations were tested and found not to interfere with the test.

	Concentration with no significant (\pm 10 %) interference
Lipemia (Intralipid®)	10 mg/mL
Bilirubin, unconjugated	0.6 mg/mL

Biotin interference

A study was conducted to evaluate biotin interference. Low and high serum control samples were spiked to final biotin concentrations of 15, 30, 60, and 600 ng/mL. The mean NSE concentration was determined for each sample and the percent recovery for each biotin concentration was calculated using the formula: Recovery (%) = 100x (Mean NSE concentration w. biotin added/Mean NSE concentration w. diluent only added).

NSE analyte level	Biotin test conc. (ng/mL)	Expected NSE conc. (μ g/L)	Observed NSE conc. (μ g/L)	Recovery (%)
Low	15	13,3	13,1	98
High 1	15	75,7	74,3	98
High 2	15	59,5	58,8	99
Low	30	13,3	13,1	98
High 1	30	75,5	73,3	97
High 2	30	59,7	57,1	96
Low	60	13,2	12,5	95
High 1	60	75,7	68,0	90
High 2	60	60,5	58,4	97
Low	600	13,3	1,75	13
High 1	600	74,8	10,3	14
High 2	600	62,8	4,97	8

Based on linear regression analysis, the lowest concentration of biotin found to influence test results (\geq 10%) was 78 ng/mL.

CLP (1272/2008) HAZARD CLASSIFICATION

The following warnings and precautions apply to

SUBS

TMB

Hazard pictograms:

Signal word: Danger

Hazard Statement: Repr. 1B: H360D May damage the unborn child.

Prevention statement: P202 Do not handle until all safety precautions have been read and understood.

Prevention: P280 Wear protective gloves / protective clothing / eye protection / face protection.

Precautionary statement response: P308+P313 IF exposed or concerned get medical advice/attention.

Precautionary statement disposal: P501 Dispose of contents / container to an approved hazardous / special waste disposal facility in accordance with local and national regulations.

Restricted to professional users.

Hazardous substances: 2- Pyrrolidone

Other hazards

None of the mixtures in the kit contains any substances considered to meet the criteria classifying them as PBT and/or vPvB.

WARRANTY

The performance data presented here were obtained using the assay procedure indicated. Any change or modification of the procedure not recommended by Fujirebio Diagnostics may affect the results, in which event Fujirebio Diagnostics disclaims all warranties expressed, implied or statutory including the implied warranty of merchantability and fitness for use.

REFERENCES

1. Paus E. and Nustad K., (1989) Immunoradiometric Assay for $\alpha\gamma$ - and $\gamma\gamma$ -Enolase (Neuron-Specific Enolase), with Use of Monoclonal Antibodies and Magnetizable Polymer Particles. *Clin. Chem.* 35: 2034-2038.
2. Dahlén U., Karlsson B., Nilsson O. and Uhl W., (1995) Development of an Enzyme Immunoassay, NSE-Enzymun Test For Determination of Neuron-Specific Enolase. XXIII International Society for Oncodevelopmental Biology and Medicine, Montréal, Québec .
3. Cooper E.H., (1994) Neuron-specific enolase. *The International Journal of Biological Markers* 9(4):205-10.
4. Cooper E.H., Pritchard J., Bailey C.C. and Ninane J., (1987) Serum neuron-specific enolase in children's cancer. *Br. J. Cancer* 56: 65–67.
5. Schneider, P. M. et al., (2002) Lung Cancer. In "Tumor markers, Physiology, Pathobiology, Technology and Clinical Applications" Eds. Diamandis E. P. et al., AACCC Press, Washington pp 287-303.
6. Bonner J. A., Sloan JA., Rowland KM., Klee GG., Kugler JW., Mailliard JA., Wiesenfeld M., Krook JE., Maksymiuk AW., Shaw EG., Marks RS and Perez EA., (2000) Significance of Neuron-specific Enolase Levels before and during Therapy for Small Cell Lung Cancer. *Clinical Cancer Research* 6: 597-601.
7. Pählman S., Esscher T., Bergvall P. and Odelstad L., (1984) Purification and characterization of human neuron-specific enolase: Radioimmunoassay development. *Tumor Biol.* 5: 127–139.
8. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
9. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).
10. Grimsey P et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. *J Pharmacokinet.* 2017; 2(4); 247-256.
11. Jenkins Colon P, Greene D.N. Biotin Interference in Clinical Immunoassays. *JALM* 2018; 2(6); 941-951.
12. Chun KY. Biotin Interference in Diagnostic Tests. *Clinical Chemistry* 2017; 63; 619-620.

CanAg[®] is a registered trademark of Fujirebio Diagnostics AB

Fujirebio Diagnostics AB
Elof Lindälvs gata 13
SE-414 58 Göteborg
Sweden
Phone + 46 31 85 70 30
Fax + 46 31 85 70 40
info@fdab.com
www.fdab.com