IVD

CE

Free 25OH Vitamin D ELISA

KAPF1991

Version : 240329

Date of issue : 29/03/2024

Revision date: 29/03/2024

History

Summary of change:

Previous Version: 230123	Current Version: 240329																																												
<p>V. REAGENTS PROVIDED</p> <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 20px;"> <tr> <td style="width: 20%;">[WASH TABLET] PBS-Tween wash buffer tablet</td> <td style="width: 15%; text-align: center;">2 tablets Lyophilized</td> <td style="width: 10%; text-align: center;">Brown</td> <td style="width: 55%; text-align: center;">Dissolve 1 tablet in 500 ml distilled water</td> </tr> </table>	[WASH TABLET] PBS-Tween wash buffer tablet	2 tablets Lyophilized	Brown	Dissolve 1 tablet in 500 ml distilled water	<p>V. REAGENTS PROVIDED</p> <p>Removal of the column relative to the color code</p> <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 20px;"> <tr> <td style="width: 20%;">[WASH SOLN CONC] Wash solution (TRIS-HCl)</td> <td style="width: 15%; text-align: center;">1 vial 10 ml</td> <td style="width: 10%; text-align: center;">Brown</td> <td style="width: 55%; text-align: center;">Dilute 200x with distilled water. Use a magnetic stirrer.</td> </tr> </table>	[WASH SOLN CONC] Wash solution (TRIS-HCl)	1 vial 10 ml	Brown	Dilute 200x with distilled water. Use a magnetic stirrer.																																				
[WASH TABLET] PBS-Tween wash buffer tablet	2 tablets Lyophilized	Brown	Dissolve 1 tablet in 500 ml distilled water																																										
[WASH SOLN CONC] Wash solution (TRIS-HCl)	1 vial 10 ml	Brown	Dilute 200x with distilled water. Use a magnetic stirrer.																																										
<p>VI. SUPPLIES NOT PROVIDED</p>	<p>VI. SUPPLIES NOT PROVIDED</p> <p>Addition of a magnetic agitator</p>																																												
<p>VII. REAGENT PREPARATION</p> <p>C. Working Wash Solution: according to the number of strips used, dissolve 1 wash buffer tablet in 500 ml distilled water, or 2 wash buffer tablets in 1L distilled water. Ensure that all the salt crystals are dissolved.</p>	<p>VII. REAGENT PREPARATION</p> <p>C. Working Wash Solution: Prepare an adequate volume of working wash solution by adding 199 volumes of distilled water to one volume of 200x concentrated Wash solution. Use a magnetic stirrer to homogenize. Dispose of unused working wash solution at the end of the day.</p>																																												
<p>VIII. STORAGE AND EXPIRATION DATING OF REAGENTS</p> <p>After reconstitution, wash buffer is stable for two weeks at room temperature</p> <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 10px;"> <thead> <tr> <th style="width: 30%;">Kit component</th> <th style="width: 70%;">In use stability</th> </tr> </thead> <tbody> <tr><td>Anti-Vitamin D coated plate</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>Sample diluent</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>100x Biot-Vit D reagent</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>Streptavidin-HRP reagent</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>Calibrator 1 - 6</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>Control 1-2</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>PBS-Tween wash buffer (reconstituted tablet)</td><td>Maximum of 2 weeks at room temperature</td></tr> <tr><td>Chromogenic solution TMB</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>Stop reagent</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>Biot-Vit D dilution buffer</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> </tbody> </table>	Kit component	In use stability	Anti-Vitamin D coated plate	Maximum of 2 weeks at 2-8 °C	Sample diluent	Maximum of 2 weeks at 2-8 °C	100x Biot-Vit D reagent	Maximum of 2 weeks at 2-8 °C	Streptavidin-HRP reagent	Maximum of 2 weeks at 2-8 °C	Calibrator 1 - 6	Maximum of 2 weeks at 2-8 °C	Control 1-2	Maximum of 2 weeks at 2-8 °C	PBS-Tween wash buffer (reconstituted tablet)	Maximum of 2 weeks at room temperature	Chromogenic solution TMB	Maximum of 2 weeks at 2-8 °C	Stop reagent	Maximum of 2 weeks at 2-8 °C	Biot-Vit D dilution buffer	Maximum of 2 weeks at 2-8 °C	<p>VIII. STORAGE AND EXPIRATION DATING OF REAGENTS</p> <p>Unused working wash solution should be disposed of at the end of the day.</p> <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 10px;"> <thead> <tr> <th style="width: 30%;">Kit component</th> <th style="width: 70%;">In use stability</th> </tr> </thead> <tbody> <tr><td>Anti-Vitamin D coated plate</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>Sample diluent</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>100x Biot-Vit D reagent</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>Streptavidin-HRP reagent</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>Calibrator 1 - 6</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>Control 1-2</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>Working wash solution</td><td>1 day at room temperature.</td></tr> <tr><td>Chromogenic solution TMB</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>Stop reagent</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> <tr><td>Biot-Vit D dilution buffer</td><td>Maximum of 2 weeks at 2-8 °C</td></tr> </tbody> </table>	Kit component	In use stability	Anti-Vitamin D coated plate	Maximum of 2 weeks at 2-8 °C	Sample diluent	Maximum of 2 weeks at 2-8 °C	100x Biot-Vit D reagent	Maximum of 2 weeks at 2-8 °C	Streptavidin-HRP reagent	Maximum of 2 weeks at 2-8 °C	Calibrator 1 - 6	Maximum of 2 weeks at 2-8 °C	Control 1-2	Maximum of 2 weeks at 2-8 °C	Working wash solution	1 day at room temperature.	Chromogenic solution TMB	Maximum of 2 weeks at 2-8 °C	Stop reagent	Maximum of 2 weeks at 2-8 °C	Biot-Vit D dilution buffer	Maximum of 2 weeks at 2-8 °C
Kit component	In use stability																																												
Anti-Vitamin D coated plate	Maximum of 2 weeks at 2-8 °C																																												
Sample diluent	Maximum of 2 weeks at 2-8 °C																																												
100x Biot-Vit D reagent	Maximum of 2 weeks at 2-8 °C																																												
Streptavidin-HRP reagent	Maximum of 2 weeks at 2-8 °C																																												
Calibrator 1 - 6	Maximum of 2 weeks at 2-8 °C																																												
Control 1-2	Maximum of 2 weeks at 2-8 °C																																												
PBS-Tween wash buffer (reconstituted tablet)	Maximum of 2 weeks at room temperature																																												
Chromogenic solution TMB	Maximum of 2 weeks at 2-8 °C																																												
Stop reagent	Maximum of 2 weeks at 2-8 °C																																												
Biot-Vit D dilution buffer	Maximum of 2 weeks at 2-8 °C																																												
Kit component	In use stability																																												
Anti-Vitamin D coated plate	Maximum of 2 weeks at 2-8 °C																																												
Sample diluent	Maximum of 2 weeks at 2-8 °C																																												
100x Biot-Vit D reagent	Maximum of 2 weeks at 2-8 °C																																												
Streptavidin-HRP reagent	Maximum of 2 weeks at 2-8 °C																																												
Calibrator 1 - 6	Maximum of 2 weeks at 2-8 °C																																												
Control 1-2	Maximum of 2 weeks at 2-8 °C																																												
Working wash solution	1 day at room temperature.																																												
Chromogenic solution TMB	Maximum of 2 weeks at 2-8 °C																																												
Stop reagent	Maximum of 2 weeks at 2-8 °C																																												
Biot-Vit D dilution buffer	Maximum of 2 weeks at 2-8 °C																																												
<p>XIII. PERFORMANCE AND LIMITATIONS</p> <p>A. Limits of detection</p> <p>The Limit of Blank (LoB) and the Limit of Detection (LoD) were determined in accordance with the CLSI guideline EP17-A2. The LoB was calculated to be 1.5 pg/ml. The LoD was calculated to be 2.4 pg/ml.</p>	<p>XIII. PERFORMANCE AND LIMITATIONS</p> <p>A. Limits of detection</p> <p>The Limit of Blank (LoB) and the Limit of Detection (LoD) were determined by converting mOD into pg/ml values with standard curve equation. The LoB was calculated to be 1.86 pg/ml. The LoD was calculated to be 2.38 pg/ml.</p>																																												

C. Precision

Intermediate Precision and repeatability was determined based on CLSI EP05-A3.

Sample	N	pg/ml		Repeatability (within run)	Intermediate Precision(Total)
Pool 1	80	5.8	SD	0.29	0.34
			CV	4.9%	5.9%
Pool 2	80	9.6	SD	0.53	0.59
			CV	5.5%	6.1%
Pool 3	80	18.4	SD	0.35	0.75
			CV	1.9%	4.0%
Pool 4	80	28.1	SD	1.26	1.76
			CV	4.5%	6.3%

SD: Standard Deviation; CV: Coefficient of Variation

D. Time delay

The time delay test between the last Calibrator and sample dispensing results is shown in the following table.

TIME DELAY				
	0 min (pg/ml)	5 min (pg/ml)	10 min (pg/ml)	15 min (pg/ml)
Sample 1	4.9	4.9	4.9	4.8
Sample 2	18.6	18.4	18.6	18.6

Assay results remain accurate even when a sample is dispensed 15 minutes after the Calibrator has been added in the coated wells.

G. Validated concentration interval

The validated concentration interval of the Free 25OH Vitamin D assay is 2.4 – 35 pg/mL Free 25OH Vitamin D.

C. Precision

Reproducibility and repeatability precisions were evaluated with two samples at medium and high concentrations. Their concentrations were measured in 10 different runs in duplicate for reproducibility and in a single run for repeatability precision.

Sample	N	pg/ml		Reproducibility Precision (Total runs)
Sample 1	20	8.22	SD	0.45
			CV	5.5%
Sample 2	20	29.42	SD	3.00
			CV	10.2%

SD: Standard Deviation; CV: Coefficient of Variation

Sample	N	pg/ml		Repeatability Precision (within a run)
Sample 1	12	9.02	SD	0.32
			CV	3.5%
Sample 2	12	32.38	SD	0.64
			CV	2.0%

SD: Standard Deviation; CV: Coefficient of Variation

D. Time delay

The time delay test between the last Calibrator and sample dispensing results is shown in the following table.

TIME DELAY					
	0 min (pg/ml)	5 min (pg/ml)	10 min (pg/ml)	15 min (pg/ml)	20 min (pg/ml)
Sample 1	5.60	5.68	5.70	5.46	5.71
Sample 2	13.17	14.70	13.07	14.27	13.38

Assay results remain accurate even when a sample is dispensed 20 minutes after the Calibrator has been added in the coated wells.

G. Validated concentration interval

The validated concentration interval of the Free 25OH Vitamin D assay is 2.38 – 34.4 pg/mL Free 25OH Vitamin D.

Read entire protocol before use.

Free 25OH Vitamin D ELISA

I. INTENDED USE

The DIAsource Free 25OH Vitamin D ELISA is an in vitro diagnostic medical device intended to be used by healthcare professionals for the quantitative measurement of free 25OH Vitamin D in human serum and can be used for screening and monitoring purposes.

II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource Free 25OH Vitamin D ELISA Kit
- B. **Catalog number :** KAPF1991 : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2 B-1348 Louvain-la Neuve, Belgium.

For technical assistance or ordering information contact :

Tel : +32 (0)10 84 99 11

Fax : +32 (0)10 84 99 90

III. CLINICAL BACKGROUND

Due to its hydrophobic nature 25OH Vitamin D, and other Vitamin D metabolites, circulate on binding proteins. About 90% of the total circulating 25OH Vitamin D is bound to the so-called VDBP or DBP. The remaining 10% are bound to albumin, the main protein of human blood plasma. Although the affinity of albumin toward 25OH Vitamin D is much lower than the affinity of VDBP, the high concentration of albumin compensates for this difference. A tiny fraction representing 0.04% of the total 25OH Vitamin D concentration circulates as the free form.

The conversion of 25OH Vitamin D into the biologically active 1,25(OH)₂ Vitamin D takes place into the cells and so requires the internalization of 25OH Vitamin D from the extracellular fluid. Different transport mechanisms are likely to be involved and some of them involve the concentration of the free ligand as one of the important parameters. In these cases, the fraction of free 25OH Vitamin D relates to the biological activity of Vitamin D, and therefore may better reflect the physiological action of Vitamin D than the total concentration of 25OH Vitamin D.

The fraction of free 25OH Vitamin D represents about 0.04% of the total concentration of 25OH Vitamin D. However, this percentage is not constant and varies according to different conditions. Although the level of albumin tends to be stable among individuals, the concentration of VDBP can fluctuate in several conditions, therefore, influencing the fraction of free 25OH Vitamin D. Pregnancy leads to an increase of the VDBP levels by about 50% while, e.g., liver failure and chronic kidney disease both result in a decrease of the VDBP concentration by also about 50%. In case of elevated binding proteins concentration, the fraction of free 25OH Vitamin D is lower, and vice-versa in the case of low-binding proteins concentration. In these conditions, and in the conditions listed later, the measurement of free 25OH Vitamin D may be a better marker of the Vitamin D activity than the classical measurement of total 25OH Vitamin D. In addition to variable concentrations VDBP also exists as different polymorphic forms. The affinity of the different VDBP forms toward 25OH Vitamin D may vary although this is still under debate. A polymorphic form with a high affinity for 25OH Vitamin D will decrease the fraction of free 25OH Vitamin D available. Conversely, a low affinity VDBP form will result in higher levels of free 25OH Vitamin D.

IV. PRINCIPLES OF THE METHOD

The Free 25OH Vitamin D ELISA is based on a two-step immunoassay procedure performed in a microtiter plate. During the first incubation step, free 25OH Vitamin D [25(OH)Vitamin D2 and D3] is bound to the anti-Vitamin D antibody coated on the well of the microtiter plate. The in vivo equilibrium between free and bound 25OH Vitamin D is minimally disturbed. After washing, a fixed amount of biotinylated 25OH Vitamin D is added to each well. The non-bound biotinylated 25OH Vitamin D is removed by washing and a streptavidin peroxidase conjugate is added. In a next step TMB chromogenic substrate is added. Finally, the reaction is stopped by addition of stop reagent and the absorbance [A450nm] is measured using a plate spectrophotometer. The concentration of free 25OH Vitamin D (pg/ml) in the sample is inversely proportional to the absorbance in each sample well.

V. REAGENTS PROVIDED

Reagents	96 Tests Kit	Reconstitution
 Anti-Vitamin D coated plate Microtiter plate coated with a mouse anti-25OH D2/D3 monoclonal antibody.	96 wells	Ready for use
[DIL SPE] Sample diluent: Specimen diluent, containing a fluorosurfactant and Proclin.	1 vial 14 ml	Ready for use
[Ag BIOT CONC] 100x Biot-Vit D reagent Biotinylated 25OH Vitamin D in preservation buffer.	1 vial 250 µl	Dilute 100 x with Biot-Vit D dilution buffer
[SAV HRP] Streptavidin-HRP reagent Peroxidase conjugated streptavidin, containing Proclin.	1 vial 14 ml	Ready for use
[CAL N] Calibrators ; N = 1 to 6 Lyophilized 25OH Vitamin D depleted human serum containing Proclin and BND. (see exact value on vial label)	6 vials Lyophilized	Add 250 µl distilled water
[CONTROL 1] Control 1 Normal human serum containing Proclin and BND. (see exact value on vial label)	1 vial Lyophilized	Add 250 µl distilled water
[CONTROL 2] Control 2 25OH Vitamin D spiked human serum containing Proclin and BND (see exact value on vial label)	1 vial Lyophilized	Add 250 µl distilled water
[WASH SOLN CONC] Wash solution (TRIS-HCl)	1 vial 10 ml	Dilute 200x with distilled water Use a magnetic stirrer
[CHROM TMB] Chromogenic solution TMB (Tetramethylbenzidine)	1 vial 12 ml	Ready for use
[STOP SOLN] Stop reagent 1M HCl.	1 vial 12 ml	Ready for use
[Ag BIOT SOLN] Biot-Vit D dilution buffer Buffer containing Proclin and BND.	1 vial 14.5 ml	Ready for use

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- Distilled water
- Pipettes calibrated for delivery of 10 µl, 35 µl - 250 µl, 3500 µl - 5000 µl (the use of accurate pipettes with disposable plastic tips is recommended).
- Vortex mixer
- Magnetic agitator
- Plate shaker incubator at 37 °C and 650 rpm
- Washer for microplates
- Microtiter plate reader capable of reading at 450 nm and 650 (bichromatic reading)

VII. REAGENT PREPARATION

Bring all reagents to room temperature (18-25°C) at least 30 minutes before use.

- Calibrators:** reconstitute the Calibrators with 250 µl distilled water and let them stand for 15 minutes at RT (18-25°C). After reconstitution, mix well and make sure that all lyophilized material has been reconstituted.
- Controls:** reconstitute the Controls with 250 µl distilled water. Let them stand for 15 minutes at RT (18-25°C). After reconstitution, mix well and make sure that all lyophilized material has been reconstituted.
- Working Wash Solution:** Prepare an adequate volume of working wash solution by adding 199 volumes of distilled water to one volume of 200x concentrated Wash solution.
Use a magnetic stirrer to homogenize. Dispose of unused working wash solution at the end of the day.
- Working Biot-Vit D solution:** prepare an adequate volume of working dilution of the Biot-Vit D reagent by mixing 100x concentrated Biot-Vit D reagent with Biot-Vit D dilution buffer according to the number of strips used, as indicated in the table below: for example, for 6 strips (48 wells): 60 µl 100x concentrated Biot-Vit D reagent is added to 6.0 ml Biot-Vit D dilution buffer.
Use an appropriate polypropylene container for preparation.
The preparation of working Biot-Vit D solution is not stable and must be discarded if not used.

No. Strips	100x Biot-Vit D reagent (µl)	Biot-Vit D dilution buffer (ml)
3	35	3.5
4	40	4.0
5	50	5.0
6	60	6.0
7	70	7.0
8	80	8.0
9	90	9.0
10	100	10.0
11	110	11.0
12	120	12.0

All other reagents provided are ready to use.

Do not use a plate seal during the incubation steps.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- Protect the Streptavidin-HRP reagent and chromogenic solution TMB from light.
- Do not store diluted Biotin-Vit D reagent, dilute only the required amount.
- Once opened, reconstitute the foil bag containing the microtiter plate with desiccant (stable for a maximum of 2 weeks at 2-8 °C).
- After reconstitution, calibrators and controls are stable for two weeks at 2-8°C.
- Unused working wash solution should be disposed of at the end of the day.

Kit component	In use stability
Anti-Vitamin D coated plate	Maximum of 2 weeks at 2-8 °C
Sample diluent	Maximum of 2 weeks at 2-8 °C
100x Biot-Vit D reagent	Maximum of 2 weeks at 2-8 °C
Streptavidin-HRP reagent	Maximum of 2 weeks at 2-8 °C
Calibrator 1 - 6	Maximum of 2 weeks at 2-8 °C
Control 1-2	Maximum of 2 weeks at 2-8 °C
Working wash solution	1 day at room temperature.
Chromogenic solution TMB	Maximum of 2 weeks at 2-8 °C
Stop reagent	Maximum of 2 weeks at 2-8 °C
Biot-Vit D dilution buffer	Maximum of 2 weeks at 2-8 °C

IX. SPECIMEN COLLECTION AND PREPARATION

- This kit is suitable for serum samples.
- Serum samples must be kept at 2 – 8°C.
- If the test is not run within 24h, sampling and storage at -20°C is recommended.
- Avoid repeated freeze-thawing cycles. Mix samples well after thawing and before testing.

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date.
Do not mix materials from different kit lots.
Bring all the reagents to room temperature (18-25°C) prior to use.
Thoroughly mix all reagents and samples by gentle agitation or swirling.
Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.
Use a clean plastic container to prepare the Working Wash Solution.
In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
For the dispensing of the Chromogenic Solution and the Stop Solution avoid pipettes with metal parts.
High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times.
To avoid drift, the time between pipetting of the first calibrator and the last sample must be limited to the time mentioned in section XIII paragraph D (Time delay).
Prepare a calibration curve for each run, do not use data from previous runs.
Dispense the Chromogenic Solution within 15 minutes following the washing of the microtiter plate.
During incubation with Chromogenic Solution, avoid direct sunlight on the microtiter plate.
Each well can only be used once.

B. Procedure

1. Select the required number of strips for the run. The unused strips should be resealed in the bag with a desiccant and stored at 2 – 8°C.
2. Secure the strips into the holding frame.
3. Pipette 90 µl of Sample diluent into all the wells.
4. Pipette 10 µl of each reconstituted Calibrator, Control or sample in duplicate into the appropriate wells (use a new pipette tip for each Calibrator, Control or sample).
5. Incubate for 90 minutes at 37 °C, shaking at 650 rpm.
6. Wash the plate 3 times with 350 µl working wash solution.
7. Pipette 100 µl of working Biot-Vit D solution into all wells.
8. Incubate for 30 minutes at 37 °C, shaking at 650 rpm.
9. Wash the plate 3 times with 350 µl working wash solution.
10. Pipette 100 µl of Streptavidin-HRP reagent into all the wells.
11. Incubate for 20 minutes at 37 °C, shaking at 650 rpm.
12. Wash the plate 3 times with 350 µl working wash solution.
13. Pipette 100 µl of the chromogenic solution into all the wells.
14. Incubate for 15 minutes at room temperature (18-25°C), stationary and protected from light.
15. Pipette 100 µl of the Stop reagent into all the well.
16. Read the absorbances at 450 nm within 5 minutes (reference filter 630nm or 650nm), and calculate the results as described in Section XI.

XI. CALCULATION OF RESULTS

1. Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).
2. Calculate the mean of duplicate determinations.
3. We recommend the use of computer assisted methods to construct the calibration curve. 4-parameter logistic function curve fitting is the preferred method. Reject obvious outliers.
4. By interpolation of the sample OD values, determine the Free 25OH Vitamin D concentrations of the samples from the calibration curve.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

Free 25OH Vitamin D ELISA		OD units
Calibrator	0.9 pg/ml	2.16
	3.1 pg/ml	1.71
	6.5 pg/ml	1.21
	11.6 pg/ml	0.83
	23.2 pg/ml	0.45
	40.3 pg/ml	0.23

Note : 1 pg/ml = 2.5 pmol/l

XIII. PERFORMANCE AND LIMITATIONS

A. Limits of detection

The Limit of Blank (LoB) and the Limit of Detection (LoD) were determined by converting mOD into pg/ml values with standard curve equation.

The LoB was calculated to be 1.86 pg/ml.

The LoD was calculated to be 2.38 pg/ml.

B. Specificity

Cross-reactivity of the antibody used in the Free 25OH Vitamin D ELISA was determined by the supplier as depicted in the following table:

Compound and Concentration	% Cross reaction
25OH-Vitamin D ₃ at 10 ng/ml	100
25OH-Vitamin D ₂ at 10 ng/ml	86
1, 25 (OH) ₂ -Vitamin D ₃ at 200 ng/ml	20
1, 25 (OH) ₂ -Vitamin D ₂ at 690 ng/ml	1.9
Vitamin D ₃ at 200 ng/ml	2.9
Vitamin D ₂ at 200 ng/ml	1.3
24, 25 (OH) ₂ -Vitamin D ₃ at 20 ng/ml	>100
25, 26 (OH) ₂ -Vitamin D ₃ at 4 ng/ml	>100
3 epi 25OH-Vitamin D ₃ at 20 µg/ml	0.1

The effect of potential interfering substances on samples tested using the DIAsource Free 25OH Vitamin D ELISA kit was evaluated according to CLSI EP07-A2. Ascorbic acid, unconjugated bilirubin, HAMA, rheumatoid factor, hemoglobin, triglycerides, biotin and cholesterol were tested in samples with different free 25OH Vitamin D concentrations. The tested substances did not affect the performance of free 25OH Vitamin D ELISA.

Substances	Concentration of interferent	Free 25OH Vitamin D (pg/mL)	Mean % Interference
Ascorbic Acid	[3mg/dl]	6.5	1%
		11.1	-1%
		16.7	0%
Unconjugated Bilirubin	[20mg/dl]	6.6	-6%
		11.2	1%
		15.7	-3%
HAMA	[2µg/ml]	6.6	-5%
		10.5	2%
		15.1	2%
Rheumatoid Factor	[600IU/ml]	6.6	-1%
		10.5	-3%
		15.1	-6%
Haemoglobin	[200mg/dl]	6.9	-8%
		11.2	-8%
		17.0	-10%
Triglycerides	[37mmol/l]	6.8	-5%
		11.1	-5%
		16.2	-5%
Biotin	[4mg/dl]	6.5	1%
		10.6	4%
		15.7	1%
Cholesterol	[13mmol/l]	6.4	0%
		10.6	-1%
		16.8	-3%

C. Precision

Reproducibility and repeatability precisions were evaluated with two samples at medium and high concentrations. Their concentrations were measured in 10 different runs in duplicate for reproducibility and in a single run for repeatability precision.

Sample	N	pg/ml		Reproducibility Precision (Total runs)
Sample 1	20	8.22	SD	0.45
			CV	5.5%
Sample 2	20	29.42	SD	3.00
			CV	10.2%

SD: Standard Deviation; CV: Coefficient of Variation

Sample	N	pg/ml		Repeatability Precision (within a run)
Sample 1	12	9.02	SD	0.32
			CV	3.5%
Sample 2	12	32.38	SD	0.64
			CV	2.0%

SD: Standard Deviation; CV: Coefficient of Variation

D. Time delay

The time delay test between the last Calibrator and sample dispensing results is shown in the following table.

TIME DELAY					
	0 min (pg/ml)	5 min (pg/ml)	10 min (pg/ml)	15 min (pg/ml)	20 min (pg/ml)
Sample 1	5.60	5.68	5.70	5.46	5.71
Sample 2	13.17	14.70	13.07	14.27	13.38

Assay results remain accurate even when a sample is dispensed 20 minutes after the Calibrator has been added in the coated wells.

E. Limitations of the test

- The test is an aid in the diagnosis and is to be used in conjunction with clinical findings.
- The performance of this assay has not been established in a pediatric population.

F. Method comparison

The correlation between rate dialysis analysis result vs Free 25OH Vitamin D ELISA is 0.9916.

G. Validated concentration interval

The validated concentration interval of the Free 25OH Vitamin D assay is 2.38 – 34.4 pg/mL Free 25OH Vitamin D.

XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls which contain azide will interfere with the enzymatic reaction and cannot be used.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises
- It is recommended that Controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.
- It is good practise to check visually the curve fit selected by the computer.

XV. EXPECTED VALUES

Dietary intake, race, season and age are known to affect the normal levels of 25OH Vitamin D.

Each laboratory should establish its own range based on their local population.

A reference range has been established based upon 109 apparently healthy individuals. The individual patient serum samples used were obtained in the

Netherlands from a local blood bank with informed consent. Samples were collected during summer. The following table provides a summary of the results:

	Free vitamin D concentration (pg/ml)
Highest concentration observed	17.1
Lowest concentration observed	<2.4
Median concentration observed	5.1
Average concentration observed	5.4

Classification of Vitamin D status has been established based on calculations performed by DIAsource, from the commonly accepted cut-off values for total 25OH Vitamin D, and from the % of free 25OH Vitamin D, and the linear correlation between both parameters, established in a large normal population.

These cut-offs are to be used by laboratories under their own responsibilities, as they have not been established by international societies yet. That said, there is also no consensus on the total 25OH Vitamin D cut-off values either.

N = 279 (normal population)

First calculation method, % Free 25OH Vitamin D = 0.020%

Second calculation method, Free 25OH Vitamin D = 0.138 Total 25OH Vitamin D + 1.4 pg/ml

Mean of the two calculation methods:

Level	Cut-off used for Total 25OH Vitamin D (ng/ml)	Proposed cut-off for Free 25OH Vitamin D (pg/ml)
Deficiency	<10	<2.4
Insufficiency	10-29	2.4-5.8
Sufficiency	30-100	>5.8

XVI. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum specimens should be in accordance with local safety procedures. All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with all reagents, Stop Solution contains HCl. In case of contact, wash thoroughly with water.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.


For more information, see Material Safety Data Sheet (MSDS)

XVII. REFERENCES

- BOUILLON R. (2016), Free 25-hydroxyvitamin D: impact of vitamin D binding protein assays on racial genotypic associations, JCEM, <http://dx.doi.org/10.1210/jc.2016-1104>.
- SOLLID S.T. (2016), Effects of vitamin D binding protein phenotypes and vitamin D supplementation on serum total 25(OH) D and directly measured free 25(OH)D, Eur. J. Endocrinol. April 1, 174:445-452.
- TANGPRICHA V. (2015), Free 25-Hydroxyvitamin D Concentrations in Cystic Fibrosis, Am. J. Med. Sci. 2015 Nov;350(5):374-9.
- ALOIA J. (2015), Free 25(OH)D and the Vitamin D Paradox in African Americans, J. Clin. Endocrinol. Metab. 2015 Jul 10:JC20152066.
- SCHWARTZ J.B. (2014), A comparison of direct and calculated free 25(OH) Vitamin D levels in clinical populations, J. Clin. Endocrinol. Metab., 99(5):1631-7.
- BIKLE D. (2013), Variability in free 25(OH) vitamin D levels in clinical populations, J. Steroid Biochem. Mol. Biol., S09600760.
- BIKLE D., Gee E, Halloran B, Haddad JG. Free 1,25-dihydroxyvitamin D levels in serum from normal subjects, pregnant subjects, and subjects with liver disease. J Clin Invest 1984;74:1966-71.

8. VAN HOOFF H.J., Swinkels LM, Ross HA, Sweep CG, Benraad TJ. Determination of non-protein-bound plasma 1,25-dihydroxyvitamin D by symmetric (rate) dialysis. Anal Biochem. 1998 May 1;258(2):176-83.

XVIII. PROTOCOL SUMMARY

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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
XVII. REFERENTIES

1. BOUILLON R. (2016), Free 25-hydroxyvitamin D: impact of vitamin D binding protein assays on racial genotypic associations, JCEM, <http://dx.doi.org/10.1210/jc.2016-1104>.
2. SOLLID S.T. (2016), Effects of vitamin D binding protein phenotypes and vitamin D supplementation on serum total 25(OH) D and directly measured free 25(OH)D, Eur. J. Endocrinol. April 1, 174:445-452.
3. TANGPRICHA V. (2015), Free 25-Hydroxyvitamin D Concentrations in Cystic Fibrosis, Am. J. Med. Sci. 2015 Nov;350(5):374-9.
4. ALOIA J. (2015), Free 25(OH)D and the Vitamin D Paradox in African Americans, J. Clin. Endocrinol. Metab. 2015 Jul 10:JC20152066.
5. SCHWARTZ J.B. (2014), A comparison of direct and calculated free 25(OH) Vitamin D levels in clinical populations, J. Clin. Endocrinol. Metab., 99(5):1631-7.
6. BIKLE D. (2013), Variability in free 25(OH) vitamin D levels in clinical populations, J. Steroid Biochem. Mol. Biol., S09600760.
7. BIKLE D., Gee E, Halloran B, Haddad JG. Free 1,25-dihydroxyvitamin D levels in serum from normal subjects, pregnant subjects, and subjects with liver disease. J Clin Invest 1984;74:1966-71.
8. VAN HOOFF H.J., Swinkels LM, Ross HA, Sweep CG, Benraad TJ. Determination of non-protein-bound plasma 1,25-dihydroxyvitamin D by symmetric (rate) dialysis. Anal Biochem. 1998 May 1;258(2):176-83.

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