
Instructions for use / Gebrauchsanweisung
3-CAT ELISA

REF

KAPL10-1600

IVD

Table of contents

| | | |
|-------|---|----|
| 1. | Introduction | 4 |
| 1.1 | Intended use and principle of the test | 4 |
| 1.2 | Clinical application | 4 |
| 2. | Procedural cautions, guidelines, warnings and limitations | 4 |
| 2.1 | Procedural cautions, guidelines and warnings | 4 |
| 2.2 | Limitations | 5 |
| 2.2.1 | Interfering substances and proper handling of specimens | 5 |
| 2.2.2 | Drug and food interferences | 5 |
| 2.2.3 | High-Dose-Hook effect | 5 |
| 3. | Storage and stability | 5 |
| 4. | Materials | 6 |
| 4.1 | Contents of the kit | 6 |
| 4.2 | Calibration and Controls | 8 |
| 4.3 | Additional materials required but not provided in the kit | 8 |
| 4.4 | Additional equipment required but not provided in the kit | 8 |
| 5. | Sample collection, handling and storage | 8 |
| 6. | Test procedure | 8 |
| 6.1 | Preparation of reagents and further notes | 9 |
| 6.2 | Sample preparation, extraction and acylation | 9 |
| 6.3 | Adrenaline ELISA | 10 |
| 6.4 | Noradrenaline ELISA | 10 |
| 6.5 | Dopamine ELISA | 11 |
| 7. | Calculation of results | 11 |
| 7.1 | Expected reference value | 12 |
| 7.2 | Typical standard curve | 12 |
| 8. | Control samples | 12 |
| 9. | Assay characteristics | 12 |
| 9.1 | Performance data | 12 |
| 10. | References/Literature | 14 |
| 11. | Changes | 14 |

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1. Introduction**1.1 Intended use and principle of the test**

Enzyme Immunoassay for the quantitative determination of adrenaline (epinephrine), noradrenaline (norepinephrine) and dopamine in plasma and urine.

Adrenaline (epinephrine), noradrenaline (norepinephrine) and dopamine are extracted by using a cis-diol-specific affinity gel, acylated and then converted enzymatically.

The subsequent competitive ELISA uses the microtiter plate format. The antigen is bound to the solid phase of the microtiter plate. The acylated standards, controls and samples compete with the solid phase bound analytes for a fixed number of antibody binding sites. After the system is in equilibrium, free antigen and free antigen-antibody complexes are removed by washing. The antibody bound to the solid phase is detected by an anti-rabbit IgG-peroxidase conjugate using TMB as a substrate resulting in a colour reaction. The reaction is monitored at a wavelength of 450 nm.

Quantification of unknown samples is achieved by comparing their absorbance with a reference curve prepared with known standard concentrations. Manual processing of the ELISA is recommended. The use of automatic laboratory equipment is the responsibility of the user. This in-vitro diagnostic is for professional use only.

1.2 Clinical application

In humans the catecholamines adrenaline (epinephrine), noradrenaline (norepinephrine) and dopamine are neurotransmitters of the sympathetic nervous system and are involved in many physiological processes. The sympathetic nervous system sets the body to a heightened state of alert, also called as the body's fight-or-flight response.

In the human body the catecholamines and their metabolites indicate the adaptation of the body to acute and chronic stress.

Next to the metanephrine/normetanephrine the catecholamines are important for the diagnosis and the follow-up of tumors of the sympathoadrenal system like the pheochromocytoma. The quantitative determination of catecholamines in urine is preferred for the diagnosis of these tumors, whereas the determination of catecholamines in plasma is medically sensible for the localization of the tumor and for function testing. Values above the cut-off can provide an indication for neuroendocrine tumors.

However, in literature various diseases like hypertension, cardiovascular diseases, schizophrenia and manic depression are described with abnormal low or high levels of catecholamines.

Therapeutic consequences should never be based on laboratory results alone, even if these results are assessed in accordance with the quality criteria of the method. Any laboratory result is only a part of the total clinical picture of the patient.

Only in cases where the laboratory results are in an acceptable agreement with the overall clinical picture of the patient, it can be used for therapeutic consequences.

2. Procedural cautions, guidelines, warnings and limitations**2.1 Procedural cautions, guidelines and warnings**

- (1) This kit is intended for professional use only. Users should have a thorough understanding of this protocol for the successful use of this kit. Only the test instruction provided with the kit is valid and must be used to run the assay. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
- (2) This assay was validated for a certain type of sample as indicated in Intended Use (please refer to Chapter 1). Any off-label use of this kit is in the responsibility of the user and the manufacturer cannot be held liable.
- (3) The principles of Good Laboratory Practice (GLP) must be followed.
- (4) In order to reduce exposure to potentially harmful substances, wear lab coats, disposable protective gloves and protective glasses where necessary.
- (5) If serious incidents should occur in connection with this product, they should be reported to the manufacturer and the competent national authorities.
- (6) All kit reagents and specimens should be brought to room temperature and mixed gently but thoroughly before use. For dilution or reconstitution purposes, use deionized, distilled, or ultra-pure water. Avoid repeated freezing and thawing of reagents and specimens.

- (7) The microplate contains snap-off strips. Unused wells must be stored at 2 – 8 °C in the sealed foil pouch with desiccant and used in the frame provided. Microtiter strips which are removed from the frame for usage should be marked accordingly to avoid any mix-up.
- (8) Duplicate determination of sample is highly recommended.
- (9) Once the test has been started, all steps should be completed without interruption. Make sure that the required reagents, materials, and devices are prepared for use at the appropriate time.
- (10) Incubation times do influence the results. All wells should be handled in the same order and time intervals.
- (11) To avoid cross-contamination of reagents, use new disposable pipette tips for dispensing each reagent, sample, standard and control.
- (12) A standard curve must be established for each run.
- (13) The controls should be included in each run and fall within established confidence limits. The confidence limits are listed in the QC-Report provided with the kit.
- (14) Do not mix kit components with different lot numbers within a test and do not use reagents beyond expiry date as shown on the kit labels.
- (15) Avoid contact with Stop Solution containing 0.25 M H₂SO₄. It may cause skin irritation and burns. In case of contact with eyes or skin, rinse off immediately with water.
- (16) TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Rinse contaminated items before reuse.
- (17) For information about hazardous substances included in the kit please refer to Safety Data Sheet (SDS). The Safety Data Sheet for this product is made available directly on the website of the manufacturer or upon request.
- (18) Kit reagents must be regarded as hazardous waste and disposed of according to national regulations.
- (19) The expected reference values reported in this test instruction are only indicative. It is recommended that each laboratory establishes its own reference intervals.
- (20) In case of any severe damage to the test kit or components, the manufacturer has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components must not be used for a test run. They must be stored properly until the manufacturer decides what to do with them. If it is decided that they are no longer suitable for measurements, they must be disposed of in accordance with national regulations.
- (21) The results obtained with this test kit should not be taken as the sole reason for any therapeutic consequence but must be correlated to other diagnostic tests and clinical observations.

2.2 Limitations

Any inappropriate handling of samples or modification of this test might influence the results.

2.2.1 Interfering substances and proper handling of specimens

Plasma

Samples containing precipitates or fibrin strands or which are hemolytic or lipemic might cause inaccurate results. Hemolytic samples (up to 4 mg/ml hemoglobin), icteric samples (up to 50 mg/dl bilirubin) and lipemic samples (up to 800 mg/dl triglycerides) have no influence on the assay results.

If the concentrations cannot be estimated and there are doubts as to whether the above limit values for hemolytic, icteric or lipemic samples are complied with, the samples should not be used in the assay.

24-hour urine

Please note the sample collection! If the percentage of the final concentration of acid is too high, this will lead to incorrect results for the urine samples.

2.2.2 Drug and food interferences

There are no known substances (drugs) which ingestion interferes with the measurement of catecholamine level in the sample.

2.2.3 High-Dose-Hook effect

No hook effect was observed in this test.

3. Storage and stability

Store kit and reagents at 2 – 8 °C until expiration date. Do not use kit and components beyond the expiry date indicated on the kit labels. Once opened, the reagents are stable for 2 months when stored at 2 – 8 °C. Once the resealable pouch of the ELISA plate has been opened, care should be taken to close it tightly again including the desiccant.

4. Materials

4.1 Contents of the kit

| | | |
|------------------|--|---|
| BA D-0090 | FOILS | Adhesive Foil – ready to use |
| Content: | Adhesive foils in a resealable pouch | |
| Number: | 3 x 4 foils | |
| BA E-0030 | WASH-CONC 50x | Wash Buffer Concentrate – concentrated 50x |
| Content: | Buffer with a non-ionic detergent and physiological pH | |
| Volume: | 3 x 20 ml/vial, purple cap | |
| BA E-0040 | CONJUGATE | Enzyme Conjugate – ready to use |
| Content: | Goat anti-rabbit immunoglobulins conjugated with peroxidase | |
| Volume: | 3 x 12 ml/vial, red cap | |
| Description: | Species is goat | |
| BA E-0055 | SUBSTRATE | Substrate – ready to use |
| Content: | Chromogenic substrate containing 3,3',5,5'-tetramethylbenzidine, substrate buffer and hydrogen peroxide | |
| Volume: | 3 x 12 ml/vial, black cap | |
| BA E-0080 | STOP-SOLN | Stop Solution – ready to use |
| Content: | 0.25 M sulfuric acid | |
| Volume: | 3 x 12 ml/vial, grey cap | |
| BA E-0131 | W ADR MN | Adrenaline Microtiter Strips – ready to use |
| Content: | 1 x 96 wells (12x8) antigen precoated microwell plate in a resealable blue pouch with desiccant | |
| BA E-0231 | W NAD NMN | Noradrenaline Microtiter Strips – ready to use |
| Content: | 1 x 96 wells (12x8) antigen precoated microwell plate in a resealable yellow pouch with desiccant | |
| BA E-0331 | W DOP | Dopamine Microtiter Strips – ready to use |
| Content: | 1 x 96 wells (12x8) antigen precoated microwell plate in a resealable green pouch with desiccant | |
| BA E-6110 | ADR-AS | Adrenaline Antiserum – ready to use |
| Content: | Rabbit anti-adrenaline antibody in buffer with proteins and non-mercury preservative, blue coloured | |
| Volume: | 1 x 6 ml/vial, blue cap | |
| Description: | Species of antibody is rabbit, species of protein in buffer is bovine | |
| BA E-6210 | NAD-AS | Noradrenaline Antiserum – ready to use |
| Content: | Rabbit anti-noradrenaline antibody in buffer with proteins and non-mercury preservative, yellow coloured | |
| Volume: | 1 x 6 ml/vial, yellow cap | |
| Description: | Species of antibody is rabbit, species of protein in buffer is bovine | |
| BA E-6310 | DOP-AS | Dopamine Antiserum – ready to use |
| Content: | Rabbit anti-dopamine antibody in buffer with proteins and non-mercury preservative, green coloured | |
| Volume: | 1 x 6 ml/vial, green cap | |
| Description: | Species of antibody is rabbit, species of protein in buffer is bovine | |
| BA E-6612 | ACYL-REAG | Acylation Reagent – ready to use |
| Content: | Acylation reagent in DMSO | |
| Volume: | 1 x 3 ml/vial, white cap | |

| | | |
|---------------------------|--|---|
| BA R-0050 | ADJUST-BUFF | Adjustment Buffer – ready to use |
| Content: | TRIS buffer | |
| Volume: | 2 x 4 ml/vial, green cap | |
| BA R-6611 | ACYL-BUFF | Acylation Buffer – ready to use |
| Content: | Buffer with light alkaline pH for the acylation | |
| Volume: | 1 x 20 ml/vial, white cap | |
| BA R-6613 | ASSAY-BUFF | Assay Buffer – ready to use |
| Content: | 1 M hydrochloric acid and a non-mercury preservative | |
| Volume: | 1 x 6 ml/vial, grey cap | |
| Hazard pictograms: |  | |
| | GHS05 | |
| Signal word: | Danger | |
| Hazard statements: | H314 Causes severe skin burns and eye damage. | |
| Precautionary statements: | P280 Wear protective gloves, protective clothing, eye protection. P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P310 Immediately call a doctor, a POISON CENTER. P501 Dispose of contents/container to an authorised waste collection point. | |
| BA R-6614 | COENZYME | Coenzyme – ready to use |
| Content: | S-adenosyl-L-methionine | |
| Volume: | 1 x 4 ml/vial, purple cap | |
| BA R-6615 | ENZYME | Enzyme – lyophilized |
| Content: | Catechol-O-methyltransferase | |
| Volume: | 6 vials, pink cap | |
| Description: | Catechol-O-methyltransferase from pig liver | |
| BA R-6617 | EXTRACT-BUFF | Extraction Buffer – ready to use |
| Content: | Buffer containing carbonate | |
| Volume: | 1 x 6 ml/vial, brown cap | |
| BA R-6618 | EXTRACT-PLATE 48 | Extraction Plate – ready to use |
| Content: | 2 x 48 well plates coated with boronate affinity gel in a resealable pouch | |
| BA R-6619 | HCL | Hydrochloric Acid – ready to use |
| Content: | 0.025 M Hydrochloric Acid, yellow coloured | |
| Volume: | 1 x 20 ml/vial, green cap | |

4.2 Calibration and Controls

Standards and Controls – ready to use

| Cat. no. | Component | Colour /Cap | Concentration [ng/ml] | | | Concentration [nmol/l] | | | Volume/ Vial |
|-----------|--------------|-------------|---|-------|-------|------------------------|-------|--------|--------------|
| | | | ADR | NAD | DOP | ADR | NAD | DOP | |
| BA E-6601 | STANDARD A | white | 0 | 0 | 0 | 0 | 0 | 0 | 4 ml |
| BA E-6602 | STANDARD B | yellow | 1 | 5 | 10 | 5.5 | 30 | 65 | 4 ml |
| BA E-6603 | STANDARD C | orange | 4 | 20 | 40 | 22 | 118 | 261 | 4 ml |
| BA E-6604 | STANDARD D | blue | 15 | 75 | 150 | 82 | 443 | 980 | 4 ml |
| BA E-6605 | STANDARD E | grey | 50 | 250 | 500 | 273 | 1,478 | 3,265 | 4 ml |
| BA E-6606 | STANDARD F | black | 200 | 1,000 | 2,000 | 1,092 | 5,910 | 13,060 | 4 ml |
| BA E-6609 | STANDARD A/B | purple | - | - | 4.5 | - | - | 29 | 4 ml |
| BA E-6651 | CONTROL 1 | green | Refer to QC-Report for expected value and acceptable range. | | | | | | 4 ml |
| BA E-6652 | CONTROL 2 | red | | | | | | | 4 ml |

Conversion: adrenaline [ng/ml] x 5.46 = adrenaline [nmol/l]
 noradrenaline [ng/ml] x 5.91 = noradrenaline [nmol/l]
 dopamine [ng/ml] x 6.53 = dopamine [nmol/l]

Content: Acidic buffer with non-mercury stabilizer, spiked with defined quantity of adrenaline, noradrenaline and dopamine

⚠ **for the determination of dopamine in plasma the additional **Standard A/B** is mandatory!*

4.3 Additional materials required but not provided in the kit

- Water (deionized, distilled, or ultra-pure)
- Absorbent material (paper towel)

4.4 Additional equipment required but not provided in the kit

- Calibrated precision pipettes to dispense volumes between 10 – 700 µl; 1 ml
- Microtiter plate washing device (manual, semi-automated or automated)
- ELISA reader capable of reading absorbance at 450 nm and if possible 620 – 650 nm
- Microtiter plate shaker (shaking amplitude 3 mm; approx. 600 rpm)
- Vortex mixer

5. Sample collection, handling and storage

Plasma

Whole blood should be collected into centrifuge tubes containing EDTA as anti-coagulant and centrifuged according to manufacturer's instructions immediately after collection.

In case of hemolytic, icteric or lipemic samples see 2.2.1.

Storage: up to 6 hours at 2 – 8 °C, for longer period (up to 6 months) at -20 °C.

Repeated freezing and thawing should be avoided.

Urine

Spontaneous urine or 24-hour urine, collected in a bottle containing 10 – 15 ml of 6 M HCl, can be used. If 24-hour urine is used please record the total volume of the collected urine.

Storage: up to 48 hours at 2 – 8 °C, up to 24 hours at room temperature, for longer periods (up to 6 months) at -20 °C. Repeated freezing and thawing should be avoided.

Avoid exposure to direct sunlight.

6. Test procedure

Allow all reagents and samples to reach room temperature and mix thoroughly by gentle inversion before use. Number the Extraction Plate and microwell plates (Microtiter Strips which are removed from the frame for usage should be marked accordingly to avoid any mix-up). Duplicate determinations are recommended.

The binding of the antisera and of the enzyme conjugate and the activity of the enzyme are temperature dependent. The higher the temperature, the higher the absorption values will be. Varying incubation times will have similar influences on the absorbance. The optimal temperature during the enzyme immunoassay is between 20 – 25 °C.

⚠️ The use of a microtiter plate shaker with the following specifications is mandatory: shaking amplitude 3 mm; approx. 600 rpm. Shaking with differing settings might influence the results.

6.1 Preparation of reagents and further notes

Wash Buffer

Dilute the 20 ml Wash Buffer Concentrate **WASH-CONC 50X** with water to a final volume of 1000 ml.
Storage: 2 months at 2 – 8 °C

Enzyme Solution

Reconstitute the content of the vial **ENZYME** with 1 ml water (deionized, distilled, or ultra-pure) and mix thoroughly. Add 0.3 ml of **COENZYME** followed by 0.7 ml of **ADJUST-BUFF**. The total volume of the Enzyme Solution is 2.0 ml.

⚠️ The Enzyme Solution has to be prepared freshly prior to the assay (not longer than 10 – 15 minutes in advance). Discard after use!

Adrenaline Microtiter Strips, Noradrenaline Microtiter Strips and Dopamine Microtiter Strips

In rare cases residues of the blocking and stabilizing reagent can be seen in the wells as small, white dots or lines. These residues do not influence the quality of the product.

Acylation Reagent

The **ACYL-REAG** (BA E-6612) has a freezing point of 18.5 °C. To ensure that it is liquid when being used, it must be ensured that the Acylation Reagent has reached room temperature and forms a homogeneous, crystal-free solution before being used.

6.2 Sample preparation, extraction and acylation

⚠️ *for the determination of dopamine in plasma the additional **Standard A/B** is mandatory!

| | | | | | | | | | |
|-------------------------------------|--|--------------------------|---------------|----------------------|--------------|-------------------------------------|--------------|--------------------------|--------------|
| 1. | Pipette 10 µl of standards, controls, urine samples and 300 µl of plasma samples into the respective wells of the EXTRACT-PLATE 48 . | | | | | | | | |
| 2. | Add 250 µl of water (deionized, distilled, or ultra-pure) to the wells with standards, controls and urine samples . | | | | | | | | |
| 3. | Pipette 50 µl of ASSAY-BUFF into all wells. | | | | | | | | |
| 4. | Pipette 50 µl of EXTRACT-BUFF into all wells. | | | | | | | | |
| 5. | Cover plate with FOILS and incubate 30 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). | | | | | | | | |
| 6. | Remove the foil. Empty plate and blot dry by tapping the inverted plate on absorbent material. | | | | | | | | |
| 7. | Pipette 1 ml of Wash Buffer into all wells. Incubate the plate for 5 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). Empty plate and blot dry by tapping the inverted plate on absorbent material. | | | | | | | | |
| 8. | Pipette another 1 ml of Wash Buffer into all wells. Incubate the plate for 5 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). Empty plate and blot dry by tapping the inverted plate on absorbent material. | | | | | | | | |
| 9. | Pipette 150 µl of ACYL-BUFF into all wells. | | | | | | | | |
| 10. | Pipette 25 µl of ACYL-REAG into all wells. | | | | | | | | |
| 11. | Incubate 15 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). | | | | | | | | |
| 12. | Empty plate and blot dry by tapping the inverted plate on absorbent material. | | | | | | | | |
| 13. | Pipette 1 ml of Wash Buffer into all wells. Incubate the plate for 10 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). Empty plate and blot dry by tapping the inverted plate on absorbent material. | | | | | | | | |
| 14. | Pipette 175 µl of HCL into all wells. | | | | | | | | |
| 15. | Cover plate with FOILS . Incubate 10 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). Remove the foil and discard. | | | | | | | | |
| ⚠️ | Do not decant the supernatant thereafter! | | | | | | | | |
| | The following volumes of the supernatant are needed for the subsequent ELISA: | | | | | | | | |
| | <table border="1"> <tr> <td>Adrenaline</td> <td>100 µl</td> <td>Noradrenaline</td> <td>20 µl</td> </tr> <tr> <td>Dopamine (standards + urine)</td> <td>25 µl</td> <td>Dopamine (plasma)</td> <td>50 µl</td> </tr> </table> | Adrenaline | 100 µl | Noradrenaline | 20 µl | Dopamine (standards + urine) | 25 µl | Dopamine (plasma) | 50 µl |
| Adrenaline | 100 µl | Noradrenaline | 20 µl | | | | | | |
| Dopamine (standards + urine) | 25 µl | Dopamine (plasma) | 50 µl | | | | | | |

6.3 Adrenaline ELISA

| |
|--|
| 1. Pipette 25 µl of the Enzyme Solution (refer to 6.1) into all wells of the Adrenaline Microtiter Strips U ADR MN . |
| 2. Pipette 100 µl of the extracted standards, controls and samples into the appropriate wells. |
| 3. Incubate for 30 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). |
| 4. Pipette 50 µl of the ADR-AS into all wells and cover plate with FOILS . |
| 5. Incubate for 2 h at RT (20 – 25 °C) on a shaker (approx. 600 rpm). |
| 6. Remove the foil. Discard or aspirate the content of the wells. Wash the plate 3 x by adding 300 µl of Wash Buffer , discarding the content and blotting dry each time by tapping the inverted plate on absorbent material. |
| 7. Pipette 100 µl of the CONJUGATE into all wells. |
| 8. Incubate for 30 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). |
| 9. Discard or aspirate the content of the wells. Wash the plate 3 x by adding 300 µl of Wash Buffer , discarding the content and blotting dry each time by tapping the inverted plate on absorbent material. |
| 10. Pipette 100 µl of the SUBSTRATE into all wells and incubate for 25 ± 5 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). ⚠ Avoid exposure to direct sunlight! |
| 11. Add 100 µl of the STOP-SOLN to all wells and shake the microtiter plate shortly. |
| 12. Read the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to 450 nm (if available a reference wavelength between 620 nm and 650 nm is recommended). |

6.4 Noradrenaline ELISA

| |
|--|
| 1. Pipette 25 µl of the Enzyme Solution (refer to 6.1) into all wells of the Noradrenaline Microtiter Strips U NAD NMN . |
| 2. Pipette 20 µl of the extracted standards, controls and samples into the appropriate wells. |
| 3. Incubate for 30 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). |
| 4. Pipette 50 µl of the NAD-AS into all wells and cover plate with FOILS . |
| 5. Incubate for 2 h at RT (20 – 25 °C) on a shaker (approx. 600 rpm). |
| 6. Remove the foil. Discard or aspirate the content of the wells. Wash the plate 3 x by adding 300 µl of Wash Buffer , discarding the content and blotting dry each time by tapping the inverted plate on absorbent material. |
| 7. Pipette 100 µl of the CONJUGATE into all wells. |
| 8. Incubate for 30 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). |
| 9. Discard or aspirate the content of the wells. Wash the plate 3 x by adding 300 µl of Wash Buffer , discarding the content and blotting dry each time by tapping the inverted plate on absorbent material. |
| 10. Pipette 100 µl of the SUBSTRATE into all wells and incubate for 25 ± 5 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). ⚠ Avoid exposure to direct sunlight! |
| 11. Add 100 µl of the STOP-SOLN to all wells and shake the microtiter plate shortly. |
| 12. Read the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to 450 nm (if available a reference wavelength between 620 nm and 650 nm is recommended). |

6.5 Dopamine ELISA

| |
|---|
| 1. Pipette 25 µl of the Enzyme Solution (refer to 6.1) into all wells of the Dopamine Microtiter Strips U DOP . |
| 2. Pipette 25 µl of the extracted standards, controls, urine samples and 50 µl of the extracted plasma samples into the appropriate wells. |
| 3. Add 25 µl of HCL to the standards, controls and urine samples . |
| 4. Incubate for 30 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). |
| 5. Pipette 50 µl of the DOP-AS into all wells and cover plate with FOILS . |
| 6. Incubate for 2 h at RT (20 – 25 °C) on a shaker (approx. 600 rpm). |
| 7. Remove the foil. Discard or aspirate the content of the wells. Wash the plate 3 x by adding 300 µl of Wash Buffer , discarding the content and blotting dry each time by tapping the inverted plate on absorbent material. |
| 8. Pipette 100 µl of the CONJUGATE into all wells. |
| 9. Incubate for 30 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). |
| 10. Discard or aspirate the content of the wells. Wash the plate 3 x by adding 300 µl of Wash Buffer , discarding the content and blotting dry each time by tapping the inverted plate on absorbent material. |
| 11. Pipette 100 µl of the SUBSTRATE into all wells and incubate for 25 ± 5 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). ⚠ Avoid exposure to direct sunlight! |
| 12. Add 100 µl of the STOP-SOLN to all wells and shake the microtiter plate shortly. |
| 13. Read the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to 450 nm (if available a reference wavelength between 620 nm and 650 nm is recommended). |

7. Calculation of results

| Measuring range | | Adrenaline | Noradrenaline | Dopamine |
|-----------------|--------|------------------|-------------------|-------------------|
| | Urine | 0.7 – 200 ng/ml | 2.5 – 1,000 ng/ml | 4.8 – 2,000 ng/ml |
| | Plasma | 18 – 6,667 pg/ml | 93 – 33,333 pg/ml | 75 – 33,333 pg/ml |

The standard curve, which can be used to determine the concentration of the unknown samples, is obtained by plotting the absorbance readings (calculate the mean absorbance) of the standards (linear, y-axis) against the corresponding standard concentrations (logarithmic, x-axis) using a concentration of 0.001 ng/ml for Standard A (this alignment is mandatory because of the logarithmic presentation of the data). Use non-linear regression for curve fitting (e. g. 4-parameter, marquardt).

⚠ *This assay is a competitive assay. This means: the OD-values are decreasing with increasing concentrations of the analyte. OD-values found below the standard curve correspond to high concentrations of the analyte in the sample and have to be reported as being positive.*

Urine samples and controls

The concentrations of the **urine samples** and the **Controls** can be read directly from the standard curve.

Calculate the 24 h excretion for each urine sample: **µg/24h = µg/l x l/24h**

Plasma samples

The read **Adrenaline and Noradrenaline** concentrations of the **plasma samples** have to be **divided by 30**.

The read **Dopamine** concentrations of the **plasma samples** have to be **divided by 60**.

Conversion:

Adrenaline [ng/ml] x 5.46 = Adrenaline [nmol/l]

Noradrenaline [ng/ml] x 5.91 = Noradrenaline [nmol/l]

Dopamine [ng/ml] x 6.53 = Dopamine [nmol/l]


7.1 Expected reference value

It is strongly recommended that each laboratory should determine its own reference values.

| | Adrenaline | Noradrenaline |
|---------------|-------------------------------|-------------------------------|
| 24-hour urine | < 20 µg/day (110 nmol/day) | < 90 µg/day (535 nmol/day) |
| Plasma | < 100 pg/ml | < 600 pg/ml |

7.2 Typical standard curve

⚠️ Examples: Do not use for calculation!

8. Control samples

It is recommended to use control samples according to national regulations. Use controls at both normal and pathological levels. Commercially obtained control samples should be treated like unknown samples. Control samples should fall within established confidence limits. The confidence limits of the kit controls are indicated on the QC-Report.

9. Assay characteristics

9.1 Performance data

| Analytical Sensitivity | | Adrenaline | Noradrenaline | Dopamine |
|-------------------------------|----------------|-------------------|----------------------|-----------------|
| Limit of Blank (LOB) | Urine [ng/ml] | 0.8 | 1.5 | 2.2 |
| | Plasma [pg/ml] | 9.3 | 32 | 43 |
| Limit of Detection (LOD) | Urine [ng/ml] | 0.9 | 1.7 | 2.5 |
| | Plasma [pg/ml] | 10 | 36 | 49 |
| Limit of Quantification (LOQ) | Urine [ng/ml] | 0.7 | 2.5 | 4.8 |
| | Plasma [pg/ml] | 18 | 93 | 75 |

| Analytical Specificity (Cross Reactivity) | | | |
|--|-----------------------------|---------------|----------|
| Substance | Cross Reactivity [%] | | |
| | Adrenaline | Noradrenaline | Dopamine |
| Derivatized Adrenaline | 100 | 0.08 | 0.02 |
| Derivatized Noradrenaline | 0.13 | 100 | 6.4 |
| Derivatized Dopamine | < 0.01 | 0.03 | 100 |
| Metanephrine | 0.18 | < 0.01 | < 0.01 |
| Normetanephrine | < 0.01 | 0.16 | 0.01 |
| 3-Methoxytyramine | < 0.01 | < 0.01 | 0.49 |
| 3-Methoxy-4-hydroxyphenylglycol | < 0.01 | < 0.01 | < 0.01 |
| Tyramine | < 0.01 | < 0.01 | 0.18 |
| Phenylalanine, Caffeinic acid, L-Dopa, Homovanillic acid, Tyrosine, 3-Methoxy-4-hydroxymandelic acid | < 0.01 | < 0.01 | < 0.01 |

| Precision | | | | | | | |
|-----------------------------------|--------|---------------|--------|------------------------------------|--------|---------------|--------|
| Intra-Assay Urine (n = 60) | | | | Intra-Assay Plasma (n = 60) | | | |
| | Sample | Range [ng/ml] | CV [%] | | Sample | Range [pg/ml] | CV [%] |
| Adrenaline | 1 | 6.2 ± 1.1 | 17.4 | Adrenaline | 1 | 64.7 ± 15.9 | 24.7 |
| | 2 | 21.4 ± 2.7 | 12.4 | | 2 | 258 ± 32.5 | 12.7 |
| | 3 | 59.4 ± 7.8 | 13.1 | | 3 | 948 ± 105 | 11.0 |
| Noradrenaline | 1 | 26.1 ± 3.6 | 13.8 | Noradrenaline | 1 | 510 ± 65 | 12.8 |
| | 2 | 97 ± 12.8 | 13.4 | | 2 | 1,358 ± 194 | 14.3 |
| | 3 | 267 ± 35 | 13.1 | | 3 | 3,363 ± 374 | 11.1 |
| Dopamine | 1 | 82 ± 16.1 | 19.7 | Dopamine | 1 | 75 ± 22 | 29.8 |
| | 2 | 253 ± 41.1 | 16.3 | | 2 | 353 ± 86 | 24.4 |
| | 3 | 714 ± 67 | 9.4 | | 3 | 1,187 ± 293 | 24.9 |
| Inter-Assay Urine (n = 33) | | | | Inter-Assay Plasma (n = 18) | | | |
| | Sample | Range [ng/ml] | CV [%] | | Sample | Range [pg/ml] | CV [%] |
| Adrenaline | 1 | 5.2 ± 0.9 | 17.9 | Adrenaline | 1 | 76.4 ± 11.1 | 14.5 |
| | 2 | 17.8 ± 2.1 | 11.7 | | 2 | 247 ± 27.5 | 11.1 |
| | 3 | 54.2 ± 6.6 | 12.1 | | 3 | 771 ± 101 | 13.1 |
| Noradrenaline | 1 | 19.5 ± 3.9 | 20.0 | Noradrenaline | 1 | 445 ± 40.9 | 9.2 |
| | 2 | 80.6 ± 10.6 | 13.2 | | 2 | 1,232 ± 134 | 10.9 |
| | 3 | 226 ± 39.5 | 17.4 | | 3 | 3,283 ± 302 | 9.2 |
| Dopamine | 1 | 79.3 ± 18.8 | 23.7 | Dopamine | 1 | 238 ± 67.0 | 28.2 |
| | 2 | 222 ± 27.0 | 12.1 | | 2 | 1,072 ± 201 | 18.8 |
| | 3 | 630 ± 69.0 | 11.0 | | 3 | 3,449 ± 491 | 14.2 |

| Recovery | | | | |
|-----------------|--------|---------------------|----------|-----------|
| | | Range | Mean [%] | Range [%] |
| Adrenaline | Urine | 4.5 – 53.5 ng/ml | 106 | 94 – 120 |
| | Plasma | 9.1 – 4,268 pg/ml | 105 | 88 – 117 |
| Noradrenaline | Urine | 58.6 – 260 ng/ml | 103 | 91 – 113 |
| | Plasma | 51 – 14,251 pg/ml | 87 | 75 – 107 |
| Dopamine | Urine | 225 – 1,306 ng/ml | 110 | 101 – 124 |
| | Plasma | 57.4 – 16,054 pg/ml | 89 | 84 – 92 |

| Linearity | | | | |
|---------------|--------|-----------------------|----------|-----------|
| | | Serial dilution up to | Mean [%] | Range [%] |
| Adrenaline | Urine | 1:512 | 108 | 92 – 123 |
| | Plasma | 1:512 | 105 | 94 – 115 |
| Noradrenaline | Urine | 1:512 | 112 | 100 – 127 |
| | Plasma | 1:512 | 112 | 102 – 125 |
| Dopamine | Urine | 1:512 | 104 | 83 – 126 |
| | Plasma | 1:512 | 106 | 85 – 132 |

10. References/Literature



- (1) Kim et al. Vitamin C prevents stress-induced damage on the heart caused by the death of cardiomyocytes, through the down-regulation of the excessive production of catecholamine, TNF- α , and ROS production in GULO(-I-) Vit C-Insufficient mice. Free Radical Biology and Medicine, 65:573-583 (2013)
- (2) Bada et al. Peripheral vasodilatation determines cardiac output in exercising humans: insight from atrial pacing. The Journal of Physiology, 590(8):2051-2060 (2012)
- (3) Parks et al. Employment and work schedule are related to telomere length in women. Occupational & Environmental Medicine 68(8):582-589 (2011)

For updated literature or any other information please contact your local supplier.

11. Changes

| Version | Release Date | Chapter | Change |
|---------|--------------|---------|---|
| 19.0 | 2023-11-28 | 4.1 | - Hazard labelling updated according to SDS |
| 19.1 | 2024-04-02 | 6.5 | - Correction antiserum symbol |

Symbols:

| | | | | | |
|---|------------------------------|---|------------------|---|-----------------------------------|
|  | Storage temperature |  | Manufacturer |  | Contains sufficient for <n> tests |
|  | Use-by date |  | Batch code |  | For in-vitro diagnostic use only! |
|  | Consult instructions for use |  | Content |  | CE marking of conformity |
|  | Caution |  | Catalogue number |  | Distributor |
|  | Date of manufacture | | | | |

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- Adrenaline ELISA ^{Fast Track}
- Noradrenaline ELISA ^{Fast Track}
- Dopamine ELISA ^{Fast Track}
- 3-CAT ELISA ^{Fast Track}

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