DHEA RIA

KIPB1138

History

Summary of change :

Previous Version : 210121	Current Version : 240412
Old DiaSource logo	New DiaSource logo on the front page
Contact address tech.support@diasource.be	Correction of the contact address products.support@diasource.be

DHEA RIA

en

Radioimmunoassay for the *in vitro* determination of dehydroepiandrosterone in human serum or plasma

KIPB1138

IN VITRO DIAGNOSTIC

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 90

1. INTENDED USE

Measurement of serum/plasma DHEA is a useful marker of adrenal androgen synthesis. His measurement contributes to the:

- Diagnosis of hyperandrogenism states in women
- Diagnosis and monitoring of CAH (congenital adrenal hyperpalsia)
- Diagnosis and differential diagnosis of premature adrenarche, diagnosis of androgen-adrenal secreting tumor.

2. PRINCIPLE OF THE ASSAY

The radioimmunoassay of dehydroepiandrosterone (DHEA) is a competition assay. Prior to the assay, samples (serum or plasma) are extracted with ethyl ether; the solvent is evaporated and the dry residues are re-dissolved in the recovery buffer of the kit. The re-dissolved extracts and calibrators are then incubated with ¹²⁵I-labeled DHEA, as tracer in antibody-coated tube. After incubation, the content of tubes is aspirated and the bound radioactivity is determined. A calibration curve is established and unknown values are determined by interpolation from the curve.

3. REAGENTS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Expiry dates printed on vial labels only apply to the long-term storage of components by the manufacturer. Do not take into account.

Anti-DHEA monoclonal antibody-coated tubes : 2 x 50 tubes (ready-to-use)

Ag	¹²⁵ I
----	------------------

¹²⁵I-labeled DHEA tracer : one 55 mL vial (ready-to-use)
The vial contains 185 kBq, at the date of manufacture, of ¹²⁵I-labeled DHEA in liquid form with proteins, sodium azide (< 0.1 %; see § Precautions) and a dye.

CAL	N
-----	---

Calibrators : six 1 mL vials (ready-to-use)
The calibrator vials contain from 0 to approximately 30 ng/mL (0 to 104 nmol/L) of DHEA in the same buffer as that employed for re-dissolving the sample extracts. The buffer contains sodium azide (<0.1%; see § Precautions). The exact concentration is indicated on each vial label. Calibrators are verified to an internal reference standard.

BUF

Recovery buffer : one 30 mL vial (ready-to-use)
This buffer is used for redissolving dry extracts. It contains bovine serum albumin and sodium azide (<0.1 %; see § Precautions).

4. MATERIAL REQUIRED BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- precision micropipets (100 µL; 300 µL).
- semi-automatic pipets (300 µL; 500 µL).
- glass pipets (5 mL or 10 mL).
- glass vials with teflon or glass stoppers (from 5 to 20 mL).
- pyrex tubes for recovery of ether phase.
- anhydrous ethyl ether (analytical grade).
- evaporator (Speed Vac type) or 37°C water bath.
- vortex-type mixer.
- aspiration system.
- gamma counter set for 125 iodine.

5. PRECAUTIONS

5.1 General remarks

- The vials with calibrators should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A calibration curve must be included with each assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once.

5.2 Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use.

Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.
- This kit contains ¹²⁵I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations

5.3 Sodium azide

Some reagents contain sodium azide as preservative. Sodium azide may react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

5.4 Material of human origin

All serum and plasma samples should be handled as if capable of transmitting hepatitis or AIDS. Waste should be discarded according to the country rules.

5.5 Ethyl ether

Ethyl ether is a volatile and highly flammable organic solvent. Extraction and evaporation must be done in a ventilated hood. Avoid any contact with a flame and do not pipet reagents by mouth.

6. SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in dry tubes or in tubes with EDTA.
- Separate serum or plasma from cells by centrifugation.
- Serum and plasma samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage, maximum 2 months, keep frozen at < -20°C after aliquoting so as to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- If samples have concentrations greater than the highest standard, they must be diluted after extraction in the recovery buffer.

Serum and EDTA plasma values for 15 samples (serum values ranging from 1.36 to 7.31 ng/mL) were compared using the RIA DHEA kit. Results are as follows:

$$[\text{EDTA-plasma}] = 0.7879[\text{serum}] + 0.1898, \\ R = 0.9562$$

7. ASSAY PROCEDURE

7.1 Extraction of samples

Note : The extraction must be done in clean glass vials or tubes, pre-rinsed with ether and fitted with teflon or glass stoppers. The organic phase must be recovered carefully so as to avoid contamination with the aqueous phase.

Samples only are extracted before assay; do not extract calibrators.

Bring samples to room temperature and mix well before starting extraction.

- Number one vial for every sample.
- Place 300 µL of each sample into corresponding vial.
- Add 3 mL of ethyl ether (see § Precautions) to each vial. Stopper carefully.
- Vortex vials vigorously.
- Keep vials at -18°C until aqueous phase freezes.
- Take off carefully organic phase without contaminating with aqueous phase and place into numbered 6 mL pyrex tubes.
- Evaporate ether phase completely with either evaporator (e.g. Speed Vac type) or by placing tubes into 37°C water bath.

Note : the tubes must be firmly attached to test tube rack, since after evaporation of the ether, they will be lighter and tend to float off.

At this stage it is possible to stopper the tubes containing the dry extracts and to store them for up to 7 days in the cold (2-8°C) prior to continuing the assay.

- Re-dissolve dry ether extracts in 300 µL of resuspension buffer. Vortex vigorously (30 sec).

7.2 Assay procedure

Bring all reagents to room temperature before pipeting.

Step 1 Additions *	Step 2 Incubation	Step 3 Counting
To antibody-coated tubes, add sequentially: - 100 µL of calibrator or extract and - 500 µL of tracer. Mix.	Cover tubes. Incubate 1 hour at 37°C in water bath.	Aspirate carefully the content of tubes (except of the 2 tubes «total cpm») Determine bound cpm (B) and total cpm (T) for 1 min.

* Add 500 µL of tracer to 2 additional tubes to obtain total cpm

8. RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of DHEA concentrations in samples measured at the same time as the calibrator.

8.1 Calibration curve

The results in the package insert were calculated using a semi-logarithmic curve fit ("spline" mode) with B/T (%) or B/B₀ (%) on vertical axis and the DHEA concentration of the calibrators on the horizontal axis (ng/mL). Other data reduction methods may give slightly different results.

Total activity: 68,628 cpm				
Calibrators	DHEA (ng/mL)	cpm (n=3)	B/T (%)	B/B ₀ (%)
0	0	45,680	66.56	100
1	0.46	41,153	59.97	90.1
2	1.0	37,426	54.53	81.9
3	3.1	26,948	39.27	59.0
4	10.0	13,921	20.28	30.5
5	31.0	6,396	9.32	14.0

(Example of calibration curve, do not use for calculation)

8.2 Samples

For each sample, locate the B/T (%) or B/B₀ (%) on the vertical axis and read off the corresponding DHEA concentration on the horizontal axis.

To convert ng/mL into nmol/L multiply results by 3,47.

9. QUALITY CONTROL

Good laboratory practices imply that control samples be used regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended that their results be analyzed using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: products.support@diasource.be

10. EXPECTED VALUES

It is suggested that each laboratory establishes its own normal values. The following values obtained with healthy subjects are indicative only.

	Age (years)	n	Median (ng/mL)	Min-Max (ng/mL)	2.5th-97.5th percentile (ng/mL)
Men	All ages (20-67)	100	4.56	1.26 - 16.18	1.44 - 13.45
	20-30	25	6.89	1.85 - 16.18	2.58 - 15.19
	31-40	25	5.04	2.11 - 12.58	2.17 - 11.75
	41-50	25	4.96	1.74 - 11.57	1.87 - 10.96
	Over 51	25	2.51	1.26 - 9.18	1.30 - 7.72
Women	All ages (19-61)	100	2.75	0.43 - 17.53	0.88 - 9.50
	19-30	25	4.06	1.59 - 17.53	1.63 - 13.59
	31-40	25	2.89	1.31 - 9.88	1.39 - 8.19
	41-50	25	2.64	1.20 - 6.13	1.62 - 6.02
	Over 51	25	1.73	0.43 - 2.99	0.48 - 2.98

CHILDREN:

Age (years)	Min-Max (ng/mL)			
	girls	n	boys	n
0 - 8	< 0.3 - 2.0	84	< 0.3 - 2.2	94
8 - 10	0.5 - 3.5	9	< 0.3 - 2.5	30
10 - 12	0.5 - 5.5	21	< 0.3 - 3.5	33
12 - 14	0.6 - 9.0	26	0.9 - 6.0	22

11. PERFORMANCE CHARACTERISTICS

(for more details, see the data sheet "appendix")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

10.1 Sensitivity

10.1.1 Analytical sensitivity: 0.17 ng/ml

10.1.2 Functional sensitivity: 0.30 ng/ml

10.2 Specificity

The antibody used in the immunoassay is highly specific for dheA. Extremely low cross reactivities were obtained against other naturally occurring steroids (aldosterone, androsterone, progesterone, etc.) Or therapeutic drugs that may be present in patient samples (danazol, etc.).

10.3 Precision

10.3.1 Intra-assay

Serum samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 6.8 %.

10.3.2 Inter-assay

Serum samples were assayed in duplicate in 10 different series. Coefficients of variation were found below or equal to 7.9 %.

10.4 Accuracy

10.4.1 Dilution test

High-concentration serum samples were serially diluted after extraction with the recovery buffer. The recovery percentage obtained were between 85.2 % and 117 %.

10.4.2 Recovery test

Low-concentration serum samples were spiked with known quantities of dheA. The recovery percentages obtained were between 94.4 % and 115 %.

10.5 Measurement range (from analytical sensitivity to highest calibrator): 0.17 to approximately 30 ng/ml.

12. LIMITATIONS OF THE PROCEDURE

The non-respect of the instructions in this package insert may affect results significantly. Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Do not use hemolyzed, lipemic or icteric samples.

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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