# **Adiponectin Elisa**

***KAPME09***

# History

---

## Summary of change :

| Previous Version :                                                              | Current Version :                                                       |
|---------------------------------------------------------------------------------|-------------------------------------------------------------------------|
| 231222                                                                          | 250703                                                                  |
| <b>5.6 Sample dilution</b><br>Suppression of alternative sample dilution method | Update of the hazard warnings in the "Warnings and Precautions" section |

# Adiponectin-ELISA

Enzyme Immunoassay for the Quantitative Determination of human Adiponectin  
in serum and plasma

KAPME09

*IN VITRO DIAGNOSTIC USE*

en

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 91

| Adiponectin ELISA KAPME09    | 96 Determinations                                                                                                                                             |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Principle of the test        | Enzyme-linked Immunoassay                                                                                                                                     |
| Duration (incubation period) | 1.75 h                                                                                                                                                        |
| Antibodies                   | monoclonal antibodies                                                                                                                                         |
| Dilution Buffer              | Ready for use                                                                                                                                                 |
| Washing buffer               | 20 fold concentrate                                                                                                                                           |
| Calibrators                  | 5 single calibrators: 2 - 100µg/L,<br>native human Adiponectin                                                                                                |
| Assay Range                  | 0.27 - 31000µg/L                                                                                                                                              |
| Control                      | 2 control sera, freeze-dried                                                                                                                                  |
| Sample                       | human serum / plasma                                                                                                                                          |
| Required sample volume       | 10 µL                                                                                                                                                         |
| Sample dilution              | 1:310                                                                                                                                                         |
| Analytical sensitivity       | ≤ 0.27 µg/L                                                                                                                                                   |
| Intra- / Interassay Variance | <10 %                                                                                                                                                         |
| Reference values             | Langkamp M, Pridzun P, Böttner A, Kiess W, Thiery J, Kratzsch J, Reference<br>Intervals for Adiponectin Levels in Human Serum. 2005 DGKL Annual Meeting Jena. |

## 1. Intended Use

The Elisa KAPME09 is intended to be used for quantitative measurement of human Adiponectin in human serum and plasma samples.

## 2. Introduction

Adiponectin is a 30kDa protein which percentage in serum proteins is 0.01%. It is mainly synthesized by adipocytes, but also muscle cells and hepatocytes have the ability to synthesize Adiponectin. It consists of a Collagen-like N-terminal and a globular C-terminal domain [1]. In vivo Adiponectin appears with different oligomers. Beside the trimer and dimer also high molecular multimers exist [1-3]. Two different receptors are known, both receptors are ubiquitarily expressed, though the distribution in the tissues varies. The Adiponectin Receptor 1 (AdipoR1) is especially in muscle- and AdipoR2 in liver tissue synthesized [4].

Several studies show that adiponectin correlates negatively with BMI and thus it could have relevance for the energy metabolism for example through the regulation of fatty acid oxidation. Beside the correlation with BMI, Adiponectin level is associated with the Insulin-Resistance [5-7] and so also linked with Type II Diabetes. Adiponectin is associated also with glucose- und lipometabolism [8, 9].

Furthermore adiponectin is involved in inflammatory processes [10-14] and therewith it is of importance for appearance of arteriosclerosis [4, 5, 15] and coronaritis [16, 17], thus the determination of Adiponectin level in plasma could serve to estimate the risk of coronary disease [18, 19]. Beside this Adiponectin influences further physiological processes as for example the angiogenesis [20, 21].

## 3. Assay principle

The DIAsource ELISA for Adiponectin KAPME09 is a so-called Sandwich-Assay using two specific and highly affine antibodies. The Adiponectin in the samples binds to the first antibody coated on the microtiter plate. In the following step the second specific anti-Adiponectin-Antibody binds in turn to the immobilised Adiponectin. The second antibody is biotinylated and will be applied in a mixture with a Streptavidin-Peroxidase-Enzyme Conjugate. In the closing substrate reaction the turn of the colour will be catalysed quantitatively depending on the Adiponectin-level of the samples

#### 4. Warnings and Precautions

##### For In Vitro Diagnostic Use only. For Professional use only.

The DIAsource kit is suitable only for in vitro diagnostics and not for internal use in humans and animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. DIAsource will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided. A Material Safety Data sheet is available on request.

Do not use obviously damaged or microbial contaminated or spilled material.

**Caution: This kit contains material of human and/or animal origin. Therefore, all components and patient's specimens should be treated as potentially infectious.**

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations.

##### Human Serum

Following components contain human serum: **Controls 1 and 2 CONTROL N, and Calibrators 1-5 CAL N.**

Source human serum for the controls provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV). No known methods can offer total security of the absence of infectious agents; therefore, all components and patient's specimens should be treated as potentially infectious.

#### 4.1 Components

All components of the product that contain hazardous substances in declarable concentrations in accordance with European Regulation 1272/2008 are listed below:

##### Dilution Buffer *DIL BUF* Antibody-Conjugate *Ab BIOT SAV HRP*

The buffers contain the hazardous ingredients listed below, based on this the following labeling of the two mixtures results:

Hazard pictograms (CLP):

GHS07

Signal word (CLP):

Warning

Hazardous components:

2-methylisothiazol-3(2H)-one (< 0,06 %)  
reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) (< 0,015 %)

Hazard statements (CLP):

H317 - May cause an allergic skin reaction.

Safety instructions (CLP):

P280 - Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.  
P272 - Contaminated work clothing should not be allowed out of the workplace.  
P261 - Avoid breathing dust/fume/gas/mist/vapours/spray.  
P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.  
P302+P352 - IF ON SKIN: Wash with plenty of Water.  
P501 - Dispose of contents/container in accordance with national disposal regulations.

##### Washing Buffer *WASH SOLN CONC*

The buffer contains a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one (< 0.015 %) as a preservative.

EUH208 Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1).

May produce an allergic reaction.

EUH210 Safety data sheet available on request.

##### Stop Solution *STOP SOLN*

The stop solution contains 0.2 M sulphuric acid.

EUH210 Safety data sheet available on request.

#### 4.2 Protective equipment and first aid measures

After contact with skin: Take off immediately all contaminated clothing. Rinse skin with water or shower. In case of skin reactions, consult a physician.

After contact with eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. In case of eye irritation consult an ophthalmologist.

After ingestion: Call a doctor if you feel unwell.

Avoid contact with the ingredients. Detailed information on personal protective equipment and first aid measures can be found in the safety data sheet for Adiponectin ELISA. This can be accessed electronically via the DIAsource website at <https://diasource-diagnostics.com/> or will be made available on request via [products.support@diasource.be](mailto:products.support@diasource.be).

## 5. SAMPLES

### 5.1 Sample type

Serum and Plasma

Serum and heparin plasma levels are comparable. In EDTA- and Citrate Plasma-samples levels were found approx. 18% lower, because of the relatively high amount of anticoagulant.

### 5.2 Specimen collection

The blood sample for serum preparation should be gained according to standardized venipuncture procedure. Hemolytic reactions have to be avoided.

### 5.3 Required sample volume: 10 µL

### 5.4 Sample stability

In firmly closable sample vials

- Storage at 20-25°C: 2 days
- Storage at -20° C: min. 2 years
- Freeze-thaw cycles max. 3

The storage of samples over a period of 2 years at -20°C showed no influence on the reading. Freezing and thawing of samples should be minimized, 3 freeze-thaw cycles showed no effect on the measured adiponectin concentration.

### 5.5 Interference

Hemoglobin, triglyceride and bilirubin in the sample do not interfere to a concentration of **5 mg/mL**, **100 mg/mL** and **100 µg/mL** respectively. However, the use of hemolytic, lipemic or icteric samples should be validated by the user.

### 5.6 Sample dilution

- Dilution: **1:310** with Dilution Buffer
- Dilute for example **300 µL** Dilution Buffer in PE-/PP-Tubes (application of a multi-stepper is recommended in larger series), add **10 µL** Serum- or Plasma (dilution: 1:31). Add **900 µL** Dilution Buffer in another PE-/PP-tube and **100 µL** of the thoroughly mixed first dilution. After mixing, use 2×100 µL from this **1:310** diluted sample in the assay.
- Sample stability after dilution of the sample: maximum 1 hour at 20-25°C.

## 6. MATERIAL

### 6.1. Material provided

The reagents listed below are sufficient for 96 wells including the calibration curve

|       |                                                                                              |      |      |                                                                                                                                                       |      |                                                                                                                                         |     |                                                                                                                                                                  |
|-------|----------------------------------------------------------------------------------------------|------|------|-------------------------------------------------------------------------------------------------------------------------------------------------------|------|-----------------------------------------------------------------------------------------------------------------------------------------|-----|------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1)    |             |      |      | <b>Microtiter plate</b> , ready for use, coated with mouse-anti-Adiponectin-antibody. Wells are separately breakable (8x12) wells                     |      |                                                                                                                                         |     |                                                                                                                                                                  |
| 2)    | CAL                                                                                          | N    |      | <b>Cal N</b> , lyophilized (native human adiponectin), concentrations are given on vial labels and on quality certificate in ng/ml. <b>5 x 750 µl</b> |      |                                                                                                                                         |     |                                                                                                                                                                  |
| 3)    | CONTROL                                                                                      |      |      | <b>Control Serum 1 &amp; 2</b> , lyophilized (human serum), concentration is given on quality certificate in ng/ml. <b>2 x 500 µl</b>                 |      |                                                                                                                                         |     |                                                                                                                                                                  |
| 4)    | <table border="1"><tr><td>Ab</td><td>BIOT</td></tr><tr><td>SAV</td><td>HRP</td></tr></table> |      |      | Ab                                                                                                                                                    | BIOT | SAV                                                                                                                                     | HRP | <b>Antibody-Conjugate, 12 ml, ready for use</b> , mouse-anti-hAdiponectin-antibody biotinylated + streptavidin horseradish peroxidase conjugat. <b>1 x 12 ml</b> |
| Ab    | BIOT                                                                                         |      |      |                                                                                                                                                       |      |                                                                                                                                         |     |                                                                                                                                                                  |
| SAV   | HRP                                                                                          |      |      |                                                                                                                                                       |      |                                                                                                                                         |     |                                                                                                                                                                  |
| 5)    | DIL                                                                                          | BUF  |      | <b>Dilution buffer</b> , ready for use. <b>1 x 125 ml</b>                                                                                             |      |                                                                                                                                         |     |                                                                                                                                                                  |
| 6)    | WASH                                                                                         | SOLN | CONC | <b>Washing Buffer</b> , 20-fold concentrated solution. <b>1 x 50 ml</b>                                                                               |      |                                                                                                                                         |     |                                                                                                                                                                  |
| 7)    | <table border="1"><tr><td>CHROM</td><td>TMB</td></tr></table>                                |      |      | CHROM                                                                                                                                                 | TMB  | <b>Chromogenic Substrate</b> , ready for use, horseradish-peroxidase-(HRP)-substrate, stabilised Tetramethylbencidine. <b>1 x 12 ml</b> |     |                                                                                                                                                                  |
| CHROM | TMB                                                                                          |      |      |                                                                                                                                                       |      |                                                                                                                                         |     |                                                                                                                                                                  |
| 8)    | <table border="1"><tr><td>STOP</td><td>SOLN</td></tr></table>                                |      |      | STOP                                                                                                                                                  | SOLN | <b>Stopping Solution</b> , ready for use, 0.2 M sulphuric acid , <b>1 x 12 ml</b>                                                       |     |                                                                                                                                                                  |
| STOP  | SOLN                                                                                         |      |      |                                                                                                                                                       |      |                                                                                                                                         |     |                                                                                                                                                                  |
| 9)    |                                                                                              |      |      | Sealing tape for covering of the microtiter plate, 2 x, adhesive.                                                                                     |      |                                                                                                                                         |     |                                                                                                                                                                  |

### 6.2. Materials Required but not provided

Distilled (Aqua destillata) or deionized water for dilution of the Washing Buffer (950 mL).

Precision pipettes and multichannel pipettes with disposable plastic tips

Polyethylene PE/Polypropylene PP tubes for dilution of samples

Vortex-mixer

Microtiter plate shaker (350 rpm)

Microtiter plate washer (recommended)

Micro plate reader ("ELISA-Reader") with filter for 450 and ≥590 nm

## 7. Technical Notes

### Storage Conditions

Store the kit at 2-8°C after receipt until its expiry date. The lyophilized reagents should be stored at –20 °C after reconstitution. Avoid repeated thawing and freezing.

### Storage Life

The shelf life of the components **after initial opening** is warranted for **4 weeks**, store the unused strips and microtiter wells **airtight** together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided. The **reconstituted components** Calibrators and Control Sera **1&2** must be stored at –20°C (max. 4 weeks). For further use, thaw quickly but gently (avoid temperature increase above room temperature and avoid excessive vortexing). Up to 3 of the freeze-thaw cycles did not influence the assay. The 1:20 diluted Washing Buffer is 4 weeks stable at 2-8°C

### Preparation of reagents

Bring all reagents to room temperature (20 - 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming. Reagents with different lot numbers cannot be mixed.

### Reconstitution

The Calibrators and Control Sera **1&2** are reconstituted with the Dilution Buffer . It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

**Dilution1&2** with the Dilution Buffer in the same ratio (1:310) as the sample.

The required volume of Washing Buffer is prepared by 1:20 dilution of the provided 20fold concentrate with Aqua dest.

### Assay Procedure

When performing the assay, Blank, Calibrators Controls **1&2** and the samples should be pipette as fast as possible (e.g. <15 minutes). To avoid distortions due to differences in incubation times, Antibody Conjugate as well as the succeeding Substrate Solution should be added to the plate in the same order and in the same time interval as the samples. Stopping Solution should be added to the plate in the same order as Substrate Solution .

All determinations (Blank, Calibrators, Controls **1&2** and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

### Incubation

**Incubation at room temperature means: Incubation at 20 - 25°C.** The Substrate Solution , stabilised H<sub>2</sub>O<sub>2</sub>-Tetramethylbencidine, is photosensitive–store and incubation in the dark.

### Shaking

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtiter plate shaker. We are recommending 350 rpm. Due to certain technical differences deviations may occur, in case the rotation frequency must be adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/ or false values, excessive shaking may result in high optical densities and/ or false values.

### Washing

Proper washing is of basic **importance** for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

The danger of handling with potentially infectious material must be taken into account.

When using an **automatic microtiter** plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

**Manual washing** is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamical swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

### 8. Assay Procedure

| Preparation of reagents                                                        |                                                                                                      | Reconstitution:                                     | Dilution                                 |
|--------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|-----------------------------------------------------|------------------------------------------|
| Calibrators                                                                    |                                                                                                      | in <b>750 µL</b> Dilution Buffer                    | -                                        |
| Control Sera                                                                   |                                                                                                      | in <b>500 µL</b> Dilution Buffer                    | <b>1:310</b> with <b>Dilution buffer</b> |
| Washing Buffer                                                                 |                                                                                                      | -                                                   | <b>1:20</b> with <b>Aqua dest.</b>       |
| Sample dilution: with Dilution Buffer <b>1:310</b>                             |                                                                                                      |                                                     |                                          |
| Before assay procedure bring all reagents to room temperature <b>20-25°C</b> . |                                                                                                      |                                                     |                                          |
| Assay Procedure in Double Determination:                                       |                                                                                                      |                                                     |                                          |
| Pipette                                                                        | Reagents                                                                                             | Position                                            |                                          |
| 100 µL                                                                         | Dilution Buffer ( <b>Blank</b> )                                                                     | A1/A2                                               |                                          |
| 100 µL                                                                         | Calibrator <b>1 (2 ng/mL)</b>                                                                        | B1/B2                                               |                                          |
| 100 µL                                                                         | Calibrator <b>2 (10 ng/mL)</b>                                                                       | C1/C2                                               |                                          |
| 100 µL                                                                         | Calibrator <b>3 (30 ng/mL)</b>                                                                       | D1/D2                                               |                                          |
| 100 µL                                                                         | Calibrator <b>4 (70 ng/mL)</b>                                                                       | E1/E2                                               |                                          |
| 100 µL                                                                         | Calibartor <b>5 (100 ng/mL)</b>                                                                      | F1/F2                                               |                                          |
| 100 µL                                                                         | Control Serum <b>1</b> (1:310 diluted)                                                               | G1/G2                                               |                                          |
| 100 µL                                                                         | Control Serum <b>2</b> (1:310 diluted)                                                               | H1/H2                                               |                                          |
| 100 µL                                                                         | Sample (1:310 diluted)                                                                               | in the rest of the wells according the requirements |                                          |
| Cover the wells with the sealing tape.                                         |                                                                                                      |                                                     |                                          |
| Sample Incubation: <b>1 h at 20-25°C, 350 rpm</b>                              |                                                                                                      |                                                     |                                          |
| 3 x 300 µL                                                                     | Aspirate the contents of the wells and <b>wash</b> 3 x with 300 µL each Washing Buffer / <b>well</b> |                                                     | In each well                             |
| 100 µL                                                                         | Antibody-POD-Conjugate                                                                               |                                                     | In each well                             |
| Cover the wells with the sealing tape.                                         |                                                                                                      |                                                     |                                          |

|                                                      |                                                                                                      |              |
|------------------------------------------------------|------------------------------------------------------------------------------------------------------|--------------|
| <b>Incubation: 30 Minutes at 20-25°C, 350 rpm</b>    |                                                                                                      |              |
| 3 x 300 µL                                           | Aspirate the contents of the wells and <b>wash</b> 3 x with 300 µL each Washing Buffer / <b>well</b> | In each well |
| 100 µL                                               | Substrate Solution                                                                                   | In each well |
| <b>Incubation: 15 Minutes in the Dark at 20-25°C</b> |                                                                                                      |              |
| 100 µL                                               | Stopping Solution                                                                                    | In each well |
|                                                      | Measure the absorbance within 30 min at <b>450 nm</b> with $\geq 590$ nm as reference wavelength.    |              |

## 9. Quality Control

Good laboratory practice requires that controls are included in each assay. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance. The test results are only valid if the test has been performed following the instructions. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state or local standards/laws. All calibrators and kit controls must be found within the acceptable ranges as stated on the QC Certificate. If the criteria are not met, the run is not valid and should be repeated. Each laboratory should use known samples as further controls.

### Quality criteria

For the evaluation of the assay it is required that the absorbance values of the blank should be below 0.25, and the absorbance of calibrator 5 should be above 1.00.

Samples, which yield higher absorbance values than calibrator 5, should be re-tested with a higher dilution.

## 10. Evaluation of results

### Establishing of the calibration curve


| Calibrator | 1 | 2  | 3  | 4  | 5   |
|------------|---|----|----|----|-----|
| ng/mL      | 2 | 10 | 30 | 70 | 100 |

- 1) Calculate the **mean absorbance** value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the blank from the mean absorbances of all other samples and calibrators
- 3) Plot the calibration concentrations on the x-axis versus the mean value of the absorbance of the calibrators on the y-axis.
- 4) Recommendation: Calculation of the calibration curve should be done by using a computer program, because the curve is in general (without respective transformation) not ideally described by linear regression. **A higher-grade polynomial, or four parametric logistic (4-PL) curve fit or non-linear regression** are usually suitable for the evaluation (as might be spline or point-to-point alignment in individual cases).
- 5) The Adiponectin concentration in ng/mL of the samples and controls **1** and **2** can be calculated by **multiplication** with the respective **dilution factor**.

### Example of a typical calibration curve

The following data is for demonstration only and cannot be used in place of data generation at the time of assay.

|                            | Blank | 1     | 2     | 3     | 4     | 5     |
|----------------------------|-------|-------|-------|-------|-------|-------|
| ng/mL                      | 0.0   | 2     | 10    | 30    | 70    | 100   |
| OD <sub>(450-620 nm)</sub> | 0.008 | 0.071 | 0.357 | 1.022 | 2.053 | 2.817 |

**Figure 1** Exemplary calibration curve

The exemplary shown calibration curve in Figure 1 **cannot** be used for calculation of your test results. You have to establish a calibration curve for each test you conduct!

#### Exemplary calculation of Adiponectin concentrations

Sample dilution: 1:310

|                                    |       |
|------------------------------------|-------|
| Measured extinction of your sample | 0.408 |
| Measured extinction of the blank   | 0.008 |

Your measurement program will calculate the Adiponectin concentration of the diluted sample automatically by using the difference of sample and blank for the calculation. You only have to determine the most suitable curve fit (here: polynomial 2<sup>nd</sup> degree).

In this exemplary case the following equation is solved by the program to calculate the adiponectin concentration in the sample:

$$\begin{aligned} 0.400 &= -6 \times 10^{-5}x^2 + 0.0344x \\ 11.89 &= x \end{aligned}$$

If the dilution factor (1:310) is taken into account the adiponectin concentration of the undiluted sample is

$$11.89 \times 310 = 3685.9 \text{ ng/mL} = 3.69 \text{ } \mu\text{g/mL}$$

The exemplary shown calibration curve in Fig.1 cannot be used for calculation of your test results. You have to establish a calibration curve for each test you conduct!

#### Interpretation of results

The test results should not be the only base for therapeutic decisions. The results should be interpreted in regard to anamnesis, further clinical observations and results of other diagnostic investigations. Further, it is recommended to establish reference and cut-off values corresponding to the relevant group of patients for each laboratory.

#### Limitation of procedure

The sensitive human Adiponectin ELISA, KAPME09 is based on monoclonal antibodies. Generally, this technique could be sensible to heterophilic antibodies or rheumatic factors in the sample. Their influence is reduced by assay design, but cannot be excluded completely.

#### 11. Reference values

The expected values for serum adiponectin, which were determined with the DIASource ELISA KAPME09 in healthy donors and analysed by Prof. Dr. J. Kratzsch, Department of Laboratory Medicine, University Hospital Leipzig, are given below (Tab. 1). These data show significant correlation between Adiponectin-Serum values and age as well as gender of the probands, in turn the correlation between the respective BMI seems to be less significant. In the samples of neonatal cord blood very high values were found. Several different statistical analyses were performed to adapt for certain individual demands. The best suited data can be chosen respectively for interpretation of the own measurements.

**Table 1a Expected values for adults**, gender specific mean as well as median, 5. and 95. percentiles are given.

| Sex    | Number [n] | Mean [µg/ml] | Median [µg/ml] | Standard Deviation | 5 <sup>th</sup> Percentile [µg/ml] | 95 <sup>th</sup> Percentile [µg/ml] |
|--------|------------|--------------|----------------|--------------------|------------------------------------|-------------------------------------|
| Female | 101        | <b>10.2</b>  | 9.1            | 4.6                | 4.0                                | 19.4                                |
| Male   | 125        | <b>6.8</b>   | 6.1            | 4.1                | 2.0                                | 13.9                                |
| total  | 226        | <b>8.3</b>   | 7.5            | 4.6                | 2.4                                | 19.3                                |

**Table 1b Expected values for children**, gender specific mean as well as median, 5. and 95. percentiles are given.

| Sex    | Number [n] | Mean [µg/ml] | Median [µg/ml] | Standard Deviation | 5 <sup>th</sup> Percentile [µg/ml] | 95 <sup>th</sup> Percentile [µg/ml] |
|--------|------------|--------------|----------------|--------------------|------------------------------------|-------------------------------------|
| Female | 131        | <b>8.71</b>  | 8.18           | 4.32               | 3.05                               | 15.6                                |
| Male   | 134        | <b>8.97</b>  | 8.12           | 5.13               | 3.36                               | 18.6                                |
| total  | 265        | <b>8.84</b>  | 8.18           | 4.74               | 3.33                               | 16.5                                |

**Table 1c Expected values for Adiponectin**, age specific mean as well as median, 5. and 95. percentiles are given.

| Age group [a] | Number [n] | Mean [µg/ml] | Median [µg/ml] | 5. Percentile [µg/ml] | 95. Percentile [µg/ml] |
|---------------|------------|--------------|----------------|-----------------------|------------------------|
| < 7.99        | 46         | <b>12.82</b> | 11.7           | 2.33                  | 26.5                   |
| 8 – 9.99      | 40         | <b>8</b>     | 8.09           | 3.96                  | 14.9                   |
| 10-11.99      | 55         | <b>8.02</b>  | 7.14           | 3.36                  | 13.8                   |
| 12 – 13.99    | 26         | <b>8.21</b>  | 7.54           | 4.5                   | 13.2                   |
| 14 – 15.99    | 59         | <b>8.12</b>  | 8.09           | 3.67                  | 13.7                   |
| 16 – 19.99    | 41         | <b>7.97</b>  | 7.79           | 2.74                  | 13.3                   |
| alle          | 267        | <b>8.88</b>  | 8.18           | 3.33                  | 16.7                   |
| 20 – 29.99    | 47         | <b>6.72</b>  | 6.38           | 2.50                  | 12.25                  |
| 30 – 39.99    | 38         | <b>7.38</b>  | 6.69           | 1.98                  | 19.29                  |
| 40 – 49.99    | 55         | <b>8.42</b>  | 8.20           | 2.41                  | 17.85                  |
| 50 – 59.99    | 47         | <b>9.61</b>  | 8.55           | 2.15                  | 19.85                  |
| > 60          | 32         | <b>9.52</b>  | 8.57           | 3.00                  | 21.10                  |
| alle          | 226        | <b>8.33</b>  | 7.5            | 2.41                  | 19.29                  |

**Table 1d Expected values for Adiponectin**, age as well as gender specific mean and median, BMI and 25. And 75. percentiles are given.

| Female             |    |              | Adiponectin (µg/ml): |          |                      |              |
|--------------------|----|--------------|----------------------|----------|----------------------|--------------|
| Age (Years):       | n: | BMI: AV ± SD | AV ± SD::            | Median : | Percentile: 25.- 75. | Min. – Max.: |
| Newborn Cord blood | 19 | --           | <b>29.80 ± 12.49</b> | 26.1     | 19.5 -35.2           | 16.9 - 61.4  |
| < 3.99             | 9  | 15.73 ± 0.79 | <b>14.43 ± 7.76</b>  | 11.2     | 8.2 - 21.8           | 2.3 - 26.7   |
| 4.0 - 7.99         | 11 | 16.01 ± 1.94 | <b>8.46 ± 4.73</b>   | 9.3      | 2.9 - 12.1           | 1.4 - 15.6   |
| 8.0 - 9.99         | 22 | 17.58 ± 3.84 | <b>7.92 ± 3.00</b>   | 8.2      | 5.2 - 10.0           | 3.6 - 15.1   |
| 10.0 - 11.99       | 33 | 17.83 ± 1.86 | <b>7.66 ± 4.59</b>   | 6.6      | 5.0 - 8.8            | 3.1 - 20.9   |
| 12.0 - 13.99       | 11 | 19.85 ± 2.31 | <b>8.22 ± 5.64</b>   | 7.5      | 6.5 - 9.2            | 4.9 - 13.2   |
| 14.0 - 15.99       | 27 | 19.91 ± 1.72 | <b>8.83 ± 9.25</b>   | 8.9      | 5.2 - 11.8           | 2.6 - 17.7   |
| 16.0 – 19.99       | 18 | 21.64 ± 2.64 | <b>9.00 ± 3.22</b>   | 8.7      | 6.9 - 11.2           | 2.7 - 14.0   |
| 20.0 - 29.99       | 24 | 23.12 ± 5.01 | <b>7.39 ± 3.35</b>   | 7.3      | 5.7 - 9.0            | 3.4 - 17.8   |
| 30.0 - 39.99       | 17 | 23.20 ± 2.86 | <b>9.19 ± 3.89</b>   | 8.6      | 7.2 - 10.4           | 3.6 - 19.3   |
| 40.0 - 49.99       | 26 | 24.50 ± 4.11 | <b>9.93 ± 3.59</b>   | 9.5      | 7.5 - 11.6           | 4.4 - 19.6   |
| 50.0 - 59.99       | 21 | 24.61 ± 3.31 | <b>11.5 ± 5.49</b>   | 10.0     | 8.0 - 15.9           | 2.0 - 23.1   |
| >60.0              | 8  | 24.63 ± 1.89 | <b>15.6 ± 4.64</b>   | 15.3     | 11.4 - 18.2          | 11.2 - 24.1  |

| Male               |    |              | Adiponectin (µg/ml): |          |                      |              |
|--------------------|----|--------------|----------------------|----------|----------------------|--------------|
| Age (Years):       | n: | BMI: AV ± SD | AV ± SD:             | Median : | Percentile: 25.- 75. | Min. – Max.: |
| Newborn Cord blood | 10 | --           | 27.80 ± 7.68         | 26.7     | 22.2 -31.0           | 15.6 - 40.6  |
| < 3.99             | 14 | 16.17 ± 1.81 | 16.57 ± 6.55         | 14.3     | 11.6 - 21.2          | 5.8 - 40.3   |
| 4.0 - 7.99         | 12 | 15.69 ± 1.05 | 11.24 ± 5.43         | 9.7      | 8.9 - 15.9           | 3.5 - 18.6   |
| 8.0 - 9.99         | 18 | 16.45 ± 1.76 | 8.11 ± 2.93          | 7.6      | 6.2 - 9.1            | 5.00 - 15.4  |
| 10.0 - 11.99       | 21 | 18.34 ± 2.18 | 8.43 ± 3.91          | 7.8      | 5.2 - 10.9           | 3.4 - 20.2   |
| 12.0 - 13.99       | 14 | 18.61 ± 2.11 | 7.59 ± 2.86          | 7.1      | 6.0 - 10.3           | 2.4 - 12.2   |
| 14.0 - 15.99       | 32 | 19.86 ± 2.00 | 7.53 ± 2.52          | 7.4      | 5.1 - 9.3            | 3.8 - 15.4   |
| 16.0 -19.99        | 23 | 22.03 ± 2.42 | 7.16 ± 3.53          | 6.9      | 4.2 - 9.6            | 2.0 - 13.9   |
| 20.0 - 29.99       | 23 | 23.43 ± 2.48 | 5.44 ± 2.29          | 5.8      | 4.0 - 6.9            | 1.3 - 10.3   |
| 30.0 - 39.99       | 21 | 23.33 ± 2.72 | 5.92 ± 4.60          | 4.4      | 2.7 - 6.7            | 1.9 - 20.6   |
| 40.0 - 49.99       | 22 | 23.79 ± 2.41 | 6.13 ± 2.92          | 5.5      | 3.8 - 8.3            | 2.1 - 11.6   |
| 50.0 - 59.99       | 23 | 26.68 ± 2.77 | 7.45 ± 4.50          | 6.7      | 5.0 - 8.8            | 1.4 - 19.6   |
| >60.0              | 24 | 25.72 ± 2.12 | 7.48 ± 3.92          | 7.6      | 4.6 - 9.2            | 3.0 - 21.1   |

n= Number of Probands AV=Average Value BMI=Body Mass Index (kg/m<sup>2</sup>) SD=Standard Deviation


## 12. Performance Characteristics

### Sensitivity

The analytical sensitivity of the ELISA KAPME09 was measured by the variability of the signal of the calibrator 0. Based on the twofold standard deviation of the blank the mean analytical sensitivity is < 0.27 ng/mL (Range 0.094 to 0.59 ng/mL).

### Specificity

Adiponectin exists in different oligomeric forms: the high, medium and low molecular weight form. Different numbers of the adiponectin monomer aggregate specifically to form a complex. In Figure 2a the five different forms of human adiponectin are shown schematically. In parallel the results of a size-exclusion chromatography of human serum measured with the DIASource KAPME09 Adiponectin ELISA are shown.

**Figure 2 Adiponectin Structure** (A) proteinaceous structure of human adiponectin including posttranslational modifications and multimeric forms (taken from Nakano et al<sup>9</sup>). (B) Results of a serum separation by size-exclusion chromatography. The sample was fractionated and the adiponectin content of each fraction was measured by DIASource KAPME09 Adiponectin, for comparison the corresponding size in kDa is indicated by black arrows.

The results shown in Figure 2b clearly demonstrate that the DIASource KAPME09 Adiponectin ELISA detects all forms of Adiponectin present in human serum: the trimer at 65 kDa, the hexamer at 150 kDa and the high molecular weight forms of >280 kDa. The DIASource KAPME09 Adiponectin ELISA therefore measures total Adiponectin.

<sup>9</sup> Nakano Y, Tahima S, Yoshimi A, Akiyama H, Tsushima M, Tanioka T, Negoro T, Tomita M, Tobe T: A novel enzyme-linked immunosorbent assay specific for high-molecular-weight adiponectin. J Lipid Res 2006, 47(7):1572-1582

## Precision

### Intra-Assay Variance & Accuracy

Intra assay variance and accuracy is exemplarily shown with two samples. The adiponectin concentration of these samples was repeatedly measured in one assay.

**Table 2 Intra-Assay Variation.** Recombinant adiponectin was diluted in dilution buffer and the adiponectin concentration of the dilution was measured repeatedly within one assay.

|          | Determinations [n] | Mean value [µg/L] | Standard deviation [µg/L] | VC [%] | Target Value [µg/L] |
|----------|--------------------|-------------------|---------------------------|--------|---------------------|
| Sample 1 | 8                  | 7.108             | 0.22                      | 3.14   | 6                   |
| Sample 2 | 8                  | 107.96            | 3.97                      | 3.67   | 100                 |

In both samples the variability is less than 5% and the deviation from the target value is <20%.

### Inter-Assay Variance

Serum samples were repeatedly measured in independent assays of different lots. On average the coefficient of variation was 7.5% (SD 1.6). The results of 5 samples are shown in table 3.

**Table 3 Inter-Assay Variation**

|              | Sample 1 | Sample 2 | Sample 3 | Sample 4 | Sample 5 |
|--------------|----------|----------|----------|----------|----------|
| Mean [µg/mL] | 4.72     | 5.25     | 8.36     | 5.45     | 22.29    |
| CV [%]       | 8.16     | 8.14     | 6.93     | 8.05     | 7.30     |
| n            | 99       | 68       | 62       | 174      | 62       |


### Linearity

Linearity of sample dilution was tested by serial dilution (1:100 – 1:4000) of human serum samples and recalculation of the adiponectin content in comparison to the mean adiponectin concentration of all dilutions (Table 4), no diluted sample showed a deviation of >30%.

**Table 4 Linearity.** Human serum sample were diluted in dilution buffer and adiponectin content was recalculated. Measurements results are shown in [mg/L]. No deviation of the mean >30% was detected.

| µg/mL    | Mean  | 1:100 | 1:200  | 1:400  | 1:800  | 1:1000 | 1:2000 | 1:4000 |
|----------|-------|-------|--------|--------|--------|--------|--------|--------|
| Sample 1 | 5.76  | 6.53  | 6.331  | 5.764  | 5.49   | 6.067  | 6.114  | 4.056  |
| Sample 2 | 11.53 | 10.93 | 12.107 | 11.395 | 11.454 | 11.567 | 12.884 | 10.362 |
| Sample 3 | 12.07 | 13.57 | 12.853 | 12.03  | 11.974 | 11.338 | 11.548 | 11.169 |
| Sample 4 | 4.89  | 4.659 | 4.886  | 4.384  | 4.425  | 5.851  | 5.13   | n/a    |

Additionally dilutions of 1:1500 to 1:96000 were evaluated with two samples. In Figure 3 the results are shown and demonstrate that in the tested samples no effect of dilution could be detected on measured adiponectin concentrations. The deviation of the target concentration of each dilution was less than 30% in all dilutions.

**Figure 3 Linearity.** Adiponectin concentration was measured in two human samples diluted 1:1500 to 1:96000.

### Recovery

Trueness and traceability of the DIAsource Adiponectin ELISA KAPME09 was evaluated by recovery of recombinant adiponectin in human serum. The recovery of recombinant Adiponectin yielded in a serum matrix on average 110%.

**Table 5 Recovery of recombinant human Adiponectin in Serum.** Recombinant Adiponectin was added in different amounts to human serum. The Adiponectin content of the so enriched samples was measured and recovery in comparison to enriched buffer calculated.

| R&D recombinant Adiponectin lot 1022911 | Dilution buffer | Serum   | Recovery   |
|-----------------------------------------|-----------------|---------|------------|
| ng/mL                                   | ng/mL           | ng/mL   | %          |
| 0                                       | 0               | 0.00778 | ---        |
| 75                                      | 87.92           | 95.71   | <b>109</b> |
| 37.5                                    | 51.84           | 59.98   | <b>116</b> |
| 18.75                                   | 26.55           | 26.7    | <b>101</b> |
| 9.375                                   | 13.35           | 15.15   | <b>113</b> |

#### Interference

Interference of bilirubin and triglycerides was tested by adding different amounts of these substances to human serum containing adiponectin. For comparison the same amount of buffer without any substance was also added to the serum. Table 6 demonstrates that neither bilirubin nor triglycerides or haemoglobin exert any influence on the measurement of adiponectin in human serum.


**Table 6** Interference. Serum samples were enriched with different amount of triglycerides, bilirubin or haemoglobin. The relative amount of adiponectin measured in comparison with native serum is shown here [%].

| Triglyceride<br>100 mg/mL | Bilirubin<br>100 µg/mL | Hemoglobin<br>1 mg/mL | Hemoglobin<br>5 mg/mL |
|---------------------------|------------------------|-----------------------|-----------------------|
| 94                        | 96                     | 90                    | 109                   |
| 90                        | 93                     | 97                    | --                    |
| 95                        | 94                     | 93                    | --                    |

### 13. COMPARISON STUDIES

DIASource Adiponectin, KAPME09 was compared with two different, commercially available test systems. In linear regression analysis both tests showed a good coefficient of determination ( $R^2=0.896$  and  $R^2=0.97$ ). Thus, the comparability of the results between the DIASource KAPME09 and other test systems is possible. In dependence on the respective system the absolute deviation of measurement results is different, but because of the excellent correlation the results are comparable after applying a factor. The results of both studies are shown in figure 4.

**Figure 4** Comparison of DIASource KAPME09 with commercially available test systems (left) radioimmunoassay (n=14) and (right) ELISA (n=84)

### 14. Instructions for use for scientific application

#### Scientific Application

In addition to serum and plasma samples Adiponectin can be determined in other human body fluids and in cell culture supernatants of human cell lines for research purposes.

#### Samples suitable for scientific application

**Serum, plasma, saliva, urine, breast milk and cell culture supernatant** of human cell lines

The recommended dilution for serum and plasma samples in Dilution Buffer : **(1:310)**.

In the other samples, the Adiponectin levels can vary considerable, the optimal dilution must be found out by the customer.

**Table 7** Results of sample matrix tests . Adiponectin was added to the respectively diluted samples. Enriched samples were measured without further dilution. Shown is the relative recovery of added Adiponectin of the concentration measured in respectively enriched assay buffer.

| Matrix Dilution | Urine | Saliva | Breast Milk | Cell Culture Medium<br>10% FCS | Cell Culture Medium |
|-----------------|-------|--------|-------------|--------------------------------|---------------------|
| 1:2             | 80    | 87     | 89          | 83                             | 95                  |
| 1:5             | 95    | 80     | 92          | 92                             | 97                  |
| 1:10            | 92    | 87     | -           | 101                            | 85                  |
| 1:20            | 94    | 99     | -           | 83                             | 91                  |

#### Species Cross-Reactivity

Serum of the cited species was diluted and used as sample in this assay system. No signal was detected in serum of the following species:

Horse, Cow, Chicken, Rabbit, Dog, Guinea pig, Sheep, Mouse, Goat, Donkey, Rat, Cat

Whether this obvious non-reactivity is species specific should be assessed individually by each customer.

#### 15. Literatur/References

- Nakano, Y., et al., *Isolation and characterization of GBP28, a novel gelatin-binding protein purified from human plasma*. J Biochem (Tokyo), 1996. **120**(4): p. 803-12.
- Pajvani, U.B., et al., *Structure-function studies of the adipocyte-secreted hormone Acrp30/adiponectin. Implications for metabolic regulation and bioactivity*. J Biol Chem, 2003. **278**(11): p. 9073-85.
- Tsao, T.S., et al., *Role of disulfide bonds in Acrp30/adiponectin structure and signaling specificity. Different oligomers activate different signal transduction pathways*. J Biol Chem, 2003. **278**(50): p. 50810-7.
- Shimada, K., T. Miyazaki, and H. Daida, *Adiponectin and atherosclerotic disease*. Clin Chim Acta, 2004. **344**(1-2): p. 1-12.
- Higashiura, K., et al., *Correlations of adiponectin level with insulin resistance and atherosclerosis in Japanese male populations*. Clin Endocrinol (Oxf), 2004. **61**(6): p. 753-9.
- Spranger, J., et al., *Adiponectin is independently associated with insulin sensitivity in women with polycystic ovary syndrome*. Clin Endocrinol (Oxf), 2004. **61**(6): p. 738-46.
- Zoico, E., et al., *Adipocytokines, fat distribution, and insulin resistance in elderly men and women*. J Gerontol A Biol Sci Med Sci, 2004. **59**(9): p. M935-9.
- Ye, J.M., et al., *PPARalpha /gamma ragaglitazar eliminates fatty liver and enhances insulin action in fat-fed rats in the absence of hepatomegaly*. Am J Physiol Endocrinol Metab, 2003. **284**(3): p. E531-40.
- Yamauchi, T., et al., *Adiponectin stimulates glucose utilization and fatty-acid oxidation by activating AMP-activated protein kinase*. Nat Med, 2002. **8**(11): p. 1288-95.
- Delaigle, A.M., et al., *Induction of adiponectin in skeletal muscle by inflammatory cytokines: in vivo and in vitro studies*. Endocrinology, 2004. **145**(12): p. 5589-97.
- Winzer, C., et al., *Plasma adiponectin, insulin sensitivity, and subclinical inflammation in women with prior gestational diabetes mellitus*. Diabetes Care, 2004. **27**(7): p. 1721-7.
- Xydakis, A.M., et al., *Adiponectin, inflammation, and the expression of the metabolic syndrome in obese individuals: the impact of rapid weight loss through caloric restriction*. J Clin Endocrinol Metab, 2004. **89**(6): p. 2697-703.
- Motoshima, H., et al., *Adiponectin suppresses proliferation and superoxide generation and enhances eNOS activity in endothelial cells treated with oxidized LDL*. Biochem Biophys Res Commun, 2004. **315**(2): p. 264-71.
- Wolf, A.M., et al., *Adiponectin induces the anti-inflammatory cytokines IL-10 and IL-1RA in human leukocytes*. Biochem Biophys Res Commun, 2004. **323**(2): p. 630-5.
- Okamoto, Y., et al., *Adiponectin reduces atherosclerosis in apolipoprotein E-deficient mice*. Circulation, 2002. **106**(22): p. 2767-70.
- Schlegel, A., *Adiponectin and risk of coronary heart disease*. Jama, 2004. **292**(1): p. 40; author reply 40.
- Choi, K.M., et al., *Inflammation, Insulin Resistance, and Glucose Intolerance in Acute Myocardial Infarction Patients without a Previous Diagnosis of Diabetes Mellitus*. J Clin Endocrinol Metab, 2004.
- Nakamura, Y., et al., *Implications of plasma concentrations of adiponectin in patients with coronary artery disease*. Heart, 2004. **90**(5): p. 528-33.
- Pischon, T., et al., *Plasma adiponectin levels and risk of myocardial infarction in men*. Jama, 2004. **291**(14): p. 1730-7.
- Shibata, R., et al., *Adiponectin stimulates angiogenesis in response to tissue ischemia through stimulation of amp-activated protein kinase signaling*. J Biol Chem, 2004. **279**(27): p. 28670-4.
- Fernandez-Real, J.M., et al., *Adiponectin is associated with vascular function independent of insulin sensitivity*. Diabetes Care, 2004. **27**(3): p. 739-45.

[Revision date : 2025-07-03]

Other translations of this Instructions for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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