# **Resistin Elisa**

***KAPME50***

# History

---

## Summary of change :

| Previous Version :       | Current Version :                                                                                                       |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------|
| 230727                   | 250703                                                                                                                  |
| Outdated hazard warnings | Update of the hazard warnings in the<br>"WARNINGS AND PRECAUTIONS" section<br><br>Spanish translation has been included |

# Resistin-ELISA

Enzyme Immunoassay for the Quantitative Determination of human Resistin  
KAPME50  
IN VITRO DIAGNOSTIC USE

en

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 91

## INTENDED USE

Measurement of human Resistin in human Serum and Plasma Sample

## INTRODUCTION

Resistin, a cysteine-rich protein of 11.3 kDa (1), was firstly found in mice (2) and constitutes together with RELM $\alpha$ , RELM $\beta$  and RELM $\gamma$  the protein family of resistin-like molecules (RELM).

In humans, resistin and RELM $\beta$  (1) but no other proteins of the RELM family were found. The human form of resistin shows a homology of 53% to the murine protein (4). It has 11 cysteine-residues, is synthesized as a propeptide of 108 amino acids and secreted as a dimer, build by a disulfide bridge of cysteine residues (22). Beside this intermolecular disulfide bridge, 5 additional intramolecular ones exist (5,6).

Appearance of multi- and oligomer formation was proved by size exclusion chromatography. Thereby it was shown, that oligomer formation is SDS-insensitive but can be inhibited by  $\beta$ -mercaptoethanol and is therefore likely to be caused by disulfide bridges (1). Further on, the resistin structure seems to be dependent on its concentration, as circular dichroism analysis shows a concentration dependent shift of  $\alpha$ -helical to  $\beta$ -sheet structure (1).

Resistin expression was demonstrated in white adipose tissue (10), pituitary (11) and pancreatic islets (12) of mice as well as in brown adipose tissue of rats. In humans, resistin expression in adipocytes can be detected but only at a very low level. But in vitro, resistin expression of non-adipocytes in fatty tissue was shown (13). Human resistin gene is also expressed in pancreatic islets (12), pre-adipocytes (14) macrophages (15) and bone marrow (39). So, resistin is of relevance for inflammation processes as well as for lipid metabolism.

Most investigation refers to the mouse model. Here, the existence of trimeric and hexameric resistin in serum was demonstrated (7). In comparison to adiponectin biology, it is highly probable that different resistin oligomers have different biologic function (8, 9).

In mice, a correlation between adiposity, insulin resistance and resistin expression was found empirically. In humans, respective study results are not clear – several studies show an association of resistin serum concentration and adiposity or insulin resistance (17, 25-31). But others failed in confirming these results (14, 16-24). Therefore, there is requirement for valid and reproducible determination of resistin serum concentration.

Relevance of resistin in other physiologic processes than energy metabolism was investigated by several different approaches. Experiments with endothelial cells gave interesting results. Here, resistin was shown to enhance expression of VCAM-1 and ICAM-1 (33, 34). By this way, resistin is potentially able to influence endothelial inflammation (35, 36) and, thereby atherosclerosis. These results were confirmed by experiments in mice, where endothelin-1 was shown to regulate resistin secretion (37, 38).

In recent research human resistin was shown to increase pre-adipocyte proliferation and lipolysis of mature adipocytes (38). By the way of modulating MAPK-signalling pathways resistin exerts crucial influence on energy metabolism.

Present research demonstrates, that Resistin exerts influence on a broad variety of physiological processes, however a clear and defined biological role of resistin remains still unexisting.

This ELISA-kit enables the user to determine the exact concentration of Resistin in human serum/plasma as well as other body fluids and thereby assists investigation of Resistin biology.

With optimized special sample buffer: reliable quantification of resistin - independent of the absolute content.

## REAGENTS PROVIDED

- 1)  **Microtiter plate**, ready for use: **Microtiter plate** with 96 wells, divided up in 12 strips with 8 wells separately breakable, coated with rabbit-anti-resistin-antibody.
- 2) 

|     |   |
|-----|---|
| CAL | N |
|-----|---|

**Calibrators 1-5**, lyophilized: contain recombinant human Resistin. Calibrator values are given on vial labels and on quality certificate. **Reconstitute in 750  $\mu$ L Calibrator Diluent.**
- 3) 

|     |     |
|-----|-----|
| DIL | CAL |
|-----|-----|

**Calibrator Diluent**, 120 ml, ready for use.
- 4) 

|         |   |
|---------|---|
| CONTROL | N |
|---------|---|

**Controls 1 and 2**, lyophilized: Contain human Serum and have to be reconstituted with **250  $\mu$ L Dilution buffer**. The Resistin target value concentration and the respective range is given on the vial label. The **dilution** of the **Control 1 and 2** in **Calibrator Diluent** should be according done to the dilution of the respected samples.
- 5) 

|      |      |      |
|------|------|------|
| BIOT | CONJ | CONC |
|------|------|------|

**Biotin Conjugate**, 120  $\mu$ L, 100-fold concentrated solution, contains rabbit biotinylated anti-Resistin antibody, please dilute before use 1:100 in Dilution buffer.
- 6) 

|     |     |      |
|-----|-----|------|
| SAV | HRP | CONC |
|-----|-----|------|

**HRP Conjugate**, 120  $\mu$ L, 100-fold concentrated solution, contains HRP (Horseradish peroxidase)-labelled Streptavidin, please dilute before use **1:100** in **Dilution Buffer**.
- 7) 

|     |     |
|-----|-----|
| DIL | BUF |
|-----|-----|

**Dilution buffer**, 25 ml, ready for use.

- 8) 

|      |      |      |
|------|------|------|
| WASH | SOLN | CONC |
|------|------|------|

**Washing Buffer**, 50 ml, 20-fold concentrated: Washing Buffer has to be diluted 1:20 with distilled or demineralised water before use.
- 9) 

|       |     |
|-------|-----|
| CHROM | TMB |
|-------|-----|

**Chromogenic Substrate**, 12 ml, ready for use, horseradish-peroxidase-(HRP) substrate, stabilised Tetramethylbencidine.
- 10) 

|      |      |
|------|------|
| STOP | SOLN |
|------|------|

**Stopping Solution**, 12 ml, ready for use, 0,2 M sulphuric acid, *Caution!*
- 11) **Sealing tape** for covering of the microtiter plate, 2 x, adhesive.

#### MATERIALS REQUIRED BUT NOT PROVIDED

Distilled (Aqua destillata) or deionized water for dilution of the **Washing Buffer (A. dest.)**, 950 mL.

Precision pipettes and multichannel pipettes with disposable plastic tips

Polyethylene PE/Polypropylene PP tubes for dilution of samples

Vortex-mixer

Microtiter plate shaker (350 rpm)

Microtiter plate washer (recommended)

Micro plate reader ("ELISA-Reader") with filter for 450 and  $\geq 590$  nm

#### WARNINGS AND PRECAUTIONS

**For in-vitro diagnostic use only. For professional use only.**

The DiaSource kit is suitable only for in vitro diagnostics and not for internal use in humans and animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. DiaSource will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided.

**Caution: This kit contains material of human and/or animal origin. Therefore, all components and patient's specimens should be treated as potentially infectious.**

Do not use obvious damaged or microbial contaminated or spilled material.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations.

#### Human Serum

Following components contain human serum: Controls 1 and 2 *CONTROL N*

Source human serum for the Controls provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV). No known methods can offer total security of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

#### COMPONENTS

All components of the product that contain hazardous substances in declarable concentrations in accordance with European Regulation 1272/2008 are listed below:

**Dilution Buffer *DIL BUF* / Calibrator Diluent *DIL CAL* / Antibody Conjugate *BIOT CONJ CONC* / Enzyme Conjugate *SAV HRP CONC***

The buffers contain the hazardous ingredients listed below, based on this the following labeling of the two mixtures results:

Hazard pictograms (CLP):

GHS07

Signal word (CLP):

Hazardous components:

Hazard statements (CLP):

Safety instructions (CLP):

Warning

2-methylisothiazol-3(2H)-one (< 0,06 %)

reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) (< 0,015 %)

H317 - May cause an allergic skin reaction.

P280 - Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

P272 - Contaminated work clothing should not be allowed out of the workplace.

P261 - Avoid breathing dust/fume/gas/mist/vapours/spray.

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention. P302+P352 - IF ON SKIN: Wash with plenty of Water.

P501 - Dispose of contents/container in accordance with national disposal regulations.

**Washing Buffer WASH SOLN CONC**

The buffer contains a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one (< 0.015 %) as a preservative.

EUH208

Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1).

May produce an allergic reaction.

EUH210

Safety data sheet available on request.

**Stop Solution STOP SOLN**

The stop solution contains 0.2 M sulphuric acid.

EUH210

Safety data sheet available on request.

**PROTECTIVE EQUIPMENT AND FIRST AID MEASURES**

After contact with skin: Take off immediately all contaminated clothing. Rinse skin with water or shower. In case of skin reactions, consult a physician.

After contact with eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. In case of eye irritation consult an ophthalmologist.

After ingestion: Call a doctor if you feel unwell.

Avoid contact with the ingredients. Detailed information on personal protective equipment and first aid measures can be found in the safety data sheet for Resistin ELISA. This can be accessed electronically via the DiaSource website at <https://diasource-diagnostics.com/> or will be made available on request via [products.support@diasource.be](mailto:products.support@diasource.be).

**METHOD**

The enzyme immunoassay for Resistin is a so-called Sandwich-Assay. It utilizes a specific high affinity polyclonal rabbit antiserum coated on the wells of a microtiter plate. The Resistin in the samples binds quantitatively to the immobilized antiserum. In the following step, the biotinylated antiserum binds in turn to Resistin. After washing, Streptavidin-Peroxidase-Enzyme conjugate will be added, which will bind highly specific to the biotin of the antiserum and will catalyse in the closing substrate reaction the turn of the colour, quantitatively depending on the Resistin level of the samples.

**SPECIMEN**

Serum and Plasma

Serum as well plasma samples are suitable (significant deviation of Resistin levels in corresponding serum-, Heparin-, EDTA-, Citrate-plasma-Samples were not found. A possible dilution of the sample by the anticoagulant must be taken into account).

Haemolytic samples appear to show falsely high Resistin levels, using such samples should be checked out critically.

By means of the special sample buffer an external sample preparation prior to the assay is not required.

**Specimen collection**

Use standard venipuncture for the blood sampling. Haemolytic reactions are to be avoided.

**Required sample volume: 15 µL****Sample stability**

In firmly closable plastic tubes.

Storage at 20-25°C max. 2 days

Storage at -20°C max. 2 years

The storage of samples, over a period of 2 years at -20°C, showed no effect on the measured value. Freeze/thaw cycles of samples should be minimized.

**Interference**

Triglyceride, bilirubin and hemolysate in the sample do not interfere to a concentration of 100 mg/mL and 100 µg/mL or 1 mg/mL. However, the use of hemolytic, lipemic or icteric samples should be validated by the user.

**Sample Preparation**

- Dilution: 1:21 with **Calibrator Diluent**
- Pipette 300 µL Calibrator Diluent in PE-/PP-Tubes (application of a multi-stepper is recommended in larger series), add 15 µL Serum- or Plasma (dilution 1:21). After mixing use 2 x 100 µL of this dilution in the assay
- According to expected Resistin levels the dilution with Calibrator Diluent can be higher or lower. The excellent linearity of this test system allows sample dilution of 1:5 to 1:400.
- Because the Calibrator Diluent is special composed for the correct determination of Resistin, the dilution should be at least 1:5!

**TECHNICAL RECOMMENDATIONS****Storage Conditions**

Store the kit at 2-8°C after receipt until its expiry date. The lyophilized reagents should be stored at -20 °C after reconstitution. Avoid repeated thawing and freezing.

#### Storage Life

The shelf life of the components **after initial opening** is warranted for **4 weeks**, store the unused strips and microtiter wells **airtight** together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided. The **reconstituted components: Calibrators 1-5 and Controls 1 and 2** must be stored at -20°C (max. 4 weeks). For further use, thaw quickly but gently (avoid temperature increase above room temperature and avoid excessive vortexing). The 1:20 diluted **Washing Buffer** is 4 weeks stable at 2-8°C.

#### Preparation of reagents

Bring all reagents to room temperature (20 - 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming. Reagents with different lot numbers cannot be mixed.

#### Reconstitution

The **Calibrators 1-5** are reconstituted with the **Calibrator Diluent**.

The **Controls 1 and 2** are reconstituted with the **Dilution Buffer**.

It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

The reconstituted calibrators and controls can be stored for 4 weeks at -20°C. Repeated freeze/thaw cycles have to be avoided.

#### Dilution

After **reconstitution**, **dilute the Controls 1 and 2** with the **Calibrator Diluent** in the same ratio (**1:21**) as the samples.

The required volume of **Washing Buffer** is prepared by 1:20 dilution of the provided 20fold concentrate with Aqua dest.

The diluted **Washing Buffer** is stable for 4 weeks at 2-8°C. It has to be at room temperature for usage!

#### Antibody and Enzyme Conjugate

Use the **Dilution Buffer** for the dilution of the **Biotin Conjugate** and **HRP Conjugate** 100 fold concentrates. The diluted solutions are only limited stable at 2-8°C and should be prepared daily fresh.

#### ASSAY PROCEDURE

When performing the assay, Blank, **Calibrators 1-5, Controls 1 and 2** and the samples should be pipetted as fast as possible (e.g. <15 minutes). To avoid distortions due to differences in incubation times, **Biotin Conjugate** and the **HRP Conjugate** as well as the succeeding **Chromogenic Substrate** should be added to the plate in the same order and in the same time interval as the samples.

**Stop Solution** should be added to the plate in the same order as **Chromogenic Substrate**.

All determinations (Blank, **Calibrators 1-5, Controls 1 and 2** and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

#### Incubation

**Incubation at room temperature means: Incubation at 20 - 25°C.** The **Chromogenic Substrate**, stabilised Tetramethylbenzidine, is photosensitive—storage and incubation in the dark.

#### Shaking

The incubation steps should be performed at mean rotation frequency of a microtiter plate shaker. We recommend 350 rpm. Depending on the design of the shaker, the shaking frequency should be adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/ or false values, excessive shaking may result in high optical densities and/ or false values.

#### Washing

Proper washing is of basic **importance** for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided **Washing Buffer** diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

When using an **automatic microtiter** plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

**Manual washing** is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamical swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

# ASSAY PROCEDURE

| Preparation of reagents |                      | Reconstitution:                      | Dilution                                                     |
|-------------------------|----------------------|--------------------------------------|--------------------------------------------------------------|
| CAL 1-5                 | Calibrators          | in 750 µL Calibrator Diluent DIL CAL | -                                                            |
| Control 1               | Control 1            | in 250 µL Dilution Buffer DIL BUF    | 1:21 with Calibrator Diluent DIL CAL                         |
| Control 2               | Control 2            | in 250 µL Dilution Buffer DIL BUF    | 1:21 with Calibrator Diluent DIL CAL                         |
| BIOT CONJ CONC          | Biotin Conjugate     | -                                    | 1:100 with Dilution Buffer DIL BUF<br>→ BLOT CONJ CONC 1:100 |
| SAV HRP CONC            | HRP Conjugate        | -                                    | 1:100 with Dilution Buffer DIL BUF<br>→ SAV HRP CONC 1:100   |
| WASH SOLN CONC          | Washing Buffer conc. | -                                    | 1:20 with Aqua dest.<br>→ WASH SOLN CONC 1:20                |

Samples + Controls 1 and 2: dilute 1:21 in Calibrator Diluent , mix immediately

Before assay procedure bring all reagents to room temperature (20-25°C).

## Assay procedure in double determination

| Pipette | Reagents                         | Position                                            |
|---------|----------------------------------|-----------------------------------------------------|
| 100 µL  | Calibrator Diluent (blank value) | A1/A2                                               |
| 100 µL  | Calibrator 1 (0.02 ng/mL)        | B1/B2                                               |
| 100 µL  | Calibrator 2 (0.1 ng/mL)         | C1/C2                                               |
| 100 µL  | Calibrator 3 (0.3 ng/mL)         | D1/D2                                               |
| 100 µL  | Calibrator 4 (0.6 ng/mL)         | E1/E2                                               |
| 100 µL  | Calibrator 5 (1.0 ng/mL)         | F1/F2                                               |
| 100 µL  | Control 1 (1:21 diluted)         | G1/G2                                               |
| 100 µL  | Control 2 (1:21 diluted)         | H1/H2                                               |
| 100 µL  | Sample SPE (1:21 diluted)        | in the rest of the wells according the requirements |

Cover the wells with the sealing tape.

## Sample Incubation: 2 h at 20-25°C, 350 rpm

|           |                                                                                                           |              |
|-----------|-----------------------------------------------------------------------------------------------------------|--------------|
| 5x 300 µL | Aspirate the contents of the wells and wash 5x with 300 µL each Washing Buffer WASH SOLN CONC 1:20 / well | In each well |
| 100 µL    | Biotin Conjugate BLOT CONJ CONC 1:100                                                                     | In each well |

Cover the wells with the sealing tape.

## Incubation: 1 h at 20-25°C, 350 rpm

|           |                                                                                                           |              |
|-----------|-----------------------------------------------------------------------------------------------------------|--------------|
| 5x 300 µL | Aspirate the contents of the wells and wash 5x with 300 µL each Washing Buffer WASH SOLN CONC 1:20 / well | In each well |
| 100 µL    | HRP Conjugate SAV HRP CONC 1:100                                                                          | In each well |

Cover the wells with the sealing tape.

## Incubation: 30 Minutes at 20-25°C, 350 rpm

|           |                                                                                                           |              |
|-----------|-----------------------------------------------------------------------------------------------------------|--------------|
| 5x 300 µL | Aspirate the contents of the wells and wash 5x with 300 µL each Washing Buffer WASH SOLN CONC 1:20 / well | In each well |
| 100 µL    | Chromogenic Substrate                                                                                     | In each well |

Measure the absorbance within 30 min at **450 nm** with  $\geq 590$  nm as reference wavelength.

Measured extinction of the blank

0.05

Your measurement program will calculate the Resistin concentration of the diluted sample automatically by using the difference of sample and blank for the calculation. You only have to determine the most suitable curve fit (here: polynomial 2<sup>nd</sup> degree).

In this exemplary case the following equation is solved by the program to calculate the Resistin concentration in the sample:

$$y = 0.0046207 + 3.3674x - 1.5709x^2$$

$$0.2686 = x$$

if the dilution factor (1:21) is taken into account the Resistin concentration of the undiluted sample is

$$0.2686 \times 21 = 5.64 \text{ ng/mL} = 0.00564 \text{ } \mu\text{g/mL}$$

## PERFORMANCE CHARACTERISTICS

### Calibrators

The calibrators are prepared from recombinant human Resistin (19.5 kDa, 2 x 92 amino acids, expressed in E. coli) in concentrations of 20, 100, 300, 600 and 1000 pg/ml (pico Gramm / ml, equal to 0.02 ng/ml-1 ng/ml).

### Sensitivity

The analytical sensitivity was determined by measuring the blank value. For this purpose, the signal intensity of blank value + 2SD was converted into a resistin concentration using the calibration curve. The analytical sensitivity of the assay is 0.012 ng / mL.

### Interference

Interference of physiological appearing substance with the Resistin measurement was investigated. Serum samples have been enriched with different concentrations of possibly interfering substances and the amount of Resistin was measured and compared with the Resistin concentration in the same sample without any enrichment. In table 1 the relative results are shown. None of the tested substances interfered significantly with Resistin measurement.

**Table 1:** Interference: Three serum samples where enriched with indicated amount of the potentially interfering substance and measured. Shown is % of Resistin of the native, non enriched serum sample

|         | <b>Triglyceride<br/>100 mg/ml</b> | <b>Bilirubin<br/>100 µg/ml</b> | <b>Haemolysate<br/>1000 µg/ml</b> |
|---------|-----------------------------------|--------------------------------|-----------------------------------|
| Serum 1 | 101                               | 93                             | 94                                |
| Serum 2 | 115                               | 99                             | 99                                |
| Serum 3 | 104                               | 103                            | 147                               |

**Table 2:** Effects of coagulation inhibitors were investigating by adding indicated amounts of inhibitors to PP enriched with 0.3 ng/ml Resistin. Relative amounts of Resistin measured in inhibitor containing samples in comparison to 0.3 ng/ml Resistin containing Sample Buffer (PP) are shown.

|              |         | <b>% of Resistin in PP</b> |           |
|--------------|---------|----------------------------|-----------|
|              |         | <b>Mean (n=3)</b>          | <b>SD</b> |
| 3.8 g/l      | Citrate | 94                         | 7.67      |
| 0.0068 mol/l | EDTA    | 93                         | 4.96      |
| 30,000 IE/l  | Heparin | 96                         | 4.89      |

### Reproducibility and Precision

The inter- and intra assay coefficients of variability are in average < 10%. Exemplary determinations are shown in table 3 and table 4.

**Table 3:** Intra-Assay-Variation

|                 | <b>Number of<br/>determinations</b> | <b>Mean value<br/>(ng/ml)</b> | <b>Standard deviation<br/>(ng/ml)</b> | <b>VC<br/>(%)</b> |
|-----------------|-------------------------------------|-------------------------------|---------------------------------------|-------------------|
| <b>Sample 1</b> | 26                                  | 2.91                          | 0.16                                  | 5.55              |
| <b>Sample 2</b> | 15                                  | 4.58                          | 0.24                                  | 5.33              |
| <b>Sample 3</b> | 17                                  | 4.60                          | 0.23                                  | 5.04              |
| <b>Sample 4</b> | 7                                   | 2.50                          | 0.09                                  | 3.37              |
| <b>Sample 5</b> | 23                                  | 4.09                          | 0.27                                  | 6.67              |

**Table 4:** Inter-Assay-Variation (results of 11 determinations, each)

|                 | Number of determinations | Mean value (ng/ml) | Standard deviation (ng/ml) | VC (%) |
|-----------------|--------------------------|--------------------|----------------------------|--------|
| <b>Sample 1</b> | 16                       | 2.81               | 0.13                       | 4.49   |
| <b>Sample 2</b> | 15                       | 4.79               | 0.24                       | 4.97   |

**Recovery and Linearity**

The DIAsource Resistin ELISA is over a very wide range dilution authentic. The linearity of serum dilutions is over a very wide range excellent (s.Tab.5).

**Table 5:** Recovery and linearity of the Sample Dilution (characteristic results of two different sera)

| Dilution | Sample 1 (native 5.5 ng/ml) |              | Sample 2 (native 2.25 ng/ml) |              |
|----------|-----------------------------|--------------|------------------------------|--------------|
|          | plus 5 ng/ml                | Recovery (%) | plus 12.25 ng/ml             | Recovery (%) |
| 1:50     | 9.71                        | 92.5         | 14.99                        | 103.4        |
| 1:100    | 10.60                       | 101.0        | 13.64                        | 94.1         |
| 1:200    | 10.44                       | 99.4         | 14.10                        | 97.2         |
| 1:400    | 10.32                       | 98.3         | 14.33                        | 98.8         |

Different human sera were spiked with recombinant human Resistin in varying concentrations (e.g. in Table 5). The recovery of Resistin yielded on average 98 % of the theoretically expected amount.

**EVALUATION OF RESULTS****Table 6:** The expected values for Resistin were determined with the DIAsource ELISA in healthy probands and analysed by Prof. Dr. J. Kratzsch. Institute for Laboratory Medicine, University of Leipzig.

| Female       |     |         |         | Resistin (ng/mL): |                      |              |
|--------------|-----|---------|---------|-------------------|----------------------|--------------|
| Age (Years): | n:  | AV Age: | AV BMI: | AV $\pm$ SD:      | 25.- 75. Percentile: | Min. – Max.: |
| 18 - 30      | 96  | 23.0    | 23.1    | 7.2 $\pm$ 2.6     | 5.4 – 8.8            | 3.1 – 14.7   |
| 31 - 40      | 63  | 36.5    | 24.3    | 8.1 $\pm$ 2.3     | 6.4 – 9.6            | 3.6 – 13.1   |
| 41 - 50      | 67  | 44.9    | 24.8    | 7.3 $\pm$ 2.5     | 5.7 – 8.1            | 4.0 – 16.1   |
| 51 - 60      | 29  | 54.7    | 25.0    | 7.2 $\pm$ 2.6     | 5.4 – 8.5            | 4.0 – 15.5   |
| 61 - 65      | 9   | 62.7    | 25.2    | 6.6 $\pm$ 1.1     | 6.0 – 6.7            | 5.4 – 9.3    |
| Male         |     |         |         | Resistin (ng/mL): |                      |              |
| Age (Years): | n:  | AV Age: | AV BMI: | AV $\pm$ SD:      | 25.- 75. Percentile: | Min. – Max.: |
| 18 - 30      | 107 | 23.9    | 24.1    | 6.4 $\pm$ 1.8     | 5.0 – 7.6            | 2.5 – 13.1   |
| 31 - 40      | 59  | 35.9    | 25.0    | 6.7 $\pm$ 3.2     | 4.8 – 7.4            | 3.8 – 26.9   |
| 41 - 50      | 66  | 45.0    | 25.2    | 6.5 $\pm$ 2.8     | 4.5 – 7.4            | 2.4 – 16.7   |
| 51 - 60      | 36  | 54.8    | 26.4    | 6.1 $\pm$ 2.1     | 4.7 – 7.2            | 3.2 – 13.3   |
| 61 - 68      | 20  | 63.2    | 25.6    | 7.2 $\pm$ 1.8     | 6.0 – 8.2            | 4.5 – 11.2   |

n=Number of Probands. AV=Average Value. BMI=Body Mass Index (kg/m<sup>2</sup>). SD=Standard Deviation

**Table 7:** Grouping of the expected values

| Sex  | Number | Mean [ng/mL] | Standard deviation | 2.5. Percentile | 9.5. Percentile |
|------|--------|--------------|--------------------|-----------------|-----------------|
| Male | 288    | 6.48         | 2.44               | 3.32            | 11.68           |

|               |     |      |      |      |       |
|---------------|-----|------|------|------|-------|
| <b>Female</b> | 264 | 7.41 | 2.47 | 3.68 | 13.60 |
| <b>Total</b>  | 552 | 6.93 | 2.49 | 3.58 | 13.12 |

## Instructions for use for scientific application

### SCIENTIFIC APPLICATION

Resistin is relevant e.g. in research of: Adiposity Insulin Resistance, Diabetes Arteriosclerosis Inflammation. Due to Assays high sensitivity (12 pg/mL) this is also suitable for measurements in cell culture supernatants and in non-serum / plasma samples (e.g. in breast milk, urine, saliva).

### SAMPLES SUITABLE FOR SCIENTIFIC APPLICATION

Serum and Plasma, common cell culture medium, saliva, breast milk and urine were found to be suitable specimens too. Resistin concentrations may be completely different in body fluids of human origin other than serum or cell culture supernatants. (See Table 8 for recovery and dilution recommendations). Due to the special Sample Buffer, external sample pre-treatment is not necessary.

**Table 8:** Samples were enriched with 0.3 ng/mL Resistin and measured in comparison to non-enriched sample. Relative recovery of added Resistin is shown.

| Matrix                   | Dilution | % Recovery |
|--------------------------|----------|------------|
| Cerebrospinal fluid      | 1:2      | 129        |
| Cerebrospinal fluid      | 1:10     | 93         |
| Cerebrospinal fluid      | 1:40     | 103        |
| Amnion fluid             | 1:10     | 85         |
| Amnion fluid             | 1:40     | 91         |
| Saliva                   | 1:10     | 99         |
| Saliva                   | 1:21     | 86         |
| Urine                    | 1:10     | 79         |
| Urine                    | 1:21     | 85         |
| Breast milk              | 1:2      | 97         |
| Breast milk              | 1:10     | 58         |
| Breast milk              | 1:21     | 63         |
| Cell culture supernatant | 1:2      | 100        |

### CROSS REACTIVITY

The cross-reactivities of several commercially available animal sera were tested in this assay. Commercially available sera from cattle, cat, chicken, dog, donkey, goat, guinea pig, horse, mouse, pig, rabbit, rat and sheep were diluted 1:10 and used in the test and the signal intensity was measured. No cross-reactivity could be detected.

### LITERATUR / LITERATURE

1. Aruna B. Ghosh S. Singh AK. Mande SC. Srinivas V. Chauhan R. Ehtesham NZ. Human recombinant resistin protein displays a tendency to aggregate by forming intermolecular disulfide linkages. *Biochemistry* 2003;42(36):10554-9.
2. Steppan CM. Lazar MA. Resistin and obesity-associated insulin resistance. *Trends Endocrinol Metab* 2002;13(1):18-23.
3. Yang RZ. Huang Q. Xu A. McLenithan JC. Eisen JA. Shuldiner AR. Alkan S. Gong DW. Eison JA. Comparative studies of resistin expression and phylogenomics in human and mouse. *Biochem Biophys Res Commun* 2003;310(3):927-35.
4. Banerjee RR. Lazar MA. Resistin: molecular history and prognosis. *J Mol Med* 2003;81(4):218-26.
5. Banerjee RR. Lazar MA. Dimerization of resistin and resistin-like molecules is determined by a single cysteine. *J Biol Chem* 2001;276(28):25970-3.
6. Raghu P. Ghosh S. Soundarya K. Haseeb A. Aruna B. Ehtesham NZ. Dimerization of human recombinant resistin involves covalent and noncovalent interactions. *Biochem Biophys Res Commun* 2004;313(3):642-6.
7. Patel SD. Rajala MW. Rossetti L. Scherer PE. Shapiro L. Disulfide-dependent multimeric assembly of resistin family hormones. *Science* 2004;304(5674):1154-8.
8. Kim KH. Zhao L. Moon Y. Kang C. Sul HS. Dominant inhibitory adipocyte-specific secretory factor (ADSF)/resistin enhances adipogenesis and improves insulin sensitivity. *Proc Natl Acad Sci U S A* 2004;101(17):6780-5.
9. Chen J. Wang L. Boeg YS. Xia B. Wang J. Differential dimerization and association among resistin family proteins with implications for functional specificity. *J Endocrinol* 2002;175(2):499-504.
10. O.Beltowski J. Adiponectin and resistin--new hormones of white adipose tissue. *Med Sci Monit* 2003;9(2):RA55-61.
11. Morash BA. Ur E. Wiesner G. Roy J. Wilkinson M. Pituitary resistin gene expression: effects of age, gender and obesity. *Neuroendocrinology* 2004;79(3):149-56.
12. Minn AH. Patterson NB. Pack S. Hoffmann SC. Gavrilova O. Vinson C. Harlan DM. Shalev A. Resistin is expressed in pancreatic islets. *Biochem Biophys Res Commun* 2003;310(2):641-5.
13. Fain JN. Cheema PS. Bahouth SW. Lloyd Hiler M. Resistin release by human adipose tissue explants in primary culture. *Biochem Biophys Res Commun* 2003;300(3):674-8.
14. Janke J. Engeli S. Gorzelniak K. Luft FC. Sharma AM. Resistin gene expression in human adipocytes is not related to insulin resistance. *Obes Res* 2002;10(1):1-5.

15. Patel L. Buckels AC. Kinghorn IJ. Murdock PR. Holbrook JD. Plumpton C. Macphee CH. Smith SA. Resistin is expressed in human macrophages and directly regulated by PPAR gamma activators. *Biochem Biophys Res Commun* 2003;300(2):472-6.
16. Youn BS. Yu KY. Park HJ. Lee NS. Min SS. Youn MY. Cho YM. Park YJ. Kim SY. Lee HK. Park KS. Plasma resistin concentrations measured by enzyme-linked immunosorbent assay using a newly developed monoclonal antibody are elevated in individuals with type 2 diabetes mellitus. *J Clin Endocrinol Metab* 2004;89(1):150-6.
17. Chen XD. Lei T. Xia T. Gan L. Yang ZQ. Increased expression of resistin and tumour necrosis factor- $\alpha$  in pig adipose tissue as well as effect of feeding treatment on resistin and cAMP pathway. *Diabetes Obes Metab* 2004;6(4):271-9.
18. Volarova de Courten B. Degawa-Yamauchi M. Considine RV. Tataranni PA. High serum resistin is associated with an increase in adiposity but not a worsening of insulin resistance in Pima Indians. *Diabetes* 2004;53(5):1279-84.
19. Shalev A. Patterson NB. Hirshberg B. Rother KI. Harlan DM. Resistin serum levels in type 1 diabetes pre- and post-islet transplantation. *Metabolism* 2004;53(4):403-4.
20. Fujinami A. Obayashi H. Ohta K. Ichimura T. Nishimura M. Matsui H. Kawahara Y. Yamazaki M. Ogata M. Hasegawa G. Nakamura N. Yoshikawa T. Nakano K. Ohta M. Enzyme-linked immunosorbent assay for circulating human resistin: resistin concentrations in normal subjects and patients with type 2 diabetes. *Clin Chim Acta* 2004;339(1-2):57-63.
21. Degawa-Yamauchi M. Bovenkerk JE. Juliar BE. Watson W. Kerr K. Jones R. Zhu Q. Considine RV. Serum resistin (FIZZ3) protein is increased in obese humans. *J Clin Endocrinol Metab* 2003;88(11):5452-5.
22. Azuma K. Katsukawa F. Oguchi S. Murata M. Yamazaki H. Shimada A. Saruta T. Correlation between serum resistin level and adiposity in obese individuals. *Obes Res* 2003;11(8):997-1001.
23. Wang H. Chu WS. Hemphill C. Elbein SC. Human resistin gene: molecular scanning and evaluation of association with insulin sensitivity and type 2 diabetes in Caucasians. *J Clin Endocrinol Metab* 2002;87(6):2520-4.
24. Brichard SM. Delparte ML. Lambert M. Adipocytokines in anorexia nervosa: a review focusing on leptin and adiponectin. *Horm Metab Res* 2003;35(6):337-42.
25. Kawanami D. Maemura K. Takeda N. Harada T. Nojiri T. Imai Y. Manabe I. Utsunomiya K. Nagai R. Direct reciprocal effects of resistin and adiponectin on vascular endothelial cells: a new insight into adipocytokine-endothelial cell interactions. *Biochem Biophys Res Commun* 2004;314(2):415-9.
26. Verma S. Li SH. Wang CH. Fedak PW. Li RK. Weisel RD. Mickle DA. Resistin promotes endothelial cell activation: further evidence of adipokine-endothelial interaction. *Circulation* 2003;108(6):736-40.
27. Kaser S. Kaser A. Sandhofer A. Ebenbichler CF. Tilg H. Patsch JR. Resistin messenger-RNA expression is increased by proinflammatory cytokines in vitro. *Biochem Biophys Res Commun* 2003;309(2):286-90.
28. Lu SC. Shieh WY. Chen CY. Hsu SC. Chen HL. Lipopolysaccharide increases resistin gene expression in vivo and in vitro. *FEBS Lett* 2002;530(1-3):158-62.
29. Zhong Q. Lin CY. Clarke KJ. Kemppainen RJ. Schwartz DD. Judd RL. Endothelin-1 inhibits resistin secretion in 3T3-L1 adipocytes. *Biochem Biophys Res Commun* 2002;296(2):383-7.
30. Barton M. Carmona R. Ortmann J. Krieger JE. Traupe T. Obesity-associated activation of angiotensin and endothelin in the cardiovascular system. *Int J Biochem Cell Biol* 2003;35(6):826-37.
31. Ort TA. A. A.; MacDougall. J.R.; Nelson. P.J.; Rothenberg. M.E.; Wu. F.; Eisen. A.; Halvorsen. Y.-D.C. Recombinant human FIZZ3/Resistin Stimulates Lipolysis in Cultured Human Adipocytes. Mouse Adipose Explants and in Normal Mice. *Endocrinology* 2005; doi:10.1210/en.2004-1421.

Other translations of this Instructions for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

**Summary of the Assay**

| Reagent preparation:                                                                                                                                                 | Reconstitution:                     | Dilution:                           |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------|-------------------------------------|
| Calibrators 1 – 5                                                                                                                                                    | in <b>750 µl</b> Calibrator Diluent |                                     |
| Control 1                                                                                                                                                            | in <b>250 µl</b> Dilution Buffer    | <b>1:21</b> with Calibrator Diluent |
| Control 2                                                                                                                                                            | in <b>250 µl</b> Dilution Buffer    | <b>1:21</b> with Calibrator Diluent |
| Biotin Conjugate                                                                                                                                                     |                                     | <b>1:100</b> with Dilution Buffer   |
| HRP Conjugate                                                                                                                                                        |                                     | <b>1:100</b> with Dilution Buffer   |
| Washing Buffer                                                                                                                                                       |                                     | <b>1:20</b> with <b>Aqua. dest.</b> |
| <b>Sample dilution: 1:21</b> (e.g. 15 µl Serum with 300 µl Calibrator Diluent).<br>Before assay procedure bring all reagents to room temperature ( <b>20-25°C</b> ). |                                     |                                     |

**Assay Procedure for Double Determination**

| Pipette                                                                                                    | Reagents                                                                      | Position        |
|------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------|-----------------|
| 100 µl                                                                                                     | Calibrator Diluent (blank value)                                              | A1/2            |
| 100 µl                                                                                                     | Calibrator <b>1</b> (0.02 ng/ml)                                              | B1/2            |
| 100 µl                                                                                                     | Calibrator <b>2</b> (0.1 ng/ml)                                               | C1/2            |
| 100 µl                                                                                                     | Calibrator <b>3</b> (0.3 ng/ml)                                               | D1/2            |
| 100 µl                                                                                                     | Calibrator <b>4</b> (0.6 ng/ml)                                               | E1/2            |
| 100 µl                                                                                                     | Calibrator <b>5</b> (1.0 ng/ml)                                               | F1/2            |
| 100 µl                                                                                                     | Control 1 (diluted)                                                           | G1/2            |
| 100 µl                                                                                                     | Control 2 (diluted)                                                           | H1/2            |
| 100 µl                                                                                                     | Sample dilution                                                               | following wells |
| Cover the wells with the sealing tape.                                                                     |                                                                               |                 |
| <b>Incubation: 2 h at 20-25°C, ≥ 350 rpm</b>                                                               |                                                                               |                 |
| 5x 300 µl                                                                                                  | Aspirate the contents of the wells and wash 3x with <b>300 µl</b> Wash Buffer | each well       |
| 100 µl                                                                                                     | <b>1:100</b> diluted Biotin Conjugate                                         | each well       |
| <b>Incubation: 1 h at 20-25°C, ≥350 rpm</b>                                                                |                                                                               |                 |
| 5x 300 µl                                                                                                  | Aspirate the contents of the wells and wash 3x with <b>300 µl</b> Wash Buffer | each well       |
| 100 µl                                                                                                     | <b>1:100</b> diluted HRP Conjugate                                            | each well       |
| <b>Incubation: 30 min at 20-25°C, ≥350 rpm</b>                                                             |                                                                               |                 |
| 5x 300 µl                                                                                                  | Aspirate the contents of the wells and wash 3x with <b>300 µl</b> Wash Buffer | each well       |
| 100 µl                                                                                                     | Chromogenic Substrate                                                         | each well       |
| <b>Incubation: 30 min in the dark at 20-25°C</b>                                                           |                                                                               |                 |
| 100 µl                                                                                                     | Stop Solution                                                                 | each well       |
| Measure the absorbance within <b>30 min</b> at <b>450 nm</b> with <b>≥ 590 nm</b> as reference wavelength. |                                                                               |                 |

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### BIBLIOGRAFÍA/LITERATURA

1. Aruna B. Ghosh S. Singh AK. Mande SC. Srinivas V. Chauhan R. Ehtesham NZ. Human recombinant resistin protein displays a tendency to aggregate by forming intermolecular disulfide linkages. *Biochemistry* 2003;42(36):10554-9.
2. Stepan CM. Lazar MA. Resistin and obesity-associated insulin resistance. *Trends Endocrinol Metab* 2002;13(1):18-23.
3. Yang RZ. Huang Q. Xu A. McLenithan JC. Eisen JA. Shuldiner AR. Alkan S. Gong DW. Eison JA. Comparative studies of resistin expression and phylogenomics in human and mouse. *Biochem Biophys Res Commun* 2003;310(3):927-35.
4. Banerjee RR. Lazar MA. Resistin: molecular history and prognosis. *J Mol Med* 2003;81(4):218-26.
5. Banerjee RR. Lazar MA. Dimerization of resistin and resistin-like molecules is determined by a single cysteine. *J Biol Chem* 2001;276(28):25970-3.
6. Raghu P. Ghosh S. Soundarya K. Haseeb A. Aruna B. Ehtesham NZ. Dimerization of human recombinant resistin involves covalent and noncovalent interactions. *Biochem Biophys Res Commun* 2004;313(3):642-6.
7. Patel SD. Rajala MW. Rossetti L. Scherer PE. Shapiro L. Disulfide-dependent multimeric assembly of resistin family hormones. *Science* 2004;304(5674):1154-8.
8. Kim KH. Zhao L. Moon Y. Kang C. Sul HS. Dominant inhibitory adipocyte-specific secretory factor (ADSF)/resistin enhances adipogenesis and improves insulin sensitivity. *Proc Natl Acad Sci U S A* 2004;101(17):6780-5.
9. Chen J. Wang L. Boeg YS. Xia B. Wang J. Differential dimerization and association among resistin family proteins with implications for functional specificity. *J Endocrinol* 2002;175(2):499-504.
10. O.Beltowski J. Adiponectin and resistin--new hormones of white adipose tissue. *Med Sci Monit* 2003;9(2):RA55-61.
11. Morash BA. Ur E. Wiesner G. Roy J. Wilkinson M. Pituitary resistin gene expression: effects of age, gender and obesity. *Neuroendocrinology* 2004;79(3):149-56.
12. Minn AH. Patterson NB. Pack S. Hoffmann SC. Gavrilova O. Vinson C. Harlan DM. Shalev A. Resistin is expressed in pancreatic islets. *Biochem Biophys Res Commun* 2003;310(2):641-5.
13. Fain JN. Cheema PS. Bahouth SW. Lloyd Hiler M. Resistin release by human adipose tissue explants in primary culture. *Biochem Biophys Res Commun* 2003;300(3):674-8.

14. Janke J. Engeli S. Gorzelniak K. Luft FC. Sharma AM. Resistin gene expression in human adipocytes is not related to insulin resistance. *Obes Res* 2002;10(1):1-5.
15. Patel L. Buckels AC. Kinghorn IJ. Murdock PR. Holbrook JD. Plumpton C. Macphee CH. Smith SA. Resistin is expressed in human macrophages and directly regulated by PPAR gamma activators. *Biochem Biophys Res Commun* 2003;300(2):472-6.
16. Youn BS. Yu KY. Park HJ. Lee NS. Min SS. Youn MY. Cho YM. Park YJ. Kim SY. Lee HK. Park KS. Plasma resistin concentrations measured by enzyme-linked immunosorbent assay using a newly developed monoclonal antibody are elevated in individuals with type 2 diabetes mellitus. *J Clin Endocrinol Metab* 2004;89(1):150-6.
17. Chen XD. Lei T. Xia T. Gan L. Yang ZQ. Increased expression of resistin and tumour necrosis factor-alpha in pig adipose tissue as well as effect of feeding treatment on resistin and cAMP pathway. *Diabetes Obes Metab* 2004;6(4):271-9.
18. Volarova de Courten B. Degawa-Yamauchi M. Considine RV. Tataranni PA. High serum resistin is associated with an increase in adiposity but not a worsening of insulin resistance in Pima Indians. *Diabetes* 2004;53(5):1279-84.
19. Shalev A. Patterson NB. Hirshberg B. Rother KI. Harlan DM. Resistin serum levels in type 1 diabetes pre- and post-islet transplantation. *Metabolism* 2004;53(4):403-4.
20. Fujinami A. Obayashi H. Ohta K. Ichimura T. Nishimura M. Matsui H. Kawahara Y. Yamazaki M. Ogata M. Hasegawa G. Nakamura N. Yoshikawa T. Nakano K. Ohta M. Enzyme-linked immunosorbent assay for circulating human resistin: resistin concentrations in normal subjects and patients with type 2 diabetes. *Clin Chim Acta* 2004;339(1-2):57-63.
21. Degawa-Yamauchi M. Bovenkerk JE. Juliar BE. Watson W. Kerr K. Jones R. Zhu Q. Considine RV. Serum resistin (FIZZ3) protein is increased in obese humans. *J Clin Endocrinol Metab* 2003;88(11):5452-5.
22. Azuma K. Katsukawa F. Oguchi S. Murata M. Yamazaki H. Shimada A. Saruta T. Correlation between serum resistin level and adiposity in obese individuals. *Obes Res* 2003;11(8):997-1001.
23. Wang H. Chu WS. Hemphill C. Elbein SC. Human resistin gene: molecular scanning and evaluation of association with insulin sensitivity and type 2 diabetes in Caucasians. *J Clin Endocrinol Metab* 2002;87(6):2520-4.
24. Brichard SM. Delporte ML. Lambert M. Adipocytokines in anorexia nervosa: a review focusing on leptin and adiponectin. *Horm Metab Res* 2003;35(6):337-42.
25. Kawanami D. Maemura K. Takeda N. Harada T. Nojiri T. Imai Y. Manabe I. Utsunomiya K. Nagai R. Direct reciprocal effects of resistin and adiponectin on vascular endothelial cells: a new insight into adipocytokine-endothelial cell interactions. *Biochem Biophys Res Commun* 2004;314(2):415-9.
26. Verma S. Li SH. Wang CH. Fedak PW. Li RK. Weisel RD. Mickle DA. Resistin promotes endothelial cell activation: further evidence of adipokine-endothelial interaction. *Circulation* 2003;108(6):736-40.
27. Kaser S. Kaser A. Sandhofer A. Ebenbichler CF. Tilg H. Patsch JR. Resistin messenger-RNA expression is increased by proinflammatory cytokines in vitro. *Biochem Biophys Res Commun* 2003;309(2):286-90.
28. Lu SC. Shieh WY. Chen CY. Hsu SC. Chen HL. Lipopolysaccharide increases resistin gene expression in vivo and in vitro. *FEBS Lett* 2002;530(1-3):158-62.
29. Zhong Q. Lin CY. Clarke KJ. Kempainen RJ. Schwartz DD. Judd RL. Endothelin-1 inhibits resistin secretion in 3T3-L1 adipocytes. *Biochem Biophys Res Commun* 2002;296(2):383-7.
30. Barton M. Carmona R. Ortmann J. Krieger JE. Traupe T. Obesity-associated activation of angiotensin and endothelin in the cardiovascular system. *Int J Biochem Cell Biol* 2003;35(6):826-37.
31. Ort TA. A. A.; MacDougall. J.R.; Nelson. P.J.; Rothenberg. M.E.; Wu. F.; Eisen. A.; Halvorsen. Y.-D.C. Recombinant human FIZZ3/Resistin Stimulates Lipolysis in Cultured Human Adipocytes. Mouse Adipose Explants and in Normal Mice. *Endocrinology* 2005; doi:10.1210/en.2004-1421.

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