TG Ab RIA

R-CI-100

History

Summary of change:

Previous Version: 230123				Current Version: 260505																																																																			
8. CALCULATION OF RESULTS “Determine the calibration curve by logarithmic plot (log / lin) of B / T (%) of each calibrator (y-axis) versus relative concentration (x-axis). Calculate the B / T (%) of each sample and determine the concentration by interpolation on the calibration curve.				Correction of a typo linked to the fixation percentage expressed, by mistake, in B/T (%) instead of B/B0 (%)																																																																			
				8. CALCULATION OF RESULTS “Determine the calibration curve by logarithmic plot (log / lin) of B/B0 (%) of each calibrator (y-axis) versus relative concentration (x-axis). Calculate the B/B0 (%) of each sample and determine the concentration by interpolation on the calibration curve.																																																																			
<table><tr><th>Description</th><th>Average cpm</th><th>B/T (%)</th><th>TG Ab (IU/ml)</th></tr><tr><td>Total activity (T)</td><td>57953</td><td>-</td><td>-</td></tr><tr><td>CAL 0</td><td>28090</td><td>100</td><td>0</td></tr><tr><td>CAL 1</td><td>24340</td><td>86.6</td><td>20</td></tr><tr><td>CAL 2</td><td>16927</td><td>60.3</td><td>60</td></tr><tr><td>CAL 3</td><td>9706</td><td>34.5</td><td>200</td></tr><tr><td>CAL 4</td><td>4100</td><td>14.6</td><td>1000</td></tr><tr><td>CAL 5</td><td>2048</td><td>7.3</td><td>2000</td></tr></table>				Description	Average cpm	B/T (%)	TG Ab (IU/ml)	Total activity (T)	57953	-	-	CAL 0	28090	100	0	CAL 1	24340	86.6	20	CAL 2	16927	60.3	60	CAL 3	9706	34.5	200	CAL 4	4100	14.6	1000	CAL 5	2048	7.3	2000	<table><tr><th>Description</th><th>Average cpm</th><th>B/B0 (%)</th><th>TG Ab (IU/ml)</th></tr><tr><td>Total activity (T)</td><td>57953</td><td>-</td><td>-</td></tr><tr><td>CAL 0</td><td>28090</td><td>100</td><td>0</td></tr><tr><td>CAL 1</td><td>24340</td><td>86.6</td><td>20</td></tr><tr><td>CAL 2</td><td>16927</td><td>60.3</td><td>60</td></tr><tr><td>CAL 3</td><td>9706</td><td>34.5</td><td>200</td></tr><tr><td>CAL 4</td><td>4100</td><td>14.6</td><td>1000</td></tr><tr><td>CAL 5</td><td>2048</td><td>7.3</td><td>2000</td></tr></table>				Description	Average cpm	B/B0 (%)	TG Ab (IU/ml)	Total activity (T)	57953	-	-	CAL 0	28090	100	0	CAL 1	24340	86.6	20	CAL 2	16927	60.3	60	CAL 3	9706	34.5	200	CAL 4	4100	14.6	1000	CAL 5	2048	7.3	2000
Description	Average cpm	B/T (%)	TG Ab (IU/ml)																																																																				
Total activity (T)	57953	-	-																																																																				
CAL 0	28090	100	0																																																																				
CAL 1	24340	86.6	20																																																																				
CAL 2	16927	60.3	60																																																																				
CAL 3	9706	34.5	200																																																																				
CAL 4	4100	14.6	1000																																																																				
CAL 5	2048	7.3	2000																																																																				
Description	Average cpm	B/B0 (%)	TG Ab (IU/ml)																																																																				
Total activity (T)	57953	-	-																																																																				
CAL 0	28090	100	0																																																																				
CAL 1	24340	86.6	20																																																																				
CAL 2	16927	60.3	60																																																																				
CAL 3	9706	34.5	200																																																																				
CAL 4	4100	14.6	1000																																																																				
CAL 5	2048	7.3	2000																																																																				

ENGLISH

RADIOIMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF ANTI-THYROGLOBULIN AUTOANTIBODIES IN HUMAN SERUM R-CI-100 - 96 Determinations

FOR IN VITRO DIAGNOSTIC USE ONLY

1. CLINICAL APPLICATIONS

Autoimmune thyroid disease is more common in women. The prevalence of antibodies in women increases with age: about 10% between the ages of 18-24 and about 30% between the ages of 55-65 for Tg Ab. For TM Ab between the ages of 18-24 around 15%, between the ages of 55-65 about 24%.

2. PRINCIPLE OF THE ASSAY

The present method is based on the Competitive Radioimmunoassay (RIA). During the incubation, the anti-Tg monoclonal antibody in the solid phase competes with the anti-Tg autoantibody of the calibrators or samples for specific binding sites of the ¹²⁵I-labeled Tg (tracer). After aspirating and washing, the radioactivity in the tubes is measured in a gamma counter. The bound radioactivity is in inverse proportion to the anti-Tg autoantibody concentration of the samples. The calibrators of the anti-TG autoantibody test were standardized against the TGAb standard MRC65 / 93.

3. REAGENTS PROVIDED WITH THE KIT

- The reagents are sufficient for 96 determinations.
- Store the kit and reagents at 2-8°C.
- The expiration date of each reagent is shown on the label.

Reagents	96 Tests Kit	Colour Code	Reconstitution
 Tubes coated with anti-human thyroglobulin monoclonal antibody	2 x 48	Green	Ready to Use
Ag ¹²⁵I TRACER : ¹²⁵I-Thyroglobulin in phosphate buffer. Preservative: NaN3 (<0.1%)	1 Vial 22ml 137 kBq	Red	Ready to Use
[CAL]0] Calibrator Zero TRIS/BSA buffer. Preservative: NaN3 (<0.1%)	1 Vial 4.2ml	Yellow	Ready to Use
[CAL]N] Calibrator – N = 5 Anti-Thyroglobulin in Human plasma. Preservative: NaN3 (<0.1%) (see exact values on vial labels)	5 Vials 0.7ml	Yellow	Ready to Use
[CONTROL N] Control - N = 2 Anti-Thyroglobulin in Human plasma. Preservative: NaN3 (<0.1%) (see exact values on vial labels)	2 Vials 0.7ml	Silver	Ready to Use
[WASH SOLN CONC] Washing Solution (50 x concentrated) Tris-HCl buffer with detergent and preservative NaN3 (<0.1%).	1 Vial 20ml	Brown	Bring to 1000 ml with distilled water

1. **Controls and Calibrators:** For the exact value, refer to the vial label.
2. **Coated Tubes:** Unused tubes must be stored at 2-8°C , protected from moisture.
3. **Washing Solution (50 x concentrated):** The diluted washing solution is stable for 2 months at 2-8°C.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Plastic test tubes.
- Test tube racks.
- Adjustable, automatic micropipettes with disposable tips.
- Vortex mixer.
- Graduated cylinder.
- Aspiration pump or automated washing device.
- Scintillation gamma counter.
- Distilled water.
- Orbital shaker adjustable at 150 rpm.

5. WARNINGS AND PRECAUTIONS FOR USERS

For in vitro diagnostic use.
Only experienced laboratory personnel should use this test and handling should be in agreement with GLP.
Radioactive Material - Not for Internal or External Use in Humans or Animals.

In order to obtain reproducible results, the following rules must be observed :

- Do not mix reagents of different lots.
- Do not use reagents beyond their expiry date.
- Use thoroughly clean glassware.
- Use distilled water, stored in clean containers.
- Avoid any contamination among samples; to this purpose disposable tips should be used for each sample and reagent.

Safety :
This kit contains ¹²⁵I (half-life: 60 days) ,emitting ionizing X (28 keV) and γ (35.5 keV) radiations. This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.
All radioactive handling should be executed in a designated area. away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.
Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.
The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.
All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.
Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up. Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

For more information, see Material Safety Data Sheet (MSDS).

6. SPECIMEN COLLECTION

Use only serum samples. Highly lipemic or hemolytic samples must be discarded. Samples can be stored at 2-8 ° C for 1 day; For longer periods it is recommended to freeze the samples aliquoted at - 20 ° C. Repeated freezing or thawing should be avoided.

7. ASSAY PROCEDURE

- Before use, all reagents and samples must be brought to room temperature (18-25°C).
 - Before use, the samples must be mixed by panning.
 - For all calibrators, a double measurement is recommended.
 - Each tube can only be used once.
1. Prepare empty tubes for total counts and coated tubes for calibrators, samples, and controls.
 2. Pipette **200µl** of each calibrator, controls and samples into the appropriate tubes.
 3. Add **200µl** of the radioactive tracer to all tubes.
 4. Mix the tubes on a vortex mixer and incubate for **90 minutes** at room temperature (18-25°C) on an orbital shaker at 150 rpm.
 5. Aspirate the contents of the tubes except for the Total Counts tubes.
 6. Add **2 ml** of the diluted wash solution to all tubes except in the Total Counts tubes. Suction the contents of all tubes thoroughly or pour gently onto absorbent paper.
 7. Count the radioactivity bound to the tubes in a gamma counter for 1 minute. We suggest measuring the background of the instrument before counting. To avoid deviations in the sensitivity of the system, the background must be reduced to a minimum or adjusted.

Note: If high anti-Tg values (higher than the last calibrator) are expected for some patient samples, the original samples must be diluted. The dilution must be carried out with the zero calibrator. **Do not use a buffer for this purpose.**

ASSAY SCHEME

Tubes Reagent	Total Activity	Calibrators	Control	Samples
Calibrators	----	20 µl	----	----
Controls	----	----	20 µl	----
Samples	----	----	----	20 µl
Tracer	200 µl	200 µl	200 µl	200 µl
- Incubate: 90 min R.T. (18-25°C), shaking (150 rpm)				
- Aspirate - washing: 1x 2ml - suction				
- Counting				

8. CALCULATION OF RESULTS

Determine the calibration curve by logarithmic plot (log / lin) of B/B0 (%) of each calibrator (y-axis) versus relative concentration (x-axis). Calculate the B/B0 (%) of each sample and determine the concentration by interpolation on the calibration curve.

EXAMPLE OF CALCULATION

The values reported below must be considered as an example and may not be used in place of experimental data.

Description	Average cpm	B/B0 (%)	TG Ab (IU/ml)
Total activity (T)	57953	-	-
CAL 0	28090	100	0
CAL 1	24340	86.6	20
CAL 2	16927	60.3	60
CAL 3	9706	34.5	200
CAL 4	4100	14.6	1000
CAL 5	2048	7.3	2000
CONTROL 1	10843	38.6	158.2
CONTROL 2	5757	20.5	604.7
P1	13998	49.8	91.5
P2	6634	23.6	464.8
P3	4270	15.2	942.9

INTERNAL QUALITY CONTROL

- If the actual values do not correspond to the nominal values indicated on the vials, the values can not be further processed without an appropriate explanation of the deviations.
- If extra controls are desired, each laboratory can create its own pool, which should be frozen in aliquots.

- Acceptance criteria for the difference between the results of the repeat tests based on the samples must be based on Good Laboratory Practice.

9. REFERENCE VALUES

It is recommended that each laboratory determines its own reference interval. Values reported below are only indicative.

Less than 30 IU/ml	Negative for Thyroglobulin autoantibodies
Between 30 - 70 IU/ml	Limit
More than 70 IU/ml	Positive for Thyroglobulin autoantibodies

10. PERFORMANCES OF THE ASSAY

SENSITIVITY

Analytical sensitivity

The sensitivity was calculated based on the calibration curve and expressed in such a way that the smallest doses give a significant difference to the zero calibrator (mean - 3 times the standard deviation).

This dosage is 6 IU/ml.

Functional sensitivity

The functional test sensitivity is the lowest value measured with a maximum precision of 20% inter-assay variance. For anti-thyroglobulin autoantibodies, this value is less than 15 IU/ml.

PRECISION

The precision was evaluated after intra-assay and inter-assay variability, at different concentrations of the analytes.

Intra-assay

Serum	Mean	± (IU/ml)	S.D.	C.V. (%)	N
1	72.6	±	6.0	8.2	20
2	345.7	±	24.7	7.1	20
3	889.9	±	51.0	5.7	20

Inter-assay

Serum	Mean	± (IU/ml)	S.D.	C.V. (%)	N
1	90.3	±	6.3	7.0	9
2	390.1	±	47.0	11.7	9
3	985.3	±	65.8	7.4	9

ACCURACY

Accuracy of the method has been checked by the recovery and parallelism tests.

Recovery Test

Samples, mixed with equal volumes of each calibrator, were tested.

Added	Expected (IU/ml)	Measured (IU/ml)	Recovery (%)
S1	-	50.4	-
S1 + CAL 0	25.2	32.0	127.0
S1 + CAL 1	35.2	39.1	111.1
S1 + CAL 2	55.2	59.7	108.2
S1 + CAL 3	125.2	111.8	89.3
S1 + CAL 4	525.2	391.4	74.5
S1 + CAL 5	1025.2	1251.1	122.0
S2	-	102.9	-
S2 + CAL 0	51.5	63.9	124.1
S2 + CAL 1	61.5	64.2	104.4
S2 + CAL 2	81.5	84.7	103.9
S2 + CAL 3	151.5	133.0	87.8
S2 + CAL 4	551.5	367.2	66.6
S2 + CAL 5	1051.5	1261.3	120.0

Parallelism Test

Sera with high analyte concentration was tested at different dilutions with the Zero Calibrator.

Dilution	Expected (IU/ml)	Measured (IU/ml)	Recovery (%)
S1 undiluted	-	1109	-
1/2	554.5	586.7	105.8
1/4	277.3	275.0	99.2
1/8	138.6	143.0	103.2
1/16	69.3	82.0	118.3

Note: Due to the heterogeneity of autoantibodies, non-linear dilution may be appropriate for some patient samples.

Other translations of this Instruction for Use can be downloaded from our website: <https://www.diasource-diagnostics.com/>

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1. Prentice L.M. et al. Geographical distribution of subclinical autoimmune thyroid disease in Britain : A study using highly sensitive direct assays for autoantibodies to thyroglobulin and thyroid peroxidase. Acta Endocrinologica , 1990, 123, 493-8.

2. Jaume J.C. et al. Thyroid peroxidase autoantibody epitopic « fingerprints » in juvenile Hashimoto's thyroiditis : evidence for conservation over time an in families. Clin. Exp. Immunol., 1996, 104, 115-123.

3. Czarnocka B. et al. Immunochemical properties of hTPO. Thyroperoxidase and Thyroid Autoimmunity, Ed. P. Carayon, J. Ruf. 1990, 207, 59-67.

4. Czarnocka B. et al. Majority of thyroid peroxidase autoantibodies in patients with autoimmune thyroid disease are directed to a single TPO domain. Autoimmunity, 1996, 23, 145-154.

5. Takamatsu J. et al. Correlation of antithyroglobulin and antithyroid-peroxidase antibody profiles with clinical and ultrasound characteristics of chronic thyroiditis. Thyroid, 1998, 8, 1101-6.

6. Smyth P. et al. Serum thyroid peroxidase autoantibodies, thyroid volume, and outcome in breast carcinoma. J. Endocrinol. Metab., 1998, 83, 2711-6.

7. Gauna A. et al. Immunological aspects of Graves' disease patients in different clinical stages. J. Endocrinol. Invest., 1989, 12, 671-7.

8. Stagnaro-Green A. et al. Detection of at-risk pregnancy by means of highly sensitive assays for thyroid autoantibodies. JAMA, 1990, 264, 1422-5.

Physikalisches Daten of ¹²⁵I
Physical characteristics of ¹²⁵I
Características físicas ¹²⁵I
φυσικά χαρακτηριστικά του ¹²⁵I
Caractéristiques physique de ¹²⁵I
Caratteristica fisica ¹²⁵I
Fysieke kenmerken van ¹²⁵I
Charakterystyka fizyczna I-¹²⁵
Características físicas ¹²⁵I
Физични характеристики на ¹²⁵I
¹²⁵I Fizikai tulajdonságai

t_{1/2} = 59.9 Tagen, days, dias, ημέρες, jours, giorni dagen, dni, días,
дена, nap

Wichtigste Emissionnen
Main emissions
Emissionnes principales
κύριες εκπομπές
Emissions principales
Emissione principale
Hoofduitzendingen
Głowne rodzaje promieniowania
Emissões principais
Основни лъчения
Főbb kibocsátási értékek

	E (MeV)	%
γ	0.035	
X	0.027	114
	0.032	25