


mAb Hu ACTH Capture, Purified Unconjugated



PRODUCT SPECIFICATION

Cat # **5300696**

Clone:	180/A2 LF3 BC8
Batch*:	<i>[Batch dependent]</i>
Size	<i>[≥ 1 mg]</i>
Format:	Purified Unconjugated, Liquid
Concentration:	<i>[≥ 1 mg/mL, Batch dependent]</i>
Buffer:	PBS 0,1M pH 7,4
Preservative:	NaN3 0,5 g/L
Purification:	Protein A affinity chromatography
Isotype:	IgG1 K
Immunogen:	Human ACTH (18-39)
Mouse / fusion partners:	Balb/c (Charles River)/ NSO
Specifications :	Recognizes Human ACTH (18-39) and Human ACTH (1-39), no cross reaction with ACTH (1-24)
Applications:	Immunoassays
Storage:	2-8°C
Expiry date:	<i>[Batch dependent]</i>

* Only the figures before the last slash should be taken into account for the batch identification.

Date

Signature